



Handy PEA

Field Portable Continuous Excitation Chlorophyll Fluorescence System

- ✓ Compact (170 x 85 x 40mm), lightweight (565gms)
- ✓ Large-scale screening capacity up to 1000 full trace data files
- ✓ High time resolution detection for discrimination of fast fluorescence induction kinetics
- ✓ Saturating high intensity focused red LED array for accurate determination of the Fm parameter
- ✓ Upload user-defined, repeatable experimental protocols for automatic field execution
- ✓ Interchangeable sensor unit cables with lengths of up to 20 metres
- ✓ Windows® data transfer & analysis software



Hansatech Instruments

Hansatech Instruments is a small, British, scientific instrument company located in the heart of rural Norfolk. For over 40 years, our efforts have been concentrated towards the design & manufacture of high quality instrumentation for teaching and research in the fields of cellular respiration and photosynthesis. Our instruments are now in use in a wide range of programs in more than 100 countries throughout the world and have gained an enviable reputation for quality, reliability and excellent price/performance.



Products

Hansatech Instruments product range covers a wide range of applications in the fields of photosynthesis and cellular respiration. We manufacture oxygen measurement systems based on Clark type polarographic oxygen sensors, chlorophyll fluorescence measurement systems for both continuous excitation and pulse-modulated measurement techniques and optical instrumentation for the measurement of sample chlorophyll content.



Support

Purchasers of Hansatech Instruments products can be assured of ongoing support and prompt and efficient attention to enquiries at all times. Customers are encouraged to register their instruments on our website which allows access to our Support Ticketing System in addition to instruments manuals and software upgrades.



Scan the code for further information.



Overview

The Handy PEA chlorophyll fluorimeter consists of a compact, light-weight control unit encapsulating sophisticated electronics providing the high time resolution essential in performing measurements of fast chlorophyll fluorescence induction kinetics.

The chlorophyll fluorescence signal received by the sensor head during recording is digitised within the Handy PEA control unit using a fast Analogue/Digital converter. The chlorophyll fluorescence signal is digitised at different rates dependent upon the different phases of the induction kinetic. Initially, data is sampled at 10µs intervals for the first 300µseconds.

Simple to configure and operate, the Handy PEA chlorophyll fluorimeter features the capacity to store up to 5 user-defined protocols for different field applications. Protocols are written using a custom Windows® software package, PEA Plus (supplied). This allows single or multiple measurement assays with optional pre-illumination periods to be defined and uploaded to the memory of Handy PEA via RS232 serial communications. The use of protocols ensures maximum reproducibility of results during field applications involving large scale screening away from a laboratory environment. A waterproof, tactile keypad allows selections and inputs to be made and a liquid crystal display module presents menu options and data.

Up to 1000 recordings of between 0.1 – 300 seconds may be saved in the memory of Handy PEA chlorophyll fluorimeter. Saved data may be viewed onscreen in numerical format with calculated parameters or transferred to the PEA Plus software where it may be viewed graphically or exported to external software packages for further statistical analysis.

The sensor unit consists of an array of 3 ultra-bright red LED's optically filtered to a peak wavelength of 650 nm, which is readily absorbed by the chloroplasts of the leaf, at a maximum intensity of up to 3500 µmol m⁻² s⁻¹ at the sample surface. The LED's are focused via lenses onto the leaf surface to provide even illumination over the area of leaf exposed by the leafclip (4mm dia). LED's have the advantage of being rugged, emitting low levels of heat, and of rising to full intensity very rapidly (typically microseconds) after being switched on. This feature eliminates the inaccuracies of Fo measurement and the constraints on speed and reliability associated with a shutter which is a necessary item in systems using filament lamps rather than LED's.

An optical feedback circuit monitors and corrects changes in the output intensity of the LED's. These changes are caused by internal heat build up in the LED's. The circuit also compensates for intensity changes caused by variation in ambient temperature. The light source is calibrated before leaving the factory but may be calibrated by the user at regular intervals using the SQS Serial Quantum Sensor. The detector is a high performance Pin photodiode and associated amplifier circuit. The optical design and filtering ensure that it responds maximally to the longer wavelength fluorescence signal and blocks the reflected shorter wavelength LED light used as the source of illumination.

Leafclips and Sample Dark Adaptation

The first step in the measurement process is to dark adapt the leaf with a small, lightweight leafclip. The clip has a small shutter plate which should be closed over the leaf when the clip is attached so that light is excluded and dark adaptation takes place.

During dark adaptation, all the reaction centres are fully oxidised and available for photochemistry and any chlorophyll fluorescence yield is quenched. This process takes a variable amount of time and depends upon plant species, light history prior to the dark transition and whether or not the plant is stressed. Typically, 15 – 20 minutes may be required to dark adapt effectively. In order to reduce waiting time before measurement, a number of leaves may be dark adapted simultaneously using several leafclips. Some users even make measurements at night, thus ensuring an adequate supply of readily dark adapted samples and zero waiting time!

Technical Specifications

Dimensions:	170 (l) x 85 (w) x 40mm (d). Weight: 565g
Communications:	RS232 bi-directional serial communications
Operating Conditions:	0 - 40°C. Non-condensing humidity
Battery:	3 x rechargeable Ni-MH 3.6V, 1.8Ahr
Battery Charger:	Integral switch mode charger 8-13.5V input (nominal 12V input)
Battery Life:	Typically 8 - 9 hours when performing 1 second recordings at 120 measurements/hr
Display:	8 line x 20 character LCD display
Illumination:	Focused array of ultra-bright red LED's with NIR short pass cut-off filters. Peak wavelength 650 nm. Spectral-line half width 22 nm.
Max. intensity:	Up to 3500 µmol m ⁻² s ⁻¹
Detector:	Fast response PIN photodiode with RG9 long pass filter
Electronics:	16 bit microprocessor, 12 bit resolution, A/D 10µsec acquisition rate, 8 bit DAC for light source control, real time clock
Record Length:	0.1 - 300 seconds
Memory:	512K battery backed RAM. Sufficient for up to 1000 one second duration recordings with full trace data

Hansatech Instruments Ltd
Narborough Road, Pentney, King's Lynn, Norfolk, UK PE32 1JL
Tel: +44 (0)1760 338877 Fax: +44 (0)1760 337303
www.hansatech-instruments.com info@hansatech-instruments.com