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
Safety instructions

- No liquid may get into the drive unit of the disperser.
- The work has to be interrupted immediately if you notice unusual noise and/ or increasing emission of liquid from the upper rinsing drilling. The slide bearing must be controlled, and replaced if necessary.
- While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles.
- Follow the health and safety regulations accident prevention regulations.
- Note the operating instructions of the disperser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.

Product information

- The shaft bearing of the dispersion tool consists of a slide bearing (PTFE).
- Content of package: shaft and dispersion head, spare part grooved piston ring bearing
- The plastic materials used are approved by FDA (Food and Drug Administration Approval).

Application instructions

- Attention! Never run the dispersion tool dry. Without cooling by the medium the slide bearing will be destroyed.
- The dispersion tool is not suitable for continuous operation.
- Dispersing causes the medium to heat up.
- Likewise the optimal dispersion duration and rotating frequency must be determined by attempts.
- The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator. Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Designation	Application	Max. permissible number of revolutions [1/min]	Working range [ml]
S 25 N-18 G	Wet crushing, suspension, emulsifying, disagglomeration, batch-operation	24000	10-1500
S 25 N-25 G	Wet crushing, suspension, emulsifying, disagglomeration, batch-operation	24000	50-2000
S 25 N-25 F	Wet crushing, suspension, emulsifying, disagglomeration, batch-operation	24000	100-2000

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs. The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Spare parts list

Item	Designation	Item	Designation	Spare parts diagram see front page. For orders of spare parts please specify the fabrication number on the shaft tube, the type of dispersion tool and the designation of the spare part.
4	Shaft tube	9	Rotor	
5	Bushing	10	Stator	
6	Friction ring	13	Shaft	
7	Grooved piston ring bearing			

Cleaning, disinfection and sterilization

Allowed procedures	Methods of sterilisation	Notice
Damp heat	Autoclave at 121°C at 2 bar positive pressure	-
Chemical procedures	By germ - killing dissolvents alcohol, phenol,formalin.... Disinfectant remainders are to be removed with germ - free water.	-
Hot air	Germ - killing by hot-air with 160 - 190 °C (approx. 30 min)	-

For pre-cleaning the dispersion tool, operate it in in a dissolvent for dissolving the residuaries of the substance. Due to the large flow rate, rotor and stator are cleaned to a large extent. The dispersion tool must be divided and cleaned immediately after working, so that the adhering residuaries of the substance don't cultivate unwanted bacterial cultures.

For cleaning purposes disassemble the dispersion tool as described below (items see front page):

- Unscrew the rotor (9) from the shaft (13) using the generator and shaft key included in the tool set. The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.
- Unscrew the stator (10) from the shaft (4) (left-hand thread). Insert flat key into the front-end slot of the shaft tube. By using the generator key the stator may now be screwed off in a clockwise direction.
- Pull the shaft (13) out of the shaft tube in a downward direction. The individual parts threaded on the shaft (5, 6 and 7) may simply be pulled off.
- The slotted bearing (7) shall be replaced in the case of wear.

The assembly of the dispersion tool takes place in reverse order.

Attention! The rotor must not be tightened too strongly, because otherwise the thread will be damaged.

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