**DE  
EN  
FR  
ES**

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# **Operation Manual** (EN)

Translation of the German original manual

Piston pump

## **LLG-uniVACUUPUMP 1**

**Lab Logistics Group GmbH**

Am Hambuch 1  
53340 Meckenheim  
Germany

T +49 2225 9211 0

F +49 2225 9211 11

[info@llg.de](mailto:info@llg.de)

[www.llg.de](http://www.llg.de)

Customer Support +49 2225 9211 0

We are constantly working on the further development of all our product models. Reprinting or reproduction of this manual, including extracts, is not allowed without the prior written permission of Co. Lab Logistics Group GmbH.

All rights under the copyright laws are expressly reserved by Co. Lab Logistics Group GmbH.

We reserve the right to make changes and amendments.

## Contents

|          |                                                       |           |
|----------|-------------------------------------------------------|-----------|
| <b>1</b> | <b>Important Information .....</b>                    | <b>4</b>  |
| 1.1      | General Information.....                              | 4         |
| 1.2      | Target Groups .....                                   | 4         |
| 1.3      | Intended Use.....                                     | 4         |
| 1.4      | Use for an Unauthorized Purpose .....                 | 4         |
| 1.5      | Safety Devices .....                                  | 5         |
| 1.6      | Meaning of the Warning notes .....                    | 5         |
| 1.7      | Product Standards, Safety Regulations .....           | 5         |
| <b>2</b> | <b>Basic Safety Instructions .....</b>                | <b>6</b>  |
| 2.1      | General Information.....                              | 6         |
| 2.2      | Electricity.....                                      | 6         |
| 2.3      | Mechanical Systems .....                              | 6         |
| 2.4      | Hazardous Substances .....                            | 7         |
| 2.5      | High Temperatures.....                                | 7         |
| <b>3</b> | <b>Description.....</b>                               | <b>8</b>  |
| 3.1      | Design .....                                          | 8         |
| 3.1.1    | Connections - Suction-/Pressure side .....            | 8         |
| 3.1.2    | Connecting to the electricity supply .....            | 8         |
| 3.2      | Protection measures against liquids in the pump ..... | 9         |
| 3.2.1    | Condensate separator (Suction side).....              | 9         |
| 3.3      | Areas of Application .....                            | 9         |
| 3.4      | Scope of Delivery .....                               | 9         |
| <b>4</b> | <b>Technical Data .....</b>                           | <b>10</b> |
| <b>5</b> | <b>Installation and Operation .....</b>               | <b>11</b> |
| 5.1      | Unpacking .....                                       | 11        |
| 5.2      | Setting up and connecting.....                        | 11        |
| 5.2.1    | Setting up .....                                      | 11        |
| 5.2.2    | Connecting.....                                       | 11        |
| 5.2.2.1  | Electrical Connection .....                           | 11        |
| 5.2.2.2  | Vacuum connection (Suction side).....                 | 11        |
| 5.2.2.3  | Exhaust connection (Pressure side).....               | 11        |
| 5.3      | Operation.....                                        | 12        |
| 5.3.1    | Start-up.....                                         | 12        |
| 5.3.2    | Decommissioning.....                                  | 12        |
| 5.4      | Storage.....                                          | 12        |
| 5.5      | Scrap Disposal .....                                  | 12        |
| <b>6</b> | <b>Maintenance and Servicing .....</b>                | <b>13</b> |
| 6.1      | Maintenance Performed by the User .....               | 13        |
| 6.2      | Maintenance by the Manufacturer.....                  | 13        |
| 6.3      | Damage Report.....                                    | 13        |
| <b>7</b> | <b>Troubleshooting .....</b>                          | <b>14</b> |

- Damage Report
- EC Declaration of Conformity

# Important Information

## 1 Important Information

### 1.1 General Information

The **Piston Pump** conforms to the following directives:

|                        |                                         |
|------------------------|-----------------------------------------|
| <b>2006 / 42 / EC</b>  | Machinery Directive                     |
| <b>2004 / 108 / EC</b> | Electromagnetic Compatibility Directive |
| <b>97 / 23 / EC</b>    | Pressure Equipment Directive            |
| <b>2014 / 29 EU</b>    | Simple Pressure Vessels Directive       |

The CE sign is located on the rating plate.

Observe the binding national and local regulations when fitting the pump into installations!

### 1.2 Target Groups

This Operating Manual is intended for the personnel planning, operating and maintaining Piston Pumps.

This group of people includes:

- Designers and fitters of vacuum apparatus
- Employees working on commercial laboratory and industrial vacuum technology applications
- Service personnel for Piston Pumps

The personnel operating and maintaining the Piston Pumps must have the technical competence required to perform the work that has to be done.

The user must authorize the operating personnel to do the work that has to be done.

The personnel must have read and understood the complete Operating Manual before using the Piston Pumps.

The Operating Manual must be kept at the place of use and be available to the personnel when required.

### 1.3 Intended Use

- The layout of the Piston Pumps must be appropriate for the conditions of use. The user bears the sole responsibility for this.
- The Piston Pumps may only be operated under the conditions stated
  - in the "Technical Data" section,
  - on the type plate, and
  - in the technical specification for the order concerned.
- Piston Pumps are approved for extracting, pumping and compressing gases and vapours. If these gases and vapours are toxic or explosive, then the user must observe the currently valid safety regulations for this application.

### 1.4 Use for an Unauthorized Purpose

It is forbidden to use the pump for applications deviating from the technical data stated on the type plate or the conditions stated in the supply contract, or to operate it with missing or defective protective devices.

## 1.5 Safety Devices

Measures such as the following are for the safety of the operating personnel:

- electrical connection with a protective conductor (operating mode S1) and an earthing plug
- Motor protection device (thermal)
- "Hot Surface" label on the pump - warning notice 

The Piston Pumps must not be operated without these elements.

## 1.6 Meaning of the Warning notes

Take note of the warning notices. They are each in the following box:

|                                                                                   |                              |
|-----------------------------------------------------------------------------------|------------------------------|
|  | <b>CAUTION ! / WARNING !</b> |
| <b>Hazard which may lead to serious injuries or material damage.</b>              |                              |

## 1.7 Product Standards, Safety Regulations

Piston Pumps meet the following product standards:

|                                                    |                                                                                                                                                                                         |
|----------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>DIN EN ISO 12100-2010</b>                       | Safety of machinery - General principles for design - Risk assessment and risk reduction                                                                                                |
| <b>DIN EN ISO 13857:2008-06</b>                    | Safety of machinery - Safety distances to prevent hazard zones being reached by upper and lower limbs                                                                                   |
| <b>DIN EN 1012-1</b><br><b>DIN EN 1012-2</b>       | Compressors and Vacuum pumps - Safety requirements - Part 1: Compressors<br>Part 2: Vacuum pumps                                                                                        |
| <b>DIN EN ISO 2151</b>                             | Acoustics - Noise test code for compressors and vacuum pumps - Engineering method (grade 2)                                                                                             |
| <b>DIN EN 60204-1</b>                              | Safety of machinery - Electrical equipment of machines - Part 1: General requirements                                                                                                   |
| <b>DIN EN 61000-6-2</b><br><b>DIN EN 61000-6-4</b> | Electromagnetic compatibility (EMC) - Part 6-2: Generic standards - Immunity for industrial environments<br>Part 6-4: Generic standards - Emission standard for industrial environments |
| <b>DIN EN 61010-1</b>                              | Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements                                                                 |
| <b>DIN EN 50110-1</b>                              | Operation of electrical installations                                                                                                                                                   |
| <b>Directive 2012/19/EC</b>                        | Electrical and electronics - old devices (WEEE)                                                                                                                                         |
| <b>Directive 2011/65/EC</b>                        | Dangerous materials in electrical and electronics devices (RoHS)                                                                                                                        |
| <b>China - RoHS</b>                                | Environment protection law - China 2007-03                                                                                                                                              |

The following additional safety regulations apply in the FR Germany:

|                                  |                                                   |
|----------------------------------|---------------------------------------------------|
| <b>BGV A3</b>                    | Electrical equipment and operating materials      |
| <b>VBG 5</b>                     | Power-driven machines                             |
| <b>BGR 120</b>                   | Guidelines for laboratories                       |
| <b>BGI 798</b>                   | Hazard assessment in the laboratory               |
| <b>BGG 919 (VBG 16)</b>          | Accident prevention regulations for "compressors" |
| <b>BGR 189 (BGR 195;192;197)</b> | Use of protective working clothes                 |

Observe the standards and regulations applying in your country when you use the Piston Pumps.

# Basic Safety Instructions

## 2 Basic Safety Instructions

### 2.1 General Information

Warning notices must be observed. Disregarding them may lead to damage to health and property.

The Piston Pumps must be operated by personnel who can detect impending dangers and take action to prevent them from materialising.

The manufacturer or authorized workshops will only service or maintain the Piston Pumps if it is accompanied by a fully completed damage report. Precise information about the contamination (also negative information if necessary) and thorough cleaning of the Piston Pumps are legally binding parts of the contract.

Contaminated Piston Pumps and their individual parts must be disposed of in accordance with the legal regulations.

The local regulations apply in foreign countries.

### 2.2 Electricity

The Piston Pump of operation mode S1 is supplied. When the location of operation mode S1 devices is changed, please note that the testing must be repeated in accordance with DIN EN 0105, DIN EN 0702 and BGV A2.

The local regulations apply in foreign countries.

Please note the following when connecting to the electrical power supply system:

- The electrical power supply system must have a protective connector according to DIN VDE 0100-410 (IEC 60364-4-41).
- The protective connector must not have any breaks.
- The connecting cable must not be damaged.

### 2.3 Mechanical Systems

Improper use can lead to injuries or material damage. Observe the following instructions:

- Only operate the pumps with hoses of the specified dimensions.
- The maximum permissible pressure of 1 bar at the suction connection must not be exceeded.
- Hazardous substances must be separated out as far as this is technically possible before they reach the pump.
- External mechanical stresses and vibrations must not be transmitted to the pump. Only use flexible laboratory hoses for connecting pumps.
- The overpressure generated at the pressure port must not exceed 3 bar.
- The pump must not be used to suck up fluids. Lay the exhaust pipe so that it slopes downwards, so allowing condensate to flow out of the pump. Collect the condensate and dispose of it in an environmentally compatible manner.
- Prevent dyes exuding.
- Maintain a space of least 20 mm between the pump and adjacent parts in order to enable the pump to cool.

|                                                                                                                                                   |                  |
|---------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                | <b>CAUTION !</b> |
| <b>Solid particles in the pumping medium impair the pumping action and can lead to damage. Prevent solid particles penetrating into the pump.</b> |                  |

## 2.4 Hazardous Substances

|                                                                                                                                                                                                                                                                      |                  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                                                                                                                                     | <b>CAUTION !</b> |
| <p><b>The operating company bears the responsibility for the use of the Piston Pump. Hazardous substances in the gases to be pumped can cause personal injuries and property damage. Pay attention to the warning notices for handling hazardous substances.</b></p> |                  |

The local regulations apply in foreign countries.

### **Combustible, aggressive and explosive Gases**

Don't pump combustible, especially aggressive or explosive gases or vapors or operate this pump in an atmosphere containing combustible or explosive gases or vapors.

Examine before switching on whether that can form gas combustible gas/air mixtures which can be promoted! Consider the regulations of the guideline 1999/92/EC.

The Piston Pump is not certified according to ATEX guidelines 94/9/EC.

The Piston Pump is not recommended for pumping acidic, basic or organic vapors.

### **Poisonous Gases**

Use a condensate separator when pumping poisonous or harmful gases. Prevent such sub-stances from leaking out of the appliance or pump. Treat these substances according to the applicable environmental protection regulations.

Test the strength and leak-tightness of the connecting lines and the connected apparatus. Prevent environmental poisons, e.g. mercury, getting into the pumps.

Fulfil the requirements, for example:

- German Hazardous Substances Regulation (GefStoffV) of 01. December 2010
- Regulation (EC) 1272/2008  
(Classification, Packaging and Labelling of hazardous substances),
- Manufacturer's safety data sheets on hazardous substances.

## 2.5 High Temperatures

The Piston Pumps may heat up as a result of the temperature of the gas being pumped and through compression heat.

Prevent the following maximum permissible temperatures from being exceeded.

- + 40 °C for the environment, and
- + 60 °C for the gas to be pumped.

The electric motor has a thermal overload protection.

## Description

### 3 Description

#### 3.1 Design

##### 3.1.1 Connections - Suction-/Pressure side

##### 3.1.2 Connecting to the electricity supply

The Piston Pump is supplied with complete electrical wiring. It is connected via a mains connection cable and a power plug. Mains connection cable and plug must comply with the requirements of the line disconnection devices (current, output). Motor is provided with a thermal overload protection ex works, protecting the vacuum pump from damage or destruction, respectively.

|                                                                                                                                                                                                                     |                  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                                                                                  | <b>WARNING !</b> |
| <p>Should the user change the electrical connection, for example for fitting into a system, then this may only be performed by a electrical specialist under observance of the accident prevention regulations.</p> |                  |

|                                                                                                                                                                               |                  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                                            | <b>WARNING !</b> |
| <p>Please keep in mind that the barometric pressure changes from day to day. This also changes the displayed pressure values because of latent air pressure fluctuations.</p> |                  |

### **3.2 Protection measures against liquids in the pump**

#### **3.2.1 Condensate separator (Suction side)**

If the possibility of penetration of larger liquid quantities we recommend the installation of a larger collection vessel (eg. Woulff bottle).

### **3.3 Areas of Application**

**The Piston Pump is intended to:**

- Pumping and compressing neutral and aggressive gases and vapours.
- Generating a vacuum down to a minimum ultimate pressure 292 mbar.
- Use in physical and chemical laboratories in trade and industry.
- Use for vacuum filtration, vacuum drying and other vacuum technology applications.

### **3.4 Scope of Delivery**

The scope of delivery is specified in the supply contract.

The LLG-uni *VACUUMPUMP* 1 comes inclusive mains connection cable with plug EU and UK and 2x Hose clamps 8 – 16 mm.

## 4 Technical Data

| Parameter                                                                                                              | Unit               | Data                                              |
|------------------------------------------------------------------------------------------------------------------------|--------------------|---------------------------------------------------|
| <b>Pumping speed</b><br>at atmospheric pressure                                                                        | m <sup>3</sup> / h | 0.55                                              |
|                                                                                                                        | l / min            | 9.2                                               |
| <b>Ultimate pressure</b>                                                                                               | mbar               | 292                                               |
| <b>Overpressure, max.</b>                                                                                              | bar (ü)            | 3.3                                               |
| <b>Suction-/ Pressure connection</b>                                                                                   | -                  | Hose nozzle DN 6<br>for hose inside diameter 6 mm |
| <b>Ambient temperature</b>                                                                                             | °C                 | + 10 to + 40                                      |
| <b>Max. Operating gas temperature</b>                                                                                  |                    | + 60                                              |
| <b>Noise level</b><br>DIN EN ISO 2151                                                                                  | dB (A)             | ≤ 45                                              |
| <b>Voltage / Frequency</b>                                                                                             | V, Hz              | 230, 50                                           |
| <b>Motor power</b>                                                                                                     | W                  | 25                                                |
| <b>Type of protection</b><br>DIN EN 60529                                                                              | -                  | IP 20                                             |
| <b>Weight</b>                                                                                                          | kg                 | 2.0                                               |
| <b>Dimensions (W/D/H)</b>                                                                                              | mm                 | 194 / 114 / 191                                   |
| <b>Order No.</b>                                                                                                       | -                  | <b>6.263 580</b>                                  |
| <b>Piston pump</b><br><i>inclusive mains connection<br/>cable with plug EU and UK and<br/>2x Hose clamps 8 – 16 mm</i> |                    | <b>LLG-uniVACUUPUMP 1</b>                         |

The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of products. It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability whatsoever in connection therewith. We does not warrant, guarantee or assume any obligation or liability in connection with this information.

## 5 Installation and Operation

### 5.1 Unpacking

Carefully unpack the **Piston Pump**.

Check the pump for:

- Transport damage,
- Conformity with the specifications of the supply contract (type, electrical supply data),
- Completeness of the delivery.

Please inform us without delay if there are discrepancies between the delivery and the contractually agreed scope of delivery, or if damage is detected.

Please take note of the general terms of business of the manufacturing firm.

**In case of a claim under warranty, the device must be returned in packaging that is suitable for protecting it during transport.**

### 5.2 Setting up and connecting

#### 5.2.1 Setting up

- Set the pump on a flat and horizontal surface.
- Remove the protective caps on the suction and pressure ports.
- Note that the cooling of the pump is guaranteed, *see chapter 2.3*.
- Note that on the installation location no moisture acts on the pump.

#### 5.2.2 Connecting

##### 5.2.2.1 Electrical Connection

Before the electrical connection of the pump review the specifications on the nameplate with the existing electrical connection conditions, *see chapter 2.2*.

**CAUTION !**

**The electric motor has a thermal overload protection.**

**After triggering the protective fuse after a certain cooling occurs, an auto restart (Reset) the motor.**

##### 5.2.2.2 Vacuum connection (Suction side)

The vacuum connection (suction port) consists of a hose nozzle DN 6.

**Suitable vacuum line:**

- PVC fabric hose 6 x 3 mm and hose clamp

Make sure that the vacuum line is kept as short as possible from the pump to the apparatus.

##### 5.2.2.3 Exhaust connection (Pressure side)

The exhaust connection (pressure port) also consists of a hose nozzle DN 6.

**Exhaust line for vacuum- and compressor application:**

- PVC fabric hose 6 x 3 mm and hose clamp

# Installation and Operation

|                                                                                   |                  |
|-----------------------------------------------------------------------------------|------------------|
|  | <b>CAUTION !</b> |
| <b>Use only for compressed air applications suitable conduit systems!</b>         |                  |

## 5.3 Operation

|                                                                                            |                  |
|--------------------------------------------------------------------------------------------|------------------|
|           | <b>CAUTION !</b> |
| <b>Observe the basic safety instructions when using the Piston Pump, <i>chapter 2</i>.</b> |                  |

### 5.3.1 Start-up

- Make sure that when you start the pump, which may arise in the connection lines condensates, cannot penetrate into the pump.  
► **Action: suction side condensate separator**
- Always try to avoid aspiration of foreign body particles!  
► **Action: suction side inline filter**
- The Piston pump is switched on the power switch.
- It is recommended to let the pump run for a few minutes before use. The warm-up improves the eligibility of the pump.

### 5.3.2 Decommissioning

- If vacuum applications with special steam load, the pump needs to run after the process with an open vacuum port about 2 minutes and rinsed with atmospheric air.  
This measure is used to cleaning the pump chamber and minimizes the corrosive attack wetted material pump parts.
- The Piston pump is switched off using the power switch.

## 5.4 Storage

The pumps are to be stored in a low-dust, interior room within the temperature range from + 5 to + 40 °C and at a relative air humidity < 90% .  
Leave the protective elements on the suction and pressure ports. Another equally good protection may be used.

## 5.5 Scrap Disposal

|                                                                                                                                                                                                            |                  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                                                                         | <b>CAUTION !</b> |
| <b>The Piston Pumps must be disposed of in accordance with the 2012/19/EC guideline and the specific national regulations.<br/>Contaminated Piston Pumps must be decontaminated according to the laws.</b> |                  |

## 6 Maintenance and Servicing

Piston pump is 100% oil-free.

The piston and cylinder of the pump are lubricated maintenance free.

All bearings are sealed and permanently lubricated.

### 6.1 Maintenance Performed by the User

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |                  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                                                                                                                                                                                                                                                                                                                                                 | <b>WARNING !</b> |
| <p><b>Only perform the work that is described here, and that which is permitted to be done by the user.</b><br/> <b>All other maintenance and service work may only be performed by the manufacturer or a dealer authorized by him.</b></p> <p><b>Beware of the pump parts being possibly contaminated by hazardous substances. Wear protective clothing if there is contamination.</b></p> <p><b>Renew defective parts, if necessary! Do not clean with compressed air!</b></p> |                  |

- Check the pump daily for unusual running noises and heat building up on the surface of the pump.
- Check the electrical and vacuum connections daily.

### 6.2 Maintenance by the Manufacturer

Repairs and maintenance going beyond the extent of the work described *in chapter 6.1* or reconditioning or modification may only be performed by the manufacturer or authorized workshops.

|                                                                                                                                                                            |                  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                                         | <b>WARNING !</b> |
| <p><b>The user shall be liable for the consequences of an incorrect damage report or a contaminated pump. The statements in the damage report are legally binding.</b></p> |                  |

### 6.3 Damage Report

You find the form of the damage report to the appendix of this manual.

|                                                                                                                                                                                                                  |                  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                                                                                                                               | <b>WARNING !</b> |
| <p><b>Incomplete or incorrectly completed damage reports may endanger the service personnel!</b></p> <p><b>Give full information in the damage report, in particular regarding a possible contaminating.</b></p> |                  |

# Troubleshooting

## 7 Troubleshooting

During the warranty period, intervention in the Piston Pumps and accessory components may only be made by manufacturing firm.

| Trouble                                                | Cause                                                  | Remedy                          |                                                                          |
|--------------------------------------------------------|--------------------------------------------------------|---------------------------------|--------------------------------------------------------------------------|
|                                                        |                                                        | by:                             | with:                                                                    |
| <b>Pump does not start</b>                             | No power supply                                        | <b>Qualified electrician</b>    | Check electrical installation                                            |
|                                                        | Motor defective                                        | <b>Service workshop</b>         | Exchange                                                                 |
| <b>Pump does not generate a vacuum or overpressure</b> | Connected apparatus and/or connecting elements leaking | <b>User or Service workshop</b> | Identify and seal the leak, replace the seals and/or hoses if necessary. |
|                                                        | Pump leaking                                           |                                 | Exchange the hoses and/or fittings if necessary.                         |
|                                                        | Pump dirty                                             |                                 | General maintenance / Cleaning                                           |
|                                                        | loose screws                                           |                                 | Tighten the screws                                                       |
| <b>Running noise</b>                                   | loose screws or connectors                             |                                 | Tighten the screws and / or connectors                                   |
| <b>Cable</b>                                           | defective and/or brittle                               | <b>Qualified electrician</b>    | Exchange of the cable                                                    |

[info@llg.de](mailto:info@llg.de)  
[www.llg.de](http://www.llg.de)

Company / Facility

Street

Post Code / Town

Order No.:

Date of Purchase:

Type No.:

Serial No.:

Type:

**1. Description of the fault:** (Which fault occurred?)

**2. Conditions of use:**

a) Were the pumps / components in use?

☐ yes

☐ no

b) Which oil was used? \*

-----

c) Field of application (load, process, use):

-----

**3. Type of contamination of pumps / components:**

a) toxic \*\* <sup>1)</sup>

☐ yes

☐ no

b) corrosive \*\* <sup>1)</sup>

☐ yes

☐ no

c) microbiological \*\* <sup>1)</sup>

☐ yes

☐ no

d) explosive <sup>1)</sup>

☐ yes

☐ no

e) radioactive \*\* <sup>1)</sup>

☐ yes

☐ no

f) other harmful substances <sup>1)</sup>

☐ yes

☐ no

\* In the case of rotary vane pumps  
and roots pumps (Important!)

\*\* Pumps / components will only be accepted if proof is provided that they have been properly cleaned.  
By arrangement, cleaning can be done by LLG GmbH.

**Type(s) of substances or process-related reaction products which came into contact with the pumps/ components:**

|    |  |  |  |  |
|----|--|--|--|--|
| 1. |  |  |  |  |
| 2. |  |  |  |  |
| 3. |  |  |  |  |
| 4. |  |  |  |  |

Are the pumps / components free of harmful substances?

☐ yes

☐ no

**4. Legally binding declaration:**

I/ We declare that the statements here are correct and complete.

The pumps / components are shipped in accordance with the legal requirements. Hereby we order the repair.

Contact:

-----

Phone:

-----

Department:

-----


Fax:

-----

Date:

-----

Legally binding signature: ----- Company Stamp:

The repair can only be made if a completely filled out damage report is received. Otherwise the work may be delayed. Our general terms of delivery and payment apply.

<sup>1)</sup> Answer all questions or mark with ☐ yes or ☐ no. If a box is not filled out, we automatically suppose that there was a contamination and the required measures are charged to you.

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