IKA

designed for scientists

Dispergierwerkzeuge / Dispersion tools / Outils dispersants

S 25 KD (- LR) - 18 G

S 25 KD (- LR) - 25 G



S 25 KD (- LR) - 25 F

S 25 KD (- LR) - 18 G - ST

S 25 KD (- LR) - 25 G - ST

S 25 KD (- LR) - 30 G

S 25 KD (- LR) - 30 F

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EN

Content

	Page
Explication of warning symbols	7
Safety instructions	7
Product information	8
Application instructions	8
Unpacking	8
Spare parts list	8
Maintenance and cleaning	8
Warranty	9
Technical Data	9

Explication of warning symbols

-  **DANGER** Indicates an (extreme) hazardous situation, which, if not avoided, will result in death, serious injury.
-  **WARNING** Indicates a hazardous situation, which, if not avoided, can result in death, serious injury.
-  **CAUTION** Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
-  **NOTICE** Indicates practices which, if not avoided, can result in equipment damage.
-  **CAUTION** Burns hazard!

Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

 **CAUTION** The dispersion tool may heat up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.

 **NOTICE** When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.

 **WARNING** Don't touch rotating parts during operation.

 **CAUTION** Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of splashing and evaporation of liquids; body parts, hair, clothing and jewelry getting caught.

- The work has to be interrupted immediately if you notice unusual noise and/or increasing emission of liquid from the upper, sideways rinsing drilling. The bearings must be controlled, and replaced if necessary.
- Note the operating instructions of the disperser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- Do not use any damaged dispersion tools.

 **NOTICE** Never run the dispersion tool dry. Without cooling by the medium the sealing will be destroyed.

- The dispersion tool is not suitable for continuous operation. Likewise the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute. Longer application times bring no improvements, increase however the sample temperature substantially.
- Dispersing causes the medium to heat up.

- Note the min./max. immersion depth in the “**Technical data**”.
- The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator (integrated in the shaft tube). Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Product information

- The shaft bearing of the dispersion tools consist of a ball bearing. The dispersing tool is sealed by means of a sealing.
- All materials are FDA (Food and Drug Administration) conform.

Application instructions

Application: Wet crushing, dispersions, emulsifying (Batch-operation).

Unpacking

Unpacking:

- Please unpack the device carefully.
- In case of any damage, a detailed report must be sent immediately (post, rail or forwarder).

Delivery scope:

- Dispersion tool according to you order
- Shaft key (see **Fig. 1**)
- Flat key (see **Fig. 1**)
- Flat key **25 F** (only with **S 25 KD (- LR)- 25 F**, see Fig. 1)
- Mesh bag (only with **S 25 KD...**)
(for collecting and packaging small parts of the dispersion tool for cleaning, e.g. washing the rotor, shaft sleeve and sealing in dishwasher.)
- Hook (only with **S 25 KD...**)
(for hook and hang the service card, see **Fig. 2**)
- Technical information
- Short instruction.

Spare parts list

S 25 KD ...

Item	Designation
1	Plastic cap
2	Shaft
3	Shaft tube
4	Ball-bearing
5	Sealing
6	Shaft sleeve
7	Stator
8	Rotor

S 25 KD - LR...

Item	Designation
1	Plastic cap
2	Shaft
2002	Shaft tube
4	Ball-bearing
5	Sealing
6	Shaft sleeve
7	Stator
8	Rotor
10	Deep-groove ball bearing
11	Coupling
12	Threaded pin

Spare parts diagram see front page.

For orders of spare parts, please specify the type of dispersion tool and the designation of the spare part.

Maintenance and cleaning

Cleaning:

The sterilization of the complete dispersion tool is **not** possible.

For pre-cleaning the dispersion tool, operate it in a detergent for dissolving the residuary of the substance. Due to the large flow rate, rotor and stator are cleaned to a large extent. The dispersion tool must be divided and cleaned immediately after working, so that the adhering residuary of the substance don't cultivate unwanted bacterial cultures.

For cleaning purposes disassemble the dispersion tool as described below (items see front page):

- The rotor (**8**) is unscrewed from the shaft (**2**) by means of the flat key and shaft key included in the tool set. Use the flat key to hold the rotor and unscrew it by turning the shaft (**2**) (counterclockwise) with the shaft key.
- Unscrew the stator (**7**) from the shaft tube (**3/2002**) (left-hand thread). Insert shaft key into the front-end slot of the shaft tube. By using the flat key the stator can now be screwed off in a clockwise direction.

- Pull the shaft (**2**) out of the shaft tube (**3/2002**) in a downward direction. The individual parts threaded on the shaft (**4**, **5** and **6**) may simply be pulled off.
- The sealing (**5**) shall be replaced in the case of wear.

Note: The sealing (**5**) mustn't be assembled in the opposite direction (see **Fig. 3**)!

The assembly of the dispersion tool takes place in reverse sequence.

The rotor must not be tightened too strongly, otherwise the thread will be damaged.

Spare parts order:

Spare parts diagram see front page.

For orders of spare parts please specify the serial number on the shaft tube, the type of dispersion tool and the designation of the spare part.

Repair:

Please send the dispersion tool for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For repair, please request the “**Decontamination Certificate**” form **IKA** or use the download printout of it from **IKA** website: **www.ika.com**.

If you require servicing, return the instrument in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs. The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Technical Data

		S 25 KD (- LR) - 18 G	S 25 KD (- LR) - 25 G	S 25 KD (- LR) - 25 F
Working range	ml	10 - 1500	50 - 2000	50 - 2000
Stator/Rotor (Ø)	mm	18/12.6	25/17	25/18
Gap between rotor and stator	mm	0.3	0.5	0.5
Max. allowable speed	rpm	25000		
Max. circumferential speed	m/s	16.5	22.2	23.6
Min./max. immersion depth	mm	40/185 (LR:40/220)		
Material in contact with medium		FFKM and PTFE (conform to the FDA, Food and Drug Administration rules), AISI 316L		
pH range		2 - 13		
Suitable for solvents		Yes		
Suitable for abrasive substances		No		
Max. working temperature	°C	150		
Sterilization methods		all meththods (except the ball bearing)		
Working range vacuum	mbar	100		
Ultimate fineness, suspensions	µm	10 - 50	15 - 50	5 - 25
Ultimate fineness, emulsions	µm	1 - 10	1 - 10	1 - 5

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IKA India Private Limited
IKA England LTD.
IKA Works (Thailand) Co. Ltd.
IKA Works Kenya Ltd.