BINDER

GAS SUPPLY SERVICE

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UK CA UKCA Declaration of Conformity

Name and address of manufacturer	BINDER GmbH Im Mittleren Ösch 5, 78532 Tuttlingen, Germany
Name and address of UK Authorised Representative	Comply Express Ltd Unit C2, Coalport House, Stafford Park 1, Telford TF3 3BD
Object of the Declaration	External Gas Cylinder Changer
Type Designation	BINDER Gas Supply Service
BINDER Art. No.	8012-2344

The Objects of the Declaration described above are in conformity with the relevant UK Regulations and UK Guidelines:

- **Electrical Equipment (Safety) Regulations 2016**
Statutory Instruments 2016 No. 1101 – Consumer Protection Health and safety
- **Electromagnetic Compatibility Regulations 2016**
Statutory Instruments 2016 No. 1091 – Electromagnetic Compatibility
- **The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012**
Statutory Instruments 2012 No. 3032 – Environmental Protection

References of standards and/or technical specifications applied for this Declaration of Conformity, or parts thereof:

S.I. 2016 No. 1101:	EN 61010-1:2010 EN 60204-1:2018
S.I. 2016 No. 1091:	EN 61326-1:2013
S.I. 2012 No. 3032:	EN IEC 63000:2018

This Declaration is issued under the sole responsibility of the manufacturer.

Tuttlingen	10.01.2024			
Place	Date	P. Wimmer Vice President	J. Bollaender Director R & D	BINDER GmbH

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BINDER Gas Supply Service (E2)

Automatic Gas Cylinder Changeover Device, external

Content

EU–Konformitätserklärung / EU Declaration of Conformity / Déclaration de conformité UE / Declaración de conformidad UE	2
UKCA-Konformitätserklärung / UKCA Declaration of Conformity / Déclaration de conformité UKCA / Declaración de conformidad UKCA.....	4
1. Safety	36
1.1. Personnel Qualification	36
1.2. Mounting and Operating Instructions	36
1.3. Legal considerations	36
1.3.1. Intellectual property.....	37
1.4. Structure of the safety instructions.....	37
1.4.1. Signal word panel	37
1.4.2. Safety alert symbol	38
1.4.3. Pictograms	38
1.4.4. Word message panel structure	38
1.5. General safety instructions on installing and operating BINDER Gas Supply Service.....	39
1.6. Precautions when working with gases	39
1.7. Precautions when handling gas cylinders.....	41
1.8. Intended use	42
1.9. Foreseeable Misuse	43
1.10. Residual Risks	44
1.11. Operating instructions by the operator	45
2. Scope of delivery, transportation, storage, and installation.....	45
2.1. Scope of delivery.....	45
2.2. Transport and storage.....	45
2.3. Installation and ambient conditions.....	45
3. Description	46
3.1. Unit overview.....	47
3.2. Technical data	48
4. Installation and connections.....	49
4.1. Connection overview.....	49
4.2. Connecting BINDER Gas Supply Service to the gas cylinders	51
4.2.1. Hose connection at the gas cylinders and setting the gas supply pressure.....	51
4.2.2. Gas cylinder connection kits (option).....	52

4.2.3.	Hose connection at BINDER Gas Supply Service.....	52
4.2.4.	Conversion table for gas inlet pressures, bar – psi	52
4.3.	Connecting BINDER Gas Supply Service to the CO ₂ incubator	53
4.3.1.	Connecting the gas hose to the CO ₂ incubator	53
4.3.2.	Mounting the Y-piece when connecting two CO ₂ incubators.....	56
4.3.3.	Hose connection at BINDER Gas Supply Service.....	56
4.4.	Leak test.....	56
4.5.	Electrical connection of BINDER Gas Supply Service.....	56
5.	Operating BINDER Gas Supply Service.....	57
5.1.	Start up.....	57
5.2.	Changeover function	57
5.3.	Operating and alert indications	58
6.	Service / contact	60
7.	Decommissioning	60
8.	Disposal.....	60
8.1.	Disposal of the transport packing.....	60
8.2.	Disposal of BINDER Gas Supply Service and the power supply unit in the Federal Republic of Germany.....	61
8.3.	Disposal of BINDER Gas Supply Service and the power supply unit in the member states of the EU except for the Federal Republic of Germany	61
8.4.	Disposal of the BINDER Gas Supply Service and the power supply unit in non-member states of the EU	62

1. Safety

1.1. Personnel Qualification

The device must only be installed, tested, and started up by personnel qualified for assembly, startup, and operation of the device. Qualified personnel are persons whose professional education, knowledge, experience and knowledge of relevant standards allow them to assess, carry out, and identify any potential hazards in the work assigned to them. They must have been trained and instructed, and be authorized, to work on the device.

The device should only be operated by laboratory personnel especially trained for this purpose and familiar with all precautionary measures required for working in a laboratory. Observe the national regulations on minimum age of laboratory personnel.

1.2. Mounting and Operating Instructions

These mounting and operating instructions are part of the components of delivery. Always keep them handy for reference in the vicinity of the device. If selling the unit, hand over the mounting and operating Instructions to the purchaser.

To avoid injuries and damage observe the safety instructions of this manual. Failure to follow instructions and safety precautions can lead to significant risks.

	 DANGER
Dangers due to failure to observe the instructions and safety precautions. Serious injuries and device damage. Risk of death. ➤ Observe the safety instructions in these mounting and operating instructions. ➤ Follow the operating procedures in these mounting and operating instructions. ➤ Carefully read the complete mounting and operating instructions of the device prior to installing and using the device. ➤ Keep the mounting and operating instructions for future reference	

Make sure that all persons who use the device and its associated work equipment have read and understood the mounting and operating instructions.

These mounting and operating instructions are supplemented and updated as needed. Always use the most recent version of the instructions. When in doubt, call the BINDER Service Hotline for information on the up-to-dateness and validity of these mounting and operating instructions.

1.3. Legal considerations

These mounting and operating instructions are for informational purposes only. They contain information for correct and safe installing, start-up, operation, and decommissioning of the product. Note: the contents and the product described are subject to change without notice.

Understanding and observing the instructions in these mounting and operating instructions are prerequisites for hazard-free use and safety during operation and maintenance. Images are to provide basic understanding. They may deviate from the actual version of the device. The actual scope of delivery can, due to optional or special design, or due to recent technical changes, deviate from the information and illustrations in these instructions. In no event shall BINDER be held liable for any damages, direct or incidental arising out of or related to the use of this manual.

These mounting and operating instructions cannot cover all conceivable applications. If you would like additional information, or if special problems arise that are not sufficiently addressed in this manual, please ask your dealer or contact us directly, e.g. by phone at the number located on page one of this manual

Furthermore, we emphasize that the contents of these mounting and operating instructions are not part of an earlier or existing agreement, description, or legal relationship, nor do they modify such a relationship. All obligations on the part of BINDER derive from the respective purchase contract, which also contains the entire and exclusively valid statement of warranty administration and the general terms and conditions, as well as the legal regulations valid at the time the contract is concluded. The statements in this manual neither augment nor restrict the contractual warranty provisions.

1.3.1. Intellectual property

These mounting and operating instructions are protected by copyright. Any unauthorized copying or disclosure to third parties is strictly prohibited. We reserve the right to take legal action and, if necessary, to assert claims for damages in the event of infringement.

Trademark Information: All BINDER trademarks relating to products or service, as well as trade names, logos and product names used on the website, products and documents of BINDER company are trademarks or registered trademarks of BINDER company (including BINDER GmbH, BINDER Inc.) in the U.S. and other countries and communities of states. This includes word marks, position marks, word/figurative marks, design configurations, figurative marks, and design patents.

Patent Information: BINDER products, categories of products, and accessories may be covered by one or more patents and/or utility models in the U.S. and other countries and communities of states. This information is provided to satisfy the virtual patent marking provisions of various jurisdictions, in particular it is intended to serve as notice under 35 U.S.C. § 287(a). Products and services listed on the BINDER website may be sold individually or as part of a combination product. Additional patent applications may also be pending in the U.S. and other countries and communities of states.

Please visit www.binder-world.com for more information.

1.4. Structure of the safety instructions

In these mounting and operating instructions, the following harmonized denominations and symbols indicate dangerous situations following the harmonization of ISO 3864-2 and ANSI Z535.6.

1.4.1. Signal word panel

Depending on the seriousness and probability of the consequences, dangers are identified with a signal word, the corresponding safety color, and if appropriate, the safety alert symbol.

 DANGER
Indicates an imminently hazardous situation that, if not avoided, will result in death or serious (irreversible) injury.

 WARNING
Indicates an imminently hazardous situation that, if not avoided, will result in death or serious (irreversible) injury.

**CAUTION**

Indicates a potentially hazardous situation which, if not avoided, may result in moderate or minor (reversible) injury.

NOTICE

Indicates a potentially hazardous situation which, if not avoided, may result in damage to the product and/or its functions or of a property in its proximity.

1.4.2. Safety alert symbol

Use of the safety alert symbol indicates **risk of injury**.

Observe all measures that are marked with the safety alert symbol in order to avoid death or injury.

1.4.3. Pictograms

Warning signs			
 Electrical hazard	 Explosive atmosphere	 Gas cylinders	 Suffocation hazard
 CO ₂ suffocation and poisoning hazard	 Explosive substances	 Fire promoting agents	
Mandatory action signs			Prohibition signs
 Mandatory regulation	 Read instructions	 Environment protection	 Do NOT spray with water
	Information to be observed in order to ensure optimum function of the product.		

1.4.4. Word message panel structure

Type / cause of hazard.

Possible consequences.

- Ø Instruction how to avoid the hazard: prohibition
- Instruction how to avoid the hazard: mandatory action

Observe all other notes and information not necessarily emphasized in the same way, in order to avoid disruptions which could result in direct or indirect injury or property damage.

1.5. General safety instructions on installing and operating BINDER Gas Supply Service

Do not install or operate BINDER Gas Supply Service in potentially explosive areas.

	 DANGER
	Danger of explosion due to introduction of flammable or explosive substances in the device. Serious injury or death from burns and / or explosion pressure. - Ø Do NOT operate the device in potentially explosive areas. ➤ KEEP explosive dust or air-solvent mixtures AWAY from the vicinity of the device.

NO liquid must pass into the inside of the gas cylinder changeover device. Make sure that the installation site is dry.

	 DANGER
	Electrical hazard by water entering the device. Deadly electric shock. Damage to the equipment. - Ø The device must NOT become wet during operation, cleaning, or maintenance. - Ø Do NOT install the device in damp areas or in puddles. ➤ Set up the device in a way that it is splash-proof.

	 CAUTION
	Risk of damage caused by opening by force or damage to the gas cylinder changeover device. Damage to BINDER Gas Supply Service. Ø Do NOT open the gas cylinder changeover device.

1.6. Precautions when working with gases

Notes on handling carbon dioxide (CO₂)

Carbon dioxide (CO₂) in high concentrations is hazardous to health. It is colorless and almost odorless and therefore practically imperceptible. Vent out any CO₂ gas that may escape via good room ventilation or a suitable connection to an exhaust system.

	The installation of a CO₂ warning system is required!
---	---

	 DANGER Danger of suffocation and poisoning by high concentration of CO₂ (> 4 Vol.-%). Death by suffocation. ⊘ Do NOT set up the device in non-ventilated recesses ➤ Ensure technical ventilation measures. ➤ Make sure that the installation area is equipped with a functioning CO₂ warning system. ➤ Observe the relevant regulations for handling CO₂. ➤ Close the CO₂ supply when decommissioning the device.
---	---

Notes on handling oxygen (O₂)

Oxygen (O₂) is colorless and almost odorless and therefore practically imperceptible. It promotes burns, which can proceed explosively. There is a fire hazard for flammable oxygenated materials, e.g. clothes and hair. O₂ is heavier than air and may accumulate in low-lying areas.

 	 DANGER Danger of fire and explosion through contact of combustible materials with O₂ with high concentration of O₂ (> 21 % O₂). Serious injury or death from burns and / or explosion pressure. ⊘ Do NOT set up units in non-ventilated recesses ➤ Ensure technical ventilation measures ➤ Observe the relevant regulations for handling O₂. ➤ Close the O₂ supply when decommissioning the device.
---	---

Take appropriate measures to prevent oxygen enrichment and fire and explosion hazards in areas where oxygen enrichment is possible.

	General information for safe handling of oxygen: • Make sure training of personnel on hazards of oxygen enrichment and necessary safety measures • Make sure adequate labeling of all oxygen equipment and facilities • Make sure gas tightness of all gas connections by checking them for leaks (e.g. with leak spray or diluted soap solution) • Close the main valve of the source of oxygen after work when not using the device • Never lubricate O₂ equipment with oil or fat. Use only materials and spare parts which are approved for use with oxygen. • Regularly inspect fire extinguishers for proper condition • Set up emergency showers where oxygen enrichment is possible • Strictest smoking ban and no ignition sources in areas where oxygen enrichment is possible • Make sure good ventilation of areas where oxygen enrichment is possible (location of the device and/or O₂ cylinders) • Persons who may have been in a possibly oxygen-enriched atmosphere must keep away from ignition sources (flames, cigarettes, etc.) and ventilate their clothes at least 15 minutes • Always keep emergency routes free
---	---

Notes on handling nitrogen (N₂)

Nitrogen (N₂) in high concentrations is hazardous to health. It is colorless and almost odorless and therefore practically imperceptible. Any N₂ gas that may escape must be safely led out via good room ventilation or a suitable connection to an exhaust system.

	 DANGER
	Risk of suffocation through high concentration of N₂. Death by suffocation. - Ø Do NOT set up units in non-ventilated recesses ➤ Ensure technical ventilation measures ➤ Observe the relevant regulations for handling N₂. ➤ Close the N₂ supply when decommissioning the device.

1.7. Precautions when handling gas cylinders

	General information for safe handling of gas cylinders: • Store and use gas cylinders only in well-ventilated locations. • Open the gas cylinder valve slowly to avoid pressure surges • Secure gas cylinders during storage and use against falling (chaining). • Transport gas cylinders with a cylinder cart, do not carry, roll, or throw them. • Always close the valve even with apparently empty cylinders; screw on the cap when not in use. Return gas cylinders with the valve closed. • Do not open gas cylinders by force. Mark them when damaged • Protect gas cylinders against fire, e.g. do not store together with flammable liquids • Observe relevant regulations for dealing with gas cylinders.
---	--

Secure the gas cylinders against falling and other mechanical damage.

	 WARNING
	Risk of injury through sudden release of the stored pressure energy when the valve safety is torn off. Injuries. ➤ Secure gas cylinders against falling (chaining). ➤ Transport gas cylinders with a cylinder cart.

The valve of the gas cylinder **always** must be closed before screwing on or unscrewing the gas hose.

	 WARNING Risk of injury through sudden release of the stored pressure energy when opening the cylinder valve of a not connected cylinder. Injuries. ➤ Close the gas cylinder valve before connecting or removing the gas hose.
---	---

1.8. Intended use

	Following the instructions in this mounting and operating instructions is part of the intended use.
---	---

Any use of the devices that does not comply with the requirements specified in these mounting and operating instructions shall be considered improper use.

Other applications than those described in this chapter are not approved.

Use

The external gas cylinder changeover device BINDER Gas Supply Service serves to automatically supply one or two CO₂ incubators with gas. It permits automatic switching to a second gas cylinder as soon as the first gas cylinder is empty.

BINDER Gas Supply Service is suitable only for non-aggressive, non-corrosive, and not inflammable gases.

In case of foreseeable use of the device there is no risk for the user through the integration of the device into systems or by special environmental or operating conditions in the sense of EN 61010-1. For this, the intended use of the device and all its connections must be observed.

The device does not dispose of any measures of explosion protection.

Contamination of the device by toxic, infectious or radioactive substances must be prevented

Medical devices

The device is not classified as a medical device as defined by Regulation (EU) No 2017/745.

Personnel Requirements

Only trained personnel with knowledge of the mounting and operating instructions can set up and install the device, start it up, operate, clean, and take it out of operation. Service and repairs call for further technical requirements (e.g. electrical know-how).

Installation site requirements

The devices are designed for setting up inside a building (indoor use).

The requirements described in the mounting and operating instructions for installation site and ambient conditions must be met.

Relevant regulations for dealing with CO₂ / O₂ / N₂ and gas cylinders must be observed.

The installation of a CO₂ warning system is required.

1.9. Foreseeable Misuse

Other applications than those described in chap. 1.8 are not approved.

This expressly includes the following misuses (the list is not exhaustive), which pose risks despite the inherently safe construction and existing technical safety equipment:

- Non-observance of the mounting and operating instructions
- Non-observance of information and warnings on the device
- Installation, startup, operation, maintenance and repair by untrained, insufficiently qualified, or unauthorized personnel
- Missed or delayed maintenance and testing
- Non-observance of traces of wear and tear
- Connection to unsuitable gases
- Failure to observe the maximum connection pressure of 2.5 bar / 36.3 *psi*
- Failure to comply with the relevant regulations for handling gas cylinders
- Failure to comply with the relevant regulations for handling CO₂ / O₂ / N₂
- Operation of the device without ventilation measures
- Operation of the device without a functioning CO₂ warning system.
- Installation, testing, service or repair in the presence of solvents
- Installation of replacement parts and use of accessories and operating resources not specified and authorized by the manufacturer
- Bypassing or changing protective systems, operation of the device without the designated protective systems
- Non-observance of messages regarding cleaning and disinfection of the device.
- Spilling water or cleaning agent on the device, water penetrating into the device during operation, cleaning or maintenance.
- Cleaning activity while device is turned on.
- Operation of the device with a damaged housing or damaged power cord.
- Continued operation of the device during an obvious malfunction
- Human error (e.g. insufficient experience, qualification, stress, exhaustion, laziness)

To prevent these and other risks from incorrect operation, it is recommended the operator issue operating instructions and standard operating procedures (SOPs).

1.10. Residual Risks

The unavoidable design features of a device, as well as its proper field of application, can also pose risks, even during correct operation. These residual risks include hazards which, despite the inherently safe design, existing technical protective equipment, safety precautions and supplementary protective measures, cannot be ruled out.

Messages on the device and in the mounting and operating instructions warn of residual risks. The consequences of these residual risks and the measures required to prevent them are listed in the mounting and operating instructions. Moreover, the operator must take measures to minimize hazards from unavoidable residual risks. This includes, in particular, issuing operating instructions.

The following list summarizes the hazards against which these mounting and operating instructions warn, and specifies protective measures at the appropriate spots:

Unpacking, Transport, Installation

- Setup of the device in unauthorized areas
- Installation of a damaged device
- Installation of a device with damaged power cord
- Inappropriate site of installation
- Missing protective conductor connection

Normal operation

- Assembly errors
- Contact with live parts in normal state

Cleaning and Decontamination

- Penetration of water into the device
- Inappropriate cleaning and decontamination agents

Malfunction and Damage

- Continued operation of the device during an obvious malfunction
- Contact with live parts during error status
- Operation of a unit with damaged power cord

Maintenance

- Maintenance work on live parts.
- Execution of maintenance work by untrained/insufficiently qualified personnel
- Electrical safety analysis during annual maintenance not performed

Trouble-shooting and Repairs

- Trouble-shooting of live parts without specified safety measures
- Absence of a plausibility check to rule out erroneous inscription of electrical components
- Performance of repair work by untrained/insufficiently qualified personnel
- Inappropriate repairs which do not meet the quality standard specified by BINDER
- Use of replacement parts other than BINDER original replacement parts
- Electrical safety analysis not performed after repairs

1.11. Operating instructions by the operator

Depending on the application and location of the device, it is recommended that the operator of the device provides the relevant information for safe operation of the device in a set of operating instructions.

Keep the operating instructions with the device at all times in a place where they are clearly visible. They must be comprehensible and written in the language of the employees.

2. Scope of delivery, transportation, storage, and installation

2.1. Scope of delivery

The external gas cylinder changeover device with monitoring system BINDER Gas Supply Service (Art. no. 8012-0408) consists of:

- Changeover device (control unit with connections)
- Power supply unit 100V-240V with socket adapters
- Y-piece
- 4 hose clamps
- Mounting instructions (Art. no. 7001-0451)

2.2. Transport and storage

- Permissible ambient temperature range during transport and storage: 0°C / 32°F up to +55°C / 131°F.
- Permissible ambient humidity: 70 % r.H. max., non-condensing

Intermediate storage in a closed and dry room.

2.3. Installation and ambient conditions

- Permissible ambient temperature range during operation: 0°C / 32°F up to +55°C / 131°F (CO₂ incubator: +18°C / 64.4°F to +30°C / 86°F.)
- Permissible ambient humidity: 70 % r.H. max., non-condensing.
- Installation height: max. 2000 m above sea level.

To completely isolate BINDER Gas Supply Service from the power supply, you must disconnect the plug of the power supply unit. Install the unit in a way that the power supply unit is easily accessible and can be easily disconnected in case of danger.

Do not install or operate the device in potentially explosive areas.

DANGER

Danger of explosion due to combustible dusts or explosive mixtures in the vicinity of the device.

Serious injury or death from burns and / or explosion pressure.

⊘ Do NOT operate the device in potentially explosive areas.

➤ KEEP explosive dust or air-solvent mixtures AWAY from the vicinity of the device.

NO liquid must pass into the inside of the gas cylinder changeover device. Make sure that the installation site is dry.

	 DANGER
Electrical hazard by water entering the device. Deadly electric shock. Damage to the equipment. - Ø Do NOT install the device in damp areas or in puddles. ➤ Set up the device in a way that it is splash-proof.	

CO₂, as well as O₂, and N₂ are harmful in high concentrations. Any gas that may escape has to be led out via good room ventilation or by connection to a suitable exhaust system.

	Observe the safety instructions in chap. 1.6 and the relevant regulations for handling CO₂, O₂, and N₂.
---	---

5,8 8,7

3. Description

The automatic gas cylinder changeover device BINDER Gas Supply Service permits safe gas supply of one or two CO₂ incubators by automatic switching to a second gas cylinder as soon as the first gas cylinder is empty. If e.g. the pressure of gas cylinder A lies below 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi, the gas cylinder changeover device triggers an alarm and automatically switches to gas cylinder B after a delay time of 2 minutes.

To operate, the gas cylinder changeover device requires a power supply, which is ensured by a power supply.

BINDER Gas Supply Service provides acoustical and optical alarms.

BINDER Gas Supply Service can be used for two CO₂ incubators max. and is suitable for cylinders of CO₂, O₂, and N₂ gas.

Functionality:

- Manual selection of the connected gas cylinder. This is only possible if the cylinder pressure (overpressure against ambient) is at least 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi.
- Automatic switching to the other gas cylinder. This takes place 2 minutes after a pressure drop below 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi and is possible only if the other gas cylinder has at least a pressure of 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi (pressure meaning overpressure against ambient).

The device is equipped with magnetic feet and can be attached to the side wall of a BINDER CO₂ incubator.

3.1. Unit overview


Figure 1: Unit front of the automatic gas cylinder changeover device

Pilot lamps A and B (ring illumination)	
	Green pilot lamp: Cylinder is in use, cylinder pressure > 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi (cylinder full)
	Pilot lamp off: Cylinder is not in use, cylinder pressure > 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi (cylinder full)
	Red pilot lamp: Cylinder is empty or closed or not connected (cylinder pressure < 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi)


Figure 2: Connections at the bottom of the gas cylinder changer control unit

The functional buttons A and B serve to manually select one of the two connected gas cylinders A and B. The button with the horn symbol serves to reset the alarm. For detailed explanation of the operating and warning indications, please see chap. 5.3.

3.2. Technical data

Dimensions: L x W x D: ca. 240 mm / 9.45 in x 165 mm / 6.50 in x 120 mm / 4.72 in

Weight Gas Supply Service: approx. 2.3 kg / 5.1 lb

Supply pressure: 0.5 bar / 7.25 psi up to 2.0 bar / 29 psi

Permissible overpressure: max. 2.5 bar / 36.3 psi

Operating voltage through power supply unit 24 V DC

Connection voltage power supply unit: 100 V-240 V

Type of protection: IP20

4. Installation and connections

4.1. Connection overview

You can connect one or two CO₂ incubators to the output (c) of BINDER Gas Supply Service. To connect a second unit, a Y-piece is delivered, providing 3 connections and to be mounted into the descending tube. 4 hose clamps are included.


Figure 3: Connection diagram for one CO₂ incubator


Figure 4: Connection diagram for two CO₂ incubators

Hose clamps ① + ② supplied with CO₂-incubator






Hose clamps ③ to ⑥ supplied with BINDER Gas Supply Service

 Gas hose with textile insert

Recommendation: gas cylinder connecting kits (chap. 4.2.2)

Overview of the connecting procedure

 Gas cylinders	Chap. 4.2.1 Hose connection at the gas cylinders and setting the gas supply pressure Chap. 4.2.2 Gas cylinder connection kits (option)
 Gas cylinders (a) (b) BINDER Gas Supply Service	Chap. 4.2.3 Hose connection at BINDER Gas Supply Service
 CO₂ incubator	Chap. 4.3.1 Connecting the gas hose to the CO₂ incubator
 CO₂ incubator CO₂ incubator Y-piece	Chap. 4.3.2 Mounting the Y-piece when connecting two CO₂ incubators
 CO₂ incubator (c) BINDER Gas Supply Service	Chap. 4.3.3 Hose connection at BINDER Gas Supply Service

4.2. Connecting BINDER Gas Supply Service to the gas cylinders

Observe the safety instructions in chap. 1.6 and the relevant regulations for handling CO₂, O₂, and N₂.

Observe the safety instructions on handling gas cylinders given in chap. 1.7.

The gas connections must be established by qualified personnel who are trained in handling the respective gases and familiar with the required safety measures.

4.2.1. Hose connection at the gas cylinders and setting the gas supply pressure

The CO₂, O₂ or N₂ gas necessary for operation must have a technical grade of 99.5 %.

The following steps are required:

Ensuring the correct output pressure

A gas supply pressure above 2.5 bar / 36 *psi* will result in unit damage to the CO₂ incubator and BINDER Gas Supply Service.

Use a pressure reducer for each cylinder and make sure to avoid any excessive outlet pressure when connecting the gas hose.

The real outlet pressure of the gas cylinders on the second manometer must **not** exceed 2.5 bar / 36 *psi*. The outlet pressure of each gas cylinder must be 2.0 bar / 29 *psi* above the ambient pressure.

NOTICE

Danger of damage by excessive outlet pressure > 2.5 bar / 36 *psi*.

Damage to the CO₂ incubator and BINDER Gas Supply Service.

- ⊘ The outlet pressure must NOT exceed the indicated value of 2.5 bar / 36 *psi*.
- Before connecting, check the outlet pressure on the pressure reducer of each cylinder.
- Adjust the outlet pressure of each cylinder to 2.0 bar / 29 *psi* above the ambient pressure.

Remark for CB CO₂ incubators before 2007: The gas connection values of the CB CO₂ incubators have been modified. If a maximum value of 2.0 bar / 29 *psi* is indicated above the gas socket of your incubator, this value will remain valid. In this case, set the gas admission pressure between 1 bar to 1.5 bar / 14 *psi* to 22 *psi* above the ambient pressure. If you connect several units with different maximum values of gas admission pressure to a common supply net, use pressure reducers for units with a lower maximum gas admission pressure.

Observe the correct outlet pressure also when replacing the gas cylinders.


Establishing the connection to the gas cylinders

Connect the gas hoses to the pressure reducers of the gas cylinders and secure the connections with hose clamps. Socket DN 6 for hose with internal diameter 6 mm / 0.24 *inches*.

Leak test

After connecting the gas cylinder, check all gas connections for leaks (e.g. with leak spray or diluted soap solution).

4.2.2. Gas cylinder connection kits (option)

Gas cylinder connection kits are available for CO₂ (Art. No. 8012-0014), O₂ (Art. No. 8012-0015) and N₂ (Art. No. 8012-0016).

The connection kit includes the following parts for connecting a gas cylinder to the CO₂ incubator:

- Pressure reducer with manometers for cylinder pressure (high pressure gauge) and outlet pressure (low pressure gauge)
- 5 m pressure hose with pre-assembled hose nozzle for quick acting closure socket
- 1 hose clamp to connect the gas hose to the pressure reducer

Figure 5: Gas cylinder connection kit

	NOTICE
	Danger of damage by excessive outlet pressure > 2.5 bar / 36 psi. Damage to the CO₂ incubator and BINDER Gas Supply Service. ⚠ The outlet pressure must NOT exceed the indicated value of 2.5 bar / 36 psi. ➤ Before connecting, check the outlet pressure on the pressure reducer of each cylinder. ➤ Adjust the outlet pressure of each cylinder to 2.0 bar / 29 psi above the ambient pressure.

	The gas connections must be established by qualified personnel who are trained in handling the respective gases and familiar with the required safety measures.
---	--

4.2.3. Hose connection at BINDER Gas Supply Service

Connect the inputs for cylinder A and for cylinder B on the gas cylinder changer control unit with both gas hoses. Socket DN 6 for hose with internal diameter 6 mm / 0.24 inches.

4.2.4. Conversion table for gas inlet pressures, bar – psi

bar	psi	bar	psi	bar	psi
1	14.5	3	43.5	5	72.5
1.5	21.7	3.5	50.7	5.5	79.7
2	29.0	4	58.0	6	87.0
2.5	36.3	4.5	65.2		

4.3. Connecting BINDER Gas Supply Service to the CO₂ incubator

Observe the safety instructions in chap. 1.6 and the relevant regulations for handling CO₂, O₂, and N₂.

The gas connections must be established by qualified personnel who are trained in handling the respective gases and familiar with the required safety measures.

4.3.1. Connecting the gas hose to the CO₂ incubator

The procedure of connecting the gas hose to the unit rear is the same for any gas connection. All quick acting closure sockets (CO₂, and N₂ and O₂ on unit with O₂ control) are degreased and supplied with a FKM gasket.

Note for units with O₂ control: The quick acting closure socket and the hose nozzle for the O₂ connection (unit with O₂ control) must be degreased.

 	 DANGER Fire and explosion hazard through contact of fat with O₂. Serious injury or death from burns and / or explosion pressure. ➤ All connection parts for the O₂ connection must be degreased.
---	--

The gas hose for connection to BINDER Gas Supply Service is already attached to the hose nozzle and secured by a hose clamp. Plug the hose nozzle into the corresponding quick acting closure socket (a) located at the rear of the CO₂ incubator.


Connect only the supplied hose nozzle to the quick acting closure socket of the CO₂ incubator.

Otherwise, the quick acting closure socket may leak, and/or it may become impossible to connect the original hose nozzle. In this case, please contact BINDER Service.


Now fit the hose nozzle in the quick acting closure socket. To remove the connection, you can pull the hose nozzle off the quick acting closure socket.

Position of the quick acting closure sockets for the gas connections in the rear control panels of BINDER CO₂ incubators


CB (E2, E3)


CB 53 without O₂ control

CB 150 / CB 210 without O₂ control


CB 53 with O₂ control

CB 150 / CB 210 with O₂ control

Figure 6: Gas connections in the rear control panel CB (E2, E3)

CB / CBF (E6, E7)

CB / CBF (E6, E7) without O₂ control

CB / CBF (E6, E7) with O₂ control

- Regular equipment: Hypoxic control range **0.2 to 20 vol. % O₂**: Only connections for CO₂, N₂
- Alternative control range **10 to 95 vol. % O₂** (option). Although the high control range is intended in particular for hyperoxic applications (> 21 vol. % O₂), it is also suitable for slightly hypoxic applications between 10 and 20 vol. % O₂. Only connections for CO₂, O₂





CB / CBF (E6, E7), without O₂ control,
regular equipmentCB / CBF (E6, E7), without O₂ control,
optional alternative control range

Figure 7: Gas connections in the rear control panel CB / CBF (E6, E7)


C (E1, E2)

C (E1, E2)

C (E3)

Figure 8: Gas connection in the rear control panel C (E1, E2, E3)

C / CB-S (E7)

C / CB-S (E7)


Figure 9: Gas connection in the rear control panel (right) C / CB-S (E7)

4.3.2. Mounting the Y-piece when connecting two CO₂ incubators

Figure 10: Connection diagram for two CO₂ incubators

① + ② + ③ Hose clamps (two supplied with CO₂ incubator, one supplied with BINDER Gas Supply Service)

 Gas hose with textile insert

4.3.3. Hose connection at BINDER Gas Supply Service

Connect the output for the incubator (c) on the gas cylinder changer control unit with the gas hose. Socket DN 6 for hose with internal diameter 6 mm / 0.24 inches.

4.4. Leak test

After establishing the gas connections, check all gas connections for leaks (e.g. with leak spray or diluted soap solution).

4.5. Electrical connection of BINDER Gas Supply Service

- Prior to connection and start-up, check the power supply voltage. Connection voltage for power supply unit: 100 V-240 V.

	NOTICE
	Danger of incorrect power supply voltage due to improper connection. Damage to the device. ➤ Check the power supply voltage before connection and start-up. ➤ Compare the power supply voltage with the connection voltage 100 V to 240 V.

- Connect the power supply unit to BINDER Gas Supply Service.
- Use only the original power supply unit according to the above specification.
- Install the socket adapter that matches the socket
- Connect the power supply unit to the power supply.
- When connecting the device, please observe the regulations specified by the local electricity supply company as well as the local and national regulations (VDE directives for Germany).
- Make sure that there is sufficient current protection in accordance with the number of devices that are to be operated. We recommend the use of a residual current circuit breaker.

- Pollution degree (acc. to IEC 61010-1): 2
- Over-voltage category (acc. to IEC 61010-1): II

See also electrical data (chap. 3.2).

To completely separate the unit from the power supply, you must disconnect the power plug. Install the unit in a way that the power plug is easily accessible and can be easily pulled in case of danger.

5. Operating BINDER Gas Supply Service

5.1. Start up

Initial state: Both gas cylinders are closed. The hose connections have been connected to the gas cylinders and to the incubator and checked for leaks. The pressure reducers are set to the correct output pressure. When CO₂ will be used, the CO₂ warning system is in operation

- Connect BINDER Gas Supply to the power supply (chap. 4.5),

The pilot lamp "A" is initially lit in green.

After 2 minutes the red pilot lamps are lit, because the supply pressure of the connected gas cylinders is still below 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi.

- Reset the alarm buzzer.
- Open the valves of both gas cylinders.

If both gas cylinders are full and open:

- The pilot lamp of the selected gas cylinder is lit in green.
- The pilot lamp of the not selected gas cylinder is off.

You can manually select gas cylinder A or B with the buttons A and B.

Gas supply is guaranteed as long as at least one of the gas cylinders delivers an output pressure higher than 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi above the ambient pressure

For details on the indications and alerts, please refer to chap. 5.3

5.2. Changeover function

- If the pressure of the selected cylinder falls below 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi for more than 2 minutes, the red pilot lamp of this empty cylinder lights up. The buzzer sounds.
- Then there is automatic switching to the other, full gas cylinder.
- You can mute the acoustical alarm by hitting the button with the horn symbol.
- If following switching to the other gas cylinder the pressure of the (probably empty) cylinder increases again to a value above 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi (e.g. when re-opening an accidentally closed gas cylinder), BINDER gas supply service will not switch back to this cylinder. But the red pilot lamp of this gas cylinder goes out.
- Having replaced the empty gas cylinder (pressure < 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi), the red pilot lamp A or B go out.
- If following alarm buzzer reset the empty gas cylinder is not replaced by a full one, then there will be a new acoustical alarm after emptying the other gas cylinder. Changeover to another cylinder is no longer possible.

	The recovery times of the gas concentrations inside the incubator following the door being opened, which are indicated in the technical data of the CO₂ incubator refer to a connection pressure of 2.0 bar / 29 <i>psi</i>. Decreasing supply pressure down to the shift point (alarm point) of the built-in gas cylinder changer of 0.5 +/- 0.1 bar / 7.25 +/- 1.5 <i>psi</i> results in longer recovery times. Check the pressure displays of your gas supply. If very short recovery times are required or the door is opened frequently, replace the gas cylinders promptly when the pressure decreases below 2.0 bar / 29 <i>psi</i>.
---	---

After using the equipment, close off the gas supply to BINDER Gas Supply Service (close the valves of the gas cylinders).

5.3. Operating and alert indications

Status during commissioning after connecting to the power supply			
A  green	B  off	Status immediately after connecting to the power supply: The power supply is present. The presence of the required gas pressure is checked	
A  red	B  red	Status 2 minutes after connecting to the power supply: - Both gas cylinders are empty or closed or not connected (cylinder pressure below 0.5 +/- 0.1 bar / 7.25 +/- 1.5 <i>psi</i>) - The audible alarm (buzzer) sounds Actions: - Acknowledge the alarm. - Make sure that both cylinders are full and connected. - Open both connected gas cylinders. Check the pressure at the pressure reducers of both cylinders: The pressure must be 1.5 bar / 21.8 <i>psi</i> to 2.0 bar / 29 <i>psi</i> during gas consumption.	  on

Normal state during operation (example)			
A  green	B  off	Status: - Gas cylinder A in use, cylinder pressure above 0.5 +/- 0.1 bar / 7.25 +/- 1.5 <i>psi</i> (full cylinder) - Gas cylinder B not in use, cylinder pressure above 0.5 +/- 0.1 bar / 7.25 +/- 1.5 <i>psi</i> (full cylinder)	

Condition with an empty gas cylinder during operation (example)

⇒ No changeover function possible!

A  green	B  red	Status: - Gas cylinder A in use, cylinder pressure above 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi (full cylinder) - Gas cylinder B is empty or closed or not connected (cylinder pressure below 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi) - The audible alarm (buzzer) sounds Actions: - Acknowledge the alarm. - Make sure that cylinder B is full, connected, and open. If necessary, replace the empty cylinder. - Check the pressure at the pressure reducer of cylinder B: The pressure must be 1.5 bar / 21.8 psi to 2.0 bar / 29 psi during gas consumption.	  on
--	--	---	--

Condition with two empty gas cylinders during operation

⇒ No changeover function and no gas supply possible!

A  red	B  red	Status: - Both gas cylinders are empty or closed or not connected (cylinder pressure below 0.5 +/- 0.1 bar / 7.25 +/- 1.5 psi) - The audible alarm (buzzer) sounds Actions: - Acknowledge the alarm. - Replace both gas cylinders as quickly as possible and carry out the commissioning again. - Make sure that both cylinders are full, connected, and open. - Check the pressure at the pressure reducers of both cylinders: The pressure must be 1.5 bar / 21.8 psi to 2.0 bar / 29 psi during gas consumption	  on
--	--	--	--

Behavior during / after power failure:

There is no changeover between the gas cylinders. Gas supply is given.

6. Service / contact

BINDER telephone hotline:	+49 (0) 7462 2005 555
BINDER fax hotline:	+49 (0) 7462 2005 93555
BINDER service hotline USA:	+1 866 885 9794 or +1 631 224 4340 x3 (toll-free in the USA)
BINDER service hotline Asia Pacific:	+852 390 705 04 or +852 390 705 03
BINDER Internet website	http://www.binder-world.com
BINDER address	BINDER GmbH, post office box 102, 78502 Tuttlingen, Germany

International customers, please contact your local BINDER distributor.

7. Decommissioning

- Turn off the CO₂ incubator at its main power switch and disconnect it from the power supply.
- Disconnect BINDER Gas Supply Service from the power supply.
- Close the valves of the gas cylinders. Remove the gas connections from BINDER Gas Supply Service.
- Remove the hose connections between the CO₂ incubator and BINDER Gas Supply Service.
- Temporal decommissioning: See indications for appropriate storage, chap. 2.2
- Final decommissioning: Dispose of the unit as described in chap. 8.

When taking into operation the system again, please observe the indications given in chap. 5 as well as the safety instructions in chap. 1.6 and 1.7.

8. Disposal

8.1. Disposal of the transport packing


Packing element	Material	Disposal
Transport box with metal clamps	Cardboard	Paper recycling
	Metal	Metal recycling
Bag for instructions manual	PE foil	Plastic recycling
Insulating air cushion foil	PE foil	Plastic recycling

If recycling is impossible, all packing parts can also be disposed of with normal waste.

8.2. Disposal of BINDER Gas Supply Service and the power supply unit in the Federal Republic of Germany

According to Annex I of Directive 2012/19/EU of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE), BINDER devices are classified as “monitoring and control instruments” (category 9) only intended for professional use“. They must not be disposed of at public collecting points.

BINDER Gas Supply Service and the power supply unit bear the symbol for the marking of electrical and electronic equipment manufactured / placed on the market in the EU after 13 August 2005 and be disposed of in separate collection according to Directive 2012/19/EU on waste electrical and electronic equipment (WEEE) and German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG). WEEE marking: crossed-out wheeled bin. A significant part of the materials must be recycled in order to protect the environment.

At the end of the device's service life, have the device disposed of according to the German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG from 20 October 2015, BGBl. I p. 1739) or contact BINDER service who will organize taking back and disposal of the device according to the German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG from 20 October 2015, BGBl. I p. 1739).


	NOTICE Danger of violation against existing law if not disposed of properly. Failure to comply with applicable law. - Ø Do NOT dispose of BINDER devices at public collecting points. ➤ Have the device disposed of professionally at a recycling company which is certified according to the German national law for electrical and electronic equipment (Elektro- und Elektronikgerätegesetz, ElektroG from 20 October 2015, BGBl. I p. 1739). <i>or</i> ➤ Instruct BINDER Service to dispose of the device. The general terms of payment and delivery of BINDER GmbH apply, which were valid at the time of purchasing the device.
--	---

Certified companies disassemble waste (used) BINDER equipment in primary substances for recycling according to Directive 2012/19/EU. The devices must be free from toxic, infectious or radioactive substances in order to eliminate any health hazards to the employees of the recycling companies.

8.3. Disposal of BINDER Gas Supply Service and the power supply unit in the member states of the EU except for the Federal Republic of Germany

According to Annex I of Directive 2012/19/EU of the European Parliament and of the Council on waste electrical and electronic equipment (WEEE), BINDER devices are classified as “monitoring and control instruments” (category 9) only intended for professional use“. They must not be disposed of at public collecting points.

BINDER Gas Supply Service and the power supply unit bear the symbol for the marking of electrical and electronic equipment manufactured / placed on the market in the EC after 13 August 2005 and be disposed of in separate collection according to the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). WEEE marking: crossed-out wheeled bin.

At the end of the device's service life, notify the distributor who sold you the device, who will take back and dispose of the device according to the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

	NOTICE Danger of violation against existing law if not disposed of properly. Failure to comply with applicable law. - Ø Do NOT dispose of BINDER devices at public collecting points. ➤ Have the device disposed of professionally at a recycling company that is certified according to conversion of the Directive 2012/19/EU into national law. <i>or</i> ➤ Instruct the distributor who sold you the device to dispose of it. The agreements apply that were agreed with the distributor when purchasing the device (e.g. his general terms of payment and delivery). ➤ If your distributor is not able to take back and dispose of the device, please contact BINDER service.
---	--

Certified companies disassemble waste (used) BINDER equipment in primary substances for recycling according to Directive 2012/19/EU. The devices must be free from toxic, infectious or radioactive substances in order to eliminate any health hazards to the employees of the recycling companies.

8.4. Disposal of the BINDER Gas Supply Service and the power supply unit in non-member states of the EU

 	NOTICE Danger of violation against existing law if not disposed of properly. Failure to comply with applicable law. Alteration of the environment. ➤ For final decommissioning and disposal of the device, please contact BINDER service. ➤ Follow the statutory regulations for appropriate, environmentally friendly disposal.
--	---

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BINDER postal BINDER GmbH, Postfach 102, 78502 Tuttlingen, Alemania