P-5 premium

EASY GTM

READ THE INSTRUCTIONS PRIOR TO PERFORMING ANY TASK!

Gas pressure and temperature measuring system Gasdruck und
Temperatur - Mess - System

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P-5 premium

EASY GTM

READ THE INSTRUCTIONS PRIOR TO PERFORMING ANY TASK!

Gas pressure and temperature measuring system

Valid for: 50.9250.00, 50.9280.00, 50.9310.00

Translation of the original operating instructions

Version 06/2025 Index 003

1 Technical data

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

The FCC ID is: WAFAMB2520T.

The evaluation software runs with WINDOWS XP / 7 / 8 © operating systems!

EASY GTM Pressure measure range	0 bar - approx. 30 bar
Allowable pressure in the bowl	0 - 20 bar
Measurement - dissolution EASY GTM	1 mbar
Measure rate	up to 4 values / s
Temperature range transmitter unit	-40°C - 110°C
Operating time with charged accu	Depending to the operating temperature approx. 1/2 year of continuous operation!
Data - transfer	USB or optional receiver USB-Stick
Radio - transfer	2,4 Ghz
Reach the radio - transfer	typically 15m inside buildings
Battery	AA 1x Energizer ® Ultimate Lithium (81.0022.00)
Change pressure-measurement	< +/- 0,5%
Temperature-dependence	< +/- 1%
Pressure hysteresis	< +/- 0,1%

2 Gas Pressure and Temperature Measuring System



The EASY GTM gas pressure and temperature measuring system is used to control the grinding process and for mechanical alloying.

NOTICE

If grinding balls with a diameter of 1 mm or smaller are used, there is a risk of the through holes in the lid getting clogged, which could falsify the pressure and temperature measurement. If these ball sizes are used nevertheless, the grinding process must not be carried out unobserved. The intermediate lid should be checked at regular intervals during the grinding process. When the MillControl software is used, the pressure and temperature are displayed. If the pressure remains constant, the cause of this may be a clogged intermediate lid.

2.1 Case contents and system design

- 1 2.6 Nm torque spanner + hexagonal bit
- 3 EASY GTM system
- 4 O-ring 88.49 x 3.53, Viton
- 5 O-ring 21 x 2.5
- 6 Coupling, gassing valve
- 7 Hexagon offset screwdriver, 2.5 mm
- 8 GTM support ring
- 9 Transmission unit

- 10 Valve insert
- 11 Lid
- 12 Bowl
- 13 Hose clamps
- 14 Compressed air hose
- 15 Cleaning rods for valves
- 16 Valve screwdriver
- 17 Hexagon screwdriver, 3 mm


2.2 Inserting / changing the battery

When inserting it for the first time and during subsequent replacement of the battery, proceed as follows:

- 1.** → Lift off the support ring.

- 2.** → Then release the 6 hexagon screws (C) for closing the system using a hexagon screwdriver (7). (See → Chapter 2.9 'Cleaning the EASY GTM system' on page 13)

- 3.** → Lift the measuring/transmission unit (9) off of the lid (11).

- 4.** → The lock ring with the 6 hexagon screws can now be removed.


Gas Pressure and Temperature Measuring System

5. → Undo the two screws (D) of the battery compartment using a hexagon offset screwdriver, size 2.5 mm (7).

6. → Lift off the battery compartment lid and place it to the side. The battery compartment can be seen.


NOTICE

Please note that the battery (2) fits exactly in place and may need to be pushed in using a certain amount of force.

The service life of the battery (2) depends on the frequency of the radio protocols (see configuration) and essentially on the temperature inside the grinding bowl.


Also pay attention to the battery's charge level. If the battery's capacity falls below 10%, it should be replaced immediately.

7. → Remove the old battery and insert the new battery (2). Battery type: 1.5 V / AA Lithium


Pay attention to the polarity when inserting the battery!

8. → Ensure that the plug contacts are in the correct position when closing the battery compartment!

2.3 Configuration of transmitter ID and data transmission frequency

The measuring/transmission unit is equipped with four configuration switches for setting the transmitter ID and data transmission frequency.

The switches illustrated are all in the **OFF** position (actuator cams set to OFF).

2.3.1 Adjustment options of the transmitter detection

In the delivery condition all switches are set to OFF. It is recognised that, in this case, bowl 1 and the transmission frequency of 1 second are the set values. The bowl with ID 1 must then also be inserted in the indicated position 1 of the mill. Similarly the EGTM system with ID 2 must be positioned in position 2 of the mill.

ID	Switch 1	Switch 2
Transmitter no. 1	OFF	OFF
Transmitter no. 2	ON	OFF

If other switch settings are set, the EGTM systems will not be recognised. Ensure that the settings shown in the table are used for your systems and that they are used in the correct position in the mill.

2.3.2 Transmission frequency setting options

Transmission every	Switch 3	Switch 4
1 second	OFF	OFF
½ second	ON	OFF

Transmission every	Switch 3	Switch 4
¼ second	OFF	ON
automatic	ON	ON

1-second transmission is the default value.

For automatic transmission, the 1 / ½ / ¼ second transmission frequency is switched over in the event of fast changes in pressure or temperature.

The data transmission frequency naturally also has an impact on the battery's service life.

1-second operation results in the longest time span.

The transmission frequency can also be checked visually: as soon as the transmission unit is switched on, the switch flashes each time data is transmitted.

2.4 Installation of the transmission unit on the EASY GTM bowl

NOTICE

The threaded holes inside the bowl, the passage holes and the Allen screws in the locking ring have to be cleaned thoroughly after each grinding to allow a safe fastening of the EASY GTM system.


The EASY GTM system must be assembled as follows before inserting the EASY GTM system in the planetary micro mill:

1. → Put on the lid and secure it in place with the 4 hexagon screws!

2. → Apply the transmitter with seal (5)!

3. → Tighten the hexagon screws crosswise with a torque spanner (1) until it disengages! To be on the safe side, re-tighten all screws again after initially tightening them. All screws are secured against becoming loose during the grinding process once they have been tightened with the torque spanner as described.

Fritsch GmbH assumes no guarantee for damage caused by grinding with loose screws!

4. → Switch on the transmitter with the On/Off button (A) (LED flashes blue) and place the support ring on the bowl. Now position the bowl with the EGTM system as described in Chapter "Clamping the grinding bowls".

NOTICE

Check the fitting of the transmission unit again after inserting the bowl!

NOTICE

You must check that the hexagon screws are sitting correctly during breaks in grinding. If screws have loosened, they must be re-tightened.

2.5 Installation of the receiver board in the P-5 premium


! DANGER

Before starting the installation, disconnect the mains plug and protect the device against being unintentionally switched back on.

Secure installation work with a warning sign.

1. → Remove the layer of rubber on the work surface and loosen the 6 screws with which the plate is attached to the housing.

2. → Insert the receiver module in the corresponding slot. (As illustrated)
3. → Position the plate for the work surface, screw it down and put the rubber mat back on top.

When the receiver board is inserted, the P-5 premium automatically recognises when the EasyGTM is switched on.

Follow the following link to view a step-by-step instructional video for installing the receiver board. ➔ [Click here.](#)

You can also view the video directly via our homepage on the product page of the device under 'Downloads - Service Videos' or on our YouTube channel 'FRITSCH Milling and Sizing'.

2.6 Entering the temperature limit

FRITSCH		
EasyGTM		
Limit values		
	Max.	Min.
Pressure	5.0	2.0 bar
Temperature	60	30 °C
Mode	Stop	

DANGER

Do not exceed a maximum temperature of 125°C. Otherwise, the EASY GTM system and the machine will be damaged.

➔ Click on the Max. and Min. fields in the 'Temperature' section and enter the maximum and minimum temperature. If the maximum temperature is reached (125° C), the device goes into cooling mode or stops until the temperature falls below the minimum temperature limit. Once the minimum temperature limit has been reached, grinding continues as normal.

If the max. temperature is reached, the device stops the grinding process and goes into cooling mode! Depending on the setting of the operating mode, see ➔ Chapter 2.8 'Selecting the operating mode' on page 13, the device stops the grinding process and opens the grinding chamber cover or reduces the speed gradually until a minimum speed is reached in order to achieve an optimal cooling effect.

2.7 Entering the pressure limit

FRITSCH		
EasyGTM		
Limit values		
	Max.	Min.
Pressure	5.0	2.0 bar
Temperature	60	30 °C
Mode	Stop	

DANGER


Do not exceed a maximum pressure of 20 bar. Otherwise, the EASY GTM system and the machine will be damaged.

➔ Click on the Max. and Min. fields in the 'Pressure' section and enter the maximum and minimum pressure. If the minimum pressure (20 bar) is reached, the device goes into cooling mode or stops until the pressure falls below the minimum pressure limit. Once the minimum pressure limit has been reached, grinding continues as normal.

If the max. pressure is reached, the device stops the grinding process and goes into cooling mode! Depending on the setting of the operating mode, see ➔ Chapter 2.8 'Selecting the operating mode' on page 13, the device stops the grinding process and opens the grinding chamber cover or reduces the speed gradually until a minimum speed is reached in order to achieve an optimal reduction of pressure.

2.8 Selecting the operating mode

You can use the 'Mode' menu item to select how the P-5 premium is to respond to the temperature or pressure limit being reached.

Slow:

With this setting, the P-5 premium gradually reduces the speed down to a minimum speed. Once the lower temperature or pressure limit has been reached, P-5 premium grinding continues as normal.

Stop:

As soon as one of the two limits is reached, the ball mill shuts down immediately. If the limit is fallen short of, restart is disabled.

Press the "<" button to confirm all entries and return to the main menu.

2.9 Cleaning the EASY GTM system

NOTICE

The threaded holes inside the bowl, the passage holes and the Allen screws in the locking ring have to be cleaned thoroughly after each grinding to allow a safe fastening of the EASY GTM system.


NOTICE

The measuring/transmission unit may not be immersed in water. It can be rubbed down with a damp cloth if necessary.

Remove the lock ring (18), incl. screws, from the measuring/transmission unit and remove all soiling.

The screws are equipped with a circlip to prevent them from falling out of the lock ring (18). Ensure that the circlips are sitting correctly! (See illustration!)

The lid (11) and the bowl (10) can be cleaned under running water as described in Chapter "Cleaning".

2.9.1 Cleaning the sensor and replacing the seal

The sensor (G) must not be subjected to any mechanical stress. It is permitted to cleaning it with a wet cloth. The O-ring (5) can be replaced if necessary.

2.10 Battery arrangement

According to the "DIRECTIVE 2006/66/EC OF THE EUROPEAN PARLIAMENT AND COUNCIL as of 6th September 2006 on batteries and rechargeable batteries and on old batteries and old rechargeable batteries and for the annulment of Directive 91/157/EEC" and article 1, §18 and article 2 of the German Act on the re-definition of the waste legislation on the product responsibility for batteries and rechargeable batteries (BattG) as of 25.06.2009, we are legally obliged as the manufacturer to inform you as the consumer about the following:

The Easy GTM product we sell contains a lithium battery, type AA, 1.5V, that can be discharged once.

When the battery is flat, it may not be disposed of along with household waste. Old batteries may contain harmful substances that can pollute the environment or damage your health. Please deliver the batteries to the regional recycling point/collection point. Please only place discharged/flat batteries in the containers intended for this purpose. Cover the terminals of lithium batteries with adhesive tape first.

All batteries and rechargeable batteries are reused. Valuable substances, such as zinc, iron and nickel can thus be recycled. Battery recycling is one of the easiest environmental protection measures. The crossed-out dustbin symbol means that (rechargeable) batteries may not be disposed of along with household waste.

You can, of course, also return discharged batteries - with sufficient postage - by post to:

Fritsch GmbH
Milling and Sizing
Industriestrasse 8
D-55743 Idar-Oberstein, Germany

For information on removing the battery, see ➔ *Chapter 2.2 'Inserting / changing the battery' on page 7.*

3 General Safety

3.1 General

This workshop manual is intended for personnel technically qualified to install, operate and maintain the above mentioned Fritsch Accessories. It contains all the necessary information for correct maintenance of the above mentioned Fritsch Accessories. The contents of this workshop manual are not contractual and cannot under any circumstance extend or restrict contract warranty clauses.

3.2 Qualification of personnel

Only qualified personnel are authorized to install, operate or maintain the above mentioned Fritsch Accessories. Any work performed by an unqualified person or failure to follow the safety instructions in this document or on the equipment may risk the safety of personnel and/or cause irreparable damage to equipment. The following personnel may be regarded as being "Qualified":

- Personnel who are trained and experienced in the adjustment and repair of Fritsch equipment (for example, installation engineers, post sales service engineers, etc).

3.3 Conformity of use

The products described in this workshop manual conform to the European Directives to which they are subjected (CE marking). However, they can only be used correctly in the context of the applications for which they are intended (described in the various documents). If all handling, transport and storage specifications are observed, and all instructions for installation, operation and maintenance are followed, the products will perform properly, with no danger to personnel or equipment.

3.4 Maintenance

It is imperative to observe the following rules when maintenance for the equipment has to be carried out.

- Safety instructions must be followed. These instructions are in the documentation or on the equipment being installed and set up.
- Always remove power before performing any operation on equipment.
- Before the equipment is powered up, check its nominal voltage to ensure that it has been adjusted to conform with the local supply voltage.
- All measures must be taken to ensure that any power return does not lead to a dangerous state which may risk personnel or the installation.
- Emergency stop and safety devices must remain effective in all the equipment's operating modes, even during abnormal operations (for example, when a wire becomes disconnected). Resetting these devices must not cause uncontrolled or improper restarts.
- Ensure that the equipment is protected against unintentional operation.
- Appropriate safety measures must be taken for the inputs and outputs (including function keys) to prevent improper states in the equipment, if no signal is received.

- Cables which carry signals must be located where they do not cause capacitive, inductive or electromagnetic interference with the equipment functions.
- Low voltage circuits should always have a ground connection in its power supply. Low voltage circuits must be protected.

3.5 Safety Checks

The BGV A3 (accident prevention regulations of the professional association) requires that electrical safety standards for electrical installations and fixed / mobile equipment are duty. Safety Checks according to these regulations were done and documented by Fritsch before delivery of the equipment.

Two standards for the examination of electrical safety (eS) after maintenance must be used:

- VDE 0701 for repair and maintenance According to the VDE 0701, after repair-, service- or maintenance work is re-quired, to do the relevant professional examinations according to BGV A3, if live parts of the unit were involved. Else it should be noted in the service report that only a basic maintenance was done and no safety examination after BGV A3 has been carried out.
- VDE 0702 for reassessment Operators are responsible to instruct appropriate tests according to the VDE 0702 by a specialist once a year. This review, with a compulsory test with certain information is indicated by a plaque on the equipment. The examination of electrical safety to VDE 0702 may be done by qualified service personnel only as an offer to the operators / users.

4 Exclusion of liability

Before using the product, be sure to have read and understood this operating manual.

The use of the product requires technical knowledge; only commercial use is permitted.

The product may be used exclusively within the scope of applications set down in this operating manual and within the framework of guidelines put forth in this operating manual and must be subject to regular maintenance. In case of non-compliance, improper use or improper maintenance, the customer assumes full liability for the functional capability of the product and for damage or injury arising from violating these obligations.

The contents of this operating manual are subject in entirety to copyright law. This operating manual and its contents may not be copied, further distributed or stored in any form, in part or in whole, without the prior written consent of Fritsch.

This operating manual has been prepared to the best of our knowledge and checked for accuracy at the time of printing. FRITSCH GMBH assumes no guarantee or liability whatsoever for the accuracy or completeness of the contents of this operating manual, including but not limited to the implied warranties of merchantability and fitness for a particular purpose, unless liability is expressly prescribed by applicable laws or jurisprudence.

FRITSCH GMBH expressly reserves the right to modify and/or update this operating manual without prior notice. The same applies to modifications and improvements to the products described in this operating manual. It is the responsibility of the user to ensure that they have the current version of this operating manual. For more information, please contact your local FRITSCH GMBH distributor or Fritsch GmbH, Industriestr. 8, D-55473 Idar-Oberstein.

Not all parts shown here are necessarily installed in the product. The buyer is not entitled to delivery of these parts. If interested, please contact your local FRITSCH GMBH distributor or Fritsch GmbH, Industriestr. 8, D-55473 Idar-Oberstein.

FRITSCH GMBH takes the greatest care to ensure that the quality, reliability and safety of your products are continuously improved and adapted to the state of the art. The supplied products as well as this operating manual conform to the current state of the art when they leave the sphere of influence of FRITSCH GMBH.

By using the product the customer agrees with this and recognizes that defects, malfunctions or errors cannot be completely excluded. To prevent risk of damage to persons or property or of other direct or indirect damage, resulting from this or other causes, the customer must implement sufficient and comprehensive safety measures for working with the product.

Fritsch GmbH excludes any liability, warranty, or other obligation to compensate for damages, regardless of whether this liability, warranty, or other obligation is explicit or implicit, contractual or arising from unlawful acts or prescribed contractually, by law, or otherwise. In no event shall the buyer be entitled to any compensation from Fritsch GmbH for any special, direct, indirect, coincidental or consequential damage, including but not limited to lost profits, lost savings, lost sales or financial loss of any kind or for compensation of third parties, for downtimes, for lost goodwill, for damage to or replacement of equipment and property, for costs or restoration of materials or goods related to the product or the use of our products, for other damage

or injury to persons (including fatal injuries) or similar. The above exclusion of liability is limited by mandatory liability as prescribed by laws or jurisprudence. Liability for negligence is excluded in all cases.

No permission is given expressly, implicitly or otherwise for the use of patents, brands or other copyrights. We also assume no liability for copyright infringements or infringements of the rights of third parties arising from the use of this product.

Neither compliance with this operating manual nor the conditions and methods used during installation, operation, use and maintenance of the product can be monitored by Fritsch GmbH. Improper execution of the installation can result in property damage and thus endanger persons. Therefore, we assume absolutely no responsibility or liability for loss, damage or costs that result from errors at installation, improper operation or improper use or improper maintenance or are in any way connected to these.

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