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
(D)

- 1 Batteriekappe
- 2 On-Off/Hold Knopf
- 3 Kalibrier Knopf
- 4 Eintauchniveau
- 5 Schutzkappe

1 ANBRINGEN DER BATTERIEKAPPE

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1 INSTALLING BATTERY CAP

The unit is shipped with the battery cap open. Close the battery cap by pressing the cap on a hard surface until the latch clicks, indicating a secure lock.

2 REPLACING THE BATTERIES

1. Lift the latch on the battery cap with a small flat object, such as a small screwdriver. Use your thumb to push the cap forwards as shown. Firmly grasp the cap and carefully separate it from the pH Meter body.

2. Ensure that the batteries are orientated correctly (correct orientation is marked inside the battery cap).

3 CALIBRATION

This tester is factory calibrated. But due to prolonged storage, the unit should be re-calibrated before use.

NOTE: Regular calibration is necessary to maintain the accuracy. Depending on usage, perform a check as regularly as possible before each series of test or if accuracy is desired. Calibrate as required. Calibration should be performed at room temperature at about 25°C or 77°F.

1 Use Essentials Buffer 7 solution for calibration.

2 Remove protective cap. Always rinse the sensor area with water and shake the tester in the same way you would use a mercury thermometer, before each and every test.

3 Place a small amount of Buffer 7 solution in a small container and place the sensor area of the pH Pen in the Buffer solution until it is fully immersed. Tap or shake a little to remove bubbles.

4 Press and hold the CAL button until it displays CAL blinking. When the buffer standard is recognised, 7.00 will be displayed in a blinking mode while it senses for an end-point.

5 When the display stops blinking and freezes, it indicates that the unit has been standardised.

6 Rinse the sensor area thoroughly with water. Calibration is completed.

7 To make measurement, press the 'ON-OFF / Hold' button to exit calibration mode. Alternative buffer standards pH 4.0 or 10.0 can be used if desired.

4 HANDLING PRECAUTIONS

Do not touch, rub or scratch the sensor. It is very delicate and might break or lose its sensitivity. Do not store unit without the protective cap. Chemical in the unit will expire faster, thus shortening its life span. Do not submerge the unit underwater. Though the unit is splash proof and water resistant, it cannot come under high pressure underwater and is beyond repair if water gets into the unit. If it is dropped into water, retrieve it immediately and wipe dry with a cloth.

Do not store the unit under high temperature or direct sunlight. This will shorten the life span of the unit. Do not clean unit with thinner or solvents. This will damage the unit. Use only a damp cloth to clean unit if needed.

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

5 TAKING A MEASUREMENT

- 1 Remove protective cap from bottom (See product layout.)
- 2 Press the 'ON-OFF / Hold' button once to switch on.
- 3 Shake the tester in the same way you would use a mercury thermometer, before each test.
- 4 Dip tester into the sample solution up to the immersion level. Shake the sensor area in solution to remove bubbles and wait for about a minute or when reading stabilizes.
- 5 Pressing the ON-OFF/Read-Hold button once will freeze the reading. Pressing again will release it for another readout.
- 6 If the sensor is dry, a slow response will result. Dip the sensor area in a cup of water or preferably Buffer 7 calibration solution for 1-2 hours before testing again.
- 7 Always rinse the sensor area with water and blot it dry before and after each test.

AUTO-LOCK MEASUREMENT

- 1 The Auto-Lock feature is user friendly, allowing the tester to take a stable reading without guesswork.
- 2 To enable auto-lock measurement, press both the ON-OFF/Read-Hold button and the CAL button and hold till a blinking small 'A' sign appears on the right side of display.
- 3 From now on, everytime you dip the tester into the solution, press the ONOFF button once, the display digit will blink while sensing for a stable reading. Once the reading is stabilised, the display digit will stop. You can now take the reading.
- 4 Press the ON-OFF button again for another reading.
- 5 To disable the Auto-Lock feature, press both the ON-OFF button and the CAL button till the small 'A' sign disappears. Now you can resume standard measurement as before.
- 6 To Switch off the tester, press the 'ON-OFF/Read-Hold' button for 3 seconds.
- 7 Replace protective cap before storing away.

6 MAINTENANCE

- When the battery symbol appears on the display, this indicates a low battery and only 2 hours of continuous use remain. Though the unit may continue to function, the accuracy of the unit will be affected beyond the 2 hours.
- To change batteries, turn to the back of the unit and lift the latch on the battery case from the top of unit, (see layout) slide the clip upward to open. Replace all four batteries according to polarity. Push cap back till latch clicks in place.
- If the unit is not used or stored for a long period of time, the sensor will become dry. This will result in a slow response to a stable reading. Soaking the sensor area in a cup of tap water or preferably Buffer 7 solution for 1-2 hours will restore sensitivity to the sensor.
- Keep in mind that all pH sensors age with time and usage. Therefore, re-calibration is necessary to maintain accurate reading.
- Note that the pH sensor has a limited life span of about 365 tests. When the unit fails to calibrate or responds very slowly, it means that the unit should be replaced. It is not possible to repair broken sensor, defective or expired unit.

7 SPECIFICATION

RANGE:	0.0 to 14.0 pH
RESOLUTION:	0.1 pH
ACCURACY:	±0.2 pH
BATTERY:	4 x 1.5V A76/LR44
BATTERY LIFE:	Approx. 150 hours
AUTO SHUT-OFF:	Approx. 15 minutes
PRODUCT LIFE:	365 tests
CASING MATERIAL:	ABS plastic
OPERATING TEMPERATURE:	0° to 50°C/32 to 122°F
SIZE (LxWxH):	170 x 32 x 15mm
WEIGHT:	Approx. 70 gm

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