[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |
[Руководство по эксплуатации](#) | [사용 지침](#) | [Bruksanvisning](#) | [使用説明書](#) |
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

Titrette®

Flaschenaufsatzbürette | Bottle-Top Burette

Impressum

BRAND GMBH + CO KG

Otto-Schott-Str. 25
97877 Wertheim (Germany)

T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

Do you need more operating manuals and translations?
Please refer to <http://www.brand.de/om> or use the following QR
Code:

The original operating manual is written in German. Other languages are translations of the original operating manual. Technical changes, errors or misprints reserved.

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The manufacturer is solely responsible for issuing this declaration of conformity with regard to meeting the essential requirements and preparing the technical documentation.

Hersteller / Manufacturer	Adresse / Address
BRAND GMBH + CO KG	Otto-Schott-Str. 25 97877 Wertheim Germany

Wertheim, 13. Mai 2022 / 13 May 2022

Patrick Ziemeck
Managing Director Technology

i.A. Siegfried Ott
Regulatory Affairs

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This document declares the accordance with the named harmonized regulations, but does not assure specific properties.

UK Declaration of Conformity

Product name: Bottle-top burette Titrette®

Product variations: 10 ml, 25 ml und 50 ml

Cat.-No. 4760141; 4760151; 4760161; 4760241; 4760251; 4760261; 4760451

The product described above meets the basic requirements that are specified in the harmonization legislation listed below:	Applied harmonized standards:
The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012 (SI 2012 No. 3032)	EN IEC 63000:2018
Electromagnetic Compatibility Regulations 2016 (SI 2016 No. 1091, as amended by SI 2019 No. 696)	EN 61326-1:2013
Other applied standards:	Scope:
EN 61010-1:2010/A1:2019/AC:2019-04, EN 61010-1:2010/A1:2019	Electrical Equipment (Safety) Regulations 2016

The manufacturer is solely responsible for issuing this declaration of conformity with regard to meeting the essential requirements and preparing the technical documentation

Manufacturer	Address
BRAND GMBH + CO KG	Otto-Schott-Str. 25 97877 Wertheim Germany

Wertheim, 13 May 2022

Patrick Ziemeck
Managing Director Technology

i.A. Siegfried Ott
Regulatory Affairs

This document declares the accordance with the named harmonized regulations, but does not assure specific properties.

Declaration of conformity – China RoHS 2

BRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in BRAND products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a “Product Conformity Assessment” (PCA) procedure was performed. As defined in GB 26572-2025 the “Maximum Concentration Value” limits (MCV) apply to these restricted substances:

+ Lead (Pb):	0.1 %	+ Hexavalent chromium (Cr(+VI)):	0.1 %
+ Mercury (Hg):	0.1 %	+ Polybrominated biphenyls (PBB):	0.1 %
+ Cadmium (Cd):	0.01 %	+ Polybrominated diphenyl ether (PBDE):	0.1 %
+ Bis(2-ethylhexyl) phthalate (DEHP):	0.1 %	+ Dibutyl phthalate (DBP):	0.1 %
+ Butyl benzyl phthalate (BBP):	0.1 %	+ Diisobutyl phthalate (DIBP):	0.1 %

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for BRAND instruments is 40 years.

Material Content Declaration for BRAND Products

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), polybrominated diphenyl ethers (PBDE), bis(2-ethylhexyl) phthalate (DEHP), butyl benzyl phthalate (BBP), dibutyl phthalate (DBP) and diisobutyl phthalate (DIBP).

Products manufactured by BRAND may enter into further devices or can be used together with other appliances. With these third party products and appliances in particular, please note the EFUP labeled on these products. BRAND will not take responsibility for the EFUP of those products and appliances.

Place, date: Wertheim, 15|07|2025

Patrick Ziemeck
 (Managing Director Technology)

i.A. Siegfried Ott
 (Regulatory Affairs)

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Table of contents

1 Introduction	58	12.3 Loosen stuck filling valve	90
1.1 Scope of delivery	58	12.4 Replacing the batteries	91
1.2 Use operating manual	58	13 Troubleshooting	92
2 Safety regulations.....	60	14 Product markings	94
2.1 General safety regulations	60	15 Ordering Information	95
2.2 Function.....	60	16 Accessories/spare parts.....	96
2.3 Limitations of Use	61	17 Repairs	101
2.4 Operating Limitations.....	61	17.1 Sending for repair.....	101
2.5 Operating exclusions	61	18 Calibration service	103
2.6 Battery specifications	61	19 Information about your laboratory in-	strument.....
2.7 Storage conditions.....	62	20 Warranty.....	104
2.8 Recommended application range ..	62	21 Disposal	105
3 Functional and operational compo-	63		
nents	63		
4 Assembly	65		
4.1 First Steps.....	65		
4.2 Priming	67		
5 Titrating.....	69		
5.1 Energy-saver mode (Auto power			
off).....	70		
5.2 Pause function	70		
6 PC interface (optional)	71		
7 Sensitive media (replacing the inspec-	72		
tion window)	72		
8 Installing a drying tube (optional)	73		
9 Error limits	74		
10 Checking the Volume (Calibration)	75		
11 Additional functions.....	76		
11.1 CAL mode (adjustment)	76		
11.2 Calibration date (GLP mode)	80		
11.3 Auto power off (APO mode)	81		
11.4 Decimal places (dP mode)	82		
12 Cleaning	83		
12.1 Standard cleaning.....	83		
12.2 Intensive cleaning	85		

1 Introduction

1.1 Scope of delivery

Bottle-top burette Titrette®, for GL 45 threaded bottles, size 10 ml, 25 ml or 50 ml, telescopic filling tube (length 170 - 330 mm), recirculation tube, 2 micro-batteries 1.5 V (AAA/UM4/LR03), 4 PP bottle-top adapters (GL 32-33, GL 38, S 40, NS 29/32), 2 tinted light-protection viewing windows, one performance certificate and this operating manual.

1.2 Use operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find updated versions of the instructions for use at www.brand.de/om.

1.2.1 Signal words and their meaning

Signal word	Meaning
 WARNING or  WARNING! ...	WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.
 CAUTION or  CAUTION! ...	CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.
NOTICE or NOTICE ...	NOTE is used to address actions that are not related to physical injury. Example: Possible property damage.

1.2.2 Symbols used in the instructions

 The sign warns of a possible hazard.

 The sign warns of hazardous media and the dangers that can arise from liquid splashes.

 The symbol indicates that the instruction manual and safety information must be observed.

 The sign suggests that a face shield be worn during work.

 The sign suggests that eye protection be worn.

 The sign suggests that you to wear protective clothing.

 The sign suggests that protective gloves be worn.

— —

1.2.3 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 General safety regulations

Please read carefully!

The laboratory device Titrette® Can be used in combination with hazardous materials, operations and equipment. However, the instructions for use cannot identify all the safety problems that may occur. It is the responsibility of the user to ensure compliance with safety and health regulations and to determine the appropriate restrictions before use.

1. Every user must read and observe these instructions for use before using the appliance.
2. Follow general hazard warnings and safety instructions, e.g. wear protective clothing, eye protection and protective gloves.
3. Observe the reagent manufacturer's instructions.
4. Do not operate the appliance in a potentially explosive atmosphere.
5. Only use the device for titrating liquids and only within the defined limits and restrictions of use.
Observe exclusions of use (see Operating exclusions, p. 61)!
If in doubt, always contact the manufacturer or dealer.
6. Always work in such a way that neither the user nor other persons are endangered. Avoid splashes. Only use suitable containers.
7. Never turn the handwheels while the titration cannula is closed with the cap.
8. Never remove the titration cannula when the glass cylinder is full.
9. Reagent may accumulate in the cap of the titration cannula. Therefore clean regularly.
10. Use a bottle holder for small bottles to avoid tipping.
11. Never carry the device mounted on the reagent bottle by the housing. Breakage or detachment of the device from the reagent bottle can lead to injuries.
12. Never use force.
13. Only use original accessories and original spare parts. Do not make any technical changes. Do not disassemble the appliance any further than described in the operating instructions!
14. Always check that the appliance is in good condition before use. Inadequately cleaned or inspected devices may come into contact with the media by the user. If you notice any malfunctions in the appliance (e.g. sluggish piston, stuck valves or leaks), stop titrating immediately and follow Troubleshooting, p. 92 . If necessary, contact the manufacturer.
15. The built-in 1.5 V micro-batteries are not rechargeable!

2.2 Function

The bottle-top burette Titrette® with electronic digital display is used for titrating aqueous and non-aqueous titration media (e.g. alcoholic KOH) up to a concentration of max. 1 mol/l. (see Recommended application range, p. 62). By using a high-precision measuring system, even the tight tolerances of class A for glass burettes can be met. The devices are labeled DE-M.

2.2.1 Handling

If handled properly, the dispensed liquid comes into contact with the following chemical-resistant materials only:

Borosilicate glass, Al_2O_3 , ETFE, PFA, FEP, PTFE, platinum-iridium, PP (screw cap). The instrument comes equipped with a recirculation valve.

2.3 Limitations of Use

This instrument is intended for titrating liquids within the following physical limitations:

- operating temperature from +15 °C to +40 °C (from 59 °F to 104 °F) of instrument and reagent
- vapor pressure up to 500 mbar
- viscosity up to 500 mm²/s
- altitude: max. 3,000 m above sea level
- relative humidity: 20 % to 90%

2.4 Operating Limitations

- Chlorinated and fluorinated hydrocarbons or chemical combinations which form deposits may make the piston difficult to move or may cause jamming.
- When working with crystallizing media, follow the cleaning instructions (see Cleaning).
- Compatibility of the instrument for the intended application (e.g., trace analysis) must be carefully checked by the user. Contact the manufacturer, if necessary.

2.5 Operating exclusions

2.5.1 Titrette

Never use the instrument for:

- liquids that attack borosilicate glass, Al_2O_3 , ETFE, PFA, FEP, PTFE or platinum-iridium (e.g., hydrofluoric acid).
- suspensions (e.g., of charcoal), as solid particles may clog or damage the instrument.
- concentrated acids and bases as well as nonpolar solvents that have a strong swelling effect on plastics (e.g., toluene, benzene).
- carbon disulfide, as this media is highly flammable.
- Autoclaving the instrument is not permitted!
- The instrument must never be exposed to an aggressive atmosphere (e.g., HCl fumes).

2.6 Battery specifications

2 micro-batteries, 1.5 V (AAA/UM4/LR03), not rechargeable!

2.7 Storage conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from -20 to 50°C (-4 to 122°F).

Relative humidity: 5 % to 95 %.

2.8 Recommended application range

The instrument can be used for the following titration media (max. concentration 1 mol/l):

Medium	Medium	Medium
Alcoholic potassium hydroxide solution	Potassium bromate solution	Oxalic acid solution
Ammonium iron (II) sulfate solution	Potassium bromide-bromate solution	Perchloric acid
Ammonium thiocyanate solution	Potassium dichromate solution	Perchloric acid in glacial acetic acid
Barium chloride solution	Potassium iodate solution	Nitric acid
Bromide bromate solution	Potassium permanganate solution*	Hydrochloric acid
Cerium (IV) sulfate solution	Potassium thiocyanate solution	Hydrochloric acid in acetone
EDTA solution	Sodium arsenite solution	Sulphuric acid
Iron (II) sulfate solution	Sodium carbonate solution	Silver nitrate solution*
Acetic acid	Sodium chloride solution	Tetra-n-butylammonium hydroxide solution
Iodine solution*	Sodium nitrite solution	Triethanolamine in acetone*
Iodine iodate solution*	Sodium thiosulfate solution	Zinc sulfate solution
Potassium hydroxide solution	Sodium hydroxide solution	


* Use light shield inspection windows (see Sensitive media (replacing the inspection window), p. 72)

This table has been carefully checked and is based on the current state of knowledge. Always follow the instructions for use of the device and the information provided by the reagent manufacturer. If you require information on chemicals that are not included in the list, please contact BRAND.

Status as of: 0219/4

3 Functional and operational components

The controls: Separate buttons for on/off, pause and CLEAR to clear the display. Handwheels with grip grooves for fast to dropwise titration.

International patents have been filed for key features of the Titrette® bottle-top burette. The installation wrench is located in the rear housing shell.

Telescoping filling tube and recirculation tube

Assembly tool

4 Assembly

⚠ WARNING


Risk of injury from hazardous media

- Wear your personal protective equipment, in particular protective clothing suitable for the medium to be used, eye protection and protective gloves!
- Before use, find out how you can use the medium safely and whether the laboratory device is suitable for it.
- Observe the chemical-specific safety data sheets and follow the operating instructions for the medium at the workplace!


4.1 First Steps

1. Inserting batteries


- a. Unscrew the air vent cap by hand or use a coin.


- b. Remove the rear housing shell.
- c. Remove the battery compartment covers and insert the batteries. Ensure that the polarity of the batteries corresponds to the symbols on the battery compartment.

- d. Close the battery compartments tightly again with the covers. Carefully press down the edges of the cover so that they fit tightly everywhere and there is no gap to the battery compartments.
- e. First hook in the rear housing shell at the top, then close it.
- f. Screw in the ventilation plug.

2. Switching the instrument on and off

- a. Briefly press the On/Off button to switch on or off.


WARNING

Follow the safety instructions

- Wear protective clothing, eye protection and protective gloves!
- Always wear protective gloves when touching the instrument or the bottle, especially when using dangerous liquids.
- Follow all safety instructions and comply with the limitations of use, see Limitations of use, p. 61.
- Comply with the operating limitations, see Operating limitations, p. 61.

3. Installing the filling tube/recirculation tube

- a. Adjust the length of the telescopic filling tube to the bottle height and attach it. Center the filling tube (side with smaller diameter) and attach it carefully to avoid damaging the olive-shaped nozzle.
- b. Insert the recirculation tube with the opening pointing outward.

NOTICE

When working with media with a high rate of crystallization, e.g. alcoholic KOH, the length of the telescopic filling tube should be adjusted to a distance of approx. 20 mm from the bottom of the bottle.

4. Mounting and aligning the instrument on a bottle

NOTICE

The included adapters are made of polypropylene (PP) and may only be used for media that do not attack PP (Accessories/spare parts, p. 96).

- a. Screw the instrument (GL 45 threads) onto the reagent bottle and then align the titrating tube with the bottle label. This is done by rotating the valve block with the titrating tube.
- b. For bottles with other thread sizes, select a suitable adapter.

The titrating tube can be adjusted by 70 mm, both horizontally and vertically.

5. Transporting the instrument

- a. When mounted to a reagent bottle, always carry the instrument as shown in the figure!

⚠ WARNING

Avoid splashing the reagent

Potential health risk, especially with hazardous media

- Never turn the hand wheels when the valve is set to 'Titrate' and the titrating tube is closed with the screw cap!
- Avoid splashing the reagent! The reagent can drip out from the titrating tube and screw cap.

4.2 Priming

NOTICE

Before first use:

After final production inspection, glycerin and ethanol residues may still be present in the instrument. To prevent media from mixing with these residues, rinse the instrument thoroughly before the first use and discard the first dispensed volumes. Avoid splattering.

b.

- a. Ensure that the screw cap of the titrating tube is screwed on tightly.
- b. Turn the valve in the direction of the arrow to 'Recirculate'.

c.

- c. First, turn the hand wheels to move the piston all the way down. To fill, rotate the piston no further than half-way up, and then empty it again.

NOTICE

If filling is not possible, see Troubleshooting, p. 92

Then rotate the hand wheel half way several times to aspirate liquid, each time using a single stroke to the lower position to empty it into the bottle. Repeat the procedure approximately 5 times until there are no more large air bubbles under the piston.

NOTICE

A few bubbles less than 1 mm in size are permitted.

e.


- d. Unscrew the screw cap of the titrating tube.
- e. Turn the valve to 'Titrate'.

f.

- f. Hold a suitable receiving vessel under the opening of the titrating tube and dispense liquid to deaerate the titrating tube until it is bubble-free. Wipe away any remaining drops from the titrating tube tip.

5 Titrating

1. Switching on the instrument


- a. Briefly press the On/Off button to switch on or off.

2. Filling the instrument

- a. Fill the instrument up to the upper stop position by gently turning the hand wheels.
b. Briefly press the CLEAR button one time to set the display value to zero.

3. Titrating

- a. Hold a suitable receiving vessel under the opening of the titrating tube.
b. Dispense liquid up to the transition point by turning the hand wheels.

NOTICE

If the fill volume is insufficient for the entire titration, gently turn the hand wheels back to the upper stop position to fill (the displayed value remains unchanged during this process). Then continue titration.

4. Filling the instrument after titration


- a. To reduce crystal deposits and evaporation, always fill the instrument all the way up to the upper stop position after titration.

5.1 Energy-saver mode (Auto power off)

If work is interrupted for more than three minutes (factory setting), the instrument automatically switches into standby mode. The display value is saved and reappears in the display after switching back on manually. The user can adjust the amount of idle time before automatically switching off (see Auto power off (APO mode), p. 81).

5.2 Pause function

If the instrument has not been completely primed, air bubbles appear in the titrating tube during titration. To prime the instrument during titration, liquid can be dispensed into a different receptacle using the Pause function without the display value changing.

- a. Press the Pause button one time.
- ⇒ The Pause signal blinks.
- b. Prime the instrument, dispense liquid, etc. See Priming, p. 67.
- c. To end the Pause function, briefly press the Pause button again.
- ⇒ The Pause signal disappears.
- d. Continue titration

6 PC interface (optional)

The device is available with an optional communication interface (RS 232) and can be connected to a USB port using a commercially available USB/RS 232 adapter and read the data from the device using the Titrette® software. This version of the device offers the following advantages:

- Transcription errors are eliminated when copying the primary data because the titration results are automatically transmitted to the PC by double-clicking on the CLEAR button. This fulfills an important GLP requirement.
- All raw data is recorded simultaneously. With each data transfer, the burette sends the titrated volume, the serial number of the instrument, the nominal volume and the adjustment value, as well as the next scheduled calibration date.

The data transmitted by the instrument is treated as keyboard input from the PC. This universal input format ensures that the instrument is compatible with all PC applications that accept keyboard inputs.

You can download and install the software free of charge at the following address:


<https://update.brand.de/Titrette/Titrette-Software-latest.exe>

The user manual for the software is available in the webshop:


<https://shop.brand.de/707538>

7 Sensitive media (replacing the inspection window)


For light-sensitive media (e.g., iodine, potassium permanganate and silver nitrate solutions), we recommend using the colored light shield inspection window.


- a.** Unscrew the air vent cap by hand or use a coin.

- b.** Remove the rear housing shell.

- c.** Unclip the rear inspection window on one side and take it out.
- d.** Insert the colored inspection window with the smaller curvature into the rear housing.

- e.** To change the front inspection window, lift up one corner of the window, e.g., with a fingernail, and take it out.
- f.** Insert the colored inspection window with the larger curvature into the front housing.
- g.** To secure the rear housing, first hang it from the top, then snap it closed and screw in the air vent cap.

8 Installing a drying tube (optional)

For moisture- or CO₂-sensitive media, it may be necessary to use a dry pipe filled with suitable absorbent (not included in the scope of delivery).

(See Accessories/spare parts)

- a. Screw in the filled drying tube in place of the air vent cap.

NOTICE

Seal the threads of the drying tube, the bottle and/or the bottle adapter with PTFE tape.

9 Error limits

Error limits in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20 °C/68 °F) of instrument, ambient environment and distilled water. The test was carried out in accordance with DIN EN ISO 8655-6 with a completely filled instrument and uniform and jerk-free dispensing.

Error limits

		Bottle-top burette Titrette®					Bottle-top burettes according to DIN EN ISO 8655-3				Glass burette Class A according to DIN EN ISO 385
Volume ml	Partial volume ml	R* ± %	± µl	VK* %	µl	R* ± %	±µl	UK* %	µl	FG** ± µl	
10'	10	0,10	10	0,05	5	0,3	30	0,1	10	20	
	5	0,20	10	0,10	5	0,6	30	0,2	10	20	
	1	1,00	10	0,50	5	3	30	1	10	20	
25'	25	0,072	18	0,024	6	0,2	50	0,1	25	30	
	12,5	0,144	18	0,048	6	0,4	50	0,2	25	30	
	2,5	0,72	18	0,24	6	2	50	1	25	30	
50'	50	0,06	30	0,02	10	0,2	100	0,1	50	50	
	25	0,12	30	0,04	10	0,4	100	0,2	50	50	
	5	0,60	30	0,20	10	2	100	1	50	50	

* R = correctness, VK = coefficient of variation ** FG = margin of error

The titration volume is displayed in 1 µl increments for device sizes 10 ml and 25 ml and in 2 µl increments for device size 50 ml. From 20 ml titration volume there is an automatic switchover to 10 µl steps.

The drop size for the 10 ml device is approx. 20 µl and for the 25 ml and 50 ml device approx. 30 µl.

NOTICE

The maximum total error for a single measurement can be calculated approximately from the sum of the error limits $FG = R + 2 \cdot VK$. This is a maximum of $25 \text{ ml} \pm 30 \text{ µl}$ for the 25 ml size and $50 \text{ ml} \pm 50 \text{ µl}$ for the 50 ml size.

This complies with the error limits for class A glass burettes in accordance with DIN EN ISO 385.

10 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. For calibration, set the 10 ml instrument to 3 decimal places. This time frame should be adjusted to correspond with individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de. You can also carry out a functional test at more frequent intervals, e.g., by titrating against a standard solution. For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. A demo version can be downloaded from www.brand.de.

Gravimetric volume testing according to DIN EN ISO 8655- 6 (for measurement conditions, see Error limits, p. 74) is performed as follows:

1. Prepare the instrument

Clean the instrument (Cleaning), fill it with test liquid (distilled water) and then prime it carefully.

2. Check the volume

- Dispense 5 drops in a separate vessel and wipe off the titration tube tip.
- Press the CLEAR button to set the display value to zero.
- 10 dispensing operations in 3 volume ranges (100 %, 50 %, 10 %) are recommended.
- Turn the hand wheels with both hands without stopping until the test volume is shown on the display. Wipe off the tip of titration tube.
- Weigh the dispensed quantity on an analytical balance. (Please refer to the user manual of the scale manufacturer.)
- Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = correction factor (e.g. 1.0029 ml/g at 20 °C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy*:

$$R\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Coefficient of variation*:

$$VK\% = \frac{100 s}{\bar{V}}$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

*) = Calculation of accuracy (A%) and coefficient of variation (CV%): A% and CV% are calculated according to the formulas of statistical quality control.

11 Additional functions

11.1 CAL mode (adjustment)


11.1.1 Adjust

After prolonged use or after replacing the dosing unit, an adjustment may be necessary to compensate for differences in accuracy up to a maximum of ± 0.999 ml. The change to the factory setting is shown on the display.

1. Calculate adjustment value


The adjustment value is the deviation of the average volume from the nominal volume (e.g: Average volume 50.024 ml, nominal volume 50 ml. Adjustment value = $50.024 \text{ ml} - 50.000 \text{ ml} = 0.024 \text{ ml}$). To calculate the average volume, see Checking the Volume (Calibration), p. 75.

2. Call up CAL mode


- a. With the appliance switched on, press and hold the CLEAR button for longer than 3 seconds until the display repeatedly shows the following modes in succession:

CAL - GLP - APO - dP

- b. As soon as CAL appears on the display, release the CLEAR button. CAL flashes and the number field is displayed.

3. Enter adjustment value

- a. The adjustment value is 0.024 ml, for example. Press the pause or on/off button until the value is reached.

4. Confirm entry

- a. Press the CLEAR button to confirm the entry of the adjustment value.
- ⇒ The change to the factory setting is made clear by the CAL symbol, which is now permanently shown on the display.

NOTICE

If the CLEAR button is not pressed within approx. 15 seconds, the initial status is retained.

11.1.2 Readjustment

The CAL symbol that is constantly shown on the display indicates that an adjustment has already been made. When the new adjustment value is entered, it is automatically added to the existing adjustment value.

1. Calculate adjustment value

The already adjusted device shows a renewed deviation of the mean volume from the nominal volume of e.g. 0.017 ml. To calculate the average volume, see Checking the Volume (Calibration), p. 75.

2. Call up CAL mode

- a.  > 3 s

- a. With the appliance switched on, press and hold the CLEAR button for longer than 3 seconds until the display repeatedly shows the following modes in succession:


CAL - GLP - APO - dP

- b.

- b. As soon as CAL appears on the display, release the CLEAR button. CAL flashes and the adjustment value of the previous adjustment appears.

3. Enter adjustment value

- a.** The adjustment value is 0.017 ml, for example. Press the pause or on/off button until the value is reached (the first time the button is pressed, the display is set to zero).

4. Confirm entry

- a.** Press the CLEAR button. The old and new adjustment values are automatically added together.
- ⇒ The change in adjustment is indicated by the CAL symbol.


NOTICE

In rare cases, the sum of the adjustment values may be zero when the new adjustment value is entered. In this case, the factory setting is restored and CAL disappears from the display.

11.1.3 Factory adjustment

The CAL symbol that is constantly shown on the display indicates that an adjustment has already been made. However, the factory setting should be restored.

1. Call up CAL mode

- a.** With the appliance switched on, press and hold the CLEAR button for longer than 3 seconds until the display repeatedly shows the following modes in succession:

CAL - GLP - APO - dP

- b.** As soon as CAL appears on the display, release the CLEAR button. CAL flashes and the adjustment value of the previous adjustment appears.

2. Restoring the factory setting


Press the On/Off button and Pause button simultaneously to restore the factory setting. The CAL symbol is hidden.

11.2 Calibration date (GLP mode)


The date for the next calibration can be saved in GLP mode (Good Laboratory Practice).

1. Entering GLP mode


- a. With the appliance switched on, press and hold the CLEAR button for longer than 3 seconds until the display repeatedly shows the following modes in succession:

CAL - GLP - APO - dP

- b. When GLP appears in the display, release the CLEAR button. The input symbol blinks and 'oFF' appears.

2. Entering the calibration date

Press and hold the Pause button until the desired date is displayed. Pressing briefly gradually extends the time until the calibration date. Pressing the On/Off button reduces the time until the calibration date. (Date input can be from 'oFF' to 12/2099)

3. Confirming input

Press the CLEAR button to confirm the input of the calibration date.


NOTICE

The saved calibration date can be called up any time the instrument is switched on. To do this, simply press and hold the On/Off button. A continuous display of GLP, the year and the month of the desired date appear. Releasing the button ends the continuous display and the instrument is switched on. (If 'oFF' is selected as the calibration date, this function is deactivated.)

11.3 Auto power off (APO mode)


In APO mode, the time for automatic power off can be set from 1 to 30 minutes. The factory default setting for automatic power off is 3 minutes. The shorter the auto power off time setting, the longer the batteries will last.

1. Entering APO mode


- a. With the appliance switched on, press and hold the CLEAR button for longer than 3 seconds until the display repeatedly shows the following modes in succession:

CAL - GLP - APO - dP

- b. When APO appears in the display, release the CLEAR button. The input symbol blinks and the factory setting is displayed.

2. Entering the time for automatic power off

- a. Press the Pause or On/Off button until the desired time entry (1 - 30 min) has been reached. 'OFF' deactivates automatic power off.

3. Confirming input

- a. Press the CLEAR button to confirm the input of desired power off time or to confirm 'OFF'.

NOTICE

When the instrument switches itself off, the last displayed value is displayed again after the instrument is switched back on. If the input value 'OFF' is confirmed, the function is deactivated and the instrument will no longer switch itself off.


11.4 Decimal places (dP mode)

In dP mode, the display can be set to show 2 or 3 decimal places (factory default setting is 2 decimal places).

NOTICE


For technical reasons, titration volumes over 20.00 ml can only be displayed with 2 decimal places.

1. Entering dP mode

- a. With the appliance switched on, press and hold the CLEAR button for longer than 3 seconds until the display repeatedly shows the following modes in succession:

CAL - GLP - APO - dP

- b. When dP appears in the display, release the CLEAR button. The input symbol blinks and the factory setting appears.

2. Changing decimal places

- a. Press the Pause button to select a display with 3 decimal places. (Pressing the button again converts the display back to 2 decimal places.)

3. Confirming input

- a. Press the CLEAR button to confirm the input of the desired decimal place format.

12 Cleaning

The instrument must be cleaned in the following situations to assure correct operation:

- Immediately when the hand wheels are more difficult to turn than usual
- prior to a reagent change
- prior to long term storage
- prior to dismantling the instrument
- regularly when using crystallizing solutions
- when liquids accumulate in the screw cap of the titrating tube


⚠ WARNING


The glass cylinder, valves, telescopic filling tube and titrating tube contain reagent! Follow the safety instructions (see Safety regulations, p. 60)!


12.1 Standard cleaning

- Set the valve to 'Recirculation' and empty the instrument completely by turning the hand wheels.
- Screw the instrument onto a bottle filled with deionized water and rinse the instrument several times by completely filling and emptying it.
- Set the valve to 'Titrate', unscrew the screw cap from the titration tube and place a suitable receiving vessel under the titrating tube. To clean the titrating tube, rinse the instrument several times by completely filling and emptying it.
- If deposits appear in the dispensing cylinder, repeat this process with a suitable cleaning agent and then rinse again with deionized water.


- e. Screw the instrument onto an empty bottle and empty the piston completely by executing several up and down piston movements in the 'Recirculation' and 'Titration' valve setting.
- f. First move the piston all the way up and then down by a half-rotation of the hand wheel.
- g. Unscrew the air vent cap by hand or use a coin.
- h. Remove the rear housing shell and take out the assembly tool.

- i. Loosen the safety ring of the dispensing unit with the assembly tool and unscrew it completely by hand.

- j. Pull the piston rod locking mechanism out to the stop.

- k. Move the top part of the instrument all the way to the top by turning the hand wheels and remove it.

- l. Remove any crystal deposits at the upper edge of the dispensing cylinder, e.g. with water and a soft bottle-brush. Then, dry it off with cellulose paper
- m. Reassemble the top part of the instrument or, if necessary, further dismantle the instrument for intensive cleaning.

NOTICE

Crystallizing media, e.g., alcoholic KOH

Depending on the frequency of use, we recommend that any crystal deposits above the piston be removed at regular intervals of approx. 8 weeks. To do this, carry out steps f-m of the standard cleaning procedure. To do this, carry out steps f-m of the standard cleaning procedure.

To reduce crystal deposits, always fill the instrument all the way up to the upper stop position after titration.

12.2 Intensive cleaning

To avoid mixing up components, do not dismantle several instruments at the same time. Calibration and, if applicable, an adjustment must always be carried out after dismantling or replacing a dispensing unit.


1. Preparing for intensive cleaning

- a. Always perform a complete standard cleaning before further dismantling the instrument.
- b. Pull out the recirculation tube and the telescopic filling tube.

2 Removing and cleaning/replacing the titration tube

- a. Set the valve to 'Recirculate' and pull the valve lever upwards to remove it (Fig. a).


- b. Hold the titrating tube as shown in the figure. To disconnect, push the discharge tube housing all the way up, then pull it forward with gentle up and down motions (Fig. b).
- c. Clean the titration tube with the integrated discharge valve (ultrasonic bath) or replace it.

3. Cleaning/replacing the dispensing unit

The dispensing unit consists of a piston and a dispensing cylinder with a valve block. If fluid is above the piston, the must be calibrated, Titrette® see Checking the Volume (Calibration), p. 75. If the instrument does not calibrate, we recommend replacing the piston and recalibrating the instrument. We, recommend replacing the complete dispensing unit.

- a. Hold the piston rod and slowly pull the piston out of the dispensing cylinder (Fig. a). Do not tilt the piston against the shaft of the piston rod.


NOTICE


If it is difficult to move the piston, open the top of the instrument, clamp the assembly tool (teeth pointing downward) between the dispensing cylinder and the top and turn the hand wheel to move the piston completely out of the dispensing cylinder (Fig. a').

- b. Use a soft cloth to Carefully clean the dispensing cylinder and the piston, or replace them. Be careful when doing this: The sealing lip of the piston is sensitive to mechanical damage

- c. To replace the piston, first slide the light gray safety ring of the piston rod **upwards** (Fig. c) and then unscrew the piston (Fig. c').
- d. Screw a new piston onto the piston rod and tighten it securely.
- e. Line up the gears of the piston and the piston rod by turning the piston back by no more than half a gear tooth.
- f. Slide the safety ring of the piston rod back **down**.

- g. Align the toothed strip (1) of the piston rod in the direction of the air vent opening (2) of the valve block. This is located opposite the titrating tube connection (3). Carefully (!) insert the piston vertically into the cleaned or new dispensing cylinder and push it about half-way in (Fig. g).


NOTICE

The sealing lip of the piston must not be damaged. Avoid contact with hard objects!

4. Installing the titrating tube

Installing the cleaned or new titrating tube:

- a. Slide the titrating tube on approx. 5 mm.
- b. Push the titrating tube housing all the way up.
- c. Then slide the titrating tube all the way in.
- d. Slide the titrating tube housing down to lock it into place.


- e. Pull up the valve lever to the 'Recirculate' position and press it in firmly.

5. Installing the top part of the housing


- a. Check that the piston rod locking mechanism has been pulled out.

- b. Attach the top part of the instrument and rotate the hand wheels to move it down while being careful that the recess of the front housing shell slides snugly over the titrating tube. Slightly rotate the top part, if necessary.




- c. Lift the safety ring of the dispensing unit and check that the tongue and groove interlock. Then, tightly screw in the safety ring by hand.


- d. Set the assembly tool on the right-hand edge of the housing and tighten it towards the left-hand edge. Then, re-insert the assembly tool into the rear housing shell.
- e. Slide the piston rod locking mechanism in all the way.
- f. First hang the rear housing shell from the top, then snap it closed and screw in the air vent cap.
- g. Carry out a function check and calibration, and make any necessary adjustments.

6. Cleaning/replacing the filling valve

- a. Remove the rear housing shell and take out the assembly tool.
- b. Pull out the telescopic filling tube and the recirculation tube.
- c. Use the assembly tool to unscrew the filling valve.
- d. If the sealing ring is contaminated or damaged, carefully remove it using a pair of angled tweezers.
- e. Clean the filling valve and sealing ring (ultrasonic bath) or replace them.
- f. Insert the cleaned or new sealing ring, if necessary.
- g. Screw in the filling valve first by hand and then tighten it with the assembly tool (1/4 turn is sufficient!).

12.3 Loosen stuck filling valve

The valve ball is firmly seated when:

- The appliance cannot be filled and if elastic resistance can be felt when the piston is turned up.
- No liquid can be dispensed even though the cylinder is filled.

In this case, loosen the valve ball, e.g. with a 200 μ l plastic pipette tip, by applying light pressure.

12.4 Replacing the batteries

If the battery symbol in the display is blinking, the battery capacity is depleted. The batteries must then be replaced.

Use only the specified battery type 1.5 V (AAA/UM4/LR03) micro batteries: The batteries are not rechargeable!

NOTICE



Always replace the battery compartment covers when replacing the batteries. These are included with the replacement batteries.

⚠ WARNING


Dispose of fully discharged batteries only in compliance with the Battery Directive. Do not discharge the battery by short-circuiting – explosion hazard!

- a. Unscrew the air vent cap by hand or use a coin.
- b. Remove the rear housing shell.

- c. Remove the battery compartment covers
- d. Remove the dead batteries using a screwdriver.
- e. Insert the new batteries and press them firmly into the sockets. Ensure that the polarity of the batteries corresponds to the symbols on the battery compartment.
- f. Close the battery compartments tightly again with the covers. Carefully press down the edges of the cover so that they fit tightly everywhere and there is no gap to the battery compartments.
- g. First hook in the rear housing shell at the top, then close it.
- h. Screw in the ventilation plug.

13 Troubleshooting

Error	Possible causes	Corrective action
Liquid is above the piston	Piston leaky	Carry out cleaning, carry out calibration, replace the dosing unit (see Cleaning, p. 83).
Piston sluggish	Dispensing cartridge contaminated or damaged by crystal deposits	Perform cleaning, replace dispensing cartridge see Cleaning, p. 83).
Filling not possible	Filling valve stuck	Clean filling valve, possibly dislodge valve ball with a 200 µl plastic pipette tip (see Loosen stuck filling valve, p. 90).
Filling not possible / liquid moves back when filling in the titration tube	Discharge valve contaminated or titration tube damaged	Clean discharge valve, or replace titration tube (see Cleaning, p. 83).
Air bubbles in the instrument	Instrument filled too quickly	Fill instrument slowly
	Filling tube is loose or damaged	Firmly attach telescopic filling tube, possibly cut the tube off approx. 1 cm from top, or replace tube.
	Filling valve loose or seal possibly not inserted	Check if the seal is inserted and tighten valve with assembly tool.
	Suction pipe does not immerse in liquid	Fill bottle, or correct the length of the telescopic filling tube.
	Recirculation tube not installed or installed incorrectly	Install the recirculation tube. The opening must face away from the bottle wall.
Titration not possible	Discharge valve stuck	Clean titration tube with integrated discharge valve, or replace titration tube (see Cleaning, p. 83). Any stuck valve ball with A 200 µl plastic pipette tip release (see Loosen stuck filling valve, p. 90)
Dispensed volume less than indicated	The instrument has not been fully primed	Prime instrument again (see Priming, p. 67).
	Seal possibly not inserted or filling valve loose	Check if the seal is inserted and tighten valve with assembly tool.

Error	Possible causes	Corrective action
	Filling valve is stuck or damaged	Clean filling valve. If necessary, replace it (see Intensive cleaning, p. 85).
The instrument does not show any function	Internal error	Restart the instrument: Remove battery, wait 1 minute and reinsert battery (see Replacing the batteries, p. 91).

14 Product markings

Symbol or number	Meaning
	General warning sign
	Note user manual
	Note eye protection
	Use hand protection
	Use protective clothing
XXZXXXXX	Serial number
	With this mark, we confirm that the product complies with the requirements set out in the EC Directives and has been subjected to the specified testing procedures.
	UKCA: United Kingdom Conformity Assessed With this mark we confirm that the product complies with the requirements specified in the UK Designated Standards.
 25	The device is marked in accordance with the German Weights and Measures Act and the Weights and Measures Ordinance. Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was added.
www.brand.de/ip	Patent information
 (here: 40 years)	China RoHS (EFUP) EFUP defines the period in years during which the hazardous substances contained in electrical and electronic equipment do not leak or change under normal operating conditions. Under normal use by the user, such electrical and electronic products do not cause severe environmental pollution, serious personal injury or damage to the user's property.
	The electrical device must not be disposed of with the household waste.

15 Ordering Information

Titrette®

	Standard	with RS 232 port
Volume	Part No.	Part No.
10 ml	4760141	4760241
25 ml	4760151	4760251
50 ml	4760161	4760261

16 Accessories/spare parts


Bottle adapter

PP, Packing unit 1 piece

External thread	for bottle thread/ground joint	Cat. No.
GL 45	GL 32- 33	704396
GL 45	GL 38	704397
GL 45	S* 40	704343
GL 45	S* 42	704349
GL 45	S* 50	704350
GL 32	NS 24/29	704424
GL 32	NS 29/32	704429

Titration tubes

With screw cap and integrated discharge and recirculation valve.

Packing unit 1 piece

For volumes	Cat. No.
10 ml	707525
25 + 50 ml	707529

Screw cap

Description	Pkg unit	Cat. No.
Screw cap with fastener.	1 piece	707528


Bottle stand

Description	Pkg unit	Cat. No.
Bottle stand. PP. Support rod 325 mm. Base plate 220 x 160 mm	1 piece	704275

Drying tube

Description	Pkg unit	Cat. No.
Drying tube Incl. PTFE seal- ing ring (with- out drying agent).	1 piece	707930

Filling valve

Description	Pkg unit	Cat. No.
Filling valve with olive- shaped nozzle and sealing ring.	1 piece	6636


Telescopic filling tube

Description	Pkg unit	Length	Cat. No.
Telescopic filling tube. FEP. Indi- vidually ad- justable lengths.	1 piece	170–330 mm	708218
	1 piece	250–480 mm	708220

Recirculation tube

Description	Pkg unit	Cat. No.
Recirculation tube	1 piece	8317

Piston

For volumes	Pkg unit	Cat. No.
10 ml	1 piece	707531
25 ml	1 piece	707530
50 ml	1 piece	707532

Dispensing cylinder with valve block

For volumes	Pkg unit	Cat. No.
10 ml	1 piece	707533
25 ml	1 piece	707535
50 ml	1 piece	707537

Viewing window

Description	Pkg unit	Cat. No.
Inspection windows, 1 set colorless and 1 set amber-colored (light shield).	1 piece	6783

Micro-batteries, 1.5 V

Description	Pkg unit	Cat. No.
Micro-batteries, 1.5 V. Not rechargeable (AAA/UM4/LR03). Incl. battery compartment covers	2 piece	7260

Venting stopper

Description	Pkg unit	Cat. No.
Venting stopper	1 piece	6659

Assembly tool

Description	Pkg unit	Cat. No.
Assembly tool	1 piece	6784

Battery compartment covers

Description	Pkg unit	Cat. No.
Battery compartment covers	2 piece	8857

RS 232 interface cable

Description	Pkg unit	Cat. No.
RS 232 interface cable. Length 2 m.	1 piece	8850

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Description	Pkg unit	Cat. No.
Titrette Software CD-ROM. German/English	1 piece	707538

Titrette® removal system for bag-in-box containers

Description	Pkg unit	Cat. No.
Titrette® extraction system for Bag-in-Box container (basic equipment)	1 piece	707550

Adapter set interface RS232 to USB for Titrette

Description	Pkg unit	Cat. No.
Adapter set interface RS232 to USB for Titrette®	1 piece	707539

17 Repairs

17.1 Sending for repair

NOTICE Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn

info@brand.co.in
www.brand.co.in

www.brand.cn.com

18 Calibration service

ISO 9001 and GLP guidelines require regular testing of your volumetric measuring devices. We recommend that you perform a volume check every 3-12 months. The cycle depends on the individual requirements of the device. In case of high frequency of use or aggressive liquids, more frequent testing should be carried out.

The complete SOP for testing can be downloaded from www.brand.de or www.brandtech.com.

BRAND also offers you the option of having your devices calibrated through our factory calibration service or through our accredited calibration laboratory. Just send us the devices to be calibrated, indicating the type of calibration you would like. You will get your devices back in a few days. A detailed calibration report (factory calibration) or an accredited calibration certificate in accordance with DIN EN ISO/IEC 17025 is enclosed with each device. More information can be obtained from your retailer or directly from BRAND. The order document is available for download at www.brand.de (Service & Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. Our service partners can forward your devices to BRAND for factory calibration, if required.

19 Information about your laboratory instrument

The online service MyProduct (<https://www.brand.de/myproduct>) offers quality certificates, accessories, and technical documentation for your laboratory instrument Titrette®. By entering the serial and article number there, you will receive information specific to your individual instrument.

On Titrette®, you will also find serialized information encoded in the quick response code. This also links to the MyProduct website, where you can access this operating manual and additional certificates for your instrument.

You will also find a data matrix code on some instruments (Transferpette® S, HandyStep® touch, and HandyStep touch® S). Scan this using a standard reader app to access the above-mentioned information via the URL <https://www.brand.de/myproduct>.

20 Warranty

We shall not be liable for the consequences of improper handling, use, maintenance, operation or unauthorized repair of the device or for the consequences of normal wear and tear, in particular of wear parts such as pistons, seals, valves or glass breakage. The same applies to non-observance of the instructions for use. In particular, we assume no liability for damage caused if the device has been disassembled further than described in the user manual or if external accessories or spare parts have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

21 Disposal

This symbol means that at the end of their service life, batteries/accumulators and electronic devices must be disposed of separately from household waste (unsorted municipal waste).

Electronic devices must be disposed of in accordance with Directive 2012/19/ EU of the European Parliament and of the Council from July 04, 2012 on waste from electrical and electronic equipment and in compliance with national disposal regulations.

Both batteries and accumulators (rechargeable batteries) contain materials that can be damaging to the environment and human health. Therefore, they must be properly disposed of in accordance with Directive 2006/66/EC of the European Parliament and of the Council from September 06, 2006 on batteries and accumulators and in compliance with national disposal regulations. Only dispose of fully discharged batteries and accumulators.

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