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
8. SPECIFICATIONS..... 32

9. ACCESSORIES..... 35

SICCO Glove Box with Transfer Chamber: Cat.No.: V 1982-08

SICCO Glove Box: Cat.No.: V 1984-08

SICCO Glove box for Exhaustion: Cat.No.: V 1988-08

SICCO Glove Box Vitrum with Transfer Chamber: Cat.No.: V 1986-08

SICCO Glove Box Antistatic with Transfer Chamber: Cat.No.: V 1992-08


SICCO Glove Box Duo 2 with Transfer Chamber: Cat.No.: V 1710-08

SICCO Glove Box Duo 2 with Transfer Chamber Antistatic Cat.No.: V 1715-08

SICCO Glove Box Duo 1 with Transfer Chambers Cat.No.: V 1730-08

SICCO Glove Box Duo 1 with Transfer Chambers Antistatic: Cat.No.: V 1735-08

SICCO Glove Box Trio with Transfer Chamber Cat.No.: V 1750-08

Please read the information listed below thoroughly. We urgently ask you to observe the instructions about security, use and maintenance.

SICCO Glove Boxes provide safe handling of toxic materials, allow production of chemical compounds and protect sensitive substances from the environment. Optionally, the chamber can be flooded with an inert gas and a special atmosphere can be created. Models with side transfer chamber/s allow to introduce products without noticeable change of the atmosphere inside. All panels are made of high-quality and transparent plastic material, except the model VITRUM with panels made of borosilicate glass 3.3. The aluminium profiles are joint together with connectors made of glass fiber-reinforced polyamide.

1. Supplied with

1 Pair of gloves respectively 1 pair of iris-ports	connecting kit for gas- filling	connecting kits for gas-filling with one transfer chamber	connecting kits for gas-filling with two transfer chambers	Operating Instructions
--	---------------------------------------	--	---	---------------------------


Glove Box with Transfer Chamber	1	0	2	0	1
Glove Box	1	1	0	0	1
Glove box for Exhaustion	1	0	0	0	1
Glove Box Vitrum with Transfer Chamber	1	0	2	0	1
Glove Box Antistatic with Transfer Chamber	1	0	2	0	1
Glove Box Duo 2 with Transfer Chamber	2	0	2	0	1
Glove Box Duo 2 with Transfer Chamber Antistatic	2	0	2	0	1
Glove Box Duo 1 with Transfer Chambers	2	3	2	3	1
Glove Box Duo 1 with Transfer Chambers Antistatic	2	3	2	3	1
Glove Box Trio with Transfer Chamber	2	1	1	1	1

2. Security Advice

- » Do not take a damaged glove box into operation.
- » Check before every use whether the glove box has been damaged.
- » Provide a safe stand for the glove box.
- » Do not work with opened swing-up door.
- » SICCO Glove Box shall be taken into operation only in rooms with sufficient ventilation.
- » Place the glove box in a reachable working height.
- » Do not expose the glove box to direct sunlight.
- » Do not store any solvents in the glove box (except for model Vitrum).
- » Do not use any sharp or pointed tools inside the glove box.
- » Do not put any explosive or inflammable materials inside the glove box.
- » Do not put any radioactive products inside the glove box.
- » Transfer the extracted air in a controlled way, e.g. by using an extractor hood.
- » Make sure that the airflow of the extractor hood does not effect or ruin the process of gas filling.
- » SICCO Glove Boxes are not suitable for overpressure or vacuum.
- » The ambient temperature range of – 20 °C to + 70 °C, respectively of -70 °C to + 150 °C for model Vitrum, shall not be exceeded.
- » Make sure to install all connecting kits for gas-filling.
- » Always use the smaller tubing as your tubing for gas-filling and the larger tubing as your aeration tubing.
- » The pressure for lead-in gas shall not exceed 0,3 bar.
- » Connect the aeration tubing at first, then the tubing for gas-filling. Disconnection has to be made in reverse order.
- » The cross section of the tubing must not be constricted.
- » Neither put the tubing for gas-filling nor the aeration tubing under tensile or compressive stress
- » Avoid unanticipated reactions which might be caused by gas residues upon opening the swing-up door or the transfer chamber.
- » Connect the ground cable (only for model "Antistatic").
- » For model for exhaustion please make sure that there is a sufficient supply air during use. Do not close any of the iris-ports, e.g. by using gloves that would seal the iris-ports.
- » These operating instructions cannot reflect all possible dangers which might occur during use. Users should check and observe all safety regulations.

3. Taking into operation

Remove the packing and the protection cover and place the glove box on a horizontal, plain bench. **Attention: Never lift the glove box at the transfer chamber. Only lift the glove box at the bottom (see picture on the right):**


Swing-up door

» Always pull both handles to open the swing-up door. Swing the door panel up to the highest point until it stops.

Please note: Due to the strong magnetic lock it is harder to open the door the longer it was closed before. This is normal and can be compensated by pulling stronger than in everyday use.

» **Attention:** The door stop defines the maximum opening of the swing-up door. Put the door on the door stop. The adjoining picture shows the maximum possible opening of the door.


Swapping with Gas

The substances and components which are stored in the glove box can be swapped with gas. The corresponding bores to lead in the gas are located on the side walls and are factory-made. For technical reasons (packing and transport), the tubing and the connecting pieces are enclosed in delivery but not yet mounted. The assembly is quite easy.

Please note: gas which is lighter than air has to be lead in from above. Gas which is heavier than air from the bottom (density of air = 1,2041 kg/m³). This information is important for further use since the tubing for gas-filling has a different cross section then the aeration tubing.

» Always use the smaller tubing (3) for gas-filling. Use the bigger tubing (4) for aeration only. This way you can prevent the accrument of overpressure and therefore any damage to your glove box.

» Open the four cable lead-ins (two cable lead-ins on the back side of the glove box and two on the transfer chamber) by unscrewing the sleeve nuts (1). Remove the sealing disks (2). Please keep the sleeve nuts and the sealing disks safe for re-use.


SICCO **Glove Box**


- » Cut all tubing with a knife to the requested length. 5 meters of each tubing are included in delivery (for models with transfer chamber 10 m, respectively 15 m, are included). The picture on the right shows the tubing for gas-filling.

- » Connect the tubing for gas-filling by sliding them on the hose connectors. Assemble the aeration tubing accordingly. The picture on the right shows the assembly of the tubing for gas-filling.
- » **Tip:** Wet the hose connectors with some water or silicone oil.

- » If the gas used for swapping is heavier than air the smaller tubing for gas-filling (3) has to be assembled on the lower cable lead-ins. The larger aeration tubing (4) has to be assembled on the upper cable lead-ins.

- » For assembling the tubing please use the included hexagon wrench. Fix the four fastening screws by turning them carefully clockwise. The picture on the right shows the aeration tubing.

- » If the gas used for swapping is lighter than air the tubing for gas-filling (3) has to be assembled on the upper cable lead-ins. Corresponding to this the aeration tubing (4) has to be assembled on the lower cable lead-ins.

Attention:

- » Transfer the extracted air in a controlled way, e.g. by using an extractor hood. Please make sure that the airflow of the extractor hood does not effect or ruin the process of gas-filling in our glove box.
- » The cross section of the tubing must not be constricted. Make sure to choose a sufficient bending radius for the tubing when installing them. Do not put any components on the tubing which could damage or crush them.
- » We recommend to purge the transfer chamber and the work space with 0,15 bar or with 0,30 bar. The following table shows you the scavenging time:

Scavenging pressure inside the inlet hose	Flow rate approx.	Scavenging time inside work space minimum	Scavenging time inside transfer chamber minimum
0,15 bar	3100 litre per h	6:00 minutes	1:00 minute
0,30 Bar	6200 litre per h	3:00 minutes	0:30 minute

Cable lead-ins

The cable lead-ins on the back side of the chamber allow to lead in any standard cables including mounted plugs or any rigid tubing. Every cable or rigid tubing with an outer diameter from 0,1 mm up to 10 mm can be assembled regardless of the cross section the cable / tubing has.


Attention:

The following steps shall not be applied to SICCO Glove box for exhaustion. The SICCO Glove box for exhaustion has iris-port cable lead-ins. Those iris-ports must never be closed during use to guarantee a sufficient supply air.

- » Unscrew the four cylinder screws. A suitable hexagon wrench is included.


- » Remove the supporting ring (5) and the supporting disk with centring collar (6).

SICCO **Glove Box**

- » Remove the silicone plug (7) and the sealing foil (8).

- » Remove the clamp ring (9).


- » Put your electrical device inside the glove box. Lead your cable through the supporting ring (5). Please pay attention to the supporting ring (5) being properly assembled with the notches in direction to yourself in order to fix the cylinder screws later.

- » Guide the plug from inside the glove box through the housing (10) of the cable lead-in.

- » Pull the cable of your electrical device until the required position is reached. Slide the clamp ring (9) on the cable.


- » Afterwards slide the silicone plug (7) on the cable as well.

SICCO **Glove Box**


- » Then slide the supporting disk with centring collar (6) on the cable with the centring collar in direction to yourself.

- » Press the silicone plug (7) into the housing (10). While doing so put pressure to the housing (10) from outside as well.

- » Shove the supporting disk with centring collar (6) along the cable. Place the clamp ring with ridge on top of the silicone plug (7). Place the supporting ring (5) on the supporting disk with centring collar (6).

- » Fix the four fastening screws by turning them carefully clockwise.

- » The assembly of a rigid tubing is identical to the assembly of a cable. Just follow the above mentioned steps.


Transfer Chamber (only for models with transfer chamber).

The transfer chamber is used to introduce further products or tools during a running process. A contamination of the working space or a loss of the filled-in gas can be prevented.

SICCO Glove Box


- » Close both doors of the transfer chamber.
- » Purge the transfer chamber with an inert gas.
- » Provide sufficient ventilation to convey the exhaust gas, e. g. with a fume hood.
- » It's the user's liability to define a sufficient period of time for gasfilling. Please pay attention to the table on page 24.

- » If the inner door of the transfer chamber is open, the interlock (11) prevents accidental opening of the outer door.
Attention: The interlock prevents accidental opening but cannot provide absolute protection!

- » Put the products into the transfer chamber. A tray facilitates the removal of the items.
- » Close the outer door.


- » It's the user's liability to define a sufficient period of time for purging the transfer chamber again with gas. Please pay attention to the table on page 24.

- » Open the inner door of the transfer chamber and remove the products or tools.

- » Close the inner door of the transfer chamber.
- » Purge the transfer chamber with inert gas.
- » Provide sufficient ventilation to convey the exhaust gas, e. g. with a fume hood.
- » It's the user's liability to define a sufficient period of time for gas-fillin Please pay attention to the table on page 24. If neces sary, the outlet tube can be connected to a metering station.


Glove Box for Exhaustion

The SICCO Glove Box for exhaustion can be connected directly to an existing exhaust system through the integrated connector. For technical reasons (packing and transport) the connector is enclosed in delivery but not yet mounted. The assembly is quite easy.

- » Open the swing-up-door and place one hand underneath the cover lid (12) for the connector.

- » Turn the screw collar (13) counter clockwise. Hold the cover lid (12) from below.
- » **Attention:** The cover lid (12) will fall into your hand.
- » Keep the collar screw (13) for further use.

- » Remove the red silicone sealing (14) from the cover lid (12).

- » If you want to use the whole cross section of your suction hose you can saw off the connector. Use a fine-toothed saw for this step. Make sure to fix the connector with a screw clamp to a solid base.

- » Remove the remaining edge with a sharp knife or use emery paper.
Attention: Risk of injury!


SICCO **Glove Box**

- » Set the red silicone sealing (14) on the connector (15).

- » Guide the connector (15) with the silicone sealing (14) through the opening from below.

- » Fix the connector (15) by screwing the collar screw (13) clockwise.
- » Do not use any tools for fixing the connector.


- » The connector is graded so that tubes with inner diameters of 100, 120, 160 and 180 mm can be assembled.

- » The iris-ports can be replaced with new ones (available as accessory, cat.-no. V 1976-07). You can remove them by slightly rotating them and simultaneously pull them off the inlet port. Assemble the new iris-ports by slightly rotating them and simultaneously slide them on the inlet port.

- » Both inlet ports can be closed with a lid (available as accessory, cat.-no. V 1790-07) to protect the interior from small particles or dust. The lids do not close gas tight.
- » **Attention:** Do not close the inlet ports during exhaustion!

The following functions are only available for versions Antistatic:

- » Connect a ground cable on the backside of the glove box and on the backside of the transfer chamber.

4. Functionality

A magnetic strip keeps the swing-up door as well as the doors of the optional transfer chamber closed tightly. The connecting kit for gas-filling allows to swap the glove box with inert gas. The natural rubber gloves, respectively EPDM gloves for versions Antistatic, allow safe working inside the glove box.

5. Service / Maintenance

- » Use only soft clothes for cleaning. Do not use any detergents which may contain solvents or which scrub.
- » Keep the magnetic door strip clean.
- » The gloves can be exchanged as follows:

- » Remove the glove from the holding ring.

- » Slide the new glove onto the holding ring.

The gloves are wearing parts and have to be exchanged regularly. Due to the diversity of possible applications, a period cannot be defined. It is each user's obligation to check regularly the state of the gloves and to exchange them if necessary.

6. Disturbances / Repair

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned glove box in its original packing to us. Then we will check whether a cost-effective repair is possible.

Objection

Reason and help

The door does not close properly.

» There are particles like dust which are on the magnetic strip and / or on the door frame. Clean the magnetic strip and the frame with a soft cloth.


The gas consumption is too high.

» Assemble the gas-saving valve V 1787-07 to the glove box and connect the outlet hose. The valve opens only to achieve a pressure compensation. This way the gas consumption can be reduced.

SICCO Gas-saving valve

- » During gas-filling the sealing ball is lifted due to the gas flow and the pass is opened.
- » When the gas-filling process is completed the sealing ball drops due to its own weight and closes the pass. An unwanted backflow of the gas or of the ambient air into the work chamber is no longer possible. Such backflow can occur e.g. through movements inside the glove box or when you slide in or off the gloves.

Attention: Since the sealing ball only drops because of its own weight the SICCO Gas-saving valve does not provide a complete protection against the loss of gas. However if moved carefully inside the glove box the gas consumption can be reduced to a minimum respectively can completely be interrupted. It's the users liability to provide a proper gas supply.

7. Waste Disposal

The glove box as well as its packing are made of recyclable materials. Please take further information on the materials used from the following specifications:

8. Specifications

Details of Transfer Chamber (for all models with Transfer Chamber)

Outer dimensions (W x H x D):	320 x 330 x 260 mm
Usable interior space (W x H x D):	220 x 275 x 320 mm
Usable volume:	ca. 24 litres
Max. usable height:	275 mm
Max. usable width:	200 mm
Distance to bottom:	40 mm
Max. all-over load:	3 kg

SICCO Glove Box

Cat.-No.:	V 1984-08
Outer dimensions (W x H x D)	890 x 700 x 600 mm
Weight:	ca. 14 kg
Usable interior space (W x H x D):	850 x 650 x 540 mm
Usable volume:	ca. 290 litres
Max. all-over load:	40 kg
Opening for gloves (O.D.):	Ø 190 mm
Distance between opening for gloves:	390 mm
Hinged lid opening (W x H):	850 x 425 mm
Opening for gas supply:	Ø 6 mm
Opening for gas outlet:	Ø 41 mm
Length of hose:	5 m tubing for gas-filling and 5 m aeration tubing
Material frame:	Aluminium with connectors made of polyamide
Material panels:	PMMA (polymethylmethacrylate)
Material base panel:	PE (polyethylene)

SICCO Glove Box with Transfer Chamber

Cat.-No.:	V 1982-08
Outer dimensions (W x H x D)	1210 x 700 x 600 mm
Weight:	ca. 20 kg
Usable interior space (W x H x D):	850 x 650 x 540
Usable volume:	ca. 290 litres
Max. all-over load:	40
Opening for gloves (O.D.):	Ø 190 mm
Distance between opening for gloves:	390 mm
Hinged lid opening (W x H):	850 x 425 mm
Opening for gas supply:	Ø 6 mm
Opening for gas outlet:	Ø 41 mm
Length of hose:	10 m tubing for gas-filling and 10 m aeration tubing
Material frame:	Aluminium with connectors made of polyamide
Material panels:	PMMA (polymethylmethacrylate)
Material base panel:	PE (polyethylene)

SICCO Glove Box with Transfer Chamber Antistatic

Cat.-No.:
 Material panels:
 Material base panel:
 Further specifications:

V 1992-08
 PC – ESD (static dissipative polycarbonate)
 PE – ESD (static dissipative polyethylene)
 see Glove Box with Transfer Chamber V 1982-08

SICCO Glove Box for Exhaustion

Cat.-No.:
 Outer dimensions (W x H x D)
 Weight:
 Usable interior space (W x H x D):
 Usable volume:
 Max. all-over load:
 Opening for Iris-Ports (O.D.):
 Distance between opening for gloves:
 Hinged lid opening (W x H):
 Connector for exhaust system (graded):
 I.D. of connector:
 Material frame:
 Material panels:
 Material base panel:

V 1988-08
 890 x 700 x 600 mm
 ca. 16 kg
 850 x 650 x 540
 ca. 290 litres
 40
 Ø 170 mm
 390 mm
 850 x 425 mm
 Ø 100, Ø 120, Ø 160, Ø 180
 Ø 90 mm
 Aluminium with connectors made of polyamide
 PMMA (Polymethylmethacrylate)
 PE (polyethylene)

SICCO Glove Box Vitrum with Transfer Chamber

Cat.-No.:
 Outer dimensions (W x H x D)
 Weight:
 Usable interior space (W x H x D):
 Usable volume:
 Max. all-over load:
 Opening for gloves (O.D.):
 Distance between opening for gloves:
 Hinged lid opening (W x H):
 Opening for gas supply:
 Opening for gas outlet:
 Length of hose:
 Material frame:
 Material panels:

 Material base panel:

V 1986-08
 1210 x 700 x 600 mm
 ca. 36 kg
 850 x 650 x 540
 ca. 290 litres
 40
 Ø 190 mm
 390 mm
 850 x 425 mm
 Ø 6 mm
 Ø 41 mm
 10 m tubing for gas-filling and 10 m aeration tubing
 Aluminium with connectors made of polyamide
 Borosilicate glass 3.3, front panel and right side panel made of aluminium
 Bottom side of base panel made of aluminium, top side of base panel made of stainless steel

SICCO Glove Box Trio with Transfer Chamber

Cat.-No.:

Outer dimensions (W x H x D)

Weight:

Usable interior space (W x H x D):

Usable volume:

Max. all-over load:

Opening for gloves (O.D.):

Distance between opening for gloves:

Hinged lid opening (W x H):

Opening for gas supply:

Opening for gas outlet:

Length of hose:

Material frame:

Material panels:

Material base panel:

V 1750-08

1590 x 690 x 590 mm

ca. 28 kg

1240 x 650 x 550 mm

ca. 440 litres

60 kg

Ø 190 mm

390 mm

1240 x 425 mm

Ø 6 mm

Ø 41 mm

10 m tubing for gas-filling and 10 m aeration tubing

Aluminium with connectors made of polyamide

PMMA (Polymethylmethacrylate)

PE (polyethylene)

SICCO Glove Box Duo 2 with Transfer Chamber Antistatic

Cat.-No.:

Weight:

Material panels:

Material base panel:

Further specifications:

V 1715-08

ca. 31,5 kg

PC – ESD (static dissipative polycarbonate)

PE – ESD (static dissipative polyethylene)

see Glove Box Duo 2 with Transfer Chamber V 1710-08

SICCO Glove Box Duo 1 with Transfer Chambers

Cat.-No.:

Outer dimensions (W x H x D):

Weight:

Usable interior space (W x H x D):

Usable volume:

Max. all-over load:

Opening for gloves (O.D.):

Distance between opening for gloves:

Hinged lid opening (W x H):

Opening for gas supply:

Opening for gas outlet:

Length of hose:

Material frame:

Material panels:

Material base panel:

V 1730-08

2345 x 690 x 590 mm

ca. 40 kg

1680 x 650 x 550 mm

ca. 600 litres

80 kg

Ø 190 mm

390 mm

1680 x 425 mm

Ø 6 mm

Ø 41 mm

15 m tubing for gas-filling and 15 m aeration tubing

Aluminium with connectors made of polyamide

PMMA (Polymethylmethacrylate)

PE (polyethylene)

SICCO Glove Box Duo 1 with Transfer Chambers Antistatic

Cat.-No.:

Weight:

Material panels:

Material base panel:

Further specifications:

V 1735-08

ca. 43,5 kg

PC – ESD (static dissipative polycarbonate)

PE – ESD (static dissipative polyethylene)

see Glove Box Duo 1 with Transfer Chambers V 1730-08

SICCO Glove Box Duo 2 with Transfer Chamber

Cat.-No.:

Outer dimensions (W x H x D):

Weight:

Usable interior space (W x H x D):

Usable volume:

Max. all-over load:

Opening for gloves (O.D.):

Distance between opening for gloves:

Hinged lid opening (W x H):

Opening for gas supply:

Opening for gas outlet:

Length of hose:

Material frame:

Material panels:

Material base panel:

V 1710-08

1200 x 710 x 890 mm

ca. 29,5 kg

850 x 665 x 845 mm

ca. 475 litres

60 kg

190 mm

390 mm

850 x 425 mm

Ø 6 mm

Ø 41 mm

10 m tubing for gas-filling and 10 m aeration tubing

Aluminium with connectors made of polyamide

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BOHLENDER GmbH