

User Manual  
English

## Portavo® 904(X) PH





## Warranty

Defects occurring within 3 years from delivery date shall be remedied free of charge at our plant (carriage and insurance paid by sender).

Sensors and accessories: 1 year

Subject to change

## Return of products under warranty

Please contact our Service Team before returning a defective device.

Ship the cleaned device to the address you have been given.

If the device has been in contact with process fluids, it must be decontaminated/disinfected before shipment. In that case, please attach a corresponding certificate, for the health and safety of our service personnel.



## Disposal

Please observe the applicable local or national regulations concerning the disposal of “waste electrical and electronic equipment”.

## Registered trademarks

The following names are registered trademarks. For practical reasons they are shown without trademark symbol in this manual.

- Calimatic®
- Memosens®
- Paraly®
- Portavo®
- Sensocheck®
- Sensoface®

|  |           |
|--|-----------|
| <b>Package Contents .....</b>                          | <b>6</b>  |
| <b>Documentation .....</b>                             | <b>7</b>  |
| <b>Overview of the Portavo 904(X) PH .....</b>         | <b>8</b>  |
| Value-Added Features .....                             | 9         |
| Protective Cover.....                                  | 10        |
| Hook.....  | 10        |
| Display.....   | 11        |
| Keypad .....   | 12        |
| <b>Start-Up .....</b>                                  | <b>13</b> |
| Inserting the Batteries .....                          | 13        |
| Batteries  |           |
| for Application in Hazardous Locations .....           | 14        |
| Connecting a Sensor .....                              | 15        |
| Switching On the Meter .....                           | 16        |
| Icons .....  | 16        |
| <b>Configuring .....</b>                               | <b>17</b> |
| <b>Calibrating .....</b>                               | <b>18</b> |
| <b>Measuring .....</b>                                 | <b>22</b> |
| Switching the Measured Value Display .....             | 22        |
| Adjusting the Temperature.....                         | 22        |
| <b>Data Logger .....</b>                               | <b>23</b> |
| Operating Modes of the Data Logger (Logger Type) ..... | 24        |
| Data Logger Menu .....                                 | 26        |
| Configuring the Data Logger .....                      | 27        |
| Starting the Data Logger using CONT .....              | 28        |
| Starting the Data Logger using START .....             | 28        |
| Displaying the Logger Data.....                        | 29        |
| Stopping the Data Logger .....                         | 30        |
| Clearing the Data Logger .....                         | 30        |
| <b>Clock .....</b>                                     | <b>31</b> |

|   |           |
|---|-----------|
| <b>Paraly SW 112 Software .....</b>         | <b>32</b> |
| <b>Error Codes and Device Messages.....</b> | <b>33</b> |
| “Sensoface” Messages.....                   | 34        |
| Error Messages .....                        | 35        |
| <b>Product Line .....</b>                   | <b>36</b> |
| Accessories .....                           | 36        |
| Sensors .....                               | 37        |
| Knick CaliMat Buffer Solutions .....        | 38        |
| <b>Specifications.....</b>                  | <b>39</b> |
| <b>Index .....</b>                          | <b>43</b> |

Check the shipment for transport damage and completeness.

The package of the Portavo 904(X) PH includes:

- The Portavo 904(X) PH incl. 4 AA batteries and premounted quiver
- Carrying strap
- Quickstart instructions in various languages
- Specific test report
- Safety instructions
- Certificates
- Data carrier with detailed user manuals and Paraly SW 112 software
- USB cable, 1.5 m

## Specific Test Report



## CD-ROM

Complete documentation:

- User manuals in different languages
- Safety Information
- Certificates
- Quickstart guides

## Safety Instructions

In official EU languages and others.

- EC Declarations of Conformity



## Certificates

- IECEx
- ATEX



## Quickstart Guides

Installation and first steps:

- Operation
- Menu structure
- Calibration
- Error messages and recommended actions



Various languages on CD-ROM and on our website:  
[www.knick.de](http://www.knick.de)



**The Portavo 904(X) PH** is a portable pH meter. A plain-text line on the high-contrast LCD screen makes operation virtually self-explanatory. The device variant 904 X PH is available for applications in hazardous locations up to Zone 0.

The meter stands out by the following features:

- Use of digital Memosens sensors
- Memosens sensors and DIN pH sensors can be used on one device.
- A detachable quiver protects the sensor and prevents it from drying out. Furthermore, it can be used for calibration.

- The rugged housing is made of a high-performance polymer. It provides high impact resistance and dimensional stability even when exposed to extreme moisture.
- Scratch-proof clear glass display, perfectly readable even after years
- Very long operating times with one set of batteries (4 x AA) or use of a Li-ion battery for reliable operation even at high or very low operating temperatures (Li-ion battery not suited for Portavo 904 X PH for application in a hazardous location)
- Data logger with 5000 values
- Micro USB port for communication with Paraly SW 112 software for data evaluation of digital sensors (Memosens)
- Sensoface icons provide single-glance information on the sensor condition (page 34)
- Calibration with “Calimatic” automatic buffer recognition (page 18)
- Manual calibration by entering individual buffer values
- Real-time clock and indication of battery charging level
- At measuring temperatures from -20 to +100 °C the temperature detector can be automatically identified.



## Value-Added Features

### Memosens

The Portavo 904 can communicate with Memosens sensors. When these digital sensors are connected to the meter, they are automatically identified and indicated by the logo shown on the right. Furthermore, Memosens allows the storage of calibration data, which will be available and can still be used when the sensor is connected to another Memosens-capable device.



### Sensoface

Sensoface provides quick information on the sensor condition. The three “smiley” faces as shown on the right represent the sensor condition during measurement and after a calibration. When the condition deteriorates, an “INFO ...” message gives a hint to the cause.



### Automatic calibration with Calimatic

Calimatic is a very convenient method for pH calibration with automatic buffer recognition. You only have to select the buffer set with the buffers used. The buffers can then be used in any order.

As delivered, this calibration method is preset. It can be adjusted or disabled in the configuration menu.



### Protective Cover

The front of the meter is protected by a cover, which can be completely flipped over and secured to the back for operation. A label on the inner side of the cover explains the control functions and device messages.



### Hook

A fold-out hook on the back allows suspending the meter. This leaves your hands free for the actual measurement. The **rating plate** is located beneath the hook.




### Protective Cover and Hook Combined

Cover and hook can be joined together to form a benchtop stand allowing comfortable and fatigue-free working at a lab bench or desk.


## Display

The meter has a three-line display for representing alphanumeric information such as measurement and calibration data, temperatures and date/time. Additional information is provided by means of icons (Sensoface, battery icon, etc.).


Some typical displays are shown here.




Measuring  
(display of measured value,  
temperature, date and time)



Calibration – step 1



End of calibration  
(display of slope)



Logger data  
(display of measured value, memory  
location, temperature, date and time)



Clock  
(display of hours and minutes,  
seconds and date).



## Keypad

The keys of the membrane keypad have a noticeable pressure point.

They have the following functions:

|               |   |
|---------------|---|
| <b>on/off</b> | Switches the meter on and displays the device and calibration data (see Start-Up)                         |
| <b>meas</b>   | Switches the meter on /<br>Activates measuring mode /<br>Stops the data logger                            |
| <b>cal</b>    | Starts calibration  |
| <b>set</b>    | Activates configuration /<br>Confirms entries   |
| <b>clock</b>  | Displays time and date, allows setting the clock using <b>set</b>   |
| <b>RCL</b>    | View stored values  |
| <b>STO</b>    | Holds and saves a measured value, allows setting and starting the logger by pressing <b>set</b> (page 23) |
| <b>▲▼</b>     | When this icon is displayed, you can use the arrow keys for navigation.                                   |

Check the shipment for transport damage and completeness (see Package Contents).



### Caution!

Do not operate the device when one of the following conditions applies:

- the device shows visible damage
- the device fails to perform the intended function
- prolonged storage at temperatures above 70 °C
- severe transport stresses

In this case, a professional routine test must be performed. This test should be carried out by the manufacturer.

### Precautions for application in hazardous locations



### Warning!

- Only open the battery compartment of the Portavo 90n X outside the hazardous location.
- Never try to open the device. If a repair should be required, return the device to our factory.
- Never use the USB port within the hazardous location.

## Inserting the Batteries







With four AA batteries, the Portavo has an operating time of over 1000 h.

Open the battery compartment on the rear of the device. Be sure to observe the correct polarity when inserting the batteries (see markings in the battery chamber). Close the battery compartment cover and screw it handtight.

A special lithium-ion battery suited to the battery compartment is available for the Portavo 904. The battery is recharged through the USB port.

**Note:** Not available for the Portavo 904 X (device variant for applications in hazardous locations).

### A battery icon in the display indicates the battery power level:

|   |                       |  |
|---|-----------------------|--|
|  | Icon fully filled     | Batteries at full capacity   |
|  | Icon partially filled | Battery capacity is sufficient   |
|  | Icon empty            | Battery capacity not sufficient;<br>calibration is possible, no logging  |
|  | Icon blinks           | Max. 10 operating hours remaining,<br>measurement is still possible<br><b>Caution!</b> It is absolutely necessary to replace<br>the batteries. |



#### Warning!

When using the Portavo 904 X (device variant for applications in hazardous locations) in a hazardous location, only the battery types listed below may be used. The batteries must be from the same manufacturer and of identical type and capacity. Never use new and used batteries together (see also Control Drawing 209.009-110).



## IECEx

## Batteries for Application in Hazardous Locations

| Batteries (4x each)     | Temp. class | Ambient temperature range                   |
|-------------------------|-------------|---|
| Duracell MN1500         | T4          | $-10\text{ °C} \leq T_a \leq +40\text{ °C}$ |
| Energizer E91           | T3          | $-10\text{ °C} \leq T_a \leq +50\text{ °C}$ |
| Power One 4106          | T3          | $-10\text{ °C} \leq T_a \leq +50\text{ °C}$ |
| Panasonic Pro Power LR6 | T3          | $-10\text{ °C} \leq T_a \leq +50\text{ °C}$ |

## Connecting a Sensor

The Portavo 904(X) PH provides several connections so that many types of sensors can be used for measurement. Note that only **one** sensor may be connected to the meter at a time.

The meter automatically recognizes a connected Memosens sensor and switches accordingly. Memosens is signaled in the display.

### Separate temperature probe

After power-on, a separate temperature probe is automatically recognized.


When you want to replace the temperature probe, you must switch off the meter and then switch it on again.



### Caution!

Always make sure that a sensor is connected to the meter before starting measurement.

Explanation: The analog pH input of the Portavo is an electrometer amplifier with an extremely high-impedance. When the sensor is not in contact with the medium or not connected to the meter, electric charges on the input can generate arbitrary, stable pH or mV values which will be shown in the display.




### Connections

- a - Micro USB port
- b - M8, 4 pins for Memosens sensors
- c - Temperature probe GND
- d - Temperature probe
- e - pH socket (DIN 19 262)

Memosens sensors have a **cable coupling**, which allows convenient replacement of sensors while the cable remains connected to the meter. The connecting cable is connected to socket **b** (M8, 4 pins for Memosens sensors).





## Switching On the Meter

When you have connected the sensor, you can switch the meter on by pressing the **on/off** or **meas** key.


When the meter is switched on with the **on/off** key, first a self test is performed and then the calibration data and settings are displayed before the meter switches to measuring mode.

When the meter is switched on with the **meas** key, it immediately switches to measuring mode.

Depending on the connected sensor and the specific measuring task, several steps for configuration and calibration must be performed as described on the following pages.

## Icons

Important information about the state of the device:







## pH Configuration

Prior to measurement, a configuration should be performed to match the connected sensor and the desired measurement performance.

Furthermore, you can select the suitable calibration method.

The following table gives you an overview. Factory settings are shown in **bold print**.

Measurement

↓ **set**

"Setup" display

|  |  |   |
|--|--|---|
| <div style="display: flex; align-items: center;"> <div style="flex: 1; border-left: 1px solid black; border-right: 1px solid black; position: relative;"> <div style="position: absolute; left: -10px; top: 0; bottom: 0; border-left: 1px solid black; border-right: 1px solid black;"></div> <div style="position: absolute; left: -10px; top: 0; bottom: 0; border-left: 1px solid black; border-right: 1px solid black;"></div> </div> <div style="flex: 1; padding: 5px;"> <div>Display 1</div> <div>Display 2</div> <div>CAL Timer</div> <div>CAL</div> <div>CAL POINTS</div> <div style="text-align: center; padding: 10px;"> <div>▲</div> <div>▼</div> </div> <div>BUFFER SET<br/>(CALIMATIC,<br/>FREE CAL)</div> <div>Auto OFF</div> <div>Temp Unit</div> <div>Time Format</div> <div>Date Format</div> <div>Default</div> </div> <div style="flex: 0.5; text-align: center; padding: 0 10px;"> <div>set</div> <div>↔</div> </div> </div> | Select using arrow keys, confirm by pressing <b>set</b> .  |   |
|  | <b>pH x.xx</b> / pH x.xxx / mV / (°C for analog pH only)   |   |
|  | <b>OFF</b> / date + time / date / time                     |   |
|  | <b>OFF</b> / 1 ... 99 days                                 |   |
|  | <b>CALIMATIC</b> /Manual/DATA INPUT/(ISFET-Zero)/FREE CAL  |   |
|  | 1 / 2 / 3 / <b>1-2-3</b> (for CALIMATIC, Manual, FREE CAL) |   |
|  | -01-   | Mettler Toledo 2.00/4.01/7.00/9.21        |
|  | <b>-02-</b>  | Knick CaliMat 2.00/4.00/7.00/9.00/12.00   |
|  | -03-   | Ciba (94) 2.06/4.00/7.00/10.00            |
|  | -04-   | NIST technical 1.68/4.00/7.00/10.01/12.46 |
|  | -05-   | NIST standard 1.679/4.006/6.865/9.180     |
|  | -06-   | HACH 4.01/7.00/10.01/12.00                |
|  | -07-   | WTW techn. buffers 2.00/4.01/7.00/10.00   |
|  | -08-   | Hamilton 2.00/4.01/7.00/10.01/12.00       |
|  | -09-   | Reagecon 2.00/4.00/7.00/9.00/12.00        |
|  | -10-   | DIN 19267 1.09/4.65/6.79/9.23/12.75       |
|  | -U1-   | loadable via Paraly SW 112 (User)         |
|  | <b>OFF</b> / 0.1h / 1h / 6h / 12h                          |   |
|  | °C / °F  |   |
|  | <b>24h</b> / 12h   |   |
|  | <b>dd.mm.yy</b> / mm.dd.yy                                 |   |
|  | <b>NO</b> / YES (reset to factory settings)                |   |
|  | <b>Note:</b> All data logger entries will be deleted.      |   |


- ▲ This icon prompts you to select a menu item using the arrow keys –
- ▼ the selection is confirmed by pressing **set**.



## CALIMATIC Calibration

### (Calibration with automatic buffer recognition)

The calibration method is selected in the configuration menu. Calibration is required to adjust the sensor to the meter. It is indispensable for achieving comparable and reproducible measurement results.




**Note:** To abort calibration, you can press **meas** at any time. This will be confirmed by the display message "CAL ABORTED". Exception: When you have selected "CAL POINTS 1-2-3" and the first calibration step has been completed, the calibration process cannot be stopped any more.



## DATA INPUT Calibration

(Calibration by entering known sensor values)

The calibration method is selected in the configuration menu.



The calibration data will be displayed successively:


Date and time

ZERO POINT

SLOPE

Then the meter switches to measuring mode.


**Note:** To abort calibration, you can press **meas** at any time.



## MANUAL Calibration

### (Manual calibration)

The calibration method is selected in the configuration menu.



**Note:** To abort calibration, you can press **meas** at any time. This will be confirmed by the display message "CAL ABORTED". Exception: When you have selected "CAL POINTS 1-2-3" and the first calibration step has been completed, the calibration process cannot be stopped any more.



## FREE CAL Calibration

(Free selection of calibration method)

FREE CAL calibration is selected in the configuration menu.

Measurement



CAL  
CALIMATIC blinks





Use ▲▼ to select the required calibration method (CALIMATIC, DATA INPUT or MANUAL).

Perform the selected calibration (see CALIMATIC, DATA INPUT or MANUAL calibration).

Once you have completed all preparations, you can start with the actual measurement.

### Keys for measurement

- 1) Connect the desired sensor to the meter. Some sensors require a special preparation. Please proceed according to the operating instructions for the sensor.
- 2) Switch the meter on using the **on/off** or **meas** key.
- 3) Depending on the measurement method and the sensor used, immerse the sensing part of the sensor in the medium to be measured.
- 4) Watch the display and wait for the reading to stabilize.
- 5) By pressing the **STO** key, you can hold and save a measured value (see data logger, page 23).



Measurement can also be controlled via the Paraly SW 112 software.

## Switching the Measured Value Display

During measurement, you can switch between pH and mV display by pressing the **meas** key.

## Adjusting the Temperature

When you connect a sensor without temperature detector, you can manually adjust the temperature for measurement or calibration:

- 1) Press **meas** to access measuring mode.  
The adjusted temperature will be displayed.
- 2) Set the desired temperature value using the ▼ or ▲ arrow.  
Holding the key depressed changes the temperature value at high speed.

## The Data Logger

The meter provides a data logger. **Prior to use**, it must be configured and then activated. You can choose from the following logger types:


- DIFF (signal-controlled logging of measured variable and temperature)
- INT (time-controlled logging at a fixed interval)
- DIFF+INT (combined time- and signal-controlled logging)
- SHOT (manual logging by pressing the **STO** key)

The data logger records up to 5000 entries and saves them in a circular buffer. Already existing entries will be overwritten.

The following data are recorded: primary value, temperature, time stamp and device status.

The Paraly SW 112 software allows convenient management of the data logger. It is always the currently selected process variable which is recorded. The “STO” icon and the memory address is displayed briefly to indicate that an entry is being saved.

### Display: Icons related to the data logger



## Operating Modes of the Data Logger (Logger Type)

### Manual logging when logger is activated (SHOT)

In this mode, a measured value is recorded when the **STO** key is pressed.

Measurement  
Logger **activated**

↓ **STO**

The measured value is saved to the address of the last recorded value + 1

### Manual logging when logger is deactivated

Measurement  
Logger **deactivated**

↓ **STO**

Measured value is maintained  
Proposed address blinks  
(address of the last recorded  
value + 1)

If desired: Select start address  
using ▲▼.


↓ **STO**

Measured value is saved to the desired address (e.g. for overwriting an incorrect measurement).

### Interval (INT)

In this mode, the measured values are cyclically recorded.


Example: Logger type "INT"






**Difference (DIFF)**

When the delta range (process variable and/or temperature) related to the last entry is exceeded, a new entry is created and the delta range is displaced upwards or downwards by the delta value. The first entry is automatically created when the data logger is started.

**Difference + Interval combined (DIFF+INT)**


When the delta range related to the last DIFF entry is exceeded, a new entry is created (example: entry **A**) and the delta range is displaced upwards or downwards by the delta value. As long as the measured value remains within the delta range, logging is performed at the preset interval. The first DIFF entry is automatically created when the data logger is started.



## Data Logger Menu

Logger display

Select using arrow keys, confirm by pressing **set**.



### Overview of data logger menu (default in bold print)

|             |          |   |   |
|-------------|----------|---|---|
| Logger type | DIFF     | Delta pH / mV                                   | OFF / pH 0.01...14.00 / <b>pH 1.00</b><br>OFF / 1... 1000 mV / <b>1 mV</b>  |
|             |          | Delta °C / °F                                   | OFF / 0.1 ... 50.0 °C / <b>1.0 °C</b><br>OFF / 0.1 ...90 °F / <b>1.0 °F</b> |
|             | INT      | Interval  | h:mm:ss<br>0:00:01 ... 9:59:59 / <b>0:01:00</b>                             |
|             | DIFF+INT | DIFF  | See logger type DIFF  |
|             |          | INT   | See logger type INT   |
|             | SHOT     | Currently selected process variable is recorded |   |

## Configuring the Data Logger

Prerequisite: The data logger is stopped (press **meas**).

Measurement

↓ **STO**

Measured value is maintained

↓ **set**

Logger: CONT blinks

↓ ▼

Logger: START blinks

↓ ▼

Logger: DEL blinks

↓ ▼

Logger: SET blinks

↓ **set**

Logger: Current logger type  
blinks


↓ **set**

Select desired logger type using ▲▼:  
DIFF, INT, DIFF+INT or SHOT.

Select the appropriate parameters using ▲▼ and confirm each selection by pressing **set**. When configuration is finished, CONT blinks. You can start the data logger by selecting START or CONT (see page 28).

## Starting the Data Logger using CONT

Prerequisite: Data logger is configured. Every time the meter has been switched off, the data logger must be restarted (exception: SHOT).




The measured value is saved to the selected start address (exception: SHOT).

“... FREE MEMORY” is displayed.

“LOGGER” and “active logger type” icons are displayed.

## Starting the Data Logger using START

Prerequisite: Data logger is configured. All existing entries are deleted. The start address for saving the values is 0001. Every time the meter has been switched off, the data logger must be restarted (exception: SHOT).



All entries will be deleted. “5000 FREE MEMORY” is displayed.

“LOGGER” and “active logger type” icons are displayed.

## Displaying the Logger Data

Pressing the **RCL** key displays all stored values. The Paraly SW 112 software allows convenient management of the data logger.

Measurement

**RCL**




The "RCL" icon and the last recorded value is displayed.

Use ▲▼ to select the desired address. Empty memory locations will also be displayed.


**RCL or meas**



Return to measurement




Example:  
Measured value stored at location  
0026



Example:  
Empty memory location 0004


## Stopping the Data Logger


You can stop the data logger at any time by pressing the **meas** key.



## Clearing the Data Logger

Selecting "DEL" deletes all data records.





Press the **clock** key to access the clock mode. Date and time will be displayed in the format as set in the configuration menu.

To set the clock, proceed as follows:

Display of  
time+date

↓ **set**

Hour display blinks  
SET HOUR



Set value.

↓ **set**

Minute display blinks  
SET MINUTE



Set value.

↓ **set**

Second display blinks and  
shows 00

**set**

Clock is started, the seconds count up.

↓ **set**

Year display blinks  
SET YEAR



Set value.

↓ **set**

Month display blinks  
SET MONTH



Set value.

↓ **set**

Day display blinks  
SET DAY



Set value.

↓ **set**

Display of  
corrected time+date

The Paraly SW 112 software supplements the Portavo series. It allows convenient management of the data that have been acquired by the meters as well as simple and clear configuration of the meters. Paraly SW 112 starts automatically when the Portavo USB port is connected to the computer.


The Paraly SW 112 software stands out by the following features:

- Intuitive Windows user interface
- Easy configuration and management of several meters
- Display of device and sensor information
- Configuration of individual buffer sets
- Convenient management and evaluation of the data logger
- Export function for Microsoft Excel
- Print function
- Updating the device software


**Note:** A detailed user manual for the Paraly SW 112 software can be found on the included data carrier.



Error messages are indicated as "ERROR ..." on the display. Information on the sensor condition is indicated by the "Sensoface" icon (friendly, neutral, sad) possibly accompanied by an info message ("INFO ...").



Example of an error message:  
ERROR 8 (identical calibration media)



Example of a "Sensoface" message:  
INFO 1 (cal timer expired)

Sensoface (the "smiley" icon) provides information on the sensor condition (maintenance request). Measurement can still be performed. After a calibration, the corresponding Sensoface icon (friendly, neutral, sad) is shown together with the calibration data. Otherwise, Sensoface is only visible in measuring mode.

The most important error messages and "Sensoface" info messages are shown on the inside of the protective cover. A complete list of messages and their meanings is provided in the following tables.



## “Sensoface” Messages

The “Sensoface” icon provides information on the sensor condition:

### Sensoface    Meaning



Sensor is okay





Calibrate the sensor soon




Calibrate or replace the sensor

The “neutral” and “sad” Sensoface icons are accompanied by an “INFO ...” message to give a hint to the cause of deterioration.

| Sensoface   | Message | Cause  |
|---|---------|--|
|   | INFO 1  | Calibration timer                            |
|   | INFO 3  | Sensocheck                                   |
|   | INFO 5  | Zero / Slope                                 |
|   | INFO 6  | Response time                                |
|   | INFO 7  | ISFET: Operating point (asymmetry potential) |
|   | INFO 8  | ISFET: Leakage current                       |
|   | INFO 9  | ORP offset                                   |

## Error Messages

The following error messages can be shown in the display.

| Message  | Cause                                       | Remedy   |
|--|---|--|
|  blinks | Battery empty                               | Replace batteries  |
| ERROR 1  | pH value out of range                       | Check whether the measurement conditions correspond to the adjusted measuring range.   |
| ERROR 2  | ORP value out of range                      |  |
| ERROR 3  | Temperature value out of range              |  |
| ERROR 4  | Sensor zero point too high/low              | Thoroughly rinse the sensor and re-calibrate. If this does not help, replace the sensor.   |
| ERROR 5  | Sensor slope too high/low                   |  |
| ERROR 8  | Calibration error:<br>Identical buffers     | Use a buffer solution with a different nominal value before starting the next calibration step.  |
| ERROR 9  | Calibration error:<br>Buffer unknown        | Make sure that you use the same buffer set as configured.  |
| ERROR 10   | Cal media interchanged                      | Repeat calibration.  |
| ERROR 11   | Measured value unstable<br>Drift too high   | Leave the sensor in the liquid until the temperature is stable. If this does not help, replace the sensor.                                 |
| ERROR 14   | Time and date invalid                       | Set time and date  |
| ERROR 18   | Configuration invalid                       | Restart, reset to factory settings (Setup: DEFAULT YES), configure and calibrate.<br>If this does not help, send in the device for repair. |
| ERROR 19   | Factory settings error                      | Device defective, send it in.  |
| ERROR 21   | Sensor error (Memosens)                     | Connect operational Memosens sensor.   |
| ERROR 25   | Buffer distance (user-defined buffer table) | Re-enter the buffer table.   |

## Accessories

| Item   | Order No.      |
|--|----------------|
| Robust field case (for meter, sensor, various small parts and user manual) | ZU 0934        |
| Adapter for BNC pH sensors to DIN socket                                   | ZU1190         |
| Replacement quiver (5 units)   | ZU 0929        |
| Memosens lab cable, M8, 4 pins   | CA/MS-001XFA-L |
| Li-ion battery   | ZU 0925        |

Please visit our website for more information on our product range: [www.knick.de](http://www.knick.de).

## Sensors

### Analog pH-sensors

|   |          |
|---|----------|
| pH/Pt1000 sensor (plastic body, length 120 mm)  | SE 101 N |
| pH/Pt1000 sensor (glass body, length 110 mm)  | SE 102 N |
| pH puncture sensor (plastic body, length 65/25 mm)  | SE 104 N |
| pH/Pt1000 sensor for measurements in hazardous areas<br>Zone 0, incl. equipotential bonding cable | ZU 6979  |
| Pt1000 temperature detector   | ZU 6959  |
| Pt 1000 temperature detector with tilted tip  | ZU 0156  |

### Digital pH sensors

|  |            |
|--|------------|
| pH/temp sensor (plastic body, length 120 mm) | SE 101 NMS |
| pH/temp sensor (glass body, length 110 mm)   | SE 102 NMS |

Memosens sensors have a **cable coupling**, which allows convenient replacement of sensors while the cable remains connected to the meter.



## Knick CaliMat Buffer Solutions

Ready-to-use quality pH buffer solutions

| pH value (20 °C) | Quantity | Order No.     |
|------------------|----------|---------------|
| 2.00 ± 0.02      | 250 ml   | CS-P0200/250  |
| 4.00 ± 0.02      | 250 ml   | CS-P0400/250  |
|                  | 1000 ml  | CS-P0400/1000 |
|                  | 3000 ml  | CS-P0400/3000 |
| 7.00 ± 0.02      | 250 ml   | CS-P0700/250  |
|                  | 1000 ml  | CS-P0700/1000 |
|                  | 3000 ml  | CS-P0700/3000 |
| 9.00 ± 0.02      | 250 ml   | CS-P0900/250  |
|                  | 1000 ml  | CS-P0900/1000 |
|                  | 3000 ml  | CS-P0900/3000 |
| 12.00 ± 0.05     | 250 ml   | CS-P1200/250  |

### Buffer sets

|                      |             |            |
|----------------------|-------------|------------|
| Set 4.00             | 3 x 250 ml  | CS-PSET4   |
| Set 7.00             | 3 x 250 ml  | CS-PSET7   |
| Set 9.00             | 3 x 250 ml  | CS-PSET9   |
| Set 4.00, 7.00, 9.00 | 250 ml each | CS-PSET479 |

|                                      |   |                         |                             |
|--------------------------------------|---|-------------------------|-----------------------------|
| pH/mV input                          | pH socket, DIN 19 262 (13/4 mm)                               |                         |                             |
| pH range                             | -2 ... 16   |                         |                             |
| Decimal places <sup>*)</sup>         | 2 or 3  |                         |                             |
|                                      | Input resistance  | 1 x 10 <sup>12</sup> Ω  | (0 ... 35 °C)               |
|                                      | Input current   | 1 x 10 <sup>-12</sup> A | (at RT, doubles every 10 K) |
| Measuring cycle                      | Approx. 1 s   |                         |                             |
| Measurement error <sup>1,2,3)</sup>  | < 0.01 pH, TC < 0.001 pH/K                                    |                         |                             |
| mV range                             | -1300 ... +1300 mV  |                         |                             |
| Measuring cycle                      | Approx. 1 s   |                         |                             |
| Measurement error <sup>1,2,3)</sup>  | < 0.1 % meas. val. + 0.3 mV, TC < 0.03 mV/K                   |                         |                             |
| <b>Temperature input</b>             | 2 x 4 mm dia. for integrated or separate temperature detector |                         |                             |
| Measuring ranges                     | NTC30 temp detector   | -20 ... +120°C          |                             |
|                                      | Pt1000 temp detector  | -40 ... +250°C          |                             |
| Measuring cycle                      | Approx. 1 s   |                         |                             |
| Measurement error <sup>1,2,3)</sup>  | < 0.2 K (Tamb = 23 °C); TC < 25 ppm/K                         |                         |                             |
| <b>Memosens pH input</b>             | M8 socket, 4 pins, for Memosens lab cable                     |                         |                             |
| Display ranges <sup>4)</sup>         | pH  | -2.00 ... +16.00        |                             |
|                                      | mV  | -2000 ... +2000 mV      |                             |
|                                      | Temperature   | -50 ... +250 °C         |                             |
| <b>Memosens pH input ISFET</b>       | M8 socket, 4 pins, for Memosens lab cable                     |                         |                             |
| Display ranges <sup>4)</sup>         | pH  | -2.00 ... +16.00        |                             |
|                                      | mV  | -2000 ... +2000 mV      |                             |
|                                      | Temperature   | -50 ... +250 °C         |                             |
| <b>Memosens ORP input</b>            | M8 socket, 4 pins, for Memosens lab cable                     |                         |                             |
| Display ranges <sup>4)</sup>         | mV  | -2000 ... +2000 mV      |                             |
|                                      | Temperature   | -50 ... +250 °C         |                             |
| Sensor standardization <sup>*)</sup> | ORP calibration (zero adjustment)                             |                         |                             |
| Permissible calibration range        | ΔmV (offset)  | -700 ... +700 mV        |                             |

<sup>\*)</sup> User-defined

1) According to EN 60746-1,  
at nominal operating conditions

2) ± 1 count

3) Plus sensor error

4) Ranges depending on Memosens sensor

|                                  |  |   |
|----------------------------------|--|---|
| <b>Sensor standardization *)</b> | pH calibration                               |   |
| Operating modes *)               | CALIMATIC                                    | Calibration with automatic buffer recognition                       |
|                                  | MANUAL                                       | Manual calibration with entry of individual buffer values           |
|                                  | DATA INPUT                                   | Data entry of zero and slope  |
| Calimatic buffer sets *)         | -01- Mettler-Toledo                          | 2.00/4.01/7.00/9.21   |
|                                  | -02- Knick CaliMat                           | 2.00/4.00/7.00/9.00/12.00   |
|                                  | -03- Ciba (94)                               | 2.06/4.00/7.00/10.00  |
|                                  | -04- NIST technical                          | 1.68/4.00/7.00/10.01/12.46  |
|                                  | -05- NIST standard                           | 1.679/4.006/6.865/9.180   |
|                                  | -06- HACH                                    | 4.01/7.00/10.01/12.00   |
|                                  | -07- WTW techn. buffers                      | 2.00/4.01/7.00/10.00  |
|                                  | -08- Hamilton                                | 2.00/4.01/7.00/10.01/12.00  |
|                                  | -09- Reagecon                                | 2.00/4.00/7.00/9.00/12.00   |
|                                  | -10- DIN 19267                               | 1.09/4.65/6.79/9.23/12.75   |
|                                  | -U1- (User)                                  | loadable via Paraly SW 112  |
| Permissible calibration range    | Zero point                                   | 6 ... 8 pH  |
|                                  | With ISFET:                                  | -750 ... +750 mV  |
|                                  | Operating point (asymmetry)                  |   |
|                                  | Slope  | approx. 74 ... 104 %<br>(possibly restricting notes from Sensoface) |
| <b>Calibration timer *)</b>      | Interval 1 ... 99 days, can be switched off  |   |
| <b>Sensoface</b>                 | Provides information on the sensor condition |   |
| Evaluation of                    | zero/slope, response, calibration interval   |   |

\*) User-defined



|                                |  |        |                          |        |                                |             |   |
|--------------------------------|--|--------|--------------------------|--------|--------------------------------|-------------|---|
| <b>Connections</b>             | 1 x pH socket, DIN 19 262<br>2 x 4-mm socket for separate temperature detector<br>1 x M8 socket, 4 pins, for Memosens lab cable<br>1 x Micro USB-B for data transmission to PC<br><b>Portavo 904 X:</b><br>Be sure to observe the safety instructions when using the USB port.   |        |                          |        |                                |             |   |
| <b>Display</b>                 | LCD STN 7-segment display with 3 lines and icons   |        |                          |        |                                |             |   |
| Sensoface                      | Status indication (friendly, neutral, sad)   |        |                          |        |                                |             |   |
| Status indicators              | Battery power level, logger  |        |                          |        |                                |             |   |
| Notices                        | Hourglass  |        |                          |        |                                |             |   |
| Keypad                         | [on/off], [cal], [meas], [set], [▲], [▼], [STO], [RCL], [clock]  |        |                          |        |                                |             |   |
| <b>Data logger</b>             | With up to 5000 memory locations   |        |                          |        |                                |             |   |
| Recording                      | Manual, interval- or event-controlled  |        |                          |        |                                |             |   |
| <b>Communication</b>           | USB 2.0  |        |                          |        |                                |             |   |
| Profile                        | HID, driverless installation   |        |                          |        |                                |             |   |
| Usage                          | Data exchange and configuration via Paraly SW 112 software   |        |                          |        |                                |             |   |
| <b>Diagnostics functions</b>   |  |        |                          |        |                                |             |   |
| Sensor data<br>(Memosens only) | Manufacturer, sensor type, serial number, operating time   |        |                          |        |                                |             |   |
| Calibration data               | Calibration date, zero, slope  |        |                          |        |                                |             |   |
| Device self-test               | Automatic memory test (FLASH, EEPROM, RAM)   |        |                          |        |                                |             |   |
| Device data                    | Device type, software version, hardware version  |        |                          |        |                                |             |   |
| <b>Data retention</b>          | Parameters, calibration data > 10 years  |        |                          |        |                                |             |   |
| <b>EMC</b>                     | EN 61326-1 (General Requirements)  |        |                          |        |                                |             |   |
| Emitted interference           | Class B (residential area)   |        |                          |        |                                |             |   |
| Immunity to interference       | Industry<br>EN 61326-2-3<br>(Particular Requirements for Transmitters)   |        |                          |        |                                |             |   |
| <b>Explosion protection</b>    | Portavo 904 X<br><table> <tr> <td>Global</td><td>IECEX Ex ia IIC T4/T3 Ga</td></tr> <tr> <td>Europe</td><td>ATEX II 1 G Ex ia IIC T4/T3 Ga</td></tr> <tr> <td>USA, Canada</td><td>               IS Class I, Division 1, Groups A,B,C,D, T4 / T3,<br/>               Ta = 40 °C / 50 °C; Entity; Type 4X<br/>               IS Class I, Zone 0, AEx ia IIC T4 / T3,<br/>               Ta = 40 °C / 50 °C; Entity; Type 4X             </td></tr> </table> For electrical parameters and further specifications, see Control Drawing No. 209.009-110 | Global | IECEX Ex ia IIC T4/T3 Ga | Europe | ATEX II 1 G Ex ia IIC T4/T3 Ga | USA, Canada | IS Class I, Division 1, Groups A,B,C,D, T4 / T3,<br>Ta = 40 °C / 50 °C; Entity; Type 4X<br>IS Class I, Zone 0, AEx ia IIC T4 / T3,<br>Ta = 40 °C / 50 °C; Entity; Type 4X |
| Global                         | IECEX Ex ia IIC T4/T3 Ga   |        |                          |        |                                |             |   |
| Europe                         | ATEX II 1 G Ex ia IIC T4/T3 Ga   |        |                          |        |                                |             |   |
| USA, Canada                    | IS Class I, Division 1, Groups A,B,C,D, T4 / T3,<br>Ta = 40 °C / 50 °C; Entity; Type 4X<br>IS Class I, Zone 0, AEx ia IIC T4 / T3,<br>Ta = 40 °C / 50 °C; Entity; Type 4X  |        |                          |        |                                |             |   |

|                                     |  |    |                         |
|-------------------------------------|--|----|-------------------------|
| <b>RoHS conformity</b>              | According to directive 2011/65/EC  |    |                         |
| <b>Power supply</b>                 |  |    |                         |
| Portavo 904                         | Batteries: 4 x AA alkaline or 4 x NiMH (rechargeable)<br>or 1 x Li-ion battery, USB chargeable |    |                         |
| Portavo 904 X                       | 4 x AA batteries<br>For battery types, see Control Drawing No. 209.009-110                     |    |                         |
| Operating time                      | Approx. 1000 h (alkaline)  |    |                         |
| <b>Nominal operating conditions</b> |  |    |                         |
| Ambient temperature                 | -10 °C ... +55 °C  |    |                         |
| Ambient temperature 904 X           | -10 °C ≤ Ta ≤ +40 °C   | T4 | Duracell MN1500         |
|                                     | -10 °C ≤ Ta ≤ +50 °C   | T3 | Energizer E91           |
|                                     | -10 °C ≤ Ta ≤ +50 °C   | T3 | Power One 4106          |
|                                     | -10 °C ≤ Ta ≤ +50 °C   | T3 | Panasonic Pro Power LR6 |
| Transport/<br>Storage temperature   | -25 ... +70 °C   |    |                         |
| Relative humidity                   | 0 ... 95 %, short-term condensing allowed  |    |                         |
| <b>Housing</b>                      |  |    |                         |
| Material                            | PA12 GF30 (silver gray RAL 7001) + TPE (black)   |    |                         |
| Protection                          | IP 66/67 with pressure compensation  |    |                         |
| Dimensions                          | Approx. (132 x 156 x 30) mm  |    |                         |
| Weight                              | Approx. 500 g  |    |                         |

0000 DELETED ("data deleted" display) 30

## A

AA batteries 13

Accessories 36

Activating the logger 28

Analog pH input 15

Arrow keys 12

ATEX 7

Automatic calibration (Calimatic) 18

## B

Batteries for application in hazardous locations 14

Battery capacity 14

Battery charge indicator 8

Battery compartment 13

Battery icon 14

Battery replacement 13

Battery types 14

Benchtop stand 10

Buffer sets 38

Buffer solutions (Knick CaliMat) 38

Buttons 12

## C

Calibration (pH), Calimatic 18

Calibration (pH), data input 19

Calibration (pH), manual 20

Calibration (pH), method selected in measuring mode 21

CaliMat buffer solutions 38

Calimatic automatic calibration 18

Calimatic, description 9

cal key 12

Carrying case (accessory) 36

CD-ROM 7

Certificates 7

Charge level of batteries 14

Clearing the datalogger 30

Clock 31

clock key 12

- Configuration (pH) 17
- Configuring the data logger 27
- Connecting a sensor 15
- Connecting cable for Memosens 15
- Connections 15
- Connection, USB 13
- Continuous recording of measured values 24
- Control elements 12
- CONT, starting the data logger 28
- Cyclic recording of measured values 24

**D**

- Data input (pH calibration) 19
- Data logger 23
- Data logger, activating 28
- Data logger, clearing 30
- Data logger configuration 27
- Data logger icons 23
- Data logger menu 26
- Data logger, stopping 30
- Data memory 23
- Data of the meter 39
- Date 31
- Deleting data logger entries 30
- Delta range (data logger) 25
- Device configuration 17
- Device messages 33
- Device properties 8, 9
- Difference (data logger mode) 25
- Difference+Interval (data logger mode) 25
- Digital sensors (pH) 37
- Display 11
- Display icons 16
- Displaying the data logger 23
- Displaying the time and date 31
- Display, switching between measured values 22
- Disposal 3
- Documentation 7
- Duracell MN1500 battery 14

## **E**

EC Declarations of Conformity 7  
Energizer E91 battery 14  
ERROR (error codes) 35  
Error messages 33  
Error messages, overview 35

## **F**

Features 8  
Field case (accessory) 36  
FREE CAL 21

## **H**

Hanging up the meter 10  
Hazardous-area application 8  
Hazardous location, batteries 14  
Holding the measured value 24  
Hook 10  
Hours, display 31

## **I**

Icons for data logger 23  
Icons in display 16  
IECEX 7  
INFO messages 34  
Inserting the batteries 13  
Interrupting the data logger 30  
Interval (data logger mode) 24  
Introduction 8

## **K**

Keypad 12  
Knick CaliMat buffer solutions 38

## **L**

Li-ion battery (accessory) 36  
Lithium-ion battery 13  
Logger 23  
Logger, activating 28  
Logger data, deleting 30  
Logger data, viewing 29  
Logger display 29  
Logger type (data logger modes) 24

**M**

Manual calibration 20  
Manual logging 24  
meas key 12  
Measured-value recording 24  
Measured values, switching the display 22  
Measuring 22  
Memory for measured values 23  
Memosens 9  
Memosens connecting cable 15  
Memosens lab cable (accessory) 36  
Memosens sensors 15  
Memosens sensors, product line 37  
Menu of data logger 26  
Menu structure of configuration 17  
Menu structure of data logger 26  
Messages 33  
Micro USB port 15  
Minutes, display 31

**O**

on/off key 12  
Operating modes of the data logger 24  
Order numbers (accessories) 36  
Overview 8  
Overview of configuration 17  
Overview of error messages 35

**P**

Package contents 6  
Panasonic Pro Power LR6 battery 14  
Paraly SW 112 software 32  
Parameter setting, data logger 27  
Parameter settings (configuration) 17  
pH buffer solutions 38  
pH configuration 17  
pH sensor 15  
pH sensors, product line 37  
pH socket, DIN 19 262 15  
Portavo 904 X 13  
Ports 15

Power-on 16  
Power One 4106 battery 14  
Product features 8  
Product presentation 8  
Protective cover 10

## Q

Quickstart guides 7

## R

Rating plate 10  
RCL, displaying the logger data 29  
RCL key 12  
Real-time clock 8  
Rechargeable battery 13  
Recorded data, display 29  
Reference numbers (accessories) 36  
Registered trademarks 3  
Replacement quiver (accessory) 36  
Replacing the batteries 14  
Return of products under warranty 3

## S

Safety instructions 7  
Saving the currently measured value 24  
Seconds, display 31  
Sensoface 9  
Sensoface messages 34  
Sensor connection 15  
Sensors 15  
Sensors, product line 37  
Sensor without temperature detector 22  
set key 12  
Setting the configuration data 17  
Setting the data logger 27  
Setting the time and date 31  
Setup (configuration) 17  
SHOT (data logger mode) 24  
Smiley face (icon) 9  
Software, Paraly SW 112 32  
Specifications 39

- Specific test report 7
- Start address (data logger) 24
- Starting the data logger using CONT 28
- Starting the data logger using START 28
- START, starting the data logger 28
- Start-up 13
- STO key 12
- STO, manual logging 24
- Stopping the data logger 30
- Structure of data logger 26
- Suspending the meter 10
- Switching on the meter 16
- Switching the measured value display 22
- Symbols in the display 16

## **T**

- T3, temperature class 14
- T4, temperature class 14
- Table of error messages 35
- Table view of configuration 17
- Technical data 39
- Temperature class 14
- Temperature, manual adjustment 22
- Temperature probe, connection 15
- Trademarks 3
- Triangle icons 12

## **U**

- USB port 15
- USB port, micro 8

## **V**

- Value-added features 9
- Viewing logger data 29
- Viewing recorded data 29

## **W**

- Warranty 3









---

**Knick**  
**Elektronische Messgeräte GmbH & Co. KG**



Beuckestr. 22  
14163 Berlin  
Germany

Phone: +49 (0)30 - 801 91 - 0  
Fax: +49 (0)30 - 801 91 - 200  
Web: [www.knick.de](http://www.knick.de)  
Email: [knick@knick.de](mailto:knick@knick.de)



085126