# Lovibond® Water Testing

Tintometer® Group

## Photometer XD 7000

### Quickstart-Guide

EN

DE

ES

FR

IT

PT

PL

NL

ZH

JA

[www.lovibond.com](http://www.lovibond.com)

ba77194y01  
71310070-1 09/2017

## XD 7000

<b>EN English</b>	Quick start guide	<b>3</b>
<b>DE Deutsch</b>	Kurzanleitung	<b>19</b>
<b>ES Español</b>	Breve manual	<b>35</b>
<b>FR Français</b>	Notice d'utilisation de base	<b>51</b>
<b>IT Italiano</b>	Manuale breve	<b>67</b>
<b>PT Português</b>	Instruções de utilização básicas	<b>83</b>
<b>PL Polski</b>	Podstawowa instrukcja obsługi	<b>99</b>
<b>NL Nederlands</b>	Basis-gebruikershandleiding	<b>115</b>
<b>ZH 中文</b>	基础版操作说明书	<b>131</b>
<b>JA 日本語</b>	取扱説明書	<b>147</b>

# 1 Safety

## 1.1 General notes

Any use other than the intended use, non-compliance with these operating instructions, deployment of insufficiently qualified personnel and unauthorized modifications on the instrument will void the manufacturer's liability and warranty for damages and subsequent damages therefrom.

The manufacturer is not liable for costs or damages incurred by the user or third party resulting from the use of this instrument, specifically in the event of unintended use of the instrument or misuse or malfunctions of the connection or the instrument.

The manufacturer is not liable for misprints.

## 1.2 Safety information

### 1.2.1 Safety information in the operating manual

This operating manual provides important information on the safe operation of the product. Read this operating manual thoroughly and make yourself familiar with the product before putting it into operation or working with it. The operating manual must be kept in the vicinity of the product so you can always find the information you need.

Important safety instructions are highlighted in this operating manual. They are indicated by the warning symbol (triangle) in the left column. The signal word (e.g. "CAUTION") indicates the level of danger:

#### **WARNING**

indicates a possibly dangerous situation that can lead to serious (irreversible) injury or death if the safety instruction is not followed.

#### **CAUTION**

indicates a possibly dangerous situation that can lead to slight (reversible) injury if the safety instruction is not followed.

#### **NOTE**

*indicates a situation where goods might be damaged if the actions mentioned are not taken.*

### 1.2.2 Safety signs on the product

Note all labels, information signs and safety symbols on the product. A warning symbol (triangle) without text refers to safety information in this operating manual.

## 1.3 Safe operation

### 1.3.1 Authorized use

The authorized use of the photometer consists exclusively of the carrying out of photometric measurements according to this operating manual. Any other use is considered **unauthorized**.

### 1.3.2 Requirements for safe operation

Note the following points for safe operation:

- The product may only be operated according to the authorized use specified above.
- The product may only be supplied with power by the energy sources mentioned in this operating manual.
- The product may only be operated under the environmental conditions mentioned in this operating manual.
- The product may not be opened.

### 1.3.3 Unauthorized use

The product must not be put into operation if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time (storing conditions, see Chapter 8 TECHNICAL DATA)

## 1.4 Handling of hazardous substances

When developing test sets, Tintometer carefully sees that the tests can be carried out as safely as possible. Some hazards by dangerous substances, however, cannot always be avoided.

If self-produced tests or solutions are used, the responsibility concerning any risks caused by those tests or solutions lies with the user (personal responsibility).

### WARNING

**Improper handling of certain reagents can cause damage to your health.**

**In any case follow the safety labels on the packing and the safety instructions of the package insert. Protective measures specified there have to be followed exactly.**

### Safety datasheets

The safety datasheets of the chemicals comprise all instructions on safe handling, occurring hazards, preventive actions and actions to take in hazardous situations. Follow these instructions in order to work safely.

## 2 Overview

### 2.1 Scope of delivery

- Spectrophotometer XD 7000
- Power pack with connection cable
- Buffer batteries 4 x AA alkaline manganese (Mignon)
- Zero cell (16 mm, round)
- Short instructions
- CD-ROM with
  - Detailed operating manual
  - Analysis instructions

### 2.2 Overview of the instrument

Front of the instrument


Figure 2-1 Front of the instrument with operating elements

### Socket field on the rear panel



Figure 2-2 Rear panel with socket field

All connections comply with SELV.

## 2.3 Keypad

### Overview

- 1 Function keys F1 to F4 (function menu-dependent)
- 2 Alphanumeric keypad
- 3 Keys with dedicated function

Figure 2-3 Keypad


**Key functions** The keys on the right side of the keypad have the following functions:

Key	Designation	Functions
	<ON/OFF>	– Switches on and off the photometer
	<HOME>	– Switches to the main menu from any operating situation. Actions that are not completed are canceled.
	<PRINT>	– Outputs the displayed measured value to an interface.
	<STORE>	– Saves a displayed measured value or spectrum
	<ZERO·BLANK>	– Starts one of the following measurements, depending on the operating situation: - Zero adjustment - Blank value measurement - Baseline measurement - User calibration
	<TIMER>	– Opens the menu, <i>Timer</i> .
	<ESC>	– Cancels the running action. Entries that have not yet been accepted are discarded. – Switches to the next higher menu level.
	<START·ENTER>	– Starts an action (e.g. measurement) – Opens a selected menu – Confirms a selection or entry
 (Arrow keys)	<▲><▼>	– Moves the selection in menus and lists one position up or down
	<◀>	– Deletes the character left of the cursor during character entries – Moves the cursor to the left in a spectrum or kinetic diagram
	<▶>	– Moves the cursor to the right in a spectrum or kinetic diagram

**Function keys** The function keys F1 to F4 have different functions depending on the operating situation. The current functions are displayed in the function key menu at the bottom edge of the display (see Section 4.2).

## 2.4 Display

### Display elements

- 1 Status line (current state, date and time)
- 2 Display range for menus and measurement results
- 3 Function keys menu

Figure 2-4 Display

### Symbols in the status line

Symbol	Designation	Function
	Save	The <STORE> key is active. You can store the displayed data with <STORE>.
	Printer	The <PRINT> key is active. You can output the displayed data to an interface with <PRINT>.

## 3 Commissioning

### 3.1 General notes on handling

Always protect the photometer XD 7000 from conditions that could damage the mechanical, optical and electronic components. Heed the following points especially:

- The temperature and humidity during operation and storage must be within the limits specified in Chapter 8 TECHNICAL DATA.
- The following influences always have to be avoided with the meter:
  - Extreme dust, moisture and wetness
  - Exposure to intensive light and heat
  - Fumes that are corrosive or contain high concentrations of solvents.
- For measuring, the meter must be placed on a flat surface.
- Spilled liquid or other material must be removed immediately.

- If a cell has broken in the cell shaft, the cell shaft must be cleaned immediately (see Section 5.2 ACTIONS IN THE CASE OF A BROKEN CELL).
- The cell shaft should always be closed when the photometer is not used.
- For the instrument to be transported the cell shaft has to be empty.
- For mobile use, we recommend using the transport case.

## 3.2 Initial commissioning

Perform the following activities:

- Insert the buffer batteries (see Section 3.2.1)
- Connect the power supply (see Section 3.2.2)
- Switch on the photometer (see Section 3.2.3)
- Set the language (see Section 3.2.3)
- Set the date and time (see Section 3.2.4)
- Carry out a zero adjustment

### Packing

This photometer is sent out in a protective transport packing.


### CAUTION

**Keep the original packing including the inner packing to protect the instrument against hard shocks if it has to be transported.**

### 3.2.1 Inserting the buffer batteries

The back-up batteries supply the integrated clock when the photometer is switched off. There are four alkaline manganese back-up batteries (type AA or Mignon) that are separately included in the delivery.

Insert the batteries as follows:

- 1 Turn the photometer upside down and place it on a soft surface.
- 2 Open the lid of the battery compartment (1).
- 3 Insert the four batteries in the battery compartment. The  $\pm$  indicators in the battery compartment must match the  $\pm$  indicators on the batteries.
- 4 Close the lid of the battery compartment.

### 3.2.2 Connecting the power supply


#### CAUTION

The line voltage of the usage location must fulfill the specifications stated on the power pack (the specifications are also given in Chapter 8 TECHNICAL DATA). Always use the supplied 12 V original power pack only.

Before plugging in the power cable check whether it is undamaged. If the power cable is damaged, the instrument must not be operated.

Connect the power pack

- 1 Connect the miniplug of the power pack to the socket (1) of the photometer.
- 2 Connect the power pack to an easily accessible power outlet.

### 3.2.3 Switch on the photometer and set the language

When the instrument is switched on during the initial commissioning, you will be guided to the setting of the instrument language automatically (<ON/OFF>).

General setup	16.01.15 9:52
German	
✓ English	
Français	
Español	
Italiano	
Bulgarian/Български	
Česko	
Chinese/ 中文	
Traditional Chinese/ 繁	
Greek/Ελληνικά	

- 1 Select a language with <▲><▼>.
- 2 Confirm the selected language with <START·ENTER>.

The language has been set.  
The display switches to the setting of the *Date* and *Time*.

### 3.2.4 Setting the date and time

During the initial commissioning, the instrument automatically guides you to the setting of the time and date after the setting of the language.

Date/Time	16.01.15 9:52
Date	16.01.2015
Time	9:52:09
	OK

The *Date/Time* menu is open.

Using <▲><▼>, select a menu item and confirm or open it with <START·ENTER>.

- 1 Select and confirm *Date*.

The input field for the current date pops up.

Date/Time	16.01.15 9:52
Date	16.01.2015
Time	9:52:09
Date	16 .01.2015
	OK

- 2 Enter the current date with <0...9> and confirm.

The input field closes.  
The date is accepted.

- 3 Select and confirm *Time*.

The input field for the current time pops up.

- 4 Set the time also.


Upon completion of the initial commissioning, you can adjust the language, data and time at any time in the *General setup* menu.

## 4 Operation

### 4.1 Switching the photometer ON/OFF

#### Switching ON


- 1 Use **<ON/OFF>** to switch the photometer ON.

The display shows

- the dialog *Self test* (with deactivated user management)
- or
- the dialog *Login* (with activated user management).

#### Self test start


- 2 Remove all cells and close the cell shaft cover.
- 3 Use **<START·ENTER>** to start the self-test.

The photometer carries out the self-test.

#### Self test

All cells must be removed and the cell shaft covers must be closed during the self-test

The self-test includes:

- the test of memory, processor, internal interfaces, filter and lamp
- a calibration for each shaft length


After the self-test is finished, the main menu will be displayed on the screen.

You can review and print the result of the self-test with the *[Info]* function key.

**Switching OFF** Press and hold the **<ON/OFF>** key until the photometer switches OFF.

## 4.2 Navigation with function keys and menus

## 5 What to do if...

### 5.1 Error causes and remedies

**Meter does not react to keystroke**

Cause	Remedy
<ul style="list-style-type: none"> <li>– Operating condition undefined or EMC load unallowed</li> </ul>	<ul style="list-style-type: none"> <li>– Processor reset: Press the <b>&lt;ON/OFF&gt;</b> and <b>&lt;ESC&gt;</b> key simultaneously</li> </ul>

**Self-test does not start.  
The instrument displays  
*Please remove cell***

Cause	Remedy
<ul style="list-style-type: none"> <li>– A cell is inserted in one of the cell shafts</li> </ul>	<ul style="list-style-type: none"> <li>– Remove the cell</li> <li>– Then press the <b>&lt;START·ENTER&gt;</b> key</li> </ul>
<ul style="list-style-type: none"> <li>– A foreign object is inserted in one of the cell shafts</li> </ul>	<ul style="list-style-type: none"> <li>– Remove foreign object</li> <li>– Then press the <b>&lt;START·ENTER&gt;</b> key</li> </ul>
<ul style="list-style-type: none"> <li>– Occasionally, the instrument carries out an automatic readjustment for the rectangular cell recognition. The informative message <i>Please remove cell</i> is displayed even when no cell is inserted.</li> </ul>	<ul style="list-style-type: none"> <li>– Press the <b>&lt;START·ENTER&gt;</b> key</li> </ul>
<ul style="list-style-type: none"> <li>– The cell shaft is contaminated</li> </ul>	<ul style="list-style-type: none"> <li>– Clean the cell shaft (see Section 5.2)</li> <li>– Restart the instrument</li> <li>– If necessary, confirm the <i>Please remove cell</i> message with <b>&lt;START·ENTER&gt;</b></li> </ul>
<ul style="list-style-type: none"> <li>– Instrument defective</li> </ul>	<ul style="list-style-type: none"> <li>– Please contact the service department.</li> </ul>

***Self test failed.***

Cause	Remedy
<ul style="list-style-type: none"> <li>– <i>System test</i>: Instrument defective</li> </ul>	<ul style="list-style-type: none"> <li>– Please contact the service department.</li> </ul>

Cause	Remedy
– <i>Filter test:</i> Instrument defective	– Please contact the service department.
– <i>Lamp test:</i> – If <u>no light</u> leaks out at the rear side of the photometer: Lamp defective – If light does leak out at the rear side of the photometer	– Replace the lamp  – Carry out a software update
– <i>Wavelength calibration:</i> – Foreign particle in the cell shaft – Lens smudged  – Instrument defective	– Remove foreign object – Clean the lens. If this happens repeatedly, check the operating conditions (see Section 8) – Please contact the service department.

## 5.2 Actions in the case of a broken cell

### WARNING

Cells can contain dangerous substances. If the contents are released, follow the safety instructions of the package insert. If necessary, take corresponding protective measures (protective goggles, protective gloves etc.).

### CAUTION

**Do not turn the photometer around to pour out the liquid!**  
When doing so, the liquid could come into contact with electronic components and damage the photometer.

The photometer has a drain device through which the contents of a broken cell can drain off without causing any damage.

### Proceeding after a cell has broken

- 1 Switch off the photometer and disconnect it from the power supply.
- 2 Let the liquid drain off into a suitable container and dispose of it properly according to the instructions of the reagent package.
- 3 Carefully remove all broken glass, e.g. with tweezers.

- 4 Carefully clean the cell shaft using a moist, lint-free cloth. If there are persistent coatings, use isopropanol for a short time. Especially clean the bottom parts of the lateral surfaces of the rectangular cell shaft where the light barriers for the automatic cell recognition are located.
- 5 Let the cell shaft dry.

After recommissioning, carry out the photometer monitoring for all measurements.

If, after recommissioning, an error occurs during the wavelength calibration, the detector lens is probably smudged. In this case, clean the detector lens.

## 6 Returns

**All instruments that are shipped back to the manufacturer must be free of measuring medium remnants and other hazardous materials. Remnants of the measuring medium on the housing may endanger humans or the environment.**

### **CAUTION**

**Keep the original packing including the inner packing to protect the instrument against hard shocks if it has to be transported.**

**The original packing is also required for the proper return of the instrument if it has to be repaired.**

**Note that damage caused by improper transport voids all warranty claims.**

## 7 Disposal

Surrender empty batteries at the appropriate collection points.

Do not dispose of the instrument in the trash bin. To dispose of the instrument, please ship it directly to us (with sufficient postage). We will dispose of the instrument in a proper and environmentally friendly manner.

## 8 Technical data

<b>Dimensions</b>	425 x 197 x 325 mm (width x height x depth)	
<b>Weight</b>	approx. 4.5 kg (without plug-in power supply)	
<b>Housing type of protection</b>	IP 30	
<b>Electrical protective class</b>	III	
<b>Test certificates</b>	CE	
<b>Allowed environmental conditions</b>	Temperature	Operation: +10 °C to + 35 °C (41 °F to 95 °F) Storage: -25 °C to + 65 °C (-13 °F to 268 °F)
	Humidity	Yearly mean: ≤ 75 % 30 days/year: 95 % Other days: 85 %
	Climatic class	2
	<b>Power supply</b>	Power pack
<b>Guidelines and norms used</b>	EMC	EC directive 2014/30/EC EN 61326-1:2013 – Interference emission: Class B – Interference immunity: IEC 61000-4-3 Tolerance extension: 0.008 E FCC Class A
	Meter safety	EC directive 2014/35/EC EN 61010-1:2010
	Climatic class	VDI/VDE 3540
<b>Communication interfaces</b>	IP protection class	EN 60529
	Ethernet	RJ45 socket
<b>Other features</b>	USB	– 1 x USB-A (for printer, USB memory devices, keyboard or bar code reader) – 1 x USB-B (for PC)
		<ul style="list-style-type: none"> <li>● Drain for spilled cell contents</li> <li>● Photometer software update and method data update possible via Internet</li> </ul>

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**Tintometer Indien Pvt. Ltd.**