

EVE™

Automatic cell counter




NanoEnTek

EVE™, Automatic cell counter

Instructions for using EVE™, Automatic cell counter described below.
For detailed instructions, refer to the manual supplied with EVE™ or
download the manual from www.nanoentek.com.


1. Load sample & insert slide

- ▶ Mix well the 10 μ L of sample and the 10 μ L of 0.4% trypan blue stain using a pipette.
- ▶ Load 10 μ L of the sample mixture on EVE™ Cell counting slide (side A) using a pipette.
- ▶ Insert the slide into the slide port of the instrument, sample side (side A) first.



2. Optimize image

- ▶ Press **Zoom** button and adjust the image using the Focus knob.
- ▶ Check the reference image for correct analysis. (next page)
- ▶ Use the grid to navigate the slide and check the image in several fields.
- ▶ Lock focusing using **Focus lock** knob if necessary



3. Count cells


- ▶ Press **Count cells** button to obtain the results.
- ▶ Record cell count.
- ▶ Press **More data** button to see data details as well as graphical representation of data.
- ▶ To read the second side of slide (side B), remove and turn around the slide, and repeat the 1 ~ 3 procedures with the side B.



The reference images for correct analysis


<Correct image>

Live cells have bright centers and dark edges.



<Incorrect image 1>

Dead cells have bright, blue centers and are counted as the live cells.



<Incorrect image 2>

Live cells have dark centers and are counted as the dead cells.



EVETM *Automatic cell counter*

Nano**EnTek**

NESMQ-EVE-001E (V.0.5)

website www.nanoentek.com

e-mail sales@nanoentek.com

NanoEnTek Inc. (HQ)

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 152-740, Korea
Tel : +82-2-6220-7940
Fax : +82-2-6220-7721

NanoEnTek USA Inc.

5627 Stoneridge Drive Suite 304, Pleasanton,
CA 94588, USA
Tel : +1-925-225-0108, +1-888-988-0108(Toll free)
Fax : +1-925-225-0109