

# SONOPULS

# **Ultrasonic Homogenizers**

Sonochemistry — Emulsifying — Homogenizing

Cell disruption — Suspending

Accelerating Reactions



### **Features**

### AMPLICHRON®-system

guarantees a constant amplitude independently from changing conditions within the sample. It ensures reproducible results for process validation. Settings within a range of 10 to 100 %. Verification of actual value at the display. Permanent control of ultrasound irradiation as well as indication of wear of the probe.

#### **Pulsation**

limits increase of temperature when processing sensitive samples. The adjustable pulsation allows cooling during rest intervals.

### **Continuous operation**

Constant sound radiation- extremely effective

### **Built-in timer**

Process duration storable. Indication of elapsed time in continuous operation or of remaining time in countdown mode.

#### Switching ON / OFF - easy to handle

either at the generator or directly at the ultrasonic converter via button or remote control.

### **Accessories**

A wide range of probes and special accessories for a vast variety of applications

### Foil keypad

easy to clean and user-friendly

### **ROHS** compliant

Devices are built lead free.

### Fail-safe during continuous operation and idling

RFI-proofed and CE-marked, also as medical device compliant to the directive for in-vitro diagnostics 98/79/EG

Features	mini20	HD 2000 series	HD 3000 series
Samples volumina	0,1 – 25 ml	1 – 1000 ml	1 – 2500 ml
Amplitude control	10 – 100 %	10 – 100 %	10 – 100 %
Power control	yes (HF power)	no	yes (HF power)
Automatic amplitude limiting	yes	no	yes, after preselection of probe
Pulsation	ON cycles 0,1-60 s OFF cycles 0,2-60 s	10-100 % - storable (duty cycle, base 1 sec)	ON cycles 0,2-600 s OFF cycles 0,3-600 s
Time modes	50 min: 59 s	99 min: 59 s continuous or timed	9 h: 59 min: 59 s continuous or timed
Safety shut down	50 min: 59 s	no	9 h: 59 min: 59 s
Display	grafic / alphanumeric liquid crystal display of amplitude, pulsation mode, time, energy	numerical seven-segment display of amplitude, pul- sation mode and time	grafic / alphanumeric liquid crystal display of amplitude, pulsation mode, time, energy and optionally temperature
Menu guided	comfortable setting of all values through "push & turn	no	comfortable setting of all values through "push & turn
Energy monitoring	in kJ	no	in kJ
Temperature monitoring and measurement	no	no	optional, 0–120 °C, temperature probe necessary, optional signal tone or switch - off
user programs	9	1	9, with software WINPULS®: 99
Remote control with PC	RS 232 (infrared)	no	RS 232 (infrared)
PC-Software, optionally available	no	no	WINPULS®
error diagnosis	yes	no	yes
Processing frequency	30 kHz	20 kHz	20 kHz
Automatic storage of the last adjusted values	yes	no	yes
Operating test	yes	no	yes
Remote control	no	foot switch	foot switch

### How to select the proper unit

Power output in watt is not the sole criterion for selecting an ultrasonic homogenizer. This value only indicates the power of the HF-generator but not the energy delivered to the sample. The amplitude at the radiating surface of the probe is the determining factor for the evaluation of the irradiation result while taking into consideration the volume of the sample.

### HF generator:

Transforming of low-frequency voltage of 50 Hz into high-frequency voltage of 20 kHz.

### **Ultrasonic converter:**

Transforming of electrical energy delivered from the generator into mechanical vibrations of 20 kHz.

### Standard and booster horns:

Horns increase the amplitude by their specially designed shape. The external thread is made for close connection of vessels.

### **Probes:**

Transmitting of ultrasonic energy into the sample. Microtips, tapered and flat tips, dia. 2, 3, 6, 13, 19 and 25 mm, for use in different volumes.



# **SONOPULS mini20**

### **SONOPULS mini20**

for volumes up to 25 ml

### Fast hand operation

Ready-to-operate for volumes from 0,5 ml to 25 ml, consisting of HF generator GM mini20, ultrasonic converter UW mini20 and microtip MS 2.5, diameter 2,5 mm.

Max. 12 Weff HF output.

Code No. 3665






Fast manual use – ideal for smallest volumes. Pulsation on pressing the bottom at the ultrasonic converter by thumb.



HF generator		GM mini20
dimensions	mm	250 × 256 × 154
weight	kg	2,0
mains supply		230 V~, 50/60 Hz (optionally 115 V~, 50/60 Hz)
converter		UW mini20
dimensions	mm	dia. ca. 50 × 160
weight	g	270
available titanium probes	dia. mm	1,5 or 2,5

# **SONOPULS HD 2000 series**

### Operating panel HD 2070 / HD 2200



### **SONOPULS HD 2070**

for volumes uo to 200 ml

### Small unit for lab routine

Ready-to-operate basic equipment for volumes from 2 ml to 50 ml consisting of generator GM 2070, ultrasonic converter UW 2070, stepped standard horn SH 70 G and microtip MS 73, 3 mm diameter. HF-output max. 70  $W_{\mbox{\footnotesize eff}}$ 

Code No. 2450



### **SONOPULS HD 2200**

for volumes uo to 1000 ml

### Standard unit for lab routine

Ready-to-operate basic equipment for volumes from 20 ml to 900 ml consisting of generator GM 2200, ultrasonic converter UW 2200, booster horn SH 213 G and titanium flat tip TT 13 of 13 mm diameter. HF-output max. 200  $W_{\mbox{\footnotesize eff}}$ 


Code No. 2530



HF generator		GM 2070	GM 2200
dimensions	mm	257 x 180 x 115	257 x 180 x 115
weight	kg	2,5	2,5
mains supply		230 V~, 50/60 Hz, optionally with	n voltage selector for 115 V~, 50/60 Hz
converter		UW 2070	UW 2200
dimensions	mm	Ø 70 x 120	dia. 70 x 120
weight	kg	1,0	1,0
available titanium probes	dia. mm	2, 3, 6, 13	2, 3, 6, 13, 19 or 25

# **SONOPULS HD 3000 series**

### Operating panel HD 3100 / HD 3200 / HD 3400



### **SONOPULS HD 3100**

for volumes uo to 200 ml

### **High-Tech for research**

Ready-to-operate for volumes from 2 ml to 50 ml, consisting of HF generator GM 3100, ultrasonic converter UW 3100, standard horn SH 70 G and microtip MS 73 diameter 3 mm. Max. 100 Weff HF output.

Code No. 3680

### **SONOPULS HD 3200**

for volumes uo to 1000 ml

### **High-Tech for reserch**

Ready-to-operate for volumes from 20 ml to 900 ml, consisting of HF generator GM 3200, ultrasonic converter UW 3200, booster horn SH 213 G and flat tip TT 13, diameter 13 mm. Max. 200  $W_{eff}$  HF output.

Code No. 3660

HF generator		GM 3100	GM 3200	
dimensions	mm	250 x 256 x 154	250 x 256 x 170	
weight	kg	2,0	2,7	
mains supply		230 V~, 50/60 Hz, optionally 115 V~, 50/60 Hz		
converter		UW 3100	UW 3200	
dimensions	mm	dia. 70 x 120	dia. 70 x 120	
weight	kg	1,0	1,0	
available titanium probes	dia. mm	2, 3, 6 or 13	2, 3, 6, 13, 19 or 25	



### **SONOPULS HD 3400**

for volumes uo to 2500 ml

### High-Tech for reserch and pilot plant stations

Ready-to-operate for volumes from 100 ml to 2500 ml, consisting of HF generator GM 3400, ultrasonic converter UW 3400, booster horn SH 3425 and extended probe VS 200 T, diameter 25 mm. Max. 400  $W_{\mbox{eff}}$  HF output.

Code No. 3690

HF generator		GM 3400
dimensions	mm	324 x 230 x 131
weight	kg	3,1
mains supply		230 V~, 50/60 Hz
converter		UW 3400
dimensions	mm	dia. 90 × 180
weight	kg	2,2
available titanium probes	dia. mm	19 or 25



## **Applications**

Ultrasonic homogenizers are used in laboratories, hospitals and in industry for scientific experiments and analysis as well as in pilot or small lot production. Here are some examples showing the vast variety of applications for ultrasonic homogenizers:

### Typical areas of application

- Disruption of cells, bacteria, virus, tissue, also mixed tissue e. q. for extraction of cell contents
- Emulsifying of hardly mixable liquids, e.g. oil and water, particle size in um range
- Deagglomeration of nanoparticles in material research (nanostructurized material) in medicine, biotechnology, automobile industry
- Acceleration of chemical reactions
- Production of dispersions



### **Analysis**

- Preparing samples for grain size determination or environmental analysis: HD 3200 or HD 2200 with tapered tip KE 76 or with extended probe VS 70 T.
- Homogenizing of cheese samples for determination of nitrates: HD 3200 or HD 2200 with MS 73

### **Biochemistry - Biology - Medicine**

- Sonication of small high-quality samples for analysis like EIA or RIA: HD 3100 and HD 2070 with microtip MS 72 or MS 73.
- Due to high amplitudes, disruption of high-resistant bacteria, cells or tissues is possible. Indirect processing of sample in cup booster BR 30 or in cup horns BB 2 G or BB 6 is recommended to avoid cross-contamination.
- Detection of prions by cyclic amplification of protein misfolding: HD 2070 with MS 73
- Simultaneous sonication of 12 samples in microplates:
- HD 3100 with MR 12-2

### **Chemistry and Sonochemistry**

Acceleration of chemical reactions or destroying of highly-molecular compounds:

HD 3200 or HD 2200 with tapered tip KE 76 and sleeve adapters NA 29 G or NA 45 G for tight fitting to a sonochemical reaction vessel.

### **Pharmacy - Cosmetic**

Production of larger volumes of long lasting emulsions, e. g. lotions and production of antigens, vaccines or liposomes:

HD 3200 or HD 2200 with flow-through cell DG 4 G





# **Applications**

Waste water samples

Aim: Homogenizing for determination of harmful substances, e.g. mineral oil, grease AOX in industrial and butcher's waste water

Quantity: 250 ml Approx.time: 5 - 10 min Unit: HD 2200/3200 with TT 13, or taller vessels with VS 70 T

Aluminium oxide suspensions

Aim: Dispersing Quantity: 100 ml Approx.time: ca. 4 min Unit: HD 3200 with KE 76

Soil samples

Aim: Extraction for determination of pH value, Mg, K, P, N – contents for recommendation of fertilizer / determination of radio nucleides to control

radioactivity in the environment (milk research)

Quantity: 50 - 100 ml / 100 - 150 ml Approx.time: a few seconds

Unit: HD 2200/3200 with KE 76 / VS 70 T

Bladder tissue Aim: Disruption Quantity: 1,5 ml

Approx.time: ca. 1,5 min

Unit: HD 2200 with MS 72, cooling necessary

Candida albicans
Aim: Disruption
Quantity: 10 ml
Approx.time: ca. 10 min
Unit: HD 2070/3100 with MS 73

ChIP (Chromatin immunoprecipitation)

Aim: DNA fragmentation Quantity: 1 ml Approx.time: ca. 2 min Unit: HD 3200 with MS 72

Large intestine tissue Aim: Disruption

Quantity: 1,5 ml Approx.time: ca. 3 min

Unit: HD 2200 with MS 72, cooling necessary

**Dispersing of solid particles Aim:** grain size analysis

Quantity: 50 - 100 ml Approx.time: ca. 2 - 5 min Unit: HD 2200/3200 with KE 76

Small intestine tissue

Aim: Disruption Quantity: 1,5 ml Approx.time: ca. 2,5 min

Unit: HD 2200 with MS 72, cooling necessary

Escherichia coli

Aim: Disruption for proteine lay off

Quantity: 10 ml

Approx.time: ca. 5 - 10 min

Unit: HD 2070/3100 with MS 73 or HD 2200 with MS 73

**Eucaryotic cells** 

Aim: Disruption for proteine lay off

Quantity: 1,5 ml Approx.time: ca. 1 min

Unit: HD 2200/3200 with BR 30 u. EH 3, cooling necessary

Meat and sausage samples

**Aim:** Homogenizing for determination of nitrates

Quantity: 100 ml Approx.time: ca. 3 min Unit: HD 2200/3200 with KE 76

Heart muscle tissue Aim: Homogenizing Quantity: 1,5 ml

Approx.time: ca. 4 min Unit: HD 2200/3200 with MS 72, cooling necessary

**Brain tissue** 

Aim: Disruption Quantity: 1,5 ml Approx.time: ca. 1 min

Unit: HD 2200 with MS 72, cooling necessary

Yeast cells

Aim: Disruption Quantity: 10 ml Approx.time: ca. 2 min Unit: HD 3200 with MS 73

Insect cells

Aim: Disruption for proteine lay off

Quantity: 20 - 50 ml

Approx.time: ca. 25 sec, pulsed

Unit: HD 2070/3100 with MS 73 and RZ 2

Liver tissue

Aim: Homogenizing for moleculargenetic tests

Quantity: 1,5 ml

Approx.time: ca. 11/2 min

Unit: HD 2200/3200 with MS 72, cooling necessary

Liposomes

Aim: Producing of small unilamellar phospholipid vesicles

Quantity: 20 ml

**Approx.time:** ca. 10 - 15 min **Unit:** HD 2070/3100 with TT 13,

cooling necessary

Lymphocytes

Aim: Disruption Quantity: 50 μl - 2 ml Approx.time: ca. 1 - 5 min

Unit: HD 2070/3100 with BR 30 and EH 3

Nano emulsions

Aim: Drop sizes within nm range

Quantity: 2 ml Approx.time: ca. 5 min Unit: HD 3100 with MS 72, cooling necessary

Nano particles

Aim: Dispersing Quantity: 100 ml Approx.time: ca. 2 min Unit: HD 3200 with KE 76

Retina

Aim: Tissue disruption Quantity: 60 ml

**Approx.time:** 15 short stokes **Unit:** mini20 with MS 2.5

Kidney tissue

Aim: Homogenizing Quantity: 1,5 ml Approx.time: ca. 40 sec Unit: HD 2200/3200 with MS 72,

cooling necessary

O/W emulsions

Aim: Finest emulsions Quantity: 10 ml Approx.time: ca. 1 min Unit: HD 3200 with KE 76,

Vessel: rosett cell

Homogenizing of aqueous ink

Aim: Dispersing of ink pigments in oil


Quantity: 200 ml Approx.time: ca. 5 min Unit: HD 2200 with VS 70 T

Carbon black dispersions

Aim: Homogenizing Quantity: 50 ml Approx.time: ca. 5 min

Unit: HD 2200 with DH 13 G, vessel: KG 3

# **Accessories and Applications**



### Saccharomyces cerevisiae

Aim: Disruption Quantity: 20 ml Approx time: 15 min

Unit: HD 2200/3200 with KE 76, addition of glass beads to


accelerate process, cooling necessary

### Staphylococcus aureus

Aim: Disruption Quantity: 10 ml

Approx time: ca. 10 min Unit: HD 2070/3100 with MS 73

# **Accessories and Applications**



### Streptococcus

Aim: Disruption Quantity: 10 ml

**Approx time:** ca. 8 - 10 min **Unit:** HD 2200/3200 with MS 73

### Detection of prions by cyclic amplification of protein misfolding

Quantity: 200 µl Approx time: ca. 1,5 min Unit: HD 2070 with MS 73

# Probes • Standard- / Booster Horns • Adapters



MS72 MS73 KE76 VS70T VS190 T VS 200 T TT13 TT19 TT25

### **Probes**

made of titanium alloy (Ti-Al6-V4) transmit mechanical longitudinal waves into the sample. They are thermo-resistant, can be treated in autoclaves and are resistant to corrosive media. Sample volume, diameter of the processing vessel and the required amplitude determine the selection of the unit and the type of probe. The higher the amplitude, the more intense the sonication.

Please note that probes are subject to wear and tear. It is advisable to order spare probes with the first homogenizer.

description			mie	crotips		tapered tip	ext	ended prol	эе
Туре		MS 1.5	MS 2.5	MS 72	MS 73	KE 76	VS 70 T	VS 190 T	VS 200 T
Code No.		3639	3652	492	529	530	494	3638	478
Diameter	mm	1,5	2,5	2	3	6	13	19	25
Length approx	mm	57	53	191	175	135	130	130	130
Standard horn for HD 2070/3100 Booster horn for HD 2200/3200 Booster horn for HD 3400		- - -	- - -	SH 70 G SH 213 G	SH 70 G SH 213 G -	SH 70 G SH 213 G -	SH 70 G SH 213 G -	- SH 219 G SH 3419	- SH 225 G SH 3425
Amplitude for HD 2070/3100 Amplitude for HD 2200/3200 Amplitude for HD 3400 Amplitude for mini20	μm <sub>SS</sub> (peak to peak)	- - - 50	- - - 70	253 / 285 282 / 286 - -	212 / 245 302 / 308 - -	165 / 191 249 / 255 - -	80 / 97 153 / 170 - -	- - 116 -	- / - 46 / 51 82 -
Volume HD 2070/3100	ml	-	-	1-25	2-50	5-100	10-200	-	-
Volume HD 2200/3200	ml	-	-	2-30	5-90	10-350	20-900	25-900	30-1000
Volume HD 3400	ml	-	-	-	-	-	-	500-1500	500-2500
Volume mini20	ml	0,1-10	0,5-25	-	-	-	-	-	-
Vessel diameter min	mm	4	6	4	6	8	17	23	29

Probe length may vary slightly due to the variations in the titanium material.

### Standard- und booster horns

(Ti-Al6-V4) are furnished with a thread for replaceable probes. With exterior (except SH 3419, SH 3425) thread to connect various vessels.

**solid standard horn** - DH 13 G - with diamond coating on the radiating area; lifetime is thirty times longer than usual.



### Flow-through standard and booster horns

material: Ti-6Al-4V, to prepare stable mixtures of non-mixable or hardly mixable liquids (oil-in-water) by direct intrusion of pre-mixed samples into the cavitation field. In combination with flow-through cell DG 4 G, the continuous treatment of 2 different media and parallel tempering is possible.



	standard horn	booster horn					diamond standard horn	flow-through standard horn	flow-through booster horn
Туре	SH 70 G	SH 213 G	SH 219 G	SH 225 G	SH 3419	SH 3425	DH 13 G	FZ 5 G	FZ 7 G
for HD	2070 / 3100	2200 / 3200		34	.00	2070 / 2200 3100 / 3200	2070 / 3100	2200 / 3200	
Code No.	486	527	3647	3634	3679	3692	403	490	452

### **Adapters**

**Sleeve adapters** made of PTFE for tight mounting on standard ground glass vessels.

NA 29 G for NS 29/32 for SH 70/213 G

NA 45 G for NS 45/40 for SH 70/213/219/225G

**Threaded adapter** made of stainless steel with external thread M 40 x 1

**GA 3 G** for SH 70/213/219/225 G



NA 29 G NA 45 G



	sle adap		threaded adapters
Type	NA 29 G	NA 45 G	GA 3 G
for HD	2070	/ 2200 / 31	100 / 3200
Code No.	540	487	473

### **Probe extensions**

for enlarging the operating depth when using flat tips.

**VS 70** between SH 70 G / 213 G and TT 13 **VS 200** between SH 225 G and TT 25



	probe extensions				
Type	VS 70	VS 200			
for HD	2070 / 2200 3100 / 3200	2200 / 3200			
Code No.	500	415			

	titanium flat tips				silica glass	probes			
TT 13	TT 19	TT 25	GS 6	GS 6 L	GS 13	GS 13 L	GS 18	GS 18 L	
497	491	532	024	048	028	050	040	054	
13	19	25	(	6	1:	3	18		
5	5	6	145	290	145	290	145	290	
SH 70 G SH 213 G -	- SH 219 G -	- SH 225 G -	SH 7	SH 70 GQ - -		SH 70 GQ - -		SH 70 GQ - -	
78 / 93 149 / 165 - -	- / - 73 / 81 - -	- / - 48 / 53 - -	- 1	13 / 13 - / - - / - - / -		13 / 13 - / - - / - - / -		13 - - -	
10-200	-	-	2-	100	25-200		25-500		
5-900	25-900	30-1000			-				
-	-	-			-				
-	-	-			-				
17	23	29	1	10 17		22			

### Silica glass probes

for connection to HD 2070/3100 with special horn SH 70 GQ.

For application in food analysis, pharmacy or environmental analysis. No intrusion of metal particles and boron compounds - ideal for trace analysis. High chemical and temperature shock resistance, no electric conductivity.





### **MULTISON®** ultrasonic probe

patent applied D 10 2004 024 214

for connection to HD 2070/3100. Composed of Multison horn MRH 12 and 12 Multison tips MRS 2, MRS 3 or MRS 2-2C . For irradiation of samples in microplates and deep well plates.

Simultaneous sonication of 12 samples. Multison tips individually replaceable.



	Multison probe c	omposed of multer 1 multison tip	tison horn with		multison tips	
Туре	MR 12-2	MR 12-2C	MR 12-3	MRS 2	MRS 3	MRS-2C
Diameter, mm	2	2	3	2	3	2
Length, mm					16	
Code No.	3626	3643	3633	3628	3629	3642

# **Direct processing**

### Processing vessel, made of stainless steel

**DG 4 G** for high-volume flow-through processing like emulsifying, dispersing or homogenizing, up to 30 l/h, The sample can be repeatedly sonicated in circulation. For connection to SH 70 G or SH 213 G with TT 13, DH 13 G.

**KG 4 G**, closed vessel with cooling jacket. Processing volume about 65 ml.

	flow-through processing vessel	cooling vessel
Type	DG 4 G	KG 4 G
for HD	2070 / 2200 3100 / 3200	2070 / 2200
Code No.	3608	3608



### Processing vessels made of glass

### Cooling vessel KG

zfor sonication of temperature-sensitve samples. The cooling jacket allows circulation of cooling liquid during sonication.

### Flow-through vessel DG

with cooling jacket for irradation of larger volumes. The cooling jacket allows circulation of cooling liquid during sonication.

### Rosett cell RZ

for homogenous and intense circulation of liquids caused by the shape of the vessel and its 3 sidearms

### Suslick cell SZ

hwith 3 sidearms for introduction of gas or measuring probes.










	cooling v	/essel	flow-throu	gh vessel	rosett cell			suslick cell			
Туре	KG 3	KG 5	DG 3	DG 5	RZ 1	RZ 2	RZ 3	RZ 4	RZ 5	SZ 3	SZ 5
for HD	2070 / 2200 3100 / 3200	2200 3200	2070 / 2200 3100 / 3200	2070 / 2200	2070 / 2	2200 / 3100	0 / 3200	2200 / 3	3200 3400	2070 3100	2200 3200
volume, ml	15	70	max. 5,6 I/H	max. 30 l/h	25	40	110	390	660	20	110
interior diameter, mm	20	35	20	53	30	42	50	75	90	20	40
height, mm	65	95	65	95	85	100	135	202	243	80	144
Code No.	536	481	538	482	3606	3607	522	3256	483	534	484

# Indirect processing

- High-intensity mini ultrasonic baths
- Processing of µl-samples without sample loss
- No cross contamination
- No contamination of samples through cavitation erosion at the probe
- Uniform sonication of several samples with simultaneous cooling
- No aerosoling when sonicating pathogenic or hazardous materials
- Application for e. g. cell disruption, preparing of liposomes







#### Microtube holder EH 6

For use in BB 6. Up to 6 samples can be treated simultaneously. Pressure plate holds tubes in place. No floating of cups.

A mixing of samples is excluded due to markings at the holder.

#### Cup horn BB 6

for indirect intense sonication. The cup horn is equipped with inlet and outlet for circulation of cooling liquid. Also useable for direct sonication.

### Cup horn BB 2 G

plastics, for indirect sonication of pathogenic material.

### Microtube holder EH 3

for use with BR 30. Up to 3 samples can be treated simultaneously. 2 exchangable discs with diameters 8,5 or 11,5 mm.

#### Inset basket BK 30

For intensive cleaning of small parts, e. g. cleaning of radioactively contaminated seeds in BR 30.

### Cup booster BR 30

For high-intensive irradiation of smallest and sensitive sample volumes, e. g. radioactive seeds or bacteria as well as for flow-through sonication of liquids like cell suspensions. During indirect processing of samples the ultrasound will be transferred by the contact liquid. The sonication will be carried out in reaction vessels or in the inset beaker BK 30. The cup booster is equipped with inlet, outlet and overflow.

A cooling as well as a flow-through process is possible.

Туре	BB 6	BB 2 G	EH 6	BR 30	BK 30	EH 3
for HD	2200 / 3200	2070 / 3100	2200 / 3200	2070 / 2200 3100 / 3200	BR 30	BR 30
Code No.	3605	552	059	082	098	078

# Stand • Sound proof boxes



Туре	HG 5	HG 10	KL 7	AT 7
für HD	2070 / 2200 3100 / 3200	3400	HG 5 HG 10	KL 7 LS 8
Code No.	459	3646	3636	3644

### **Stand**

#### Stainless steel stand

with lab clamp and non-slip mat to hold processing vessels securely in place

Clamping device KL 7 (DE 20 2006 005 654.98) for HG 5 / HG 10 with rod, swivelling clamp for reaction vessels dia. 15 mm to dia. 100 mm

**Suporting table AT 7** suitable for KL 7 or in LS 8 with non-slip mat to hold sample vessels securely in place





LS 11

### Sound proof boxes

reduce the noise level considerably. Precut holes at the backside allow connections for gas supply and flow-through processing. Acrylic door permits process monitoring.

#### 194

Plastics coated walls, 10 dB-AU damping.

#### 1 5 8

made of stainless steel, with damping material. 20 dB-AU damping. The damping material is water resistant for easy cleaning.

With rod, swivelling clamp and quick clamp for height adjustment of sample vessels.

Clamping belt for safe fixing of sample vessels with different sizes.

Also applicable for sonication of samples in glass vessels with round bottoms or with inlets from below.

Special support UG 6 is available for inverted position of the box during indirect sonication with horn BB 6 or cup booster BR 30. Ultrasonic converter is fixed safely through a special clamp.

### **LS 11**


stainless steel with damping material, 20 dB-AU damping, The damping material is water resistant for easy cleaning

**Distance tube** for direct processing with long probes

**AH 6**: for MS 72/73, KE 76, VS 70 with TT 13, VS 200 mit TT 25 / VS 200 T, VS 70 T, GS ...

BD 8: damping material for sound proof boxes

Туре	LS 4	LS 8	LS 11	UG 6	AH 6	BD 8
for HD	2070 / 2200 3100 / 3200	2070 / 2200 3100 / 3200	3400	2070 / 2200 3100 / 3200	LS 7	LS 8 / LS 11
dB-AU damping	10	20	20	-	-	-
Code No.	416	3653	3663	3616	3619	3661




### WINPULS® remote control

For process control with PC for operation systems MICROSOFT® WINDOWS® 2000 and MICROSOFT® WINDOWS® XP.

With different additional functions like test logging and comfortable data storage (up to 99 storages).

Set composed of WINPULS® software and infrared adapter IR 1 for interface RS 232

Туре	WINPULS® software with ifraded adapter IR 1			
for HD	3100 / 3200 / 3400			
Code No.	3625			



WINPULS® up PC

### Foot switch remote control

for easy switching ON/OFF of the HF generator. With 3 m cable.

Туре	TS 8				
for HD	2070 / 2200 3100 / 3200 / 3400				
Code No.	531				



### **Temperature sensor**

for measuring the sample temperature from 0 up to 120  $^{\circ}\text{C}$  Diameter sensor: 4 mm

Туре	TM 100				
for HD	3100 / 3200 / 3400				
Code No.	3622				



# **Ultrasound in Laboratory**

# SONOREX TECHNIK Vortex reactor

### **Applications**

- Intensifying of industrial, biotechnological and chemical processes
- Degassing
- Support of disinfection (bacterial elimination)
- Disinfection of liquids
- Producing of finest polishing pastes for wafer industry
- Homogenizing

Vortex reactor consisting of: Vortex reactorbloc WB and HF generator LG 2002 T



# Tube reactor SONOBLOC® Applications

- Ultrasonic intensive treatment of fibrous and bandshaped products
- Support of industrial and biotechnological processes
- Wire cleaning
- Degassing
- Support of disinfection (bacterial elimination)
- Acceleration of disintegration

SONOREX TECHNIK

Dispersing of solid particles in liquids

### Tube reactor SONOBLIOC®

consisting of: Tube reactorbloc RB and HF generator LG 1001 T



SB 8-1002.01

Technical Data	Vort	ex reactorbloc - WB	Tube reactorbloc - RB		
Туре	WB 4-1402	WB 4-1503	WB 4-1604	RB 8-1002	RB 8-1004
Flow-through rate		1 - 50 l/min	1 - 100 l/min		
Internal pressure, max.		10 bar		10 bar	
Solid particles	< 5 mm			-	
Power density, max.	480 W/I 520 W/I 550 W/I 500 W/I			N/I	
Power, max.	1400 W 1500 W 1600 W 1000 W		W		
Frequency	25 kHz	25 und 40 kHz	40 kHz	25 kHz 40 kHz	
Tube material / dimensions		ISI 316 Ti / dia. 139.7 × dia. 104 × 2 mm	Stainless steel AISI 316 Ti / dia. 60.3 × 3.6 mm		
Housing dimensions (I × w × h)	29	0 × 290 × 642 mm	260 × 150 × 990 mm		
Weight, net	approx. 50 kg			approx.35 kg	
HF-Generator (separat)	LG 1510 T LG 2002 T			LG 1001 T	

Units are equipped with standard victaulic connection. Further connection variants on request.

Separate documentation on request.

### SONOREX ultrasonic baths

- standard baths from 1 200 litres
- compact construction
- easy to handle
- optionally with heating
- integrated timer
- from 3 ltr. upwards with drain



- cost-saving reduced cleaning time
- intensive, gentle, thorough
- ⇒ large range of accessories

Further documentation on request.



### BANDELIN electronic

specialized in manufacturing of ultrasonic units, maintain a quality management system complying with the requirements of EN ISO 9001 / 12.2000 and EN ISO 13485:2003 + AC:2007

5732 e/2009-04

All units are CE marked.
Illustrations exemplarily, not to scale

Subject to technical alterations without notice. Decoration products are not included in delivery.

The general delivery terms apply.



www.bandelin.com www.sonopuls.info info@bandelin.com 60 years of experience in ultrasound technology