

Cole-Parmer® Narrow-Mouth Bottles,

HDPE, LDPE, PP

Keep samples contained with leakproof design

- Leakproof PP closure for safer transport
- HDPE bottles are FDA-compliant
- LDPE bottles feature transparent, food-grade design
- PP bottles offer superior chemical resistance and are autoclavable



“I will continue to buy this product.”



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk	Qty/ cs	Catalog number	Price/ cs
HDPE							
30 (1)	22 mm	12	GY-06038-61		72	GY-06038-62	
60 (2)	22 mm	12	GY-06038-63		72	GY-06038-64	
125 (4)	28 mm	12	GY-06038-65		72	GY-06038-66	
250 (8)	28 mm	12	GY-06038-67		72	GY-06038-68	
500 (16)	28 mm	12	GY-06038-69		48	GY-06038-70	
1000 (32)	28 mm	6	GY-06038-71		24	GY-06038-72	
2000 (64)	45 mm	1	GY-06038-73		6	GY-06038-74	
LDPE							
60 (2)	22 mm	12	GY-06038-21		72	GY-06038-22	
125 (4)	28 mm	12	GY-06038-23		72	GY-06038-24	
250 (8)	28 mm	12	GY-06038-25		72	GY-06038-26	
500 (16)	28 mm	12	GY-06038-27		48	GY-06038-28	
1000 (32)	28 mm	6	GY-06038-29		24	GY-06038-30	
PP Autoclavable							
30 (1)	22 mm	12	GY-06038-41		72	GY-06038-42	
60 (2)	22 mm	12	GY-06038-43		72	GY-06038-44	
125 (4)	28 mm	12	GY-06038-45		72	GY-06038-46	
250 (8)	28 mm	12	GY-06038-47		72	GY-06038-48	
500 (16)	28 mm	12	GY-06038-49		48	GY-06038-50	
1000 (32)	28 mm	6	GY-06038-51		24	GY-06038-52	
2000 (64)	45 mm	1	GY-06038-53		6	GY-06038-54	

1396 Series Square Bottles

Maximize storage space

- Square design provides easier storage with no wasted space between bottles
- Drip-free pouring
- Autoclavable bottle and PP cap



Capacity mL (oz)	Thread size	Qty/ cs	Catalog number	Price/ cs
100 (3)	GL32	10	GY-34594-22	
250 (8)	GL45	10	GY-34594-24	
500 (16)	GL45	10	GY-34594-26	
1000 (32)	GL45	10	GY-34594-28	

Narrow-Mouth Bottles, LDPE

Won't break or leak

- See samples more easily with translucent material
- Very low metal content, ideal for trace metal analysis
- Include leakproof, PP caps



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
8 (1/4)	20 mm	12	GY-06033-03	
15 (1/2)			GY-06033-05	
30 (1)			GY-06033-10	
60 (2)			GY-06033-20	
125 (4)	24 mm	12	GY-06033-30	
250 (8)	24 mm		GY-06033-50	
500 (16)	28 mm		GY-06033-60	
500 (16)	38-430		GY-06033-80	
1000 (32)	38-430	6	GY-06033-70	
Replacement caps, PP				
—	24 mm	12	GY-63000-04	
—	28 mm		GY-63000-06	
—	38-430		GY-63000-18	

Cole-Parmer® Round Jugs, HDPE

Easier carrying with molded-in handle

- Lightweight design won't add unnecessary weight
- Leak-free polypropylene (PP) cap



NEW



Capacity L (gal.)	Qty/ pk	Catalog number	Price/ pk
2 (1/2)	24	GY-06089-51	
4 (1)	12	GY-06089-53	

DON'T RISK IT!

Our chemical compatibility database provides important information on chemical resistance—check it out at ColeParmer.com/ChemChart

