

SAFETY DATA SHEET

Issuing Date 24-Apr-2014 Revision Date 06-Nov-2015 Revision Number 1

Section 1. Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product Name N26, N44, N75, N140, N250, N415, N750, N1000, N1000(KU), N1400, C200, G200,

S200(KU), S200, C350, G350, N350, N350(KU), C600, G600, N600, S600(KU), S600,

SV350, SV1400, VP190, VP960

1.2. Relevant identified uses of the substance or mixture and uses advised against

Recommended Use Viscometer and/or density measurement equipment calibration and performance

verification reference standard

Uses advised against No information available

1.3. Details of the supplier of the safety data sheet

Supplier

Cannon Instrument Company

2139 High Tech Rd.

State College, PA 16803-1733

TEL: (814) 353-8000; (800) 676-6232

For further information, please contact

E-mail Address No information available.

1.4. Emergency telephone number

Emergency Telephone (800) 255-3924 Domestic CHEM-TEL Inc.

Number +1 (813) 248-0585 Overseas CHEM-TEL Inc. (Please Call Collect)

Europe 112

Section 2. Hazards identification

2.1. - Classification of the substance or mixture

REGULATION (EC) No 1272/2008

Not classified

Physical Hazards

None

2.2. Label Elements

Not classified

2.3. Other information

May be irritating to eyes, respiratory system and skin.

Section 3. Composition/information on ingredients

3.1. Substances

3.2. Mixtures

Chemical Name	EC-No	CAS-No	Weight %	EU - GHS Substance	REACH No.
				Classification	

1-Decene, homopolymer,	-	68037-01-4	100	No data available
hydrogenated				

For the full text of the H-Statements mentioned in this Section, see Section 16

Section 4. First aid measures

4.1. Description of first-aid measures

Eye Contact Rinse thoroughly with plenty of water, also under the eyelids. Get medical attention if

irritation persists.

Skin Contact Wash off with warm water and soap. Get medical attention if symptoms occur.

Ingestion Do NOT induce vomiting. Drink plenty of water. Never give anything by mouth to an

unconscious person. Get medical attention.

Inhalation IF INHALED: If breathing is difficult, remove to fresh air and keep at rest in a position

comfortable for breathing. Get medical attention.

4.2. Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects Irritation.

4.3. Indication of immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

Section 5. Fire-fighting measures

5.1. Extinguishing media

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Extinguishing media which must not be used for safety reasons

No information available.

5.2. Special hazards arising from the substance or mixture

Special exposure hazards arising from the substance or preparation itself, combustion products, resulting gases None in particular.

5.3. Advice for firefighters

Special protective equipment for fire-fighters

As in any fire, wear self-contained breathing apparatus and full protective gear.

Section 6. Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing. Refer to Section 8 for personal protective equipment.

6.2. Environmental precautions

Prevent further leakage or spillage if safe to do so.

6.3. Methods and materials for containment and cleaning up

Prevent further leakage or spillage if safe to do so.

Dam up. Contain and collect spillage with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local/national regulations (see Section 13). Clean contaminated surface thoroughly.

6.4. Reference to other sections

See Section 12 for additional information.

Section 7. Handling and storage

7.1. Precautions for Safe Handling

Handling

Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes and clothing. Remove and wash contaminated clothing before re-use. Wear personal protective equipment.

Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice.

7.2. Conditions for safe storage, including any incompatibilities

Keep containers tightly closed in a dry, cool and well-ventilated place.

7.3. Specific end use(s)

Exposure Scenario

No information available.

Other Guidelines

No information available.

Section 8. Exposure controls/personal protection

8.1. Control parameters

This product, as supplied, does not contain any hazardous materials with occupational **Exposure Limits**

exposure limits established by the region specific regulatory bodies.

Derived No Effect Level Predicted No Effect Concentration (PNEC)

No information available. No information available.

8.2. Exposure controls

Engineering Measures Personal protective equipment

Eye Protection

Ensure adequate ventilation, especially in confined areas.

Safety glasses with side-shields.

Skin and Body Protection Hand ProtectionLong sleeved clothing.
Protective gloves.

Respiratory Protection When workers are facing concentrations above the exposure limit they must use

appropriate certified respirators.

Environmental Exposure Controls No information available.

Section 9. Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical State Liquid Appearance Colorless

Odor Faint hydrocarbon

<u>Property</u> <u>Values</u> <u>Remarks/ - Method</u>

рΗ No data available None known **Melting Point/Range** No data available None known >346 °C **Boiling Point/Boiling Range** None known > 224 °C Flash Point Closed cup **Evaporation rate** No data available None known Flammability (solid, gas) No data available None known

Vapor Pressure No data available None known Vapor Density No data available None known **Relative Density** at 15 °C 0.82-0.85 Water Solubility Insoluble in water. None known Solubility in other solvents Soluble in solvents. None known Partition coefficient: n-octanol/waterNo data available None known **Autoignition Temperature** No data available None known **Decomposition Temperature** No data available None known 25-2500 cSt @ 40°C **Viscosity** None known

Explosive Properties No data available Oxidizing Properties No data available

9.2. Other information

VOC Content (%)

Flammability Limits in Air

No information available

No data available

Section 10. Stability and reactivity

10.1. Reactivity

Not reactive under normal conditions.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

None under normal processing.

10.4. Conditions to avoid

Keep away from open flames, hot surfaces and sources of ignition.

10.5. Incompatible materials

Strong oxidizing agents.

10.6. Hazardous decomposition products

Carbon oxides.

Section 11. Toxicological information

11.1. Information on toxicological effects

Acute Toxicity

Product Information

Inhalation Inhalation of mist may cause irritation to the respiratory system.

Eye Contact Contact with eyes may cause irritation.

Skin Contact May cause skin irritation and/or dermatitis. Repeated or prolonged contact may cause

localized dermal effects including contact dermatitis, dry skin, or rash.

Ingestion Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea.

No information available. Sensitization No information available. **Mutagenic Effects Carcinogenic Effects** No information available. No information available. **Reproductive Toxicity Developmental Toxicity** No information available. STOT - single exposure No information available. STOT - repeated exposure No information available. **Aspiration Hazard** No information available.

Section 12. Ecological information

12.1. Toxicity

Ecotoxicity Effects

Contains no substances known to be hazardous to the environment at concentrations that would be significant.

Chemical Name	Toxicity to Algae	Toxicity to Fish	Toxicity to Microorganisms	Daphnia Magna (Water Flea)
1-Decene, homopolymer, hydrogenated		LC50 Rainbow trout: >1000 mg/L (96-hr)		EC50: 190 mg/l (48-hr)

12.2. Persistence and degradability

No information available.

12.3. Bioaccumulative potential

No information available.

12.4. Mobility in soil

Adsorbs on soil.

12.5. Results of PBT and vPvB assessment

No information available.

12.6. Other adverse effects

This product does not contain any known or suspected endocrine disruptors.

Section 13. Disposal considerations

13.1. Waste treatment methods

Waste from Residues / Unused Products

Dispose of in accordance with local regulations.

Contaminated Packaging

Empty containers should be taken to an approved waste handling site for recycling or disposal.

Section 14. Transport information

Note: The information provided below may not apply to all shipping situations. Consult

appropriate Dangerous Goods Regulations for additional requirements and mode-specific,

material-specific, or quantity-specific shipping requirements.

IMDG/IMO

14.1. UN-NumberNot regulated.14.2. Proper Shipping NameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.

14.5. Marine Pollutant None.14.6. Special Provisions None.

14.7. Transport in bulk according No information available.

to Annex II of MARPOL 73/78 and

the IBC Code

RID

14.1. UN-NumberNot regulated.14.2. Proper Shipping NameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.

14.5. Environmental hazard None. **14.6. Special Provisions** None.

ADR

14.1. UN-NumberNot regulated.14.2. Proper Shipping NameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.

14.5. Environmental hazard None. **14.6. Special Provisions** None.

ICAO

14.1. UN-NumberNot regulated.14.2. Proper shipping nameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.14.5. Environmental hazardNone.

14.5. Environmental hazard None. **14.6. Special Provisions** None.

IATA

14.1. UN-NumberNot regulated.14.2. Proper Shipping NameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.

14.5. Environmental hazard None. **14.6. Special Provisions** None.

Section 15. Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

International Inventories

TSCA Complies **EINECS/ELINCS** Complies DSL/NDSL Complies **PICCS** Complies Complies **ENCS** Complies **IECSC** Complies **AICS KECL** Complies

Legend

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

PICCS - Philippines Inventory of Chemicals and Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

AICS - Australian Inventory of Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

15.2. Chemical Safety Assessment

No information available

Section 16. Other information

Key literature references and sources for data

www.ChemADVISOR.com/

Issuing Date 24-Apr-2014

Revision Date 06-Nov-2015

Revision Note Name change.

This safety data sheet complies with the requirements of Commission Regulation (EU) No 453/2010 of 20 May 2010 amending Regulation (EC) No. 1907/2006

Revision Date 06-Nov-2015

WPS-CAN-014 - N26, N44, N75, N140, N250, N415, N750, N1000, N1000(KU), N1400, C200, G200, S200(KU), S200, C350, G350, N350, N350(KU), C600, G600, N600, S600(KU), S600, SV350, SV1400, VP190, VP960

General Disclaimer

The information provided on this SDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of Safety Data Sheet
