


Replacement Parts for L/S® Pump Heads

L/S® Standard Pump Head Parts



**Wait! There's More...**  
Mounting hardware for the Standard pump heads, see page 688. Easy Load® pump heads, see page 690.

**A. End Bell Assemblies.** Two required for a complete head assembly.

| For pump heads with CRS rollers |                             |       | For pump heads with SS rollers |                             |       |
|---------------------------------|-----------------------------|-------|--------------------------------|-----------------------------|-------|
| For pump heads                  | Catalog number              | Price | For pump heads                 | Catalog number              | Price |
| 07013-20                        | <a href="#">GY-07013-81</a> |       | 07013-21                       | <a href="#">GY-07013-91</a> |       |
| 07014-20                        | <a href="#">GY-07014-81</a> |       | 07014-21                       | <a href="#">GY-07014-91</a> |       |
| 07015-20                        | <a href="#">GY-07015-81</a> |       | 07015-21                       | <a href="#">GY-07015-91</a> |       |
| 07016-20                        | <a href="#">GY-07016-81</a> |       | 07016-21                       | <a href="#">GY-07016-91</a> |       |
| 07017-20                        | <a href="#">GY-07017-81</a> |       | 07017-21                       | <a href="#">GY-07017-91</a> |       |
| 07018-20                        | <a href="#">GY-07018-81</a> |       | 07018-21                       | <a href="#">GY-07018-91</a> |       |
| 07024-20                        | <a href="#">GY-07024-81</a> |       | 07024-21                       | <a href="#">GY-07024-91</a> |       |
| —                               | —                           | —     | 07035-21                       | <a href="#">GY-07035-91</a> |       |


**B. Rotor Assemblies.** Cold-rolled steel (CRS) or stainless steel (SS).

| For pump heads with CRS rollers |                             |       | For pump heads with SS rollers |                             |       |
|---------------------------------|-----------------------------|-------|--------------------------------|-----------------------------|-------|
| For models                      | Catalog number              | Price | For models                     | Catalog number              | Price |
| 07013-20, -42                   | <a href="#">GY-07013-80</a> |       | 07013-21, -52                  | <a href="#">GY-07013-95</a> |       |
| 07014-20, -42                   |                             |       | 07014-21, -52                  |                             |       |
| 07016-20, -42                   |                             |       | 07016-21, -52                  |                             |       |
| 07017-20, -42                   |                             |       | 07017-21, -52                  |                             |       |
| 07018-20, -42                   |                             |       | 07018-21, -52                  |                             |       |
| 07015-20, -42                   | <a href="#">GY-07015-80</a> |       | 07015-21, -52                  | <a href="#">GY-07015-90</a> |       |
| 07024-20, -42                   |                             |       | 07024-21, -52                  |                             |       |
| 07035-20                        |                             |       | 07035-21                       |                             |       |

| Description    | Catalog number              | Price/pk of 10 |
|----------------|-----------------------------|----------------|
| Thrust washers | <a href="#">GY-07021-04</a> |                |

[GY-07013-90](#) L/S tubing loading key for Standard pump heads


L/S® Easy-Load® Pump Head Parts



**D-E Main Body and Rotor/Bearing Assemblies**

| For pump heads  | Body/bearing material | Catalog number              | Price |
|---|-----------------------|-----------------------------|-------|
| <b>D Main body assemblies; occlusion bed with rear bearing and screws</b> |                       |                             |       |
| 07516-00, -02, -20, -22;<br>07518-00, -02, -20, -22                       | White PPS/CRS         | <a href="#">GY-07516-07</a> |       |
| 07516-10, -12, -30, -32;<br>07518-10, -12, -30, -32                       | White PPS/SS          | <a href="#">GY-07516-17</a> |       |
| 07518-60, -62, -80, -82   | Black PPS/SS          | <a href="#">GY-07518-67</a> |       |
| <b>E Rotor/bearing support assemblies</b>                                 |                       |                             |       |
| 07516-00, 07518-00  | White PPS/CRS         | <a href="#">GY-07516-08</a> |       |
| 07516-02, 07518-02  | White PPS/CRS         | <a href="#">GY-07516-09</a> |       |
| 07516-10, 07518-10  | White PPS/SS          | <a href="#">GY-07516-18</a> |       |
| 07516-12, 07518-12  | White PPS/SS          | <a href="#">GY-07516-19</a> |       |
| 07518-60  | Black PPS/SS          | <a href="#">GY-07518-68</a> |       |
| 07518-62  | Black PPS/SS          | <a href="#">GY-07518-69</a> |       |

L/S Easy-Load® II Pump Head Parts



**F. Rotor Assemblies**

| For pump head                          | Body/rotor material† | Catalog number   | Price |
|--|----------------------|--|-------|
| 77200-50, 77200-52                     | PPS/CRS              | <a href="#">GY-77200-58</a><br><a href="#">GY-77200-59</a> |       |
| 77200-60, 77201-60, 77200-62, 77201-62 | PPS/SS               | <a href="#">GY-77200-68</a><br><a href="#">GY-77200-69</a> |       |

†PPS—Polyphenylene sulfide PSF—Polysulfone PC—Polycarbonate  
CRS—Cold-rolled steel SS—Stainless steel

Tubing Retainer Kits (not shown)

| For pump head                | Catalog number              | Price |
|------------------------------|-----------------------------|-------|
| 77200-50, 77200-60, 77201-60 | <a href="#">GY-77200-90</a> |       |
| 77200-52, 77200-62, 77201-62 | <a href="#">GY-77200-92</a> |       |
| 77202-50, 77202-60           | <a href="#">GY-77200-94</a> |       |

L/S High-Performance 77250-62 Pump Head Parts

- G. [GY-77250-69](#) Rotor assembly with clip**
- H. [GY-77250-66](#) Tubing retainer kit**
- I. [GY-77250-67](#) Occlusion bed assembly with clip**
- [GY-77250-01](#) Mounting screws set (not shown)**

