

SAFETY DATA SHEET

Issuing Date 24-Apr-2014 Revision Date 24-Apr-2014 Revision Number 0

Section 1. Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product Name S3, C6, G6, S6, SV3 Contains Petroleum distillates, hydrotreated light naphthenic

1.2. Relevant identified uses of the substance or mixture and uses advised against

Recommended Use Viscometer and/or density measurement equipment calibration and performance

verification reference standard

Uses advised against No information available

1.3. Details of the supplier of the safety data sheet

Supplier

Cannon Instrument Company

2139 High Tech Rd.

State College, PA 16803-1733

TEL: (814) 353-8000; (800) 676-6232

For further information, please contact

E-mail Address No information available.

1.4. Emergency telephone number

Emergency Telephone (800) 255-3924 Domestic CHEM-TEL Inc.

Number +1 (813) 248-0585 Overseas CHEM-TEL Inc. (Please Call Collect)

Europe 112

Section 2. Hazards identification

2.1. - Classification of the substance or mixture

REGULATION (EC) No 1272/2008

Aspiration Toxicity Category 1

Physical Hazards

None

Classification according to EU Directives 67/548/EEC or 1999/45/EC

For the full text of the R-phrases mentioned in this Section, see Section 16

Symbol(s) Xn - Harmful R-code(s) Xn;R65

2.2. Label Elements



Signal Word Danger

Hazard Statements

H304 - May be fatal if swallowed and enters airways

Precautionary Statements - EU (§28, 1272/2008)

P308 + P313 - IF exposed or concerned: Get medical advice/ attention

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician

P331 - Do NOT induce vomiting

2.3. Other information

May cause central nervous system depression with nausea, headache, dizziness, vomiting, and incoordination Repeated or prolonged skin contact may cause allergic reactions with susceptible persons

Prolonged skin contact may defat the skin and produce dermatitis

Section 3. Composition/information on ingredients

3.1. Substances

Chemical Name	EC-No	CAS-No	Weight %	Classification	EU - GHS Substance Classification	REACH No.
Petroleum distillates, hydrotreated light naphthenic	Present	64742-53-6	100	Carc.Cat.2; R45 Xn;R65	Asp. Tox. 1 (H304) Carc. 1B (H350)	No data available

For the full text of the R-phrases mentioned in this Section, see Section 16
For the full text of the H-Statements mentioned in this Section, see Section 16

Section 4. First aid measures

4.1. Description of first-aid measures

Eye Contact Rinse thoroughly with plenty of water, also under the eyelids. Obtain medical attention if

irritation persists.

Skin Contact Wash off with warm water and soap. Remove and wash contaminated clothing before

re-use. Get medical attention if symptoms occur.

Ingestion Do NOT induce vomiting. Rinse mouth. Never give anything by mouth to an unconscious

person. Drink plenty of water. Get medical attention.

Inhalation IF INHALED: If breathing is difficult, remove to fresh air and keep at rest in a position

comfortable for breathing. Get medical attention

4.2. Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects Drowsiness. Dizziness. Irritation.

4.3. Indication of immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

Section 5. Fire-fighting measures

5.1. Extinguishing media

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Extinguishing media which must not be used for safety reasons

No information available.

5.2. Special hazards arising from the substance or mixture

Special exposure hazards arising from the substance or preparation itself, combustion products, resulting gases None in particular.

5.3. Advice for firefighters

Special protective equipment for fire-fighters

As in any fire, wear self-contained breathing apparatus and full protective gear.

Section 6. Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing. Use personal protective equipment.

6.2. Environmental precautions

Prevent further leakage or spillage if safe to do so.

6.3. Methods and materials for containment and cleaning up

Contain and collect spillage with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see Section 13).

Clean contaminated surface thoroughly.

6.4. Reference to other sections

See Section 12 for additional information.

Section 7. Handling and storage

7.1. Precautions for Safe Handling

Handling

Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes and clothing. Wear personal protective equipment.

Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice.

7.2. Conditions for safe storage, including any incompatibilities

Keep containers tightly closed in a dry, cool and well-ventilated place. Keep away from open flames, hot surfaces and sources of ignition.

7.3. Specific end use(s)

Exposure Scenario

No information available.

Other Guidelines

No information available.

Section 8. Exposure controls/personal protection

8.1. Control parameters

Exposure Limits

This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

Derived No Effect Level

Predicted No Effect Concentration

(PNEC)

No information available No information available.

8.2. Exposure controls

Engineering Measures

Ensure adequate ventilation, especially in confined areas.

Personal protective equipment

Eye Protection Safety glasses with side-shields.

Skin and Body Protection Hand Protection

Long sleeved clothing. Protective gloves.

Respiratory Protection

When workers are facing concentrations above the exposure limit they must use

appropriate certified respirators.

Environmental Exposure Controls

No information available.

Section 9. Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical State Liquid. Appearance Pale yellow.

Odor Hydrocarbon-like.

Property Values Remarks/ - Method

No data available None known Ha Melting Point/Range No data available None known **Boiling Point/Boiling Range** >193 °C None known > 96 °C **Flash Point** Open cup No data available None known **Evaporation rate** Flammability (solid, gas) No data available None known

Vapor Pressure No data available None known **Vapor Density** No data available None known **Relative Density** 0.87-0.88 at 15 °C Water Solubility Insoluble in water. None known Solubility in other solvents Soluble in solvents. None known Partition coefficient: n-octanol/waterNo data available None known **Autoignition Temperature** No data available None known **Decomposition Temperature** No data available None known 3-6 cSt @ 40°C **Viscosity** None known

Explosive Properties No information available Oxidizing Properties No information available

9.2. Other information

VOC Content (%)

Flammability Limits in Air

No information available.

No information available.

Section 10. Stability and reactivity

10.1. Reactivity

Not reactive under normal conditions.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

None under normal processing.

10.4. Conditions to avoid

Keep away from open flames, hot surfaces and sources of ignition.

10.5. Incompatible materials

Strong oxidizing agents.

10.6. Hazardous decomposition products

None under normal use.

Section 11. Toxicological information

11.1.

Acute Toxicity

Ingestion

Product Information

Inhalation Inhalation of mist may cause irritation to the respiratory system.

Eye Contact Contact with eyes may cause irritation.

Skin Contact No known hazard in contact with skin. Prolonged skin contact may defat the skin and

produce dermatitis. Prolonged or repeated contact may dry skin and cause irritation. Potential for aspiration if swallowed. Aspiration may cause pulmonary edema and

pneumonitis. Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea.

May be fatal if swallowed and enters airways.

SensitizationNo information available.Mutagenic EffectsNo information available.

Carcinogenic Effects The classification listed for the petroleum distillates in this product pertains to those that

contain more than 3% DMSO extract as measured by IP 346. The petroleum distillates in

this product do not meet that criteria to be classified as carcinogens.

Reproductive Toxicity

Developmental Toxicity

STOT - single exposure

STOT - repeated exposure

No information available.

No information available.

No information available.

Target Organ Effects Skin. Central nervous system (CNS). Eyes. Respiratory system.

Aspiration Hazard No information available.

Section 12. Ecological information

12.1. Toxicity

Ecotoxicity Effects

The environmental impact of this product has not been fully investigated.

Chemical Name	Toxicity to Algae	Toxicity to Fish	Toxicity to Microorganisms	Daphnia Magna (Water Flea)
Petroleum distillates, hydrotreated light		LC50 96 h: > 5000 mg/L (Oncorhynchus mykiss)		EC50 48 h: > 1000 mg/L (Daphnia magna)
naphthenic		(, , , , , , , , , , , , , , , , , , ,		(1

12.2. Persistence and degradability

No information available.

12.3. Bioaccumulative potential.

No information available.

12.4. Mobility in soil

Adsorbs on soil.

12.5. Results of PBT and vPvB assessment

No information available.

12.6. Other adverse effects

This product does not contain any known or suspected endocrine disruptors.

Section 13. Disposal considerations

13.1. Waste treatment methods

Waste from Residues / Unused

Products

Dispose of in accordance with local regulations.

Contaminated Packaging Empty containers should be taken to an approved waste handling site for recycling or

disposal.

Section 14. Transport information

Note: The information provided below may not apply to all shipping situations. Consult

appropriate Dangerous Goods Regulations for additional requirements and mode-specific,

material-specific, or quantity-specific shipping requirements.

IMDG/IMO

14.1. UN-NumberNot regulated.14.2. Proper Shipping NameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.

14.5. Marine Pollutant None.14.6. Special Provisions None.

14.7. Transport in bulk according No information available.

to Annex II of MARPOL 73/78 and

the IBC Code

RID

14.1. UN-NumberNot regulated.14.2. Proper Shipping NameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.

14.5. Environmental hazard None. **14.6. Special Provisions** None.

ADR

14.1.UN-NumberNot regulated.14.2.Proper Shipping NameNot regulated.14.3.Hazard ClassNot regulated.14.4.Packing GroupNot regulated.

Description Not applicable.

14.5. Environmental hazard None.14.6. Special Provisions None.

ICAO

14.1. UN-NumberNot regulated.14.2. Proper shipping nameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.

14.5. Environmental hazard None. **14.6. Special Provisions** None.

IATA

14.1. UN-NumberNot regulated.14.2. Proper Shipping NameNot regulated.14.3. Hazard ClassNot regulated.14.4. Packing GroupNot regulated.DescriptionNot applicable.

14.5. Environmental hazard None.14.6. Special Provisions None.

Section 15. Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

International Inventories

TSCA -

EINECS/ELINCS

DSL/NDSL

PICCS

Complies

ENCS

Complies

Legend

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

PICCS - Philippines Inventory of Chemicals and Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

AICS - Australian Inventory of Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

15.2. Chemical Safety Assessment

No information available

Section 16. Other information

Full text of R-phrases referred to under Sections 2 and 3

R45 - May cause cancer

R65 - Harmful: may cause lung damage if swallowed

Full text of H-Statements referred to under sections 2 and 3

H304 - May be fatal if swallowed and enters airways

H350 - May cause cancer if swallowed

Key literature references and sources for data

www.ChemADVISOR.com/

Issuing Date 24-Apr-2014

Revision Date 24-Apr-2014

Revision Note Initial Release.

This safety data sheet complies with the requirements of Commission Regulation (EU) No 453/2010 of 20 May 2010 amending Regulation (EC) No. 1907/2006

General Disclaimer

The information provided on this SDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of Safety Data Sheet