

Jars, PFA

Samples won't stick to surface

- Nonstick surface is easy to clean
- Resist corrosive liquids, solvents, acids, sensitive reagents, and staining solutions
- Use 2200-mL size for EPA TCLP soil sampling (40 CFR, Part 261 et. atl.)
- Autoclavable



Capacity mL (oz)	Dimensions (H x OD)		Neck size	Catalog number	Price/ea
	in.	cm			
60 (2)	1 7/8 x 2	4.7 x 5.1	51 mm	GY-34607-58	
120 (4)	2 x 2 1/2	5.1 x 6.4	66 mm	GY-34607-60	
240 (8)	4 1/4 x 3 1/4	10.8 x 8.3	84 mm	GY-34607-62	
360 (12)	3 1/2 x 3 1/4	8.8 x 8.3	84 mm	GY-34607-64	
480 (16)	3 3/4 x 2 1/2	9.4 x 6.4	65 mm	GY-34607-66	
1000 (32)	6 x 4 1/4	15.2 x 10.7	107 mm	GY-34607-68	
2000 (64)	10 1/2 x 4	26.9 x 10.4	105 mm	GY-34607-70	
2200 (70)	11 1/2 x 4 1/4	29.3 x 10.7	107 mm	GY-34607-72	

Wide-Mouth Sample Containers, PS

Won't break like glass



- Get glass-like clarity with more rigidity
- Included screw caps keep samples secure
- Containers are polystyrene (PS); lids are polypropylene (PP) with polyethylene foam lines

Capacity mL (oz)	Dimensions (H x OD)		Cap size	Catalog number	Qty/ pk	Price/pk
	in.	mm				
60 (2)	2 x 1 3/4	53 x 46	53 mm	GY-62263-20	48	
120 (4)	2 3/4 x 2	68 x 53	58 mm	GY-62263-30	36	
240 (8)	3 1/2 x 2 1/2	89 x 65	89 mm	GY-62263-40	36	
480 (16)	3 3/4 x 3 1/2	96 x 89	89 mm	GY-62263-50	24	
960 (32)	4 3/4 x 3 3/8	120 x 98	120 mm	GY-62263-60	24	

Cole-Parmer® Cylindrical Jars, HDPE

Get easy sample access

- Wide neck for easy powder sampling
- Include screw caps to keep samples intact
- Won't break—made of high-density polyethylene for sturdy storage



Capacity mL (oz)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	mm			
70 (2.4)	2 5/16 x 2	60 x 50	GY-63100-10	10	
120 (4)	2 3/4 x 2 1/4	71 x 56	GY-63100-12		
250 (8)	3 3/4 x 2 3/4	94 x 69	GY-63100-14		
500 (16)	4 1/4 x 3 3/8	107 x 86	GY-63100-16		
1000 (32)	5 x 4 3/8	128 x 111	GY-63100-18	10	
1500 (51)	7 7/16 x 4 3/8	182 x 111	GY-63100-20		
2000 (64)	9 1/4 x 4 3/8	235 x 111	GY-63100-22		

Collapsible Containers, LDPE

Easily expand when filled to save space

- Stack neatly when empty, just unscrew cap and fill to expand
- Include polypropylene caps



Capacity L (gal.)	Cap size	Catalog number	Qty/ pk	Price/ pk
1 (1/4)	38 mm	GY-06100-10	12	
4 (1)		GY-06100-20		
10 (2 1/2)		GY-06100-30		
19 (5)		GY-06100-40		