

Centrifuge Bottles

Get better separation of large samples

Conical Bottles

- Measure more precisely with molded-in graduations
- Include leakproof, PP cap

Flat-Bottom Bottles

- Use model 06107-20 at low speeds, filled completely, and in refrigerated centrifuge to maintain leak-free spinning
- Bottles 06107-20 and 06106-10 require sealing cap assembly 06327-57 for spinning at high speeds
- 1000 mL bottles require sealing cap 06327-69 for leak-free spinning at maximum speed



Capacity (mL)	Material	Dimensions (H x OD)	Max RCF (x g)	Autoclavable	Bottle			Adapter required			
					Catalog number	Qty/pk	Price/pk	Catalog number	Qty/pk	Price/pk	
Conical bottles											
175	PPCO	139.2 x 61.5 mm	27,500	Yes	GY-06107-28	4		GY-06106-78	4		
	PC	140.0 x 61.4 mm			GY-06106-08						
Flat-bottom bottles											
250	PPCO	127.7 x 61.8 mm	13,200	Yes	GY-06105-20	4		GY-17404-06†	6		
	HDPE	131.9 x 61.2 mm	8000	No	GY-06107-10	6					
	FEP	128.8 x 60.0 mm	4000	Yes	GY-06107-20	1					
250	PC	127.6 x 61.8 mm	27,500	Yes	GY-06106-10	4		GY-17404-06†	6		
		135.1 x 61.5 mm			GY-06106-12						
		133.4 x 60.7 mm			GY-06106-14						
450	PPCO	160.0 x 69.5 mm	13,700	Yes	GY-06105-55	4		—	—	—	
	PC				GY-06106-35						
500	PPCO	169.8 x 73.8 mm	4800	Yes	GY-06105-50	4		—	—	—	
	PPCO	170.2 x 69.5 mm	13,700		GY-06106-25						
	PC	169.6 x 69.5 mm	13,700		GY-06106-30						
1000	PPCO	184.5 x 97.5 mm	7100	Yes	GY-06105-60	4		—	—	—	
	PC	188.4 x 97.6 mm			GY-06106-40						
	PPCO	179.0 x 97.7 mm			GY-06107-75						
	PC	180.2 x 98.1 mm			GY-06107-80						

†These bottles require the noted adapter to fit the rotor.

[GY-06327-57](#) Sealing cap assembly, 38-mm. Pack of 2

[GY-06327-69](#) Replacement PP sealing cap for 450, 500, and 1000 mL bottles. Pack of 2

BOD Bottles and Bottle Racks

Easily identify samples

- Clear glass BOD bottles come with glass stoppers and a white write-on marking space

Use a bottle rack for storing and transporting your BOD bottles. Racks are corrosion resistant and plastisol coated. For our full range of benchtop dissolved oxygen meters including BOD, see your local representative.

Specifications

Capacity: 300 mL

Description	Number sequence	Qty/ cs	Catalog number	Price/ cs
Penny head stopper BOD bottles, clear	Unnumbered	24	GY-03660-14	
Robotic stopper BOD bottles, clear	Unnumbered	24	GY-03660-15	
Robotic stopper BOD bottles, clear	1–24	24	GY-03660-16	
	25–48		GY-03660-17	
	49–72		GY-03660-18	
	73–96		GY-03660-19	



Accessories

[GY-03660-09](#) Cap and insert for use with 03660-15 BOD bottles. Pack of 50

[GY-03660-12](#) Rubber bottle leash holds stopper to bottle. Pack of 25

[GY-03660-20](#) BOD bottle rack, holds twelve 300-mL bottles