

Racks for Chest Freezers

TENAK has 3 different kinds of racks for chest freezers, which fit into any freezer brand- and model on the world market.

Eco AluRack

- Low weight
- Made of aluminium

An inexpensive racking option for basic storage. It has no safety system. This rack is assembled without netting or welding, which is only available from TENAK.

Classic Rack

- With locking rod
- Made of stainless steel

This is the most commonly used type of rack for mechanical chest freezers. It is supplied with a locking rod to prevent the boxes from slipping out.




Comfort Rack

- With TENAK safety spring lock
- Made of stainless steel

The Comfort Rack is for the critical user for application where safety and/or ergonomics matter. This type of rack is highly recommended for use in LN2 freezers. The rack has the unique TENAK Safety spring lock design.

Max box size 137 x 137 x 54 mm.

Usability



Racks for Chest Freezers

Eco AluRack	Classic Rack	Comfort Rack			
Article no.	Article no.	Article no.	For boxes	HxDxW	No. of boxes
TE11458	TE21350	TE24458	32 mm	288x140x140	8
TE11460	TE21352	TE24460	32 mm	360x140x140	10
TE11462	TE21354	TE24462	32 mm	432x140x140	12
TE11464	TE21356	TE24464	32 mm	504x140x140	14
TE11465	TE21357	TE24465	32 mm	540x140x140	15
TE11466	TE21358	TE24466	32 mm	576x140x140	16
TE11468	TE21360	TE24468	32 mm	648x140x140	18
TE11470	TE21362	TE24470	32 mm	720x140x140	20
TE11500	TE21448	TE24500	50 mm	170x140x140	3
TE11502	TE21449	TE24502	50 mm	280x140x140	5
TE11504	TE21450	TE24504	50 mm	335x140x140	6
TE11506	TE21451	TE24506	50 mm	390x140x140	7
TE11508	TE21440	TE24508	50 mm	445x140x140	8
TE11510	TE21441	TE24510	50 mm	500x140x140	9
TE11512	TE21452	TE24512	50 mm	555x140x140	10
TE11514	TE21454	TE24514	50 mm	610x140x140	11
TE11516	TE21456	TE24516	50 mm	665x140x140	12
TE11518	TE21478	TE24518	50 mm	720x140x140	13
TE11520	TE21430	TE24520	50 mm	775x140x140	14
TE11522	TE21432	TE24522	50 mm	832x140x140	15
TE11530	TE21455	TE24530	75 mm	166x140x140	2
TE11532	TE21457	TE24532	75 mm	246x140x140	3
TE11534	TE21458	TE24534	75 mm	326x140x140	4
TE11536	TE21460	TE24536	75 mm	406x140x140	5
TE11538	TE21462	TE24538	75 mm	486x140x140	6
TE11540	TE21464	TE24540	75 mm	551x140x140	7
TE11542	TE21466	TE24542	75 mm	629x140x140	8
TE11544	TE21467	TE24544	75 mm	726x140x140	9
TE11550	TE21465	TE24550	85 mm	607x140x140	7
TE11552	TE21715	TE24552	85 mm	693x140x140	8
TE11560	TE21468	TE24560	100 mm	320x140x140	3
TE11562	TE21470	TE24562	100 mm	425x140x140	4
TE11564	TE21472	TE24564	100 mm	530x140x140	5
TE11566	TE21474	TE24566	100 mm	635x140x140	6
TE11570	TE25100	TE24570	130 mm	270x140x140	2
TE11572	TE25102	TE24572	130 mm	404x140x140	3
TE11574	TE25104	TE24574	130 mm	536x140x140	4
TE11576	TE25106	TE24576	130 mm	669x140x140	5