


3-way stopcocks, with solid glass plug



3-way stopcocks, with solid glass plug

Acc. to DIN 12554 and in addition to DIN standard, T-bore, with solid glass plug and screwthread retaining device, material: borosilicate glass 3.3, leakage rate to DIN 12540-1

Stopcock(NS)	Bore (mm)	Acc. to DIN	Pack (qty.)	Cat. No.
14,5	1,5	DIN 12554	1	2 5014 01
14,5	2,5	-	1	2 5014 02
18,8	2,5	DIN 12554	1	2 5018 02
18,8	4,0	-	1	2 5018 04
24	4,0	DIN 12554	1	2 5024 04
29,2	6,0	DIN 12554	1	2 5029 06

