

KERN 313-05

Weight set, E2, 1 mg - 500 g, Knob, stainless steel polished, in wooden case



Article codes, product groups & Status information

Model article-code	313-05
Model-Serie	313
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main article
EAN code (Model)	4045761003132
Planned start date	2004-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 mg - 500 g
Case wood	313-050-100
Weight set case	In wooden case

Technical data - Construction

Material	Stainless steel polished
Is weight set	Yes
Number of weights in set	24
Type of case	In wooden case

Technical data - Measuring system

Nominal value	1 mg - 500 g
OIML class	E2
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×180×110 mm
Gross weight	3.2 kg
Net weight	3.2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-304
-------------------	---------

Functions