

Injection Systems

Needles and Syringes

BD Discardit™ II – 2-Piece Syringes, Luer-Slip Tip – 2, 5, 10, 20 ml

- Transparent barrel for good visualization of syringe contents
- Smooth plunger movement - reduced sliding force
- BD Discardit™ syringes are tested during manufacture in compliance with ISO 7886-1 to make them leak-tight
- Accurate graduation suited to user's needs
- With a retaining ring to prevent accidental plunger rod withdrawal
- Design of finger grip flanges ensures stability and comfort during injection
- Ergonomically designed plunger for single-handed operation
- Sterilization method: ethylene oxide
- Sterile, single use
- CE marked



Capacity	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Discardit™ Syringes with Luer-Slip Tip				
2 ml	concentric	0,1 ml	100 / 3000	300928
5 ml	eccentric	0,2 ml	100 / 1800	309050
10 ml	eccentric	0,5 ml	100 / 1200	309110
20 ml	eccentric	1 ml	80 / 960	300296

Capacity	Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Discardit™ Syringes with Luer-Slip Tip and detached BD Microlance™ Needle							
2 ml	21G	0,8	40	1 1/2"	green	100 / 3000	300318
2 ml	22G	0,7	30	1 1/4"	black	100 / 3000	300320
2 ml	22G	0,7	40	1 1/2"	black	100 / 3000	300321
2 ml	23G	0,6	25	1"	blue	100 / 3000	300322
2 ml	23G	0,6	30	1 1/4"	blue	100 / 3000	300323
5 ml	22G	0,7	40	1 1/2"	black	100 / 1800	300506
5 ml	21G	0,8	40	1 1/2"	green	100 / 1800	300512
5 ml	22G	0,7	30	1 1/4"	black	100 / 1800	300513
5 ml	23G	0,6	30	1 1/4"	blue	100 / 1800	300516
10 ml	22G	0,7	30	1 1/4"	black	100 / 1200	300626
10 ml	21G	0,8	40	1 1/2"	green	100 / 1200	301001
10 ml	22G	0,7	40	1 1/2"	black	100 / 1200	301504
20 ml	21G	0,8	40	1 1/2"	green	80 / 960	300330