

BR 010 900

ClassicLine safety shower test cart, capacity 30 litres, mobile



Technical description

ClassicLine safety shower test cart, capacity 30 litres, mobile
With the safety shower test cart the prescribed monthly function test in accordance with laboratory regulations can be carried out quickly and simply. The trolley can be dismantled into three parts, thus saving space during storage. The assembled test cart is 2.10 metre in height and is positioned under the shower for testing. The collecting device holds approx. 30 litre and can be emptied with ease over a floor drain by means of the drain plug.

- mobile collecting device made of red plastic fitted with 4x castors
- two recessed grips integrated in the container
- drain plug for draining
- funnel made of plastic, black, diameter 440 mm, height 160 mm
- total height 2100 mm
- pipe material aluminum 47 mm, pluggable
- tank capacity 30 litre
- total weight 6 kg
- detachable in 3 parts for simple storage

Manufacturer: B-SAFETY or equal

Article-No.: BR 010 900

Technical specifications

Tank capacity:	30 litre
Total weight:	6 kg
Dimensions (H x W x D):	2100 x 440 x 440 mm

Material-Overview	Material	Material description
Collecting device	plastic, red	-
Funnel:	plastic, black	-
Pipe:	Aluminum	-

Product overview

BR 010 900	ClassicLine safety shower test cart, capacity 30 litres, mobile
-------------------	---

