

Chemware PTFE Coated Stainless Steel Tongs



Description • Opposing grip handle for easy operation

- Coated tension screw will not corrode
- Contoured jaw enables handling of crucibles and other objects
- Tapered tips for enhanced grip; rounded for safety
- Use for handling chemically active materials
- Nickel-plated stainless steel tongs are covered with extra-heavy, non-stick PTFE coating that eliminates surface build-up and loss of valuable samples
- Will not scratch or damage delicate surfaces

Image(s) are representative of the product group and not necessarily the individual product.

Order	02-Mar-2026 20:45: PM PST	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
142-51987-EA	9.4 in. (24 cm) PTFE-Coated Stainless Steel Tongs - 1 ea	Sign In	Saint Gobain Performance Plastics Corporation	D1069290			EA
142-51987-EA	9.4 in. (24 cm) PTFE-Coated Stainless Steel Tongs - 1 ea		Saint Gobain Performance Plastics Corporation	D1069290			CS(6)