


In addition to the ground stopper, a glass cap with ground joint is supplied. This provides an improved seal and protection against acid vapours.

Cat. No.	Description	Capacity (mL)	d (mm)	h (mm)	Neck	Pack Unit
212752408		100	55	145	19/17	10
212753601		250	75	180	19/17	10
212754409		500	82	220	24/20	10
212755405		1 000	109	260	29/32	10
Components						
212732406	Bottle	100	55	104	19/17	10
212733608	Bottle	250	75	133	19/17	10
212734407	Bottle	500	87	166	24/20	10
212735403	Bottle	1 000	108	208	29/32	10
212742407	Cap for Bottle 100 mL		48	73		10
212743609	Cap for bottle 250 mL		55	75		10
212744408	Cap for Bottle 500 mL		66	87		10
212745404	Cap for Bottle 1 000 mL		75	103		10

DURAN® Acid Bottle

with standard ground "pennyhead" stopper; conical shoulders, interchangeable glass cap




In addition to the ground stopper, a glass cap with ground joint is supplied. This provides an improved seal and protection against acid vapours.

Cat. No.	Description	Capacity (mL)	d (mm)	h (mm)	Neck	Pack Unit
212752465		100	55	145	19/17	10
212753667		250	75	180	19/17	10
212754466		500	82	220	24/20	10
212755462		1 000	109	260	29/32	10
Components						
212732463	Bottle	100	55	104	19/17	10
212733665	Bottle	250	76	133	19/17	10
212734464	Bottle	500	87	166	24/20	10
212735469	Bottle	1 000	108	208	29/32	10
212742464	Cap for Bottle 100 mL		48	73		10
212743666	Cap for Bottle 250 mL		55	75		10
212744465	Cap for Bottle 500 mL		66	87		10
212745461	Cap for Bottle 1 000 mL		75	103		10

DURAN® Acid Bottle Amber

with standard ground „pennyhead“ stopper; conical shoulders, interchangeable glass cap




For dosing use with the dropping pipette. Spare pipettes, clear glass, Cat. No. 232711709 and 232712405 (Quantity 10); Rubber teats, transparent, Cat. No. 292000102 (Quantity 100).

Cat. No.	Capacity (mL)	h (mm)	Neck	Pack Unit
232701708	50	79	14/15	10
232702404	100	105	14/15	10


Dropping Bottle from soda-lime glass

with interchangeable clear glass pipette standard ground joint, complete with rubber teats



Dropping Bottle Amber from soda-lime glass

with interchangeable clear glass pipette standard ground joint, complete with rubber teats




For dosing use the dropping pipette. Spare pipettes, clear glass, Cat. No. 232711709 and 232712405 (Quantity 10): rubber teats, transparent, Cat. No. 292000102 (Quantity 100).

Cat. No.	Capacity (mL)	h (mm)	Neck	Remark	Pack Unit
232701765	50	79	14/15		10
232702461	100	105	14/15	from borosilicate 3.3	10

Rubber Teat transparent


from natural rubber



Cat. No.	d ₁ (mm)	h (mm)	Pack Unit
292000102	15	35	100

DURAN® Aspirator (levelling) Bottle

with plain neck and bottom sidearm



Useful as a delivery or storage container for solutions. Outlet facilitates attachment of flexible tubing.

Cat. No.	Capacity (mL)	d (mm)	d ₁ (mm)	d ₂ (mm)	h (mm)	Pack Unit
247083603	250	73	11	5	130	10
247084402	500	89	11	5	164	10
247085407	1 000	111	11	5	200	1




Dosing of liquids is possible via an outlet.

Cat. No.	Capacity (mL)	d (mm)	h (mm)	Tubulature (NS)	Pack Unit
247014404	500	86	164	19/26	10
247015409	1 000	107	200	19/26	10
247016302	2 000	134	249	19/26	1
247017307	5 000	182	320	29/32	1
247018603	10 000	228	398	29/32	1
247019102	20 000	289	492	29/32	1

DURAN® Aspirator (levelling) Bottle

tubulated with standard ground joint, without stoppers, neck unground



ISO
4796-3

Retrace
Code


A
121 °C

Dosing of liquids is possible via a stopcock.

Cat. No.	Capacity (mL)	d (mm)	h (mm)	Neck	Tubulature (NS)	Pack Unit
247024405	500	86	164	24/29	19/26	10
247025401	1 000	107	200	29/32	19/26	10
247026303	2 000	134	249	29/32	19/26	1
247027308	5 000	182	320	45/40	29/32	1
247028604	10 000	228	398	60/46	29/32	1
247029103	20 000	289	492	60/46	29/32	1

DURAN® Aspirator (levelling) Bottle

tubulated with standard ground joint, complete with standard ground stopcock and standard ground stopper



ISO
4796-3

Retrace
Code


A
121 °C

Complete with screw connection cap, silicone seal and stopcock with PTFE spindle. Dosing of liquids is possible via a stopcock.

Cat. No.	Capacity (mL)	d (mm)	h (mm)	DIN Thread (GL)	Pack Unit
247035402	1 000	101	225	45	1
247036304	2 000	136	260	45	1
247037309	5 000	182	330	45	1
247038605	10 000	230	410	45	1

DURAN® Aspirator (levelling) Bottle

neck with DIN thread GL 45, tabulator with GL 32




ISO
4796-3

Retrace
Code

A
121 °C

DURAN® Aspirator (levelling) Bottle


neck with DIN thread GLS 80®, tabulator with GL 32



Complete with screw connection cap, silicone seal and stopcock with PTFE spindle. Dosing of liquids is possible via a stopcock.

Cat. No.	Capacity (mL)	d (mm)	h (mm)	DIN Thread (GL)	Pack Unit
247047301	5 000	182	330	80	1

Stopcock with Standard Ground Joint for Aspirator Bottle




Spare part for aspirator bottle.

Cat. No.	Capacity (mL)	Neck	Pack Unit
241480307	500 – 2 000	19/26	1
241480401	5 000 – 20 000	29/32	1

Stopcock for Aspirator Bottle

with PTFE spindle, for GL 32 screw thread



Cat. No.	Capacity (mL)	Hole (mm)	Remark	Pack Unit
241470306	1 000 + 2 000	6	suitable silicone sealing ring: Cat. No. 292361 004	1
241470409	5 000 + 10 000	8	suitable silicone sealing ring: Cat. No. 292361 201	1