



designed for scientists



H 30 Oil bath attachment, 1.5L

/// Data Sheet

The oil bath attachment H 30 can be used as oil baths together with an IKAMAG® magnetic stirrer with heating or with an IKATHERM® heating plate having a diameter of 135 mm.

- Positioning border prevents sliding on the heating plate
- Safety grip protects from bath attachment burns
- The bath attachment is made of aluminum. This ensures good heat transfer and quick heating-up of the tempering medium
- Easy cleaning
- The bath attachment can only be used as an oil bath
- Inner diameter 136-190 mm



designed for scientists

- Height 110 mm
- Volume 1,5 Liter



designed for scientists

Technical Data

Inner diameter [mm]	190
Bottom shape	hollow cylindrical
Useable volume max. [ml]	1500
Working temperature max. [°C]	300
Material in contact with medium	Aluminum cast
Dimensions (W x H x D) [mm]	310 x 110 x 134
Weight [kg]	0.948