




designed for scientists



IKAFLON 140 slide (1 pc.)

/// Data Sheet

Magnetic stirring bar made of PTFE, with centric PEEK ring for reduction of friction and better concentricity with uneven vessel bottoms. For large vessels, turbulent mixing at low speeds.

Length: 140 mm, diameter: 27 mm, diameter with ring: 35 mm.

All materials are FDA compliant.



designed for scientists

Technical Data

Length [mm]	140
Diameter [mm]	35
Shape	round
Coating	PTFE
Material in contact with medium	PEEK + PTFE
Material: Magnet	SmCo
Magnetic attraction force	strong
Working temperature max. [°C]	280
Weight [kg]	0.22