



designed for scientists



PT 1000.90 Temperature sensor, coated

/// Data Sheet

For use with RCT basic, RET basic and C-MAG HS digital series.

Immersion depth: 230 mm

Diameter: 3 mm

Coating material: Non-stick coating based on perfluoroalkoxy and polytetrafluoroethylene



designed for scientists

Technical Data

Probe type	Immersion probe
Sensor type	PT1000 (DIN IEC 751 Class A)
Accuracy of temperature measurement [K]	$\pm (0,15 + 0,2\% \text{ o.Mv.})$
Response time $t_{9/10}$ [sec]	10
Temperature measuring range [°C]	-10 - +260
Material of probe tube	stainless steel 1.4571
Coating of probe tube	SafeCoat
Diameter of probe tube [mm]	3
Length of probe tube [mm]	230
Immersion depth min. [mm]	25
Cable length [mm]	1000
Probe connection	DIN 45322 (MAS6)
Weight [kg]	0.026