




designed for scientists



R 1375 Paddle stirrer

/// Data Sheet

Tangential flow, minimum turbulence, good heat exchange, gentle treatment of product. Used at low to medium speeds.

Stirrer Ø: 70 mm

Shaft Ø: 8 mm

Shaft length: 550 mm



designed for scientists

Technical Data

Stirrer diameter [mm]	70
Shaft diameter [mm]	8
Shaft length [mm]	550
Speed max. [rpm]	800
Material in contact with medium	AISI 316L
Weight [kg]	0.266