




designed for scientists



TRIKA 35 (5 pcs.)

/// Data Sheet

Dimensions (L x Ø): 35 x 9 mm

Suitable for mixing difficult to dissolve samples or to prevent sedimentation.

Provides strong turbulence at low speeds. The scraper-like shape is perfect for the dissolution of solids and stirring up precipitates.

The PTFE coating is chemically inert and resistant to aggressive chemicals.



designed for scientists

Technical Data

Length [mm]	35
Diameter [mm]	9
Coating	PTFE
Material in contact with medium	PTFE
Material: Magnet	Alnico V
Magnetic attraction force	normal
Working temperature max. [°C]	280
Weight [kg]	0.024