



Tubing

[back](#)[\[Laboratory Tubing\]](#)[\[Silicone Tubing\]](#)[\[Tubing made of Viton®\]](#)[\[Butyl Rubber Tubing\]](#)[\[Latex Tubing\]](#)[\[Perbunan Rubber Tubing\]](#)[\[PVC Tubing\]](#)[\[PVC Tubing with reinforcing webbing\]](#)[\[Polyethylene \(PE\) Tubing\]](#)[\[PTFE Tubing\]](#)

Silicone Tubing



Quality

Silicone Rubber (VMQ), meets the following standards:

BFR XV, FDA CFR 21§ 177.2600, Europ.Pharmacopoeia

Test for biocompatibility according to 10993-1

transparent, Durometer 60 Shore A

temperature range -60°C to +200°C,

short periods up to +260°C

Inner-Ø mm	Outer-Ø mm	Wall Thickness mm	Roll Length m	Catalog No
0,5	1,3	0,4	25	310 0504
0,5	2,5	1	25	310 0525
0,8	1,6	0,4	25	310 0804
1	1,8	0,4	25	310 0118
1	3	1	25	310 0103
1	4	1,5	25	310 0104
1	6	2,5	25	310 0106
1,5	2,3	0,4	25	310 1523
1,5	3	0,75	25	310 1503
1,5	3,5	1	25	310 1535
1,5	4	1,25	25	310 1504
2	2,6	0,3	25	310 0226

2	2,8	0,4	25	310 0228
2	4	1	25	310 0204
2	5	1,5	25	310 0205
2	6	2	25	310 0206
2,5	3,3	0,4	25	310 2533
2,5	4,5	1	25	310 2545
3	3,8	0,4	25	310 0338
3	5	1	25	310 0305
3	6	1,5	25	310 0306
Inner-Ø mm	Outer-Ø mm	Wall Thickness mm	Roll Length m	Catalog No
3	7	2	25	310 0307
4	5	0,5	25	310 0405
4	6	1	25	310 0406
4	7	1,5	25	310 0407
4	8	2	25	310 0408
4	9	2,5	25	310 0409
5	6	0,5	25	310 0506
5	7	1	25	310 0507
5	8	1,5	25	310 0508
5	9	2	25	310 0509
5	10	2,5	25	310 0510
6	8	1	25	310 0608
6	9	1,5	25	310 0609
6	10	2	25	310 0610
7	9	1	25	310 0709
7	10	1,5	25	310 0710
7	11	2	25	310 0711
8	10	1	25	310 0810
8	11	1,5	25	310 0811
Inner-Ø mm	Outer-Ø mm	Wall Thickness mm	Roll Length m	Catalog No
8	12	2	25	310 0812
8	14	3	25	310 0814
8	15	3,5	25	310 0815
9	12	1,5	25	310 0912
9	13	2	25	310 0913
10	13	1,5	25	310 1013
10	14	2	25	310 1014

10	15	2,5	25	310 1015
10	16	3	25	310 1016
10	18	4	25	310 1018
11	15	2	25	310 1115
12	16	2	25	310 1216
12	18	3	25	310 1218
13	19	3	25	310 1319
14	20	3	25	310 1420
15	21	3	20	310 1521
16	21	2,5	20	310 1621
18	24	3	20	310 1824
20	26	3	20	310 2026
25	33	4	15	310 2534

Vacuum Quality				
Inner-Ø mm	Outer-Ø mm	Wall Thickness mm	Roll Length m	Catalog No
3	9	3	25	311 0309
4	12	4	25	311 0412
5	15	5	15	311 0515
6	12	3	25	311 0612
7	13	3	25	311 0713
7	15	4	25	311 0715
7	17	5	15	311 0717
8	16	4	25	311 0816
8	18	5	15	311 0818
10	20	5	15	311 1020
12	22	5	15	311 1222