

Multi-position mantles

- Available in three or six positions, heating only or heating and stirring models
- Maximum element temperature of 450°C
- Polypropylene top case and power coated aluminium base provide excellent chemical resistance
- Grounded stainless steel earth screen over every heating position provides added safety
- Thermally insulated heater cartridge for maximum heat transfer to the flask
- Unique airflow within the mantle prevents the outer casing from getting hot



Three-position Heating Mantles EME Controlled Mantles

- Multi-position heating only mantles
- Individual temperature control for each position





Three-position Heating and Stirring EMEA Controlled Mantles

- Allow the solution to be stirred and heated simultaneously
- Provides stirring speeds of between 50 and 1000rpm
- One stirring control operates each set of three recesses

Part Number	Capacity	Dimensions, mm	Weight
EME30100/CEB	100ml	270 x 640 x 110	6.2 kg
EME30250/CEB	250ml	270 x 640 x 110	6.2 kg
EME30500/CEB	500ml	285 x 755 x 130	7.4 kg
EME31000/CEB	1000ml	285 x 755 x 130	7.4 kg


Part Number	Capacity	Dimensions, mm	Weight
EMEA30100/CE	100ml	270 x 640 x 110	7 kg
EMEA30250/CE	250ml	270 x 640 x 110	7 kg
EMEA30500/CE	500ml	285 x 755 x 130	8.4 kg
EMEA31000/CE	1000ml	285 x 755 x 130	8.4 kg



Six-position Heating Mantles EME Controlled Mantles

- Multi-position heating only mantles
- Individual temperature control for each position

Part Number	Capacity	Dimensions, mm	Weight
EME60100/CEB	100ml	270 x 1210 x 110	10.1 kg
EME60250/CEB	250ml	270 x 1210 x 110	10.1 kg
EME60500/CEB	500ml	285 x 1450 x 130	12.5 kg
EME61000/CEB	1000ml	285 x 1450 x 130	12.5 kg



Six-position Heating and Stirring EMEA Controlled Mantles

- Allow the solution to be stirred and heated simultaneously
- Provides stirring speeds of between 50 and 1000rpm
- One stirring control operates each set of six recesses

Part Number	Capacity	Dimensions, mm	Weight
EMEA60100/CE	100ml	270 x 1210 x 110	10.1 kg
EMEA60250/CE	250ml	270 x 1210 x 110	10.1 kg
EMEA60500/CE	500ml	285 x 1450 x 130	12.5 kg
EMEA61000/CE	1000ml	285 x 1450 x 130	12.5 kg

Quality. Reliability. **Safety.** Electrothermal.