


E 144

Insert

For the optimum loading of up to 18 beakers.



- Suitable for 18 glass beakers
- equipped with 18 x 3 stay bars
- For use in upper and lower basket


EAN: 4002511115996 / Material number: 03808710 / Old Material Number: 69514401D

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Area of application	
Reprocessing of laboratory glassware	•
Product properties	
Material	Stainless steel, Wire, rilsan-coated
Colour	Stainless Steel, White
Capacity	
Beakers [quantity]	18
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	131

E 144

Insert

For the optimum loading of up to 18 beakers.



EAN: 4002511115996 / Material number: 03808710 / Old Material Number: 69514401D

External dimensions, net width in mm	200
External dimensions, net depth in mm	445
External dimensions, gross height in mm	180
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	490
Net weight in kg	1,37
Gross weight in kg	1,38