


A 2

Cover net

For fixing delicate wash items.




- Suitable for 1/2 inserts
- Plastic coated metal frame with plastic netting
- Size: 216 x 456 mm

EAN: 4002511115590 / Material number: 03830460 / Old Material Number: 69500201D

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Cover	•
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Product properties	
Material	Wire, rilsan-coated, Plastic
Colour	Green, White, Black
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net width in mm	216
External dimensions, net depth in mm	456
External dimensions, gross height in mm	20
External dimensions, gross width in mm	230

A 2
Cover net
For fixing delicate wash items.



EAN: 4002511115590 / Material number: 03830460 / Old Material Number: 69500201D

External dimensions, gross depth in mm	465
Net weight in kg	0,39
Gross weight in kg	0,49