


E 105/1

Insert

For the optimum loading of test tubes up to 12 x 165 mm.



- Suitable for about 200 Test tubes (up to 12 x 165 mm)
- Capacity including lid A 13
- The mesh width is 9 x 9 mm

EAN: 4002516853978 / Material number: 12667380 / Old Material Number: 69510503D

Product category

Laboratory washers	•
Large-chamber laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•

Machine affiliation

G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•

Product category

Insert/module for test tubes	•
------------------------------	---

Area of application

Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•

Product properties

Material	Stainless steel
Colour	Stainless Steel

Capacity

Test tubes [number]	200
---------------------	-----

Dimensions and weight


External dimensions, net height in mm	212
External dimensions, net width in mm	200
External dimensions, net depth in mm	210
External dimensions, gross height in mm	235
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	250
Net weight in kg	0.88

E 105/1

Insert

For the optimum loading of test tubes up to 12 x 165 mm.

EAN: 4002516853978 / Material number: 12667380 / Old Material Number: 69510503D



Gross weight in kg

1.16

Companion pack

A 13 lid [number]

1