

Petri Dishes and Contact Dishes

In different sizes and formats

Greiner Bio-One offers petri dishes with diameters of 35, 60, 94, 100 and 145 mm. Manufactured from highly transparent polystyrene (PS), they are heat-resistant up to 60 °C for use with hot agar. Depending on the application vented and non-vented design can be selected.


Besides the standard round design, a square-profile dish (120 mm x 120 mm x 17 mm) and the new OneWell Plate™ in SBS format are available.

Furthermore, Greiner Bio-One offers contact dishes and graduated germ count dishes which are used with hygiene monitoring including detection of contaminations as well as testing the effectiveness of cleaning and disinfection on plane surfaces by means of contact cultures.

The product range is completed by special models with two or three compartments which enable parallel experiments with different sampels or parameters in one dish.


Key Facts Petri Dishes

- Available in different sizes
- Easy stacking
- With or without vents
- Manufactured from crystal clear PS
- Sterilised dishes are also available




Key Facts Contact Dishes/OneWell Plate™

- Available in different sizes
- With or without vents
- Manufactured from crystal clear PS
- Contact dish with graduation to enable quick and easy analyses
- If applicable, sterilised dishes are also available



Key Facts Special Models

- Easy stacking
- With vents
- Manufactured from crystal clear PS
- Sterilised dishes are also available



Ordering Information			
Cat.-No.	Product Description	Qty/ Bag	Qty/ Case
Petri Dishes			
627 102	Petri dish, Ø 35 x 10 mm, with vents	10	740
628 102	Petri dish, Ø 60 x 150 mm, with vents	20	600
632 102	Petri dish, Ø 94 x 16 mm, without vents, heavy design	20	480
632 180	Petri dish, Ø 94 x 16 mm, without vents	20	480
633 102	Petri dish, Ø 94 x 16 mm, with vents, heavy design	20	480
633 180	Petri dish, Ø 94 x 16 mm, with vents	20	480
633 185	Petri dish, Ø 94 x 15 mm, with vents	20	480
663 102	Petri dish, Ø 100 x 15 mm, with vents, heavy design	20	420
664 102	Petri dish, Ø 100 x 20 mm, with vents, heavy design	15	360
639 102	Petri dish, Ø 145 x 20 mm, with vents	15	120

Contact Dishes / Square Petri Dish / OneWell Plate™			
629 161	Contact dish, Ø 65 x 15 mm, sterile, graduation with 10 mm squares, without vents	20	600
629 180	Contact dish, Ø 65 x 15 mm, sterile, graduation with 10 mm squares, with vents	20	600
670 102	OneWell Plate™, microplate format 127.8 x 85.5 x 14.4 mm, with vents	8	32
688 102	Petri dish, square, 120 x 120 x 17 mm, with vents	10	240

Petri Dishes - Special Models			
633 175	Germ count dish, Ø 94 x 16 mm, graduated, with vents	20	480
635 102	Petri dish, Ø 94 x 15 mm, 2 compartments, 20 ml volume/well, with vents	20	480
637 102	Petri dish, Ø 94 x 15 mm, 3 compartments, 12 ml volumen/well, with vents	20	480
657 102	Multiwell plate with lid, 6 well, 16 ml volume/well, with vents, numeric coding	2	100