

LABORHITS

Immersion Circulator CORIO™

Precise temperature control

The CORIO™ C Immersion Circulator is the basic model of the CORIO™ circulator portfolio. The bath attachment clamp is included in delivery and facilitates mounting of the circulator on any bath tank up to 30 liters. With a few simple steps you can install a pump set onto the CORIO™ CD Immersion Circulator. In no time at all, your circulator is ready for temperature control of an external application.

- The jet nozzle permits continuous adjustment of the pump stream in your bath
- Bright, white, easy to read display
- Whisper quiet
- Easy operation



Please contact us! We'll be pleased to support you.

	CORIO™ C	CORIO™ CD
Description	For bath tank up to 30 liters	For bath tank up to 50 liters
	Class I (NFL) according to DIN 12876-1	Class III (FL) according to DIN 12876-1
	Circulation capacity 6 l / min 0.1 bar	Circulation capacity 15 l / min 0.35 bar
	for internal circulation	Pump change-over between internal and external circulation
	Immersion depth 7.5 ... 16.5 cm	Immersion depth 7.5 ... 16.5 cm
		USB interface
Working Temp. range °C	+20 ... +100	+20 ... +150
Temp. stability ± °C	0,03	0,03
Heating capacity kW	2	2
PK	1	1
Cat. No.	7.673 349	4.658 066

Delivery without bath please order separately.

To adapt the CORIO™ Immersion Circulator to your individual application we offer a comprehensive range of accessories (racks, bath lids, and more) or the Pump set for CORIO™ CD with connectors M16x1, Art. No. 4.658 027