

Test weights CLASS E2

Weight sets, compact shape, polished stainless steel



- **Test weight material:** Polished stainless steel
- **Case material:** Lined plastic

Weight set		+	DKD certificate		=	Package price	
KERN			KERN			KERN	
312-024	1 g - 50 g		962-315			312-024-600	
312-034	1 g - 100 g		962-316			312-034-600	
312-044	1 g - 200 g		962-317			312-044-600	
312-054	1 g - 500 g		962-318			312-054-600	
312-064	1 g - 1 kg		962-319			312-064-600	
312-074	1 g - 2 kg		962-320			312-074-600	
312-084	1 g - 5 kg		962-321			312-084-600	

Weight sets, cylindrical shape, polished stainless steel



- **Test weight material:** Individual weights - polished stainless steel, milligram weights - aluminium/German silver
- **Case material:** Lined plastic. Milligram weights 1 mg - 500 mg in removable plastic box

Weight set		+	DKD certificate		=	Package price	
KERN			KERN			KERN	
318-22	1 mg - 500 mg		962-350			318-221-600	
313-024	1 mg - 50 g		962-301			313-024-600	
313-034	1 mg - 100 g		962-302			313-034-600	
313-044	1 mg - 200 g		962-303			313-044-600	
313-054	1 mg - 500 g		962-304			313-054-600	
313-064	1 mg - 1 kg		962-305			313-064-600	
313-074	1 mg - 2 kg		962-306			313-074-600	
313-084	1 mg - 5 kg		962-307			313-084-600	
314-024	1 g - 50 g		962-315			314-024-600	
314-034	1 g - 100 g		962-316			314-034-600	
314-044	1 g - 200 g		962-317			314-044-600	
314-054	1 g - 500 g		962-318			314-054-600	
314-064	1 g - 1 kg		962-319			314-064-600	
314-074	1 g - 2 kg		962-320			314-074-600	
314-084	1 g - 5 kg		962-321			314-084-600	

Weight sets, cylindrical shape, polished stainless steel



- **Test weight material:** Individual weights - polished stainless steel, milligram weights - aluminium/German silver
- **Box material:** Lined wood. Milligram weights 1 mg - 500 mg in removable plastic box

Weight set		+	DKD certificate		=	Package price	
KERN			KERN			KERN	
318-22	1 mg - 500 mg		962-350			318-221-600	
313-02	1 mg - 50 g		962-301			313-021-600	
313-03	1 mg - 100 g		962-302			313-031-600	
313-04	1 mg - 200 g		962-303			313-041-600	
313-05	1 mg - 500 g		962-304			313-051-600	
313-06	1 mg - 1 kg		962-305			313-061-600	
313-07	1 mg - 2 kg		962-306			313-071-600	
313-08	1 mg - 5 kg		962-307			313-081-600	
313-09	1 mg - 10 kg		962-308			313-091-600	
314-02	1 g - 50 g		962-315			314-021-600	
314-03	1 g - 100 g		962-316			314-031-600	
314-04	1 g - 200 g		962-317			314-041-600	
314-05	1 g - 500 g		962-318			314-051-600	
314-06	1 g - 1 kg		962-319			314-061-600	
314-07	1 g - 2 kg		962-320			314-071-600	
314-08	1 g - 5 kg		962-321			314-081-600	
314-09	1 g - 10 kg		962-322			314-091-600	