



Figure 1 — General design of glass filter funnel
(plain or analytical)

Table 1 — Dimensions for plain filter funnels

Dimensions in millimetres

Bowl		Stem		
Internal diameter	Wall thickness min.	Nominal		Wall thickness min.
		External diameter	Length	
35 ± 5	1	6	35	0,8
55 ± 5	1	8	55	1,3
75 ± 5	1	8	75	1,3
100 ± 5	1,2	10	100	1,3
125 ± 5	1,2	16	125	1,3
150 ± 10	1,2	16	150	1,6
200 ± 10	1,5	24	150	1,6