

LABORHITS

Universal drying oven LLG-uniOVEN 42 and LLG-uniOVEN 110



NEW

42 liter or 110 liter volume

**Quick and intuitive
operation**

Features

- Microprocessor controlled universal drying oven with forced convection
- Digital display for time and temperature
- Timer: 1 minute to 99.9 h and continuous
- Delayed start timer: 1 minute to 99.9 h
- Acoustical and visible information for end of program
- Acoustical and optical alarm at temperature differences
- Outlet for vapour exhaustion
- Interior made of stainless steel
- Adjustable over temperature protection up to 300 °C

**WARRANTY
3
YEARS**



Specifications

Internal volume: 42 l or 110 l
Working temperature: +5 °C up to +40 °C
Temperature range: RT +7 °C up to 250 °C
Temperature accuracy: ± 1 °C (at 70 °C)
Temperature uniformity: ± 2 °C (at 70 °C)
Maximum humidity: 80 %, non-condensing
Power supply: 230 V, 50/60 Hz

Scope of supply: LLG-uniOVEN 42 or 110 incl. 2 shelves

Type	Power consumption W	Internal dimensions mm (W x D x H)	External dimensions mm (W x D x H)	Weight kg	PK	Cat. No.
LLG-uniOVEN 42	850	420 x 280 x 360	705 x 505 x 540	39	1	6.263 670
LLG-uniOVEN 110	1550	500 x 445 x 500	790 x 665 x 680	66	1	6.263 675

LABORHITS

Universal drying LLG-uniOVEN 150



LLG-uniOVEN 115

150 liter volume

Quick and intuitive operation



Specifications

Internal volume: 150 l

Working temperature: +5 °C up to +35 °C

Temperature range: RT +10 °C up to 300 °C

Temperature accuracy: ± 0,5 °C (at 70 °C)

Temperature uniformity: ± 2 °C (at 70 °C)

Maximum humidity: 80 %, non-condensing

Power supply: 230 V, 50/60 Hz

Features

- Microprocessor controlled universal drying oven with forced convection
- Integrated auto-diagnostic system
- 5 different programs with 5 steps can be programmed by the user
- Controller can be protected by password
- Lockable door handle
- LED display
- Interior made of stainless steel (AISI-304)
- Fan speed adjustable
- Hot indicator for temperatures above 50 °C
- Adjustable over temperature protection

Type	Power consumption W	Internal dimensions mm (W x D x H)	External dimensions mm (W x D x H)	Weight kg	PK	Cat. No.
LLG-uniOVEN 150, with EU plug	3000	495 x 440 x 705	700 x 665 x 1010	85	1	6.263 680
LLG-uniOVEN 150, with UK plug	3000	495 x 440 x 705	700 x 665 x 1010	85	1	6.263 682