

Data Sheet

Page 1/1



CHROMABOND columns C18 Hydra volume: 3 mL, content of sorbent: 200 mg material: glass, with glass fibre filter element pack of 50

Part Number: MN-730296G
Manufacturer: Macherey-Nagel
Quantity: 1 pc

Description

CHROMABOND columns C18 Hydra volume: 3 mL, content of sorbent: 200 mg material: glass, with glass fibre filter element pack of 50

Technical Information

Brand:	CHROMABOND	pH Minimum:	2.0
Product Type:	SPE Cartridges	pH Maximum:	8.0
Particle Size:	45 µm	Mode of Separation:	Reversed Phase Chromatography (RP)
Chemistry:	C18 (RP-18, ODS, Octadecyl)	Matrix:	Silica
Chemistry Producers Term:	C18 Hydra	Particle Morphology:	fully porous
Pore Size:	60 Å		
Carbon Load:	15.0 %		
Bet Surface Area [m²/g]:	500.0		
Endcapping:	non endcapped		
Country of Origin:	Germany		
Volume:	3 ml		
Bed Mass [mg]:	200 mg		