

Product Range BRAND

Liquid Handling · Life Science Consumables · Pipetting Robots
Volumetric Instruments · General Lab Products



BRAND. For lab. For life.®





accu-jet®, BIO-CERT®, BLAUBRAND®, BRAND®, BRAND. For lab. For life.®, BRANDplates®, cell-culture™, cellGrade™, Dispensette®, EASYCAL™, HandyStep®, hydroGrade™, immunoGrade™, inertGrade™, lipoGrade™, nano-cap™, pureGrade™, QuikSip™, SafetyPrime™, seripettor®, Titrette®, Transferpette®, as well as the the figurative marks depicted here and the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany.

All other trademarks mentioned or depicted here are the property of the respective owners.

Contents

	General Information Dispensette® S · Titrette® · HandyStep® touch · Transferpette® S PCR products	8
	Liquid Handling Dispensing · Titrating · Pipetting	14
	Pipetting Robots Liquid Handling Station	48
	Life Science Consumables Pipette Tips · Micro Tubes · PCR · Microplates · Cuvettes, etc.	51
	Volumetric Instruments Bulb Pipettes · Graduated Pipettes · Volumetric Flasks · Burettes, etc.	77
	General Lab Products Beakers · Funnels · Connectors · Bottles · Detergents, etc.	113

Dispensette® S

Innovative ideas with trusted technology – the new bottle-top dispenser Dispensette® S.

- New discharge tube with or without recirculation valve
- New valve system, no sealing rings necessary
- Faster priming due to improved flow technology
- Less force needed during dispensing especially for instruments with large volumes
- Volume selection with interior scalloped track for analog devices, enhances setting reproducibility
- New 1 ml size digital and analog

The bottle-top dispenser Dispensette® S has all the features that make dispensing safer and convenient.
Innovative ideas – trusted technology.

Easy Calibration Technique

Calibration and adjustments according to ISO 9001 and GLP are done within seconds.

NEW!

Discharge valve with safety ball

closes when discharge tube is not mounted

NEW!

Discharge tube

without recirculation valve

NEW!

Large

viewport

NEW!

Hinged screw cap

doesn't get in the way when dispensing

Dispensette® S Organic Digital

Calibration mechanism

Dispensette® S Analog-adjustable

NEW!

Volume selection

using interior scalloped track

NEW!

Discharge tube

with recirculation valve

360° rotating

valve block with GL 45 thread

NEW!

Filling valve olive-shaped

for firmer attachment of the filling tube

Recirculation tube

for discharge tubes with recirculation valve

Telescoping

filling tube

Dispenser Selection Chart

Reagent	Dispensette® S	Dispensette® S Organic
Acetaldehyde	+	+
Acetic acid (glacial), 100%	+	+
Acetic acid, ≤ 96%	+	+
Acetic anhydride		+
Acetone	+	+
Acetonitrile	+	+
Acetophenone		+
Acetyl chloride		+
Acetylacetone	+	+
Acrylic acid	+	+
Acrylonitrile	+	+
Adipic acid	+	
Allyl alcohol	+	+
Aluminium chloride	+	
Amino acids	+	
Ammonia, ≤ 20%	+	+
Ammonia, 20-30%		+
Ammonium chloride	+	
Ammonium fluoride	+	
Ammonium sulfate	+	
n-Amyl acetate	+	+
Amyl alcohol (Pentanol)	+	+
Amyl chloride (Chloropentane)		+
Aniline	+	+
Barium chloride	+	
Benzaldehyde	+	+
Benzene (Benzol)	+	+
Benzine (Petroleum benzin), bp 70-180 °C		+
Benzoyl chloride	+	+
Benzyl alcohol	+	+
Benzylamine	+	+
Benzylchloride	+	+
Boric acid, ≤ 10%	+	+
Bromobenzene	+	+
Bromonaphthalene	+	+
Butanediol	+	+
1-Butanol	+	+
n-Butyl acetate	+	+
Butyl methyl ether	+	+
Butylamine	+	+
Butyric acid	+	+
Calcium carbonate	+	
Calcium chloride	+	
Calcium hydroxide	+	
Calcium hypochlorite	+	
Carbon tetrachloride		+
Chloro naphthalene	+	+
Chloroacetaldehyde, ≤ 45%	+	+
Chloroacetic acid	+	+
Chloroacetone	+	+
Chlorobenzene	+	+
Chlorobutane	+	+
Chloroform		+
Chlorosulfonic acid		+
Chromic acid, ≤ 50%	+	+
Chromosulfuric acid	+	
Copper sulfate	+	
Cresol		+
Cumene (Isopropyl benzene)	+	+

Reagent	Dispensette® S	Dispensette® S Organic
Cyclohexane		+
Cyclohexanone	+	+
Cyclopentane		+
Decane	+	+
1-Decanol	+	+
Dibenzyl ether	+	+
Dichloroacetic acid		+
Dichlorobenzene	+	+
Dichloroethane		+
Dichloroethylene		+
Dichloromethane		+
Diesel oil (Heating oil), bp 250-350 °C		+
Diethanolamine	+	+
Diethyl ether		+
Diethylamine	+	+
1,2 Diethylbenzene	+	+
Diethylene glycol	+	+
Dimethyl sulfoxide (DMSO)	+	+
Dimethylaniline	+	
Dimethylformamide (DMF)	+	+
1,4 Dioxane		+
Diphenyl ether	+	+
Essential oil		+
Ethanol	+	+
Ethanolamine	+	+
Ethyl acetate	+	+
Ethylbenzene		+
Ethylene chloride		+
Fluoroacetic acid		+
Formaldehyde, ≤ 40%	+	
Formamide	+	+
Formic acid, ≤ 100%		+
Glycerol	+	+
Glycol (Ethylene glycol)	+	+
Glycolic acid, ≤ 50%	+	
Heating oil (Diesel oil), bp 250-350 °C		+
Heptane		+
Hexane		+
Hexanoic acid	+	+
Hexanol	+	+
Hydriodic acid, ≤ 57% **	+	+
Hydrobromic acid		+
Hydrochloric acid, ≤ 20%	+	+
Hydrochloric acid, 20-37% **		+
Hydrogen peroxide, ≤ 35%		+
Isoamyl alcohol	+	+
Isobutanol	+	+
Isooctane		+
Isopropanol (2-Propanol)	+	+
Isopropyl ether	+	+
Lactic acid	+	
Methanol	+	+
Methoxybenzene	+	+
Methyl benzoate	+	+
Methyl butyl ether	+	+
Methyl ethyl ketone	+	+
Methyl formate	+	+
Methyl propyl ketone	+	+

Reagent	Dispensette® S	Dispensette® S Organic
Methylene chloride		+
Mineral oil (Engine oil)	+	+
Monochloroacetic acid	+	+
Nitric acid, ≤ 30%	+	+
Nitric acid, 30-70% */**		+
Nitrobenzene	+	+
Oleic acid	+	+
Oxalic acid	+	
n-Pentane		+
Peracetic acid		+
Perchloric acid	+	+
Perchloroethylene		+
Petroleum, bp 180-220 °C		+
Petroleum ether, bp 40-70 °C		+
Phenol	+	+
Phenylethanol	+	+
Phenylhydrazine	+	+
Phosphoric acid, ≤ 85%	+	+
Phosphoric acid, 85% + Sulfuric acid, 98%, 1:1	+	+
Piperidine	+	+
Potassium chloride	+	
Potassium dichromate	+	
Potassium hydroxide	+	
Potassium permanganate	+	
Propionic acid	+	+
Propylene glycol (Propanediol)	+	+
Pyridine	+	+
Pyruvic acid	+	+
Salicylaldehyde	+	+
Scintillation fluid	+	+
Silver acetate	+	
Silver nitrate	+	
Sodium acetate	+	
Sodium chloride	+	
Sodium dichromate	+	
Sodium fluoride	+	
Sodium hydroxide, ≤ 30%	+	
Sodium hypochlorite	+	
Sulfuric acid, ≤ 98%	+	+
Tartaric acid	+	
Tetrachloroethylene		+
Tetrahydrofuran (THF) */**		+
Tetramethylammonium hydroxide	+	
Toluene		+
Trichloroacetic acid		+
Trichlorobenzene		+
Trichloroethane		+
Trichloroethylene		+
Trichlorotrifluoro ethane		+
Triethanolamine	+	+
Triethylene glycol	+	+
Trifluoro ethane		+
Trifluoroacetic acid (TFA)		+
Turpentine		+
Urea	+	
Xylene		+
Zinc chloride, ≤ 10%	+	
Zinc sulfate, ≤ 10%	+	

The above recommendations reflect testing completed prior to publication. Always follow instructions in the operating manual of the instrument as well as the reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g., biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact BRAND. Status as of: 1116/13

* use ETFE/PTFE bottle adapter


** use PTFE seal for valve block

Note! For dispensing HF, we recommend the use of the Dispensette® S Trace Analysis bottle-top dispenser with platinum-iridium valve spring.



Titrette®

The control elements of the Titrette® bottle-top burette have an intuitive layout. Separate buttons for On/Off and Pause; Clear button allows user to reset the display and select functions. Easy-grip hand wheel action and smooth precision gears for fast or drop-wise titration make handling extremely simple and reliable.



Comparison of error limits

		Titrette® bottle-top burette				Bottle-top burettes acc. to DIN EN ISO 8655-3				Glass burettes Class A acc. to DIN EN ISO 385 and ASTM	
Volume ml	Partial volume, ml	A* ≤ ± %	μl	CV* ≤ %	μl	A* ≤ ± %	μl	CV* ≤ %	μl	EL** ± μl	
10	10	0,10	10	0,05	5	0,3	30	0,1	10	20	
	5	0,20	10	0,10	5	0,6	30	0,2	10	20	
	1	1,00	10	0,50	5	3	30	1	10	20	
25	25	0.07	18	0.025	6	0.2	50	0.1	25	30	
	12.5	0.14	18	0.05	6	0.4	50	0.2	25	30	
	2.5	0.70	18	0.25	6	2	50	1	25	30	
50	50	0.06	30	0.02	10	0.2	100	0.1	50	50	
	25	0.12	30	0.04	10	0.4	100	0.2	50	50	
	5	0.60	30	0.20	10	2	100	1	50	50	

* Error limits related to the nominal capacity (= maximum volume) indicated on the instrument, obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation.

** Error limit: EL = A + 2CV, according to DIN EN ISO 8655-6 Annex B
(A = accuracy, CV = coefficient of variation, EL = error limit)

The titration volume is displayed in steps of 1 μl at instruments with 10 ml and 25 ml size and in steps of 2 μl for 50 ml size instruments. For titration volumes above 20 ml the display will automatically switch to steps of 10 μl.

Note:

If you need an official certification which confirms the error limits that are much stricter than those of DIN EN ISO 8655-3, we recommend a calibration certificate from an accredited calibration laboratory (e.g., the DAKS laboratory at BRAND).

Material and reagents

The instrument can be used for the following titration media (maximum concentration 1 mol/l):

Acetic acid	Iron (II) sulfate solution	Potassium thiocyanate solution
Alcoholic potassium hydroxide solution	Nitric acid	Silver nitrate solution*
Ammonium iron (II) sulfate solution	Oxalic acid solution	Sodium arsenite solution
Ammonium thiocyanate solution	Perchloric acid	Sodium carbonate solution
Barium chloride solution	Perchloric acid in glacial acetic acid	Sodium chloride solution
Bromide bromate solution	Potassium bromate solution	Sodium hydroxide solution
Cerium (IV) sulfate solution	Potassium bromide bromate solution	Sodium nitrite solution
EDTA solution	Potassium dichromate solution	Sodium thiosulfate solution
Hydrochloric acid	Potassium hydroxide solution	Sulfuric acid
Hydrochloric acid in Acetone	Potassium iodate solution	Tetra-n-butylammonium hydroxide solution
Iodide Iodate solution*	Potassium permanganate solution*	Triethanolamine in Acetone*
Iodine solution*		Zinc sulfate solution

* Use light shield inspection window

When the instrument is properly handled, dispensed liquid will only come into contact with the following chemically resistant materials: borosilicate glass, Al₂O₃, ETFE, PFA, FEP, PTFE, platinum-iridium; PP (screw cap).

Limitations of use

Chlorinated and fluorinated hydrocarbons or chemical combinations which form deposits may make the piston difficult to move or may cause jamming.

Compatibility of the instrument for a special application (e.g., trace material analysis) must be checked by the user. For additional information, please contact the manufacturer.

The instrument is not autoclavable!

Operating limits

This instrument is designed for titrating liquids, observing the following physical limits:

- +15 °C to +40 °C (59 °F to 104 °F) of instrument and reagent
- Vapor pressure up to 500 mbar
- Viscosity up to 500 mm²/s
- Altitude: maximum 3000 m above sea level
- Relative humidity: 20% to 90%



12.796

HandyStep® touch · HandyStep® touch S

Overview

- Touchscreen operation with intuitive menu navigation and integrated help function
- Two models available: HandyStep® touch or HandyStep® touch S with expanded functionality
- Menu languages: German, English, Spanish, French and Chinese
- Variable volume adjustment from 1.0 µl to 50 ml
- Inductive charging with optional charging stand

Versatile

HandyStep® touch with the standard modes Multi-Dispensing, Auto-Dispensing, Pipetting and Favorites; HandyStep® touch S with the additional functions Sequential Dispensing, Multi-Aspiration and Titration

Long-lasting

Many hours of continuous operation with a single battery charge

Robust

High-quality materials provide long service life, and smooth surfaces make cleaning/disinfection easy

Easy

Automatic size detection of BRAND PD-Tips //

Efficient

In addition to BRAND PD-Tips, compatible third party tips can also be used.

Simple

Automatic tip ejection by the press of a button prevents contamination

USB-C connection
For quick charging by cable

STEP button
Ergonomically arranged STEP button for effortless volume filling and delivery

Convenient

Inductive charging with the optional charging stand



Intuitive operation

Be flexible and have the right solution at hand for every application. Even complex procedures are easy to carry out – no problem with the expanded functionality of the HandyStep® touch S.

Mode	Function	HandyStep® touch	HandyStep® touch S
Multi Dispensing (MULTI-DISP)	Successive dispensing of a set volume	✓	✓
Auto Dispensing (AUTO-DISP)	Automatic dispensing by the instrument	✓	✓
Pipetting (PIP)	Pipetting a defined volume	✓	✓
Sequential Dispensing (SEQ-DISP)	Successive dispensing of different volumes	—	✓
Multi Aspirating (MULTI-ASP)	Aspirating different volumes	—	✓
Titration	Manual or successive delivery of liquids for titration	—	✓
Favorites	Saving favorites for quick access	✓	✓
Settings	Individual customization of the device (e.g. brightness or device name)	✓	✓



Expanded functionality on the HandyStep® touch S:



The clearly structured menus make fast work with intuitive operation possible, while the most important functions and parameters are always in view.



Easy to understand

Main menu with all modes for quick selection



Clearly structured

Working mode with all important information



Organized

All options visible at a glance



Helpful

Internal help right where it is needed



The perfect team

Automatic size detection with PD-Tips // from BRAND

The HandyStep® touch saves time and prevents errors through automatic tip size recognition of the PD-Tips // from BRAND, with patented coding on the pistons. After inserting the tip, the size is automatically recognized and displayed, making it easy to select the volume to be dispensed. When a new PD-Tip of the same size is inserted, all instrument settings are maintained.




Transferpette® S

Whether you are using the single-channel pipette or the 8- and 12-channel models, the Transferpette® S is always the first choice for simple and efficient pipetting in demanding applications.

The single-channel pipettes cover a wide volume range from 0.1 µl to 10 ml. Pipettes with small nominal volumes provide high precision for small-volume applications such as those found in molecular biology, e.g. for pipetting enzymes.

Large nominal volumes are used for water and environmental analysis, for example, where durability is critical.

For standardized applications with specific volumes, we offer economical fixed volume models.



Easy Calibration technology



Volume-change protection



Integrated shaft coupling



Label on the finger rest

Transferpette® S -8 and S -12


Work in microplates and PCR plates is efficient with the multi-channel models of the Transferpette® S.

These 8- and 12-channel pipettes cover the volume range

from 0.5 µl to 300 µl. They provide all the advantages of the single-channel models, plus additional features that make work and maintenance easier.



Transferpette® S -12



Transferpette® S -8

Adapts to you

Pipetting unit can be rotated 360° for individualized comfort

Low forces

Easy mounting and ejection of tips due to O-rings

Easy cleaning and maintenance

Individual nose cones and seals can be easily removed* and for easy cleaning or replacement

* patented



Narrow centrifuge tubes



Sample preparation



Working with PCR plates



Easy maintenance

PCR plates Thermal Cycler Compatibility

One glance at the table will tell you which BRAND PCR plates are compatible with your thermal cycler. We are continuously updating the table with information from manufacturers and feedback from our customers. Ask for a free sample of our PCR plates to check compatibility with your thermal cycler (www.brand.de), with no obligation.

	non-skirted										semi-skirted						r.s.-s.*		skirted									
	7814 11	7814 12	7814 15	7814 16	7813 66	7813 67	7813 68	7813 69	7813 50	7813 54	7814 00	7813 57	7813 64	7813 71	7813 72	7813 75	7813 76	7813 73	7813 74	7813 77	7813 78	7813 45	7813 47	7813 58				
White color		x		x		x		x		x		x	x		x		x		x		x			x				
Number of wells	24		48		96		96		96		96	96	96	96	96	96	96	96	96	96	384	384	384					
Standard well	x		x				x		x		x					x												
Low profile well					x				H12				x		x			x		x								
Cut corner					H12		A12		H12		H12	H12	A12	A12	A12	A12	A1	A1	H1	H1	A24 + P24	A24	A24 + P24					
Agilent Technologies																												
AriaMx																												
Analytik Jena																												
Flex Cycler² 96 Series																												
Flex Cycler² Twin 48																												
Speed Cycler² SPR																												
qTower 2.0 und 2.2																												
APPLIED BIOSYSTEMS®																												
2700/2720	•		•				•		•		•						•				•		•		•			
3100	•		•				•		•								•						•					
3130	•		•				•		•								•						•					
3500							•		•								•						•					
3700	•		•				•		•								•						•					
3730/ 3730x	•		•				•		•								•						•					
9600	•		•				•		•		•						•						•					
9700	•		•				•		•		•						•						•		•			
9800 Fast																			•				•		•			
5700	•		•				•		•								•											
7000	•		•				•		•						•		•											
7300	•		•				•		•				•				•											
7500	•		•		•		•		•								•											
7700	•		•				•		•						•		•											
7900 HT							•								•		•						•					
7500 Fast																	•											
Quant Studio™ 6																	•											
Quant Studio™ 96-well																	•				•							
Quant Studio™ 96-well fast					•												•											
Step One Plus																	•											
Veriti 0.2 ml															•		•											
Veriti 0.1 ml																	•											
ViiA™ 7					•												•											
																	•											
AMERSHAM® Bioscience																												
MegaBace® 500																							•					
MegaBace® 1000																							•					
MegaBace® 4000																							•					
BECKMAN®																												
CEQ																							•					
MegaBace® 4000																							•					
BIOMETRA® Analytik Jena AG																												
Uno	•		•		•		•		•		•						•				•		•		•			
Uno II	•		•		•		•		•		•						•				•		•		•			
T1 Thermal Cycler	•		•		•		•		•		•						•				•		•		•			
Tgradient	•		•		•		•		•		•						•				•		•		•			
Trobot	•		•		•		•		•		•						•				•		•		•			
BIORAD®/MJ RESEARCH®																												
CFX 96 Touch™											—										—		•					
CFX 96 Connect™											—										—		•					
CFX Automated System II																							•					
T100	•		•				•		•		•						•				•		•		•			
Genecycler																							•					
C1000/ C1000 touch	•		•		•		•		•		•						•				•		•		•			
S1000	•		•		•		•		•		•						•				•		•		•			
PTC-200 DNA-Engine	•		•		•		•		•		•						•				•		•		•			
PTC-220 DNA-Eng. Dyad	•		•		•		•		•		•						•				•		•		•			
PTC-221 Dyad-Disciple	•		•		•		•		•		•						•				•		•		•			
PTC-225 DNA-Tetrad	•		•		•		•		•		•						•				•		•		•			
PTC-240 DNA-Tetrad 2																							•					
Base Station																							•					
iCycler	•		•				•		•								•				•							
MyCycler	•		•				•		•								•				•							
MyiQ	•		•				•		•								•				•				•			
iQ4	•		•				•		•								•				•				•			
iQ5	•		•				•		•								•				•				•			
CFX96																							•					
Chromo 4																							•					
Opticon																							•					
Opticon 2																							•					
Opticon 4																							•					
Corbett Research																												
Palm Cycler 96																							•					
Palm Cycler 384																							•					
Eppendorf®																												
Mastercycler® Gradient	•		•				•		•		•						•				•		•		•			
Mastercycler® Nexus	•		•				•		•		•						•				•		•		•			
Mastercycler® ep	•		•		•		•		•								•				•		•		•			
Mastercycler®					•		•		•								•				•		•		•			
Mastercycler® M384																							•					
Mastercycler® ep Realplex																							•					

* r. s.-s. = raised skirt semi-skirted

	non-skirted										semi-skirted						r.s.-s.*		skirted							
	7814 11	7814 12	7814 15	7814 16	7813 66	7813 67	7813 68	7813 69	7813 50	7813 54	7814 00	7813 57	7813 64	7813 71	7813 72	7813 75	7813 76	7813 73	7813 74	7813 77	7813 78	7813 45	7813 47	7813 58		
White color		x		x		x		x		x		x		x		x		x		x					x	
Number of wells	24		48		96		96		96		96	96	96	96	96	96	96	96	96	96	384	384	384	384		
Standard well	x		x				x		x		x					x										
Low profile well					x								x	x				x		x						
Cut corner					H12		A12		H12		H12	H12	A12	A12	A12	A12	A12	A1		H1	A24 + P24	A24	A24	A24 + P24		
Ericomp																										
Single Block	●		●		●		●		●								●									
Twin Block	●		●		●		●		●								●									
Power Block					●		●		●						●		●									
Delta Cyclor	●		●		●		●		●						●		●									
Esco																										
Swift							●																		●	
Flexi																										
Gene							●														●				●	
Genius							●														●					
G-Storm																										
GS1/ GS4/ GSX					●		●								●		●									
Labnet																										
MultiGene™ OptiMax							●										●									
Thermo Hybaid																										
Multi Block System	●		●		●		●		●		●				●		●		●		●		●		●	
Omnigene	●		●		●		●		●		●				●		●		●		●		●		●	
Omn-E	●		●		●		●		●		●				●		●		●		●		●		●	
PCR-Express	●		●		●		●		●		●				●		●		●		●		●		●	
PCR-Sprint	-		-		●		●		-		-				●		●		●		●		●		●	
pxe	●		●		●		●		●						●		●		●		●		●		●	
px2	●		●		●		●		●						●		●		●		●		●		●	
Touchdown	●		●		●		●		●		●				●		●		●		●		●		●	
MWG®																										
Primus 96	●		●		●		●		●		●				●		●		●		●		●		●	
Primus 384																					●		●		●	
The Q LifeCycler					●		●														●		●		●	
PCRmax																										
Alpha Cycler 1					●		●										●		-		●					
Alpha Cycler 2					●		●										●		-		●					
Alpha Cycler 4					●		●										●		-		●					
PEQLAB®																										
peqStar 96	●		●		●		●		●		●		●		●		●		-		●					
peqSar 96X Universal	●		●		●		●		●		●		●		●		●		●		●					
peqSar 96X Univ. Gradient	●		●		●		●		●		●		●		●		●		●		●					
peqStar 96X HPL	●		●		●		●		●		●		●		●		●		●		●					
peqStar 96X HPLGradient	●		●		●		●		●		●		●		●		●		●		●					
peqStar 384X HPL																							●		●	
Roche®																										
LightCycler® 480													●												●	
STRATAGENE®																										
Robocycler	●		●				●		●		-				●		●				●		●			
Robocycler Gradient					●		●						●		●						●					
Mx 3000	●		●						●																	
Mx 3000 P					●		●								●		●									
Mx 3005 P					●		●								●		●									
Mx 4000	●		●		●		●		●																	
TaKaRa®																										
TP2400																					●		●		●	
TP3000	●		●		●		●		●						●		●				●		●		●	
TECHNE®																										
TC-412	●		●		●		●		●		●				●		●				●		●		●	
Flexigene	●		●		●		●		●		●				●		●				●		●		●	
Genius	●		●						●		●				●		●				●		●		●	
Genius Quad					●		●						●		●		●				●		●		●	
TC-512	●		●		●		●		●		●				●		●				●		●		●	
Touchgene Gradient									●		●				●		●				●		●		●	
Touchgene X											-												-		-	
TC-3000X	●		●						●		-												-		-	
Quantica																										
Cyclogene					●		●								●		●				●		●		●	
TRANSGENOMIC®																										
Wave System																					●				●	
VWR																										
XT96	●		●		●		●		●		●		●		●		●		●		-					
UNO ⁹⁶	●		●		●		●		●		●		●		●		●		●		●					
UNO ⁹⁶ Gradient	●		●		●		●		●		●		●		●		●		●		●					
UNO ⁹⁶ HPL	●		●		●		●		●		●		●		●		●		●		●					
UNO ⁹⁶ HPL Gradient	●		●		●		●		●		●		●		●		●		●		●					
UNO ³⁸⁴ HPL																							●		●	

* r. s.-s. = raised skirt semi-skirted

• = compatible

- = not compatible

□ = no information

Stand: 0519

The following brands and trademarks as referenced in the current table are property of their respective third-party owners as listed below:

Referenced brand/trademark	Owner	Referenced brand/trademark	Owner
Agilent®	Agilent Technologies, Inc., USA	MegaBace®	GE Healthcare Bio-Sciences AB, GB
AMERSHAM®	GE HEALTHCARE LIMITED, GB	MJ RESEARCH®	BIO-RAD LABORATORIES, INC., USA
APPLIED BIOSYSTEMS®	Applied Biosystems LLC, USA	MWG®	MWG BIOTECH AG, Deutschland
AriaMx®	Agilent Technologies, Inc., USA	PEQLAB®	PEQLAB Biotechnologie GmbH, Deutschland
BECKMAN®	Beckman Coulter, Inc., USA	Roche®	F. Hoffmann-La Roche AG, Schweiz
BIOMETRA®	Biometra Biomedizinische Analytik GmbH, Deutschland	STRATAGENE®	Agilent Technologies, Inc., USA
BIO-RAD®	BIO-RAD LABORATORIES GMBH, Deutschland	TaKaRa®	TAKARA SHUZO CO., LTD., Japan
Eppendorf®	Eppendorf AG, Deutschland	TECHNE®	TECHNE INCORPORATED, USA
LightCycler®	Roche Diagnostics GmbH, Deutschland	TRANSGENOMIC®	Transgenomic, Inc., USA
Mastecycler®	Eppendorf AG, Deutschland		

Bottle-top dispenser Dispensette® S, Digital, DE-M



Cat. No.	Description	From*	Unit
4600310	Dispensette® S, Digital, DE-M, 0,1-1ml, without recirculation valve	1	piece(s)
4600311	Dispensette® S, Digital, DE-M, 0,1-1ml, with recirculation valve	1	piece(s)
4600320	Dispensette® S, Digital, DE-M, 0,2-2ml, without recirculation valve	1	piece(s)
4600321	Dispensette® S, Digital, DE-M, 0,2-2ml, with recirculation valve	1	piece(s)
4600330	Dispensette® S, Digital, DE-M, 0,5-5ml, without recirculation valve	1	piece(s)
4600331	Dispensette® S, Digital, DE-M, 0,5-5ml, with recirculation valve	1	piece(s)
4600340	Dispensette® S, Digital, DE-M, 1-10ml, without recirculation valve	1	piece(s)
4600341	Dispensette® S, Digital, DE-M, 1-10ml, with recirculation valve	1	piece(s)
4600350	Dispensette® S, Digital, DE-M, 2,5-25ml, without recirculation valve	1	piece(s)
4600351	Dispensette® S, Digital, DE-M, 2,5-25ml, with recirculation valve	1	piece(s)
4600360	Dispensette® S, Digital, DE-M, 5-50ml, without recirculation valve	1	piece(s)
4600361	Dispensette® S, Digital, DE-M, 5-50ml, with recirculation valve	1	piece(s)

Bottle-top dispenser Dispensette® S, Analog-adjustable, DE-M



Cat. No.	Description	From*	Unit
4600100	Dispensette® S, Analog-adj., DE-M, 0,1-1ml, without recirculation valve	1	piece(s)
4600101	Dispensette® S, Analog-adj., DE-M, 0,1-1ml, with recirculation valve	1	piece(s)
4600120	Dispensette® S, Analog-adj., DE-M, 0,2-2ml, without recirculation valve	1	piece(s)
4600121	Dispensette® S, Analog-adj., DE-M, 0,2-2ml, with recirculation valve	1	piece(s)
4600130	Dispensette® S, Analog-adj., DE-M, 0,5-5ml, without recirculation valve	1	piece(s)
4600131	Dispensette® S, Analog-adj., DE-M, 0,5-5ml, with recirculation valve	1	piece(s)
4600140	Dispensette® S, Analog-adj., DE-M, 1-10ml, without recirculation valve	1	piece(s)
4600141	Dispensette® S, Analog-adj., DE-M, 1-10ml, with recirculation valve	1	piece(s)
4600150	Dispensette® S, Analog-adj., DE-M, 2,5-25ml, without recirculation valve	1	piece(s)
4600151	Dispensette® S, Analog-adj., DE-M, 2,5-25ml, with recirculation valve	1	piece(s)
4600160	Dispensette® S, Analog-adj., DE-M, 5-50ml, without recirculation valve	1	piece(s)
4600161	Dispensette® S, Analog-adj., DE-M, 5-50ml, with recirculation valve	1	piece(s)
4600170	Dispensette® S, Analog-adj., DE-M, 10-100ml, without recirculation valve	1	piece(s)
4600171	Dispensette® S, Analog-adj., DE-M, 10-100ml, with recirculation valve	1	piece(s)

Bottle-top dispenser Dispensette® S, Fixed-volume, DE-M

Cat. No.	Description	From*	Unit
4600210	Dispensette® S, Fixed-vol., DE-M, 1ml, without recirculation valve	1	piece(s)
4600211	Dispensette® S, Fixed-vol., DE-M, 1ml with recirculation valve	1	piece(s)
4600220	Dispensette® S, Fixed-vol., DE-M, 2ml, without recirculation valve	1	piece(s)
4600221	Dispensette® S, Fixed-vol., DE-M, 2ml, with recirculation valve	1	piece(s)
4600230	Dispensette® S, Fixed-vol., DE-M, 5ml, without recirculation valve	1	piece(s)
4600231	Dispensette® S, Fixed-vol., DE-M, 5ml, with recirculation valve	1	piece(s)
4600240	Dispensette® S, Fixed-vol., DE-M, 10ml, without recirculation valve	1	piece(s)
4600241	Dispensette® S, Fixed-vol., DE-M, 10ml, with recirculation valve	1	piece(s)
4600290	Dispensette® S, Fixed-vol., DE-M, your choice of vol., without recircul. valve	1	piece(s)
4600291	Dispensette® S, Fixed-vol., DE-M, your choice of vol., with recircul. valve	1	piece(s)



Bottle-top dispenser Dispensette® S Organic, Digital, DE-M

Cat. No.	Description	From*	Unit
4630330	Dispensette® S Organic, Digital, DE-M, 0,5-5ml, without recircul. valve	1	piece(s)
4630331	Dispensette® S Organic, Digital, DE-M, 0,5-5ml, with recircul. valve	1	piece(s)
4630340	Dispensette® S Organic, Digital, DE-M, 1-10ml, without recircul. valve	1	piece(s)
4630341	Dispensette® S Organic, Digital, DE-M, 1-10ml, with recircul. valve	1	piece(s)
4630350	Dispensette® S Organic, Digital, DE-M, 2,5-25ml, without recircul. valve	1	piece(s)
4630351	Dispensette® S Organic, Digital, DE-M, 2,5-25ml, with recircul. valve	1	piece(s)
4630360	Dispensette® S Organic, Digital, DE-M, 5-50ml, without recircul. valve	1	piece(s)
4630361	Dispensette® S Organic, Digital, DE-M, 5-50ml, with recircul. valve	1	piece(s)



Bottle-top dispenser Dispensette® S Organic, Analog-adjustable, DE-M

Cat. No.	Description	From*	Unit
4630130	Dispensette® S Organic, Analog-adj., DE-M, 0,5-5ml, w/o recircul. valve	1	piece(s)
4630131	Dispensette® S Organic, Analog-adj., DE-M, 0,5-5ml, with recircul. valve	1	piece(s)
4630140	Dispensette® S Organic, Analog-adj., DE-M, 1-10ml, w/o recircul. valve	1	piece(s)
4630141	Dispensette® S Organic, Analog-adj., DE-M, 1-10ml, with recircul. valve	1	piece(s)
4630150	Dispensette® S Organic, Analog-adj., DE-M, 2,5-25ml, w/o recircul. valve	1	piece(s)
4630151	Dispensette® S Organic, Analog-adj., DE-M, 2,5-25ml, with recircul. valve	1	piece(s)
4630160	Dispensette® S Organic, Analog-adj., DE-M, 5-50ml, w/o recircul. valve	1	piece(s)
4630161	Dispensette® S Organic, Analog-adj., DE-M, 5-50ml, with recircul. valve	1	piece(s)
4630170	Dispensette® S Organic, Analog-adj., DE-M, 10-100ml, w/o recircul. valve	1	piece(s)
4630171	Dispensette® S Organic, Analog-adj., DE-M, 10-100ml, with recircul. valve	1	piece(s)



Bottle-top dispenser Dispensette® S Organic, Fixed-volume, DE-M



Cat. No.	Description	From*	Unit
4630230	Dispensette® S Organic, Fixed-vol., DE-M, 5ml, without recircul. valve	1	piece(s)
4630231	Dispensette® S Organic, Fixed-vol., DE-M, 5ml, with recircul. valve	1	piece(s)
4630240	Dispensette® S Organic, Fixed-vol., DE-M, 10ml, without recircul. valve	1	piece(s)
4630241	Dispensette® S Organic, Fixed-vol., DE-M, 10ml, with recircul. valve	1	piece(s)
4630290	Dispensette® S Organic, Fixed-v., DE-M, vol. of choice, w/o recircul. valve	1	piece(s)
4630291	Dispensette® S Organic, Fixed-v., DE-M, vol. of choice, with recircul. valve	1	piece(s)

Bottle-top dispenser Dispensette® S Trace Analysis, Analog-adjustable, DE-M



Cat. No.	Description	From*	Unit
4640040	Dispensette® S Trace Analysis, DE-M, PT/IR, 1-10 ml, w/o. recirculation valve	1	piece(s)
4640041	Dispensette® S Trace Analysis, DE-M, PT/IR, 1-10 ml, w. recirculation valve	1	piece(s)
4640240	Dispensette® S Trace Analysis, DE-M, Tantalum, 1-10 ml, w/o. recirculation valve	1	piece(s)
4640241	Dispensette® S Trace Analysis, DE-M, Tantalum, 1-10 ml, w. recirculation valve	1	piece(s)

Bottle-top dispenser seripettor®



Cat. No.	Description	From*	Unit
4720120	seripettor®, 0,2-2 ml	1	piece(s)
4720140	seripettor®, 1-10 ml	1	piece(s)
4720150	seripettor®, 2,5-25 ml	1	piece(s)

Bottle-top dispenser seripettor® pro



Cat. No.	Description	From*	Unit
4720420	seripettor® pro, 0,2-2 ml	1	piece(s)
4720440	seripettor® pro, 1-10 ml	1	piece(s)
4720450	seripettor® pro, 2,5-25 ml	1	piece(s)

Closure caps for seripettor®



Cat. No.	Description	From*	Unit
704552	seripettor® cap, PP, for 2+10 ml, cap for valve block	1	piece(s)
704554	seripettor® cap, PP, for 25 ml, cap for valve block	1	piece(s)

Adapter for seripettor® pro

Cat. No.	Description	From*	Unit
6208	Adapter for discharge tube seripettor® pro, ETFE, with seal	1	piece(s)
6707	Adapter for filling valve seripettor® pro, PP, with seal	1	piece(s)



Pump assembly for seripettor®

Cat. No.	Description	From*	Unit
704541	Pump assembly seripettor®, 2 ml, PC / stainless steel	1	piece(s)
704542	Pump assembly seripettor®, 10 ml, PC / stainless steel	1	piece(s)
704544	Pump assembly seripettor®, 25 ml, PC / stainless steel	1	piece(s)
704548	Pump assembly seripettor® pro, 10 ml, PPO /PEI (UV protection)	1	piece(s)
704549	Pump assembly seripettor® pro, 25 ml, PPO /PEI (UV protection)	1	piece(s)
704551	Pump assembly seripettor® pro, 2 ml, PPO /PEI (UV protection)	1	piece(s)



Filling valves for seripettor®

Cat. No.	Description	From*	Unit
6790	Valve set for seripettor (1 filling valve, 1 discharge valve and 2 seals)	1	pack



Filling tubes for seripettor®

Cat. No.	Description	From*	Unit
704532	Filling tube, PP, for seripettor®, 250 mm	1	pack
704534	Filling tube, PP, for seripettor®, 500 mm	1	pack
704536	Filling tube, PP, for seripettor®, for sterile applications, length 250 mm, O-ring	1	piece(s)
704538	Filling tube, PP, for seripettor®, for sterile applications, length 500 mm, O-ring	1	piece(s)



Flexible discharge tubing for seripettor®

Cat. No.	Description	From*	Unit
704522	Discharge tube seripettor® /seripettor® pro for 2 ml and 10 ml, PTFE, flexible, coiled, length 800 mm	1	piece(s)
704523	Discharge tube seripettor® /seripettor® pro for 25 ml, PTFE, flexible, coiled length 800 mm	1	piece(s)



Dispensing cartridge for seripettor® and QuikSip™



Cat. No.	Description		From*	Unit
704500	Dispenser cartridge for seripettor® /seripettor® pro/QuikSip, 2 ml, PE/PP	1 bag = 3 piece(s)	1	bag
704502	Dispenser cartridge for seripettor® /seripettor® pro/QuikSip, 10 ml	1 bag = 3 piece(s)	1	bag
704504	Dispenser cartridge for seripettor® /seripettor® pro/QuikSip, 5 - 25 ml	1 bag = 3 piece(s)	1	bag
704506	Dispenser cartridge for seripettor® /seripettor® pro/QuikSip 10 ml, sterile	1 pack = 7 piece(s)	1	pack
704507	Dispenser cartridge for seripettor® /seripettor® pro/QuikSip 2 ml, sterile, PE/PP	1 pack = 7 piece(s)	1	pack
704508	Dispenser cartridge for seripettor® /seripettor® pro/QuikSip 25 ml, sterile	1 pack = 5 piece(s)	1	pack

Discharge tubes for Dispensette® and seripettor®



Cat. No.	Description		From*	Unit
704518	Discharge tube, for seripettor® for 2 ml, fine tip		1	piece(s)
704520	Discharge tube, for seripettor® for 1 - 10 ml and 2,5 - 25 ml		1	piece(s)
707915	Discharge tube with integ. valve for Dispensette®, seripettor® pro, up to 10 ml, fine tip, 90 mm		1	piece(s)
707916	Discharge tube with integ. valve for Dispensette®, seripettor® pro, 5, 10 ml, 90 mm		1	piece(s)
707918	Discharge tube with integ. valve for Dispensette®, seripettor® pro, up to 100 ml, fine tip, 120 mm		1	piece(s)
707928	Discharge tube with Luer-Lock attachment for micro filter, for seripettor® pro		1	piece(s)

Bottle stand for Dispensette® and Titrette®



Cat. No.	Description		From*	Unit
704275	Bottle stand PP Dispensette®, Titrette® support rod 325 mm base plate 220 x 160 mm		1	piece(s)

Discharge tubes for Dispensette® S



Cat. No.	Description		From*	Unit
708002	Discharge tube for Dispensette® S 1, 2, 5 and 10 ml, fine tip, 105mm		1	piece(s)
708005	Discharge tube for Dispensette® S 5 and 10 ml, standard, 105 mm		1	piece(s)
708006	Discharge tube for Dispensette® S 25, 50 and 100 ml, fine tip, 135mm		1	piece(s)
708008	Discharge tube for Dispensette® S 25, 50 and 100 ml, standard, 135mm		1	piece(s)
708102	Discharge tube with recirculation valve for Dispensette® S 1, 2, 5 and 10 ml, fine tip, 108mm		1	piece(s)
708104	Discharge tube with recirculation valve for Dispensette® S 5 and 10 ml, standard, 108 mm		1	piece(s)
708106	Discharge tube with recirculation valve for Dispensette® S 25, 50 and 100 ml, fine tip, 135mm		1	piece(s)
708109	Discharge tube with recirculation valve for Dispensette® S 25, 50 and 100 ml, standard, 135mm		1	piece(s)

Discharge tubes for Dispensette® S Organic

Cat. No.	Description	From*	Unit
708012	Discharge tube for Dispensette® S Organic, 5 and 10 ml, fine tip, 108mm	1	piece(s)
708014	Discharge tube for Dispensette® S Organic 5 and 10 ml, standard, 108mm	1	piece(s)
708016	Discharge tube for Dispensette® S Organic 25, 50 and 100ml, fine tip, 135mm	1	piece(s)
708019	Discharge tube for Dispensette® S Organic 25, 50 and 100ml, standard, 135mm	1	piece(s)
708112	Discharge tube with recirculation valve for Dispensette® S Organic, 5 and 10 ml, fine tip, 108mm	1	piece(s)
708114	Discharge tube with recirculation valve for Dispensette® S Organic 5 and 10 ml, standard, 108mm	1	piece(s)
708116	Discharge tube with recirculation valve for Dispensette® S Organic 25, 50 and 100ml, fine tip, 135mm	1	piece(s)
708119	Discharge tube with recirculation valve for Dispensette® S Organic 25, 50 and 100ml, standard, 135mm	1	piece(s)



Discharge tubes for Dispensette® S Trace Analysis

Cat. No.	Description	From*	Unit
708022	Discharge tube f.Dispensette® S Trace Anal.,PT/IR valve spring, w/o.recirc.valve	1	piece(s)
708024	Discharge tube f.Dispensette® S Trace Anal., TA valve spring,w/o.recirc.valve	1	piece(s)
708122	Discharge tube f.Dispensette® S Trace Anal.,PT/IR valve spring, w.recirc.valve	1	piece(s)
708124	Discharge tube f.Dispensette® S Trace Anal., TA valve spring, w.recirc.valve	1	piece(s)



Flexible discharge tubing for Dispensette® S / S Organic / S Trace Analysis

Cat. No.	Description	From*	Unit
708132	Flexible discharge tubing with recircul. valve f.Dispensette® S / S Or. 2, 5 and 10ml, PTFE, 80cm	1	piece(s)
708134	Flexible discharge tubing with recircul. valve f.Dispensette® S / S Org. 25, 50, 100ml, PTFE, 80cm	1	piece(s)



Bottle adapter for Dispensette®, seripettor®, Titrette® and QuikSip™



Cat. No.	Description	From*	Unit
704322	Thread adapter, PP, outer thread GL 32, for bottle thread GL 22	1	piece(s)
704325	Thread adapter, PP, outer thread GL 32, for bottle thread GL 24-25	1	piece(s)
704328	Thread adapter, PP, outer thread GL 32, for bottle thread GL 28	1	piece(s)
704330	Thread adapter, PP, outer thread GL 32, for bottle thread GL 30	1	piece(s)
704338	Thread adapter, PP, outer thread GL 32, for bottle thread GL 38	1	piece(s)
704340	Thread adapter, PP, outer thread GL 32, for bottle thread S 40	1	piece(s)
704343	Thread adapter, PP, outer thread GL 45, for bottle thread S 40	1	piece(s)
704345	Thread adapter, PP, outer thread GL 32, for bottle thread GL 45	1	piece(s)
704348	Thread adapter, PP, outer thread GL 45, for bottle thread S 60	1	piece(s)
704350	Thread adapter, PP, outer thread GL 45, for bottle thread S 50	1	piece(s)
704375	Thread adapter, ETFE, outer thread GL 32, for bottle thread GL 24-25	1	piece(s)
704378	Thread adapter, ETFE, outer thread GL 32, for bottle thread GL 28	1	piece(s)
704388	Thread adapter, ETFE, outer thread GL 32, for bottle thread GL 38	1	piece(s)
704390	Thread adapter, ETFE, outer thread GL 32, for bottle thread S 40	1	piece(s)
704391	Thread adapter, PTFE, outer thread GL 45, for bottle thread S 40	1	piece(s)
704395	Thread adapter, ETFE, outer thread GL 32, for bottle thread GL 45	1	piece(s)
704396	Thread adapter, PP, outer thread GL 45, for bottle thread GL 32-33	1	piece(s)
704397	Thread adapter, PP, outer thread GL 45, for bottle thread GL 38	1	piece(s)
704398	Thread adapter, ETFE, outer thread GL 45, for bottle thread GL 32-33	1	piece(s)
704399	Thread adapter, ETFE, outer thread GL 45, for bottle thread GL 38	1	piece(s)
704419	Thread adapter, PP, outer thread GL 32, for ground joint NS 19/26	1	piece(s)
704424	Thread adapter, PP, outer thread GL 32, for ground joint NS 24/29	1	piece(s)
704429	Thread adapter, PP, outer thread GL 32, for ground joint NS 29/32	1	piece(s)
704430	Thread adapter, PP, outer thread GL 45, for bottle thread S 54	1	piece(s)
704431	Thread adapter, PP, outer thread GL 45, for bottle thread GL 35	1	piece(s)

Filling valves for Dispensette® S / S Organic



Cat. No.	Description	From*	Unit
6734	Filling valve with olive for Dispensette® S 1, 2, 5 + 10 ml, PFA/glass/ceramic	1	piece(s)
6735	Filling valve with olive for Dispensette® S 25, 50 and 100 ml, PFA/glass/ceramic	1	piece(s)

Filling valves for Dispensette®



Cat. No.	Description	From*	Unit
6697	Filling valve with sealing washer, Dispensette® III/seripettor® pro/Digital Burette, 0.5, 1, 2, 5, 10 ml	1	piece(s)
6698	Filling valve with sealing washer, Dispensette® III/Organic/seripettor® pro, 25, 50, 100 ml	1	piece(s)


Valves for Dispensette® S / S Organic / S Trace Analysis

Cat. No.	Description	From*	Unit
6727	Discharge valve for Dispensette® S 5 and 10 ml, PFA/glass/ceramic/platinum-iridium, no valve marking	1	piece(s)
6728	Discharge valve for Dispensette® S 25, 50 and 100 ml, PFA/glass/ceramic/platinum-iridium, no valve marking	1	piece(s)
6729	Discharge valve for Dispensette® S Organic, 5 and 10ml, PFA/glass/ceramic/tantalum, valve marked "ORG"	1	piece(s)
6730	Discharge valve for Dispensette® S Organic 25, 50 and 100ml, PFA/glass/ceramic/tantalum, valve marked "ORG"	1	piece(s)
6732	Discharge valve for Dispensette® S Trace Analysis, 10 ml, with Pt/Ir-spring	1	piece(s)
6733	Discharge valve for Dispensette® S Trace Analysis, 10 ml, with Tantalum-spring	1	piece(s)
6749	Discharge valve for Dispensette® S, 1 and 2 ml, valve marked with "1+2"	1	piece(s)



Telescoping filling tubes for Dispensette® S / S Organic / S Trace Analysis / seripettor pro / Dispensette® III / TA

Cat. No.	Description	From*	Unit
708210	Telescoping filling tube for Dispensette® S / S Organic / S Trace Analysis / seripettor® pro with nominal volume 1, 2, 5 and 10 ml, FEP, 70-140 mm	1	piece(s)
708212	Telescoping filling tube for Dispensette® S / S Organic / S Trace Analysis / seripettor® pro with nominal volume 1, 2, 5 and 10 ml, FEP, 125-240 mm	1	piece(s)
708214	Telescoping filling tube for Dispensette® S / S Organic / S Trace Analysis / seripettor® pro with nominal volume 1, 2, 5 and 10 ml, FEP, 195-350 mm	1	piece(s)
708216	Telescoping filling tube for Dispensette® S / S Organic / S Trace Analysis / seripettor® pro with nominal volume 1, 2, 5 and 10 ml, FEP, 250-480 mm	1	piece(s)
708218	Telescoping filling tube for Dispensette® S / S Organic / seripettor® pro with nominal volume 25, 50, 100 ml as well as Titrette®, FEP, 170-330 mm	1	piece(s)
708220	Telescoping filling tube for Dispensette® S / S Organic / seripettor® pro with nominal volume 25, 50, 100 ml as well as Titrette®, FEP, 250-480 mm	1	piece(s)



Seals

Cat. No.	Description	From*	Unit
704486	Seal for valve block, PTFE, for highly volatile reagents	1	piece(s)



Remote dispensing system for drum dispensing Dispensette®



Cat. No.	Description	From*	Unit
704209	Filling tubes, FEP, for Dispensette® 0,5 - 10 ml	1	piece(s)
704261	Remote disp. system for Dispensette®, standard supply (without Dispensette)	1	piece(s)
704267	Filling tube for remote dispensing system Dispensette®, 10m, outer diameter 7.6 mm	1	piece(s)
704268	Support rod connector for wall mounting unit for remote dispensing system Dispensette®, stain. steel	1	piece(s)
704269	Filling tube for remote dispensing system Dispensette®, 1m, outer diameter 6.9 mm	1	piece(s)
704270	Thread adapter for remote dispensing system Dispensette®, stain. steel, outer diameter 2"	1	piece(s)
704272	Table, shelf clamp for wall mounting unit, for remote dispensing system Dispensette®	1	piece(s)
704281	Thread adapter for remote dispensing system Dispensette®, PTFE, outer thread 3/4", GL 32	1	piece(s)
704282	Thread adapter for remote dispensing system Dispensette®, PTFE, inner thread 3/4", GL 32	1	piece(s)

Dispensing cartridge for Dispensette® S Trace Analysis



Cat. No.	Description	From*	Unit
708035	Dispensing cartridge with safety ring for Dispensette® S Trace Analysis	1	piece(s)

Drying tube for Dispensette®, seripettor®, Titrette®



Cat. No.	Description	From*	Unit
707930	Drying tube, PE-LD, for Dispensette®, Titrette® and Digital Burette III	1	piece(s)

Venting stopper, Micro filter with Luer-cone



Cat. No.	Description	From*	Unit
704495	Air vent cap for micro filter with Luer-cone, PP, for Dispensette®	1	piece(s)

Bottle-top burette Titrette®

Cat. No.	Description	From*	Unit
4760141	Titrette®, 10 ml, DE-M, with titration/recirc. valve	1	piece(s)
4760151	Titrette®, 25 ml, DE-M, with titration/recirc. Valve	1	piece(s)
4760161	Titrette®, 50 ml, DE-M, with titration/recirc. Valve	1	piece(s)
4760241	Titrette®, 10 ml, DE-M, with interf. RS 232	1	piece(s)
4760251	Titrette®, 25ml, DE-M, with interf. RS 232	1	piece(s)
4760261	Titrette®, 50 ml, DE-M, with interf. RS 232 with titr./recirc. Valve	1	piece(s)
4760451	Titrette® SH, DE-M, according to Soxhlet-Henkel (4 °SH = 1 ml)	1	piece(s)



Titration tube for Titrette®

Cat. No.	Description	From*	Unit
707525	Titration tube with screw cap for Titrette® 10ml, w.integ.discha./recirc.valve	1	piece(s)
707526	Titration tube with screw cap for Titrette® 25ml, w.integ.discha./recirc.valve	1	piece(s)
707529	Titration tube with screw cap for Titrette® 25 ml+50 ml, w. integ. discha./recirc. valve from SN 01K	1	piece(s)



Dispensing cylinder with valve block for Titrette®

Cat. No.	Description	From*	Unit
707533	Dispensing cylinder with valve block for Titrette® 10 ml	1	piece(s)
707535	Dispensing cylinder with valve block for Titrette® 25 ml, from serial-number 01K	1	piece(s)
707537	Dispensing cylinder with valve block for Titrette® 50 ml, from serial-number 01K	1	piece(s)



Piston head for Titrette®

Cat. No.	Description	From*	Unit
707530	Piston head for Titrette® 25 ml	1	piece(s)
707531	Piston head for Titrette® 10 ml	1	piece(s)
707532	Piston head for Titrette® 50 ml	1	piece(s)



Inspection window for Titrette®

Cat. No.	Description	From*	Unit
6783	Inspection windows for Titrette® (colorless / amber colored)	1 pack = 2 piece(s)	pack



Filling valve for Titrette®

Cat. No.	Description	From*	Unit
6636	Filling valve with olive-shaped nozzle and sealing ring, for Titrette®	1	

Remote system Titrette® for Titripac® packing system

Cat. No.	Description	From*	Unit
707550	Remote system Titrette® for Titripac® packing system	1	piece(s)

Bottle-top aspirator QuikSip™ BT-Aspirator

Cat. No.	Description	From*	Unit
4723150	QuikSip BT-Aspirator incl. cell-culture unit	1	piece(s)

cell-culture™ unit

Cat. No.	Description	From*	Unit
25950	cell-culture™ (aspiration controller)	1	piece(s)

Filling tube for QuikSip™ BT-Aspirator

Cat. No.	Description	From*	Unit
704575	Filling tube for QuikSip BT-Aspirator PP, with filling valve of PP/EPDM	1	piece(s)

Discharge valve for QuikSip™ BT-Aspirator

Cat. No.	Description	From*	Unit
704580	Discharge valve for QuikSip BT-Aspirator PP/EPDM, with sealing of EPDM	1	piece(s)

8-channel manifold for QuikSip™ BT-Aspirator

Cat. No.	Description	From*	Unit
704526	8-channel manifold, PP, for QuikSip, BT-Aspirator, autoclavable	1	piece(s)



Adapter for cell-culture™-unit

Cat. No.	Description	From*	Unit
25933	Adapter, PVC, for capillaries	1	piece(s)
25960	Adapter, SI/PVC, for pasteur pipettes	1	piece(s)
25961	Adapter, PP, for pipette tips 100 µl	1	piece(s)



Seal for QuikSip™ BT-Aspirator

Cat. No.	Description	From*	Unit	
6788	Seals for QuikSip™ (Set of 5 pieces)	1 pack = 5 piece(s)	1	pack



Tubing for cell-culture™-unit

Cat. No.	Description	From*	Unit
25962	Suction tube for cell-culture™, SI, 2 metres long	1	piece(s)



Single channel micropipette Transferpette® S, adjustable, CE-IVD, DE-M

Cat. No.	Description	From*	Unit
705868	Transferpette® S, 0,1-1 µl, adjustable volume, D-1, DE-M	1	piece(s)
705869	Transferpette® S, 0,1-2,5 µl, adjustable volume, D-2,5, DE-M	1	piece(s)
705870	Transferpette® S, 0,5-10 µl, adjustable volume, D-10, DE-M	1	piece(s)
705872	Transferpette® S, 2-20 µl, adjustable volume, D-20, DE-M	1	piece(s)
705873	Transferpette® S, 5-50 µl, adjustable volume, D-50, DE-M	1	piece(s)
705874	Transferpette® S, 10-100 µl, adjustable volume, D-100, DE-M	1	piece(s)
705878	Transferpette® S, 20-200 µl, adjustable volume, D-200, DE-M	1	piece(s)
705880	Transferpette® S, 100-1000 µl, adjustable volume, D-1000, DE-M	1	piece(s)
705882	Transferpette® S, 500-5000 µl, adjustable volume, D-5000, DE-M	1	piece(s)
705884	Transferpette® S, 1000-10000 µl, adjustable volume, D-10000, DE-M	1	piece(s)



Single channel micropipette Transferpette® S, Fixed-volume, CE-IVD, DE-M



Cat. No.	Description	From*	Unit
705808	Transferpette® S, 10 µl, fixed volume, F-10, DE-M	1	piece(s)
705816	Transferpette® S, 20 µl, fixed volume, F-20, DE-M	1	piece(s)
705820	Transferpette® S, 25 µl, fixed volume, F-25, DE-M	1	piece(s)
705828	Transferpette® S, 50 µl, fixed volume, F-50, DE-M	1	piece(s)
705838	Transferpette® S, 100 µl, fixed volume, F-100, DE-M	1	piece(s)
705844	Transferpette® S, 200 µl, fixed volume, F-200, DE-M	1	piece(s)
705854	Transferpette® S, 500 µl, fixed volume, F-500, DE-M	1	piece(s)
705862	Transferpette® S, 1000 µl, fixed volume, F-1000, DE-M	1	piece(s)
705864	Transferpette® S, 2000 µl, fixed volume, F-2000, DE-M	1	piece(s)

Pipette sets Transferpette® S, adjustable, CE-IVD, DE-M



Cat. No.	Description	From*	Unit
705890	Transferpette® S Starter-Kit Micro: 3 x pipettes (0,1-2,5 µl, 0,5-10 µl, 10-100 µl), 3 x TipBox pipette tips (2 x 0,5-20 µl, 2-200 µl), 3 shelf mounts	1	piece(s)
705891	Transferpette® S Starter-Kit Midi: 3 x pipettes (2-20 µl, 20-200 µl, 100-1000 µl), 3 x TipBox pipette tips (2 x 2-200 µl, 50-1000 µl), 3 shelf mount	1	piece(s)
705892	Transferpette® S Starter-Kit Macro: 3 x pipettes (100-1000 µl, 500-5000 µl, 1000-10000 µl), 3 x TipBox pipette tips (50-1000 µl, 0,5-5 ml, 1-10 ml), 3 shelf mounts	1	piece(s)
705893	Transferpette® S Starter-Kit Standard: 3 x pipettes (0,5-10 µl, 10-100 µl, 100-1000 µl), 3 x TipBox pipette tips (0,5-20 µl, 2-200 µl, 50-1000 µl), 3 shelf mounts	1	piece(s)
705896	Transferpette® S package 1: 5 x pipettes (0,1-2,5 µl, 0,5-10 µl, 2-20 µl, 20-200 µl, 100-1000 µl), 3 x TipBox pipette tips (0,5-20 µl, 2-200 µl, 50-1000 µl), benchtop rack f. 6 pipettes	1	piece(s)
705897	Transferpette® S package 2: 5 x pipettes (0,5-10 µl, 10-100 µl, 20-200 µl, 50-1000 µl, 500-5000 µl), 4 x TipBox pipette tips (0,5-20 µl, 2-200 µl, 50-1000 µl, 0,5-5 ml), benchtop rack f. 6 pipettes	1	piece(s)

Multichannel micropipette Transferpette® S-8, adjustable, CE-IVD, DE-M



Cat. No.	Description	From*	Unit
705900	Transferpette® S-8, adjustable volume, DE-M, M8-10, 0,5-10 µl	1	piece(s)
705906	Transferpette® S-8, adjustable volume, DE-M, M8-50, 5-50 µl	1	piece(s)
705908	Transferpette® S-8, adjustable volume, DE-M, M8-100, 10-100 µl	1	piece(s)
705910	Transferpette® S-8, adjustable volume, DE-M, M8-200, 20-200 µl	1	piece(s)
705912	Transferpette® S-8, adjustable volume, DE-M, M8-300, 30-300 µl	1	piece(s)

Multichannel micropipette Transferpette® S -12-12, adjustable, CE-IVD, DE-M

Cat. No.	Description	From*	Unit
705920	Transferpette® S-12, adjustable volume, DE-M, M12-10, 0.5-10 µl	1	piece(s)
705926	Transferpette® S-12, adjustable volume, DE-M, M12-50, 5-50 µl	1	piece(s)
705928	Transferpette® S-12, adjustable volume, DE-M, M12-100, 10-100 µl	1	piece(s)
705930	Transferpette® S-12, adjustable volume, DE-M, M12-200, 20-200 µl	1	piece(s)
705932	Transferpette® S-12, adjustable volume, DE-M, M12-300, 30-300 µl	1	piece(s)


Single channel micropipette Transferpette® electronic, CE-IVD, DE-M

Cat. No.	Description	From*	Unit
705299	Transferpette® electronic with AC charger (Europe), 0.5-10 µl, DE-M	1	piece(s)
705300	Transferpette® electronic with AC adapter (Europe), 2-20 µl, DE-M	1	piece(s)
705303	Transferpette® electronic with AC adapter (Europe), 10-200 µl, DE-M	1	piece(s)
705306	Transferpette® electronic with AC adapter (Europe), 50-1000 µl, DE-M	1	piece(s)
705307	Transferpette® electronic with AC adapter (Europe), 250-5000 µl, DE-M	1	piece(s)
705309	Transferpette® electronic with AC adapter (UK/Ireland), 0.5-10 µl, DE-M	1	piece(s)
705310	Transferpette® electronic with AC adapter (UK/Ireland), 2-20 µl, DE-M	1	piece(s)
705313	Transferpette® electronic with AC adapter (UK/Ireland), 10-200 µl, DE-M	1	piece(s)
705316	Transferpette® electronic with AC adapter (UK/Ireland), 50-1000 µl, DE-M	1	piece(s)
705317	Transferpette® electronic with AC adapter (UK/Ireland), 250-5000 µl, DE-M	1	piece(s)
705319	Transferpette® electronic with AC adapter (USA/Japan), 0.5-10 µl, DE-M	1	piece(s)
705320	Transferpette® electronic with AC adapter (USA/Japan), 2-20 µl, DE-M	1	piece(s)
705323	Transferpette® electronic with AC adapter (USA/Japan), 10-200 µl, DE-M	1	piece(s)
705326	Transferpette® electronic with AC adapter (USA/Japan), 50-1000 µl, DE-M	1	piece(s)
705327	Transferpette® electronic with AC adapter (USA/Japan), 250-5000 µl, DE-M	1	piece(s)
705329	Transferpette® electronic with AC adapter (Australia), 0.5-10 µl, DE-M	1	piece(s)
705330	Transferpette® electronic with AC adapter (Australia), 2-20 µl, DE-M	1	piece(s)
705333	Transferpette® electronic with AC adapter (Australia), 10-200 µl, DE-M	1	piece(s)
705336	Transferpette® electronic with AC adapter (Australia), 50-1000 µl, DE-M	1	piece(s)
705337	Transferpette® electronic with AC adapter (Australia), 250-5000 µl, DE-M	1	piece(s)
705339	Transferpette® electronic without AC adapter, 0.5-10 µl, DE-M	1	piece(s)
705340	Transferpette® electronic without AC adapter, 2-20 µl, DE-M	1	piece(s)
705343	Transferpette® electronic without AC adapter, 10-200 µl, DE-M	1	piece(s)
705346	Transferpette® electronic without AC adapter, 50-1000 µl, DE-M	1	piece(s)
705347	Transferpette® electronic without AC adapter, 250-5000 µl, DE-M	1	piece(s)



Multi channel micropipette Transferpette® -8 electronic, CE-IVD, DE-M



Cat. No.	Description	From*	Unit
705399	Transferpette® -8 electronic with AC adapter (Europe), 0,5-10 µl, DE-M	1	piece(s)
705400	Transferpette® -8 electronic with AC adapter (Europe), 1-20 µl, DE-M	1	piece(s)
705403	Transferpette® -8 electronic with AC adapter (Europe), 5-100 µl, DE-M	1	piece(s)
705404	Transferpette® -8 electronic with AC adapter (Europe), 10-200 µl, DE-M	1	piece(s)
705406	Transferpette® -8 electronic with AC adapter (Europe), 15-300 µl, DE-M	1	piece(s)
705409	Transferpette® -8 electronic with AC adapter (UK/Ireland), 0,5-10 µl, DE-M	1	piece(s)
705410	Transferpette® -8 electronic with AC adapter (UK/Ireland), 1-20 µl, DE-M	1	piece(s)
705413	Transferpette® -8 electronic with AC adapter (UK/Ireland), 5-100 µl, DE-M	1	piece(s)
705414	Transferpette® -8 electronic with AC adapter (UK/Ireland), 10-200 µl, DE-M	1	piece(s)
705416	Transferpette® -8 electronic with AC adapter (UK/Ireland), 15-300 µl, DE-M	1	piece(s)
705419	Transferpette® -8 electronic with AC adapter (USA/Japan), 0,5-10 µl, DE-M	1	piece(s)
705420	Transferpette® -8 electronic with AC adapter (USA/Japan), 1-20 µl, DE-M	1	piece(s)
705423	Transferpette® -8 electronic with AC adapter (USA/Japan), 5-100 µl, DE-M	1	piece(s)
705424	Transferpette® -8 electronic with AC adapter (USA/Japan), 10-200 µl, DE-M	1	piece(s)
705426	Transferpette® -8 electronic with AC adapter (USA/Japan), 15-300 µl, DE-M	1	piece(s)
705429	Transferpette® -8 electronic with AC adapter (Australia), 0,5-10 µl, DE-M	1	piece(s)
705430	Transferpette® -8 electronic with AC adapter (Australia), 1-20 µl, DE-M	1	piece(s)
705433	Transferpette® -8 electronic with AC adapter (Australia), 5-100 µl, DE-M	1	piece(s)
705434	Transferpette® -8 electronic with AC adapter (Australia), 10-200 µl, DE-M	1	piece(s)
705436	Transferpette® -8 electronic with AC adapter (Australia), 15-300 µl, DE-M	1	piece(s)

Multi channel micropipette Transferpette® -12 electronic, CE-IVD, DE-M

Cat. No.	Description	From*	Unit
705449	Transferpette® -12 electronic with AC adapter (Europe), 0,5-10 µl, DE-M	1	piece(s)
705450	Transferpette® -12 electronic with AC adapter (Europe), 1-20 µl, DE-M	1	piece(s)
705453	Transferpette® -12 electronic with AC adapter (Europe), 5-100 µl, DE-M	1	piece(s)
705454	Transferpette® -12 electronic with AC adapter (Europe), 10-200 µl, DE-M	1	piece(s)
705456	Transferpette® -12 electronic with AC adapter (Europe), 15-300 µl, DE-M	1	piece(s)
705459	Transferpette® -12 electronic with AC adapter (UK/Ireland), 0,5-10 µl, DE-M	1	piece(s)
705460	Transferpette® -12 electronic with AC adapter (UK/Ireland), 1-20 µl, DE-M	1	piece(s)
705463	Transferpette® -12 electronic with AC adapter (UK/Ireland), 5-100 µl, DE-M	1	piece(s)
705464	Transferpette® -12 electronic with AC adapter (UK/Ireland), 10-200 µl, DE-M	1	piece(s)
705466	Transferpette® -12 electronic with AC adapter (UK/Ireland), 15-300 µl, DE-M	1	piece(s)
705469	Transferpette® -12 electronic with AC adapter (USA/Japan), 0,5-10 µl, DE-M	1	piece(s)
705470	Transferpette® -12 electronic with AC adapter (USA/Japan), 1-20 µl, DE-M	1	piece(s)
705473	Transferpette® -12 electronic with AC adapter (USA/Japan), 5-100 µl, DE-M	1	piece(s)
705474	Transferpette® -12 electronic with AC adapter (USA/Japan), 10-200 µl, DE-M	1	piece(s)
705476	Transferpette® -12 electronic with AC adapter (USA/Japan), 15-300 µl, DE-M	1	piece(s)
705479	Transferpette® -12 electronic with AC adapter (Australia), 0,5-10 µl, DE-M	1	piece(s)
705480	Transferpette® -12 electronic with AC adapter (Australia), 1-20 µl, DE-M	1	piece(s)
705483	Transferpette® -12 electronic with AC adapter (Australia), 5-100 µl, DE-M	1	piece(s)
705484	Transferpette® -12 electronic with AC adapter (Australia), 10-200 µl, DE-M	1	piece(s)
705486	Transferpette® -12 electronic with AC adapter (Australia), 15-300 µl, DE-M	1	piece(s)



Stands and mounts for Transferpette® S / Transferpette® electronic

Cat. No.	Description	From*	Unit
703440	Individual stand for Transferpette®-8/12 and Transferpette® S-8/12 for 1 Transferpette®	1	piece(s)
704807	Benchtop rack (round) for pipettes Transferpette® S and S -8/-12	1	piece(s)
704811	Rack mount for pipettes Transferpette® S and S -8/-12	1	piece(s)
704812	Wall mount for pipettes Transferpette® S and S -8/-12	1	piece(s)
705385	Individual stand for pipette Transferpette® electronic, except 0,5-5 ml	1	piece(s)
705386	Individual stand for pipette Transferpette® electronic 0,25-5ml or Transferpette® S 2ml or 0,5-5ml	1	piece(s)
705390	Charging stand for 3 Transferpette® electronic, AC adapter Europe (continent), except 0,5-5 ml	1	piece(s)
705391	Charging stand for 3 Transferpette® electronic, AC adapter UK/Ireland, except 0,5-5 ml	1	piece(s)
705392	Charging stand for 3 Transferpette® electronic, AC adapter USA/Japan, except 0,5-5 ml	1	piece(s)
705393	Charging stand for 3 Transferpette® electronic, AC adapter Australia, except 0,5-5 ml	1	piece(s)



Filter for Transferpette® S / Transferpette® electronic



Cat. No.	Description		From*	Unit
704652	Filter for Transferpette® S / electronic 0,5-5 ml, Fix 2 ml	1 pack = 25 piece(s)	1	pack
704653	Filter for Transferpette® 10000 µl, for Transferpette® S	1 pack = 25 piece(s)	1	pack

Reagent reservoirs for multichannel micropipette



Cat. No.	Description		From*	Unit
703409	Reagent reservoir PP for multichannel pip. 60 ml, sterile, translucent, autoclavable	1 pack = 200 piece(s)	1	pack
703411	Reagent reservoir PP for multichannel pip. 60 ml, sterile, translucent, autoclavable	1 pack = 100 piece(s)	1	pack
703459	Reagent reservoir f. multi-channel pipettes, 60 ml, w. lid, PP	1 pack = 10 piece(s)	1	pack

Battery for Transferpette® electronic



Cat. No.	Description		From*	Unit
705500	Battery for Transferpette® electronic		1	piece(s)

Pipette tips, 0,1 - 20 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732002	Pipette tips, bulk, 0,1- 20 µl, non-sterile, PP, IVD, colorless		2000	piece(s)
732022	Pipette tips, bulk, 0,1-20 µl, non-sterile, PP, XXL, IVD, colorless		10000	piece(s)
732102	Pipette tips, racked, TipRack, 0,1-20 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732122	Pipette tips, racked, TipRack, 0,1-20 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732202	Pipette tips, racked, TipBox, 0,1-20 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1	pack
732222	Pipette tips, racked, TipBox, 0,1- 20 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack



Pipette tips, 0,5 - 20 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732004	Pipette tips, bulk, 0,5- 20 µl, non-sterile, PP, IVD, colorless		2000	piece(s)
732024	Pipette tips, bulk, 0,5-20 µl, non-sterile, PP, XXL, IVD, colorless		10000	piece(s)
732104	Pipette tips, racked, TipRack, 0,5-20 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732124	Pipette tips, racked, TipRack, 0,5-20 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732204	Pipette tips, racked, TipBox, 0,5-20 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1	pack
732224	Pipette tips, racked, TipBox, 0,5- 20 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732244	Pipette tips, racked, TipStack™, 0,5- 20 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732264	Pipette tips, racked, TipStack™, 0,5-20 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack



Pipette tips, 1 - 50 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732006	Pipette tips, bulk, 1-50 µl, non-sterile, PP, IVD, colorless		2000	piece(s)
732026	Pipette tips, bulk, 1-50 µl, non-sterile, PP, XXL, IVD, colorless		10000	piece(s)
732106	Pipette tips, racked, TipRack, 1-50 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732126	Pipette tips, racked, TipRack, 1-50 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732206	Pipette tips, racked, TipBox, 1-50 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1	pack
732226	Pipette tips, racked, TipBox, 1-50 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack





Pipette tips, 2 - 200 µl, PP, CE-IVD

Cat. No.	Description	From*	Unit
732008	Pipette tips, bulk, 2-200 µl, non-sterile, PP, IVD, yellow	1000	piece(s)
732028	Pipette tips, bulk, 2-200 µl, non-sterile, PP, XXL, IVD, yellow	10000	piece(s)
732108	Pipette tips, racked, TipRack, 2-200 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1
732128	Pipette tips, racked, TipRack, 2-200 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1
732208	Pipette tips, racked, TipBox, 2-200 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1
732228	Pipette tips, racked, TipBox, 2-200 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1
732248	Pipette tips, racked, TipStack™, 2-200 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1
732268	Pipette tips, racked, TipStack™, 2-200 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1



Pipette tips, 5 - 300 µl, PP, CE-IVD

Cat. No.	Description	From*	Unit
732010	Pipette tips, bulk, 5-300 µl, non-sterile, PP, IVD, colorless	1000	piece(s)
732030	Pipette tips, bulk, 5-300 µl, non-sterile, PP, XXL, IVD, colorless	10000	piece(s)
732110	Pipette tips, racked, TipRack, 5-300 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1
732130	Pipette tips, racked, TipRack, 5-300 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1
732210	Pipette tips, racked, TipBox, 5-300 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1
732230	Pipette tips, racked, TipBox, 5-300 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1



Pipette tips, 50 - 1000 µl, PP, CE-IVD

Cat. No.	Description	From*	Unit
732012	Pipette tips, bulk, 50-1000 µl, non-sterile, PP, IVD, blue	1000	piece(s)
732032	Pipette tips, bulk, 50-1000 µl, non-sterile, PP, XXL, IVD, blue	5000	piece(s)
732112	Pipette tips, racked, TipRack, 50-1000 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1
732132	Pipette tips, racked, TipRack, 50-1000 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1
732212	Pipette tips, racked, TipBox, 50-1000 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1
732232	Pipette tips, racked, TipBox, 50-1000 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1
732252	Pipette tips, racked, TipStack™, 50-1000 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1
732272	Pipette tips, racked, TipStack™, 50-1000 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1



Pipette tips, 500 - 5000 µl, PP, CE-IVD

Cat. No.	Description	From*	Unit
702595	Pipette tips bulk, PP, 0,5-5 ml, colorless, non-sterile, IVD	200	piece(s)
702600	Pipette tips bulk, PP, 0,5-5 ml, colorless, non-sterile, XXL, IVD	1000	piece(s)
702605	Pipette tips racked, PP, TipBox 5 ml, round, with 28 pipette tips 0,5-5 ml	1 pack = 28 piece(s)	1

Pipette tips, 1 - 10 ml, PP, CE-IVD

Cat. No.	Description	From*	Unit
702603	Pipette tips bulk, PP, 1-10 ml, colorless, non-sterile, IVD	200	piece(s)
702604	Pipette tips bulk, PP, 1-10 ml, colorless, non-sterile, XXL, IVD	1000	piece(s)
702608	Pipette tips racked, PP, TipBox 10 ml, round, with 18 pipette tips 1-10 ml	1 pack = 18 piece(s)	1



Ultra Low Retention pipette tips, 0,1 - 20 µl, PP, CE-IVD

Cat. No.	Description	From*	Unit
732302	Pipette tips Ultra Low Retention, racked, TipBox, 0,1- 20 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1
732322	Pipette tips Ultra Low Retention, racked, TipBox, 0,1-20 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1



Ultra Low Retention pipette tips, 0,5 - 20 µl, PP, CE-IVD

Cat. No.	Description	From*	Unit
732304	Pipette tips Ultra Low Retention, racked, TipBox, 0,5- 20 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1
732324	Pipette tips Ultra Low Retention, racked, TipBox, 0,5-20 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1
732344	Pipette tips Ultra Low Retention, racked, TipStack™, 0,5-20 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1
732364	Pipette tips Ultra Low Retention, racked, TipStack™, 0,5-20 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1



Ultra Low Retention pipette tips, 1 - 50 µl, PP, CE-IVD

Cat. No.	Description	From*	Unit
732306	Pipette tips Ultra Low Retention, racked, TipBox, 1-50 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1
732326	Pipette tips Ultra Low Retention, racked, TipBox, 1-50 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1





Ultra Low Retention pipette tips, 2 - 200 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732308	Pipette tips Ultra Low Retention, racked, TipBox, 2-200 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1	pack
732328	Pipette tips Ultra Low Retention, racked, TipBox, 2-200 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732348	Pipette tips Ultra Low Retention, racked, TipStack™, 2-200 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732368	Pipette tips Ultra Low Retention, racked, TipStack™, 2-200 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack



Ultra Low Retention pipette tips, 5 - 300 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732310	Pipette tips Ultra Low Retention, racked, TipBox, 5-300 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1	pack
732330	Pipette tips Ultra Low Retention, racked, TipBox, 5-300 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack



Ultra Low Retention pipette tips, 50 - 1000 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732312	Pipette tips Ultra Low Retention, racked, TipBox, 50-1000 µl, non-sterile, PP, IVD, colorless	1 pack = 480 piece(s)	1	pack
732332	Pipette tips Ultra Low Retention, racked, TipBox, 50-1000 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732352	Pipette tips Ultra Low Retention, racked, TipStack™, 50-1000 µl, non-sterile, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack
732372	Pipette tips Ultra Low Retention, racked, TipStack™, 50-1000 µl, BIO-CERT®, PP, IVD, colorless	1 pack = 960 piece(s)	1	pack



Filter tips, 0,1 - 1 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732502	Filter tips, bulk, 0,1-1 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732602	Filter tips, racked, TipRack, 0,1-1 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732622	Filter tips, racked, TipRack, 0,1-1 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732702	Filter tips, racked, TipBox, 0,1-1 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732722	Filter tips, racked, TipBox, 0,1-1 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack

Filter tips, 0,5 - 10 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732504	Filter tips, bulk, 0,5-10 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732604	Filter tips, racked, TipRack, 0,5-10 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732624	Filter tips, racked, TipRack, 0,5-10 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732704	Filter tips, racked, TipBox, 0,5-10 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732724	Filter tips, racked, TipBox, 0,5-10 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Filter tips, 1 - 20 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732506	Filter tips, bulk, 1-20 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732606	Filter tips, racked, TipRack, 1-20 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732626	Filter tips, racked, TipRack, 1-20 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732706	Filter tips, racked, TipBox, 1-20 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732726	Filter tips, racked, TipBox, 1-20 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Filter tips, 2 - 20 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732508	Filter tips, bulk, 2-20 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732608	Filter tips, racked, TipRack, 2-20 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732628	Filter tips, racked, TipRack, 2-20 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732708	Filter tips, racked, TipBox, 2-20 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732728	Filter tips, racked, TipBox, 2-20 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Filter tips, 5 - 100 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732510	Filter tips, bulk, 5-100 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732610	Filter tips, racked, TipRack, 5-100 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732630	Filter tips, racked, TipRack, 5-100 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732710	Filter tips, racked, TipBox, 5-100 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732730	Filter tips, racked, TipBox, 5-100 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack





Filter tips, 5 - 200 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732512	Filter tips, bulk, 5-200 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732612	Filter tips, racked, TipRack, 5-200 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732632	Filter tips, racked, TipRack, 5-200 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732712	Filter tips, racked, TipBox, 5-200 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732732	Filter tips, racked, TipBox, 5-200 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Filter tips, 50 - 1000 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732514	Filter tips, bulk, 50-1000 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732614	Filter tips, racked, TipRack, 50-1000 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732634	Filter tips, racked, TipRack, 50-1000 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack
732714	Filter tips, racked, TipBox, 50-1000 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732734	Filter tips, racked, TipBox, 50-1000 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Ultra Low Retention filter tips, 0,1 - 1 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732802	Filter tips Ultra Low Retention, racked, TipBox, 0,1- 1 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732822	Filter tips Ultra Low Retention, racked, TipBox, 0,1- 1 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Ultra Low Retention filter tips, 0,5 - 10 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732804	Filter tips Ultra Low Retention, racked, TipBox, 0,5-10 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732824	Filter tips Ultra Low Retention, racked, TipBox, 0,5-10 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Ultra Low Retention filter tips, 1 - 20 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732806	Filter tips Ultra Low Retention, racked, TipBox, 1-20 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732826	Filter tips Ultra Low Retention, racked, TipBox, 1-20 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack

Ultra Low Retention filter tips, 2 - 20 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732808	Filter tips Ultra Low Retention, racked, TipBox, 2-20 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732828	Filter tips Ultra Low Retention, racked, TipBox, 2-20 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Ultra Low Retention filter tips, 5 - 100 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732810	Filter tips Ultra Low Retention, racked, TipBox, 5-100 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732830	Filter tips Ultra Low Retention, racked, TipBox, 5-100 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Ultra Low Retention filter tips, 5 - 200 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732812	Filter tips Ultra Low Retention, racked, TipBox, 5-200 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732832	Filter tips Ultra Low Retention, racked, TipBox, 5-200 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



Ultra Low Retention filter tips, 50 - 1000 µl, PP, CE-IVD

Cat. No.	Description		From*	Unit
732814	Filter tips Ultra Low Retention, racked, TipBox, 50-1000 µl, non-sterile, PP/PE-filter, IVD, colorless	1 pack = 480 piece(s)	1	pack
732834	Filter tips Ultra Low Retention, racked, TipBox, 50-1000 µl, BIO-CERT®, PP/PE-filter, IVD, colorless	1 pack = 960 piece(s)	1	pack



TipBox, empty, with tip tray

Cat. No.	Description		From*	Unit
732990	TipBox, empty, with gray tip-tray for tips up to 50 µl		1	piece(s)
732992	TipBox, empty, with yellow tip-tray for 200 µl tips		1	piece(s)
732994	TipBox, empty, with green tip-tray for 300 µl tips		1	piece(s)
732996	TipBox, empty, with blue tip-tray for 1000 µl tips		1	piece(s)



Pipette leak testing unit PLT unit



Cat. No.	Description	From*	Unit
703970	Pipette leak testing unit PLT unit, with AC adapter 100-240V/50-60Hz	1	piece(s)

AC adapter for PLT unit



Cat. No.	Description	From*	Unit
703979	Universal AC adapter for PLT unit	1	piece(s)

Pipette adapter for PLT unit



Cat. No.	Description	From*	Unit
703975	Pipette adapter for PLT unit single-channel, for testing with tips	1	piece(s)
703976	Pipette adapter for PLT unit single-channel, for testing w/o tips	1	piece(s)
703977	Pipette adapter for PLT unit multichannel, for testing with and w/o tips	1	piece(s)

Filter for PLT unit



Cat. No.	Description	From*	Unit
703978	Filter for pipette adapter PLT unit 1 pack = 10 piece(s)	1	pack

PLTconnect software, for PLT unit



Cat. No.	Description	From*	Unit
703980	Software for PLT unit	1	piece(s)

Positive displacement pipette Transferpettor, Digital, DE-M

Cat. No.	Description	From*	Unit
701807	Positive displacement pipette Transferpettor Digital type, DE-M, 2,5- 10 µl, with glass capillaries	1	piece(s)
701812	Positive displacement pipette Transferpettor Digital type, DE-M, 5 - 25 µl, with glass capillaries	1	piece(s)
701817	Positive displacement pipette Transferpettor Digital type, DE-M, 10 - 50 µl, with glass capillaries	1	piece(s)
701822	Positive displacement pipette Transferpettor Digital type, DE-M, 20 - 100 µl, with glass capillaries	1	piece(s)
702804	Positive displacement pipette Transferpettor Digital type, DE-M, 100 - 500 µl, with caps made of PP	1	piece(s)
702806	Positive displacement pipette Transferpettor Digital type, DE-M, 200 - 1000 µl, with caps made of PP	1	piece(s)
702810	Positive displacement pipette Transferpettor Digital type, DE-M, 1,0- 5,0 ml, with caps made of PP	1	piece(s)
702812	Positive displacement pipette Transferpettor Digital type, DE-M, 2,0-10,0 ml, with caps made of PP	1	piece(s)



Positive displacement pipette Transferpettor, Fixed-volume, DE-M

Cat. No.	Description	From*	Unit
701842	Positive displacement pipette Transferpettor Fix type, DE-M, 1 µl, with glass capillaries	1	piece(s)
701844	Positive displacement pipette Transferpettor Fix type, DE-M, 2 µl, with glass capillaries	1	piece(s)
701853	Positive displacement pipette Transferpettor Fix type, DE-M, 5 µl, with glass capillaries	1	piece(s)
701858	Positive displacement pipette Transferpettor Fix type, DE-M, 10 µl, with glass capillaries	1	piece(s)
701863	Positive displacement pipette Transferpettor Fix type, DE-M, 20 µl, with glass capillaries	1	piece(s)
701864	Positive displacement pipette Transferpettor Fix type, DE-M, 25 µl, with glass capillaries	1	piece(s)
701868	Positive displacement pipette Transferpettor Fix type, DE-M, 50 µl, with glass capillaries	1	piece(s)
701873	Positive displacement pipette Transferpettor Fix type, DE-M, 100 µl, with glass capillaries	1	piece(s)
701878	Positive displacement pipette Transferpettor Fix type, DE-M, 200 µl, with glass capillaries	1	piece(s)



Station for Transferpettor

Cat. No.	Description	From*	Unit
701960	Transferpettor station for 4 instruments up to 200 µl	1	piece(s)
702890	Transferpettor station for 2 instruments, 0,5 - 10 ml	1	piece(s)



Caps for Transferpettor



Cat. No.	Description	From*	Unit
701900	Transferpettor caps, 5 µl, white glass cap. with color code	100	piece(s)
701902	Transferpettor caps, 10 µl, orange glass cap. with color code	100	piece(s)
701904	Transferpettor caps, 20 µl, black glass cap. with color code	100	piece(s)
701906	Transferpettor caps, 25 µl, 2xwhite glass cap. with color code	100	piece(s)
701908	Transferpettor caps, 50 µl, green glass cap. with color code	100	piece(s)
701910	Transferpettor caps, 100+200 µl, blue/red glass cap. color code	50	piece(s)
702852	Transferpettor-caps PP, for 100 - 500 µl, color code green	10	piece(s)
702854	Transferpettor-caps PP, for 200 -1000 µl, color code yellow	10	piece(s)
702858	Transferpettor-caps PP, for 1,0 - 5,0 ml, color code red	10	piece(s)
702860	Transferpettor-caps PP, for 2,0 -10,0 ml, color code orange	10	piece(s)

Seals for Transferpettor



Cat. No.	Description	From*	Unit
701920	Seals PTFE for Transferpettor, 20/25 µl 1 pack = 3 piece(s)	1	pack
701922	Seals PTFE for Transferpettor, 50 µl 1 pack = 3 piece(s)	1	pack
701924	Seals PTFE for Transferpettor, 100/200 µl 1 pack = 3 piece(s)	1	pack
702864	Seals, PE, for Transferpettor, for 100 - 500 µl (10 pcs.) 1 pack = 10 piece(s)	1	pack
702866	Seals, PE, for Transferpettor, for 200 -1000 µl (10 pcs.) 1 pack = 10 piece(s)	1	pack
702870	Seals, PE, for Transferpettor, for 1,0 - 5,0 ml (10 pcs.) 1 pack = 10 piece(s)	1	pack
702872	Seals, PE, for Transferpettor, for 2,0 -10,0 ml (10 pcs.) 1 pack = 10 piece(s)	1	pack

Combi-pack caps and seal for Transferpettor

Cat. No.	Description	From*	Unit
702883	Combi-Pack Caps (PP, 2 pieces) and Seal (PE, 1 piece), 100-500 µl	1	pack
702884	Combi-Pack Caps (PP, 2 pieces) and Seal (PE, 1 piece), 200-1000 µl	1	pack
702885	Combi-Pack Caps (PP, 2 pieces) and Seal (PE, 1 piece), 1000-5000 µl	1	pack
702886	Combi-Pack Caps (PP, 2 pieces) and Seal (PE, 1 piece), 2000-10000 µl	1	pack

Repair set for Transferpettor



Cat. No.	Description	From*	Unit
701964	Repair set for Transferpettor, 1/2/3/4/5 µl	1	pack
701965	Repair set for Transferpettor, 10 µl	1	pack
701966	Repair set for Transferpettor, 20/25 µl	1	pack
701967	Repair set for Transferpettor, 50 µl	1	pack
701968	Repair set for Transferpettor, 100/200 µl	1	pack

Piston rods for Transferpettor

Cat. No.	Description		From*	Unit
701928	Piston rod for Transferpettor, Fix type 1/2/3/4/5 µl, stainless steel	1 pack = 3 piece(s)	1	pack
701930	Piston rod for Transferpettor, Fix+Dig. type 10 µl, stainless steel	1 pack = 3 piece(s)	1	pack
701932	Piston rod for Transferpettor, Fix type 10 µl/Dig.type 25 µl stainless steel	1 pack = 3 piece(s)	1	pack
701934	Piston rod for Transferpettor, Fix+Dig. type 50 µl, stainless steel	1 pack = 3 piece(s)	1	pack
701936	Piston rod for Transferpettor, Fix+Dig. type 100 µl, stainless steel	1 pack = 3 piece(s)	1	pack
701938	Piston rod for Transferpettor, Fix type 200 µl, stainless steel	1 pack = 3 piece(s)	1	pack



Repetitive pipette HandyStep® touch, DE-M

Cat. No.	Description		From*	Unit
705200	HandyStep® touch, DE-M, inclusive universal AC adapter and shelf/rack mount		1	piece(s)
705201	Package HandyStep® touch, DE-M, inclusive universal AC adapter and charging stand		1	piece(s)
705210	HandyStep® touch S, DE-M, inclusive universal AC adapter and shelf/rack mount		1	piece(s)
705211	Package HandyStep® touch S, DE-M, inclusive universal AC adapter and charging stand		1	piece(s)



Charging stand for HandyStep® touch

Cat. No.	Description		From*	Unit
705220	Charging stand for HandyStep® touch and HandyStep® touch S, w/o universal AC adapter		1	piece(s)



Support stand for HandyStep® touch

Cat. No.	Description		From*	Unit
705230	Support stand for HandyStep® touch und HandyStep® touch S		1	piece(s)



Universal AC adapter for HandyStep® touch

Cat. No.	Description		From*	Unit
705250	Universal AC adapter for HandyStep® touch and HandyStep® touch S		1	piece(s)



Universal holder for HandyStep® touch



Cat. No.	Description	From*	Unit
705235	Universal holder for HandyStep® touch and HandyStep® touch S	1	piece(s)

Lithium ion battery for HandyStep® touch



Cat. No.	Description	From*	Unit
705225	Lithium ion battery for HandyStep® touch and HandyStep® touch S	1	piece(s)

Repetitive pipette HandyStep® S, DE-M



Cat. No.	Description	From*	Unit
705110	Repetitive pipette HandyStep® S	1	piece(s)

Shelf/rack mount for HandyStep® S



Cat. No.	Description	From*	Unit
705130	Shelf/Rack mount for HandyStep® S	1	piece(s)

Dispenser tips PD-Tips II, non-sterile, bulk



Cat. No.	Description	From*	Unit
705700	PD-Tips II, 0,1 ml, non-sterile, bulk, cylinder PP / piston LCP, type encoded	100	piece(s)
705702	PD-Tips II, 0,5 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	100	piece(s)
705704	PD-Tips II, 1 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	100	piece(s)
705706	PD-Tips II, 1,25 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	100	piece(s)
705708	PD-Tips II, 2,5 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	100	piece(s)
705710	PD-Tips II, 5 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	100	piece(s)
705712	PD-Tips II, 10 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	100	piece(s)
705714	PD-Tips II, 12,5 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	100	piece(s)
705716	PD-Tips II, 25 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	50	piece(s)
705718	PD-Tips II, 50 ml, non-sterile, bulk, cylinder PP / piston PE-HD, type encoded	25	piece(s)

Dispenser tips PD-Tips II, BIO-CERT®, individually wrapped

Cat. No.	Description	From*	Unit
705730	PD-Tips II, 0,1 ml, BIO-CERT®, individually wrapped, PP/LCP, type encoded	100	piece(s)
705732	PD-Tips II, 0,5 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	100	piece(s)
705734	PD-Tips II, 1 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	100	piece(s)
705736	PD-Tips II, 1,25 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	100	piece(s)
705738	PD-Tips II, 2,5 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	100	piece(s)
705740	PD-Tips II, 5 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	100	piece(s)
705742	PD-Tips II, 10 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	100	piece(s)
705744	PD-Tips II, 12,5 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	100	piece(s)
705746	PD-Tips II, 25 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	25	piece(s)
705748	PD-Tips II, 50 ml, BIO-CERT®, individually wrapped, PP/PE-HD, type encoded	25	piece(s)



Dispenser tips PD-Tips II, Set non-sterile

Cat. No.	Description	From*	Unit
705720	Set PD-Tips II, 20 pcs. of each size 0,5, 1, 1,25, 2,5, 5, 10 und 12,5 ml, non-sterile, bulk, type encoded	1 pack = 140 piece(s)	1



fAdapter or PD-Tips II 25 / 50 ml

Cat. No.	Description	From*	Unit
702398	Adapter for PD-Tips II, 25-50ml, PP, non-sterile	1 pack = 10 piece(s)	1
702399	Adapter for PD-Tips II, 25-50ml, PP, BIO-CERT®	5	piece(s)



Pipetting aid accu-jet® pro



Cat. No.	Description	From*	Unit
26300	Pipette controller accu-jet® pro, dark blue, AC adapter Europe 100-240 V/50-60 Hz	1	piece(s)
26301	Pipette controller accu-jet® pro, magenta, AC adapter Europe 100-240 V/50-60 Hz	1	piece(s)
26302	Pipette controller accu-jet® pro, green, AC adapter Europe 100-240 V/50-60 Hz	1	piece(s)
26303	Pipette controller accu-jet® pro, royal blue, AC adapter Europe 100-240 V/50-60 Hz	1	piece(s)
26304	Pipette controller accu-jet® pro, dark blue, w/o AC adapter	1	piece(s)
26310	Pipette controller accu-jet® pro, dark blue, AC adapter UK/Ireland 100-240 V/50-60 Hz	1	piece(s)
26311	Pipette controller accu-jet® pro, magenta, AC adapter UK/Ireland 100-240 V/50-60 Hz	1	piece(s)
26312	Pipette controller accu-jet® pro, green, AC adapter UK/Ireland 100-240 V/50-60 Hz	1	piece(s)
26313	Pipette controller accu-jet® pro, royal blue, AC adapter UK/Ireland 100-240 V/50-60 Hz	1	piece(s)
26320	Pipette controller accu-jet® pro, dark blue, AC adapter Australia 100-240 V/50-60 Hz	1	piece(s)
26321	Pipette controller accu-jet® pro, magenta, AC adapter Australia 100-240 V/50-60 Hz	1	piece(s)
26322	Pipette controller accu-jet® pro, green, AC adapter Australia 100-240 V/50-60 Hz	1	piece(s)
26323	Pipette controller accu-jet® pro, royal blue, AC adapter Australia 100-240 V/50-60 Hz	1	piece(s)
26330	Pipette controller accu-jet® pro, dark blue, AC adapter USA 100-240 V/50-60 Hz	1	piece(s)
26331	Pipette controller accu-jet® pro, magenta, AC adapter USA 100-240 V/50-60 Hz	1	piece(s)
26332	Pipette controller accu-jet® pro, green, AC adapter USA 100-240 V/50-60 Hz	1	piece(s)
26333	Pipette controller accu-jet® pro, royal blue, AC adapter USA 100-240 V/50-60 Hz	1	piece(s)
26340	Pipette controller accu-jet® pro, dark blue, AC adapter Japan 100-240 V/50-60 Hz	1	piece(s)
26341	Pipette controller accu-jet® pro, magenta, AC adapter Japan 100-240 V/50-60 Hz	1	piece(s)
26342	Pipette controller accu-jet® pro, green, AC adapter Japan 100-240 V/50-60 Hz	1	piece(s)
26343	Pipette controller accu-jet® pro, royal blue, AC adapter Japan 100-240 V/50-60 Hz	1	piece(s)

Pipetting aid macro



Cat. No.	Description	From*	Unit
26200	macro pipette controller for pipettes 0.1-200 ml, grey, with spare membrane filter	1	piece(s)
26201	macro pipette controller for pipettes 0.1-200 ml, green, with spare membrane filter	1	piece(s)
26202	macro pipette controller for pipettes 0.1-200 ml, blue, with spare membrane filter	1	piece(s)
26203	macro pipette controller for pipettes 0.1-200 ml, magenta, with spare membrane filter	1	piece(s)

Pipetting aid micro



Cat. No.	Description	From*	Unit
25800	Micro pipette controller for pipette up to 1 ml and disposable micro pipette with ringmark	1	piece(s)
25900	Pipette controller micro-classic for small-vol. pip. with 2 spare suction tubes	1	piece(s)

Pipette filler

Cat. No.	Description	From*	Unit
25300	Pipette filler, NR, with 3 valves, standard model, for pip. up to, 10 ml	1	piece(s)
25315	Pipette filler, NR, with 3 valves, universal model, for pip. up to 100 ml	1	piece(s)
25400	Pipette filler, NR, with 2 valves, Flip model, for pip. up to 100 ml	1	piece(s)



Membrane filter for accu-jet® pro/ macro

Cat. No.	Description	From*	Unit
26052	Macro pipette controller/membrane filt. 3 µm hydrophobic non-sterile, blister pack	1	piece(s)
26056	Macro pipette controller/membrane filt. 3 µm hydrophobic non-sterile	1	pack
26530	Membrane filter 0,2 µm, for accu-jet® pro pipette controller, ETO-sterile SAL 10EXP-6	1	piece(s)
26535	Membrane filter 0,2 µm, for accu-jet® pro pipette controller, non sterile	1	pack



AC adapters for accu-jet® pro

Cat. No.	Description	From*	Unit
26601	Recharger for accu-jet® pro, Europe 100-240 V/50-60 Hz	1	piece(s)
26602	Recharger for accu-jet® pro, UK/Ireland 100-240 V/50-60 Hz	1	piece(s)
26603	Recharger for accu-jet® pro, USA, 100-240 V/50-60 Hz	1	piece(s)
26604	Recharger for accu-jet® pro, Japan 100-240 V/50-60 Hz	1	piece(s)
26605	Recharger for accu-jet® pro, Australia 100-240 V/50-60 Hz	1	piece(s)




Adapter support for accu-jet® pro/ macro

Cat. No.	Description	From*	Unit
26220	Macro pipette controller/adappter support, total length 49 mm, PP, grey	1	piece(s)
26652	Adapter support, PP, for accu-jet® pro, dark blue	1	piece(s)
26653	Adapter support, PP, for accu-jet® pro, magenta	1	piece(s)
26654	Adapter support, PP, for accu-jet® pro, green	1	piece(s)
26655	Adapter support, PP, for accu-jet® pro, royal blue	1	piece(s)



Silicone adapter for micro-classic

Cat. No.	Description	From*	Unit
25931	Spare adapter for micro-classic, with suction tube	1	pack



Wall support for accu-jet® pro



Cat. No.	Description	From*	Unit
26540	Wall support for accu-jet® pro pipette controller, single	1	piece(s)

Snap tapes for accu-jet® pro



Cat. No.	Description	From*	Unit
26542	Snap tapes, 12 X 65 mm, 2 sets for accu-jet® pro pipette controller	1	pack

NiMH accu pack for accu-jet® pro



Cat. No.	Description	From*	Unit
26630	NiMH accu pack for accu-jet pro; 2.4 V / 730 mAh	1	piece(s)

Valve system for macro pipette controller



Cat. No.	Description	From*	Unit
26128	Valve system for macro pipette controller, PP, PTFE, silicone	1	piece(s)

Suction bellows for macro pipette controller



Cat. No.	Description	From*	Unit
26225	Macro pipette controller/suction bellows, silicone, with screw ring, PP	1	piece(s)

Adapter for accu-jet® pro/ macro



Cat. No.	Description	From*	Unit
26146	Adapter silicone for macro pipette controller, 44 mm	1	piece(s)
26508	Pipette adapter with non-return valve, silicone, for accu-jet® pro	1	piece(s)

Replacement suction system for micro pipette controller



Cat. No.	Description	From*	Unit
25805	Suction system for micro pipette controller, silicone	1	pack


EASYCAL™ test tubes

Cat. No.	Description		From*	Unit
708462	EASYCAL™ testing tubes for pipettes < 50 µl	1 pack = 250 piece(s)	1	pack




Lid for micro-weighing container

Cat. No.	Description		From*	Unit
708472	Set of lids for micro weighing vessel, 3 spare lids	1 pack = 3 piece(s)	1	pack



Tubes for micro-weighing container

Cat. No.	Description		From*	Unit
708470	Micro-weighing container incl. 10 filters and 3 cover caps		1	pack



Filter pack for micro-weighing container

Cat. No.	Description		From*	Unit
708471	Filter pack for micro weighing vessel, 20 spare filters	1 pack = 20 piece(s)	1	pack

Pipetting robot Liquid Handling Station



Cat. No.	Description	From*	Unit
709400	Liquid Handling Station pipetting robot	1	piece(s)
709402	Liquid Handling Station flow pipetting robot with FlowBox	1	piece(s)

Holder for Liquid Ends



Cat. No.	Description	From*	Unit
709463	Holder for 3 Liquid Ends, POM, single/multichannel	1	piece(s)
709465	Holder for 5 Liquid Ends, POM, single/multichannel	1	piece(s)

Liquid Ends for Liquid Handling Station



Cat. No.	Description	From*	Unit
709410	Liquid End für Liquid Handling Station, 1-Channel, 1 - 50 µl	1	piece(s)
709413	Liquid End für Liquid Handling Station, 1-Channel, 10 - 200 µl	1	piece(s)
709416	Liquid End für Liquid Handling Station, 1-Channel, 10 - 200 µl	1	piece(s)
709420	Liquid End für Liquid Handling Station, 8-Channel, 1 - 50 µl	1	piece(s)
709423	Liquid End für Liquid Handling Station, 8-Channel, 20 - 300 µl	1	piece(s)

Adapter for Liquid Handling Station



Cat. No.	Description	From*	Unit
709430	Height Adapter for Liquid Handling Station, 60 mm, Aluminum	1	piece(s)
709432	Height Adapter for Liquid Handling Station, 30 mm, Aluminum	1	piece(s)
709434	Tip Adapter for Liquid Handling Station, Aluminum	1	piece(s)
709436	Cuvette Rack for Liquid Handling Station suitable for 16 (BRAND) cuvettes, POM	1	piece(s)
709443	Reservoir Rack for Liquid Handling Station for 4 Reservoirs à 40 ml (701460/62), Aluminum	1	piece(s)
709446	PCR Adapter for Liquid Handling Station, 96-well, POM	1	piece(s)
709448	PCR Adapter for Liquid Handling Station, 384-well, optimized for BRAND plates (781345, 781347 und 781358), POM	1	piece(s)
709449	PCR Adapter for Liquid Handling Station, 384-well, universal, POM	1	piece(s)
709450	Microtube Rack for Liquid Handling Station suitable for 24 tubes of 1,5 ml each, POM	1	piece(s)
709452	Microtube Rack for Liquid Handling Station suitable for 24 tubes of 0,5 ml each, POM	1	piece(s)
709453	Microtube Rack for Liquid Handling Station suitable for 9 tubes of 5 ml each, POM	1	piece(s)
709454	Microtube Rack for Liquid Handling Station suitable for 24 tubes of 0,5 ml each, POM	1	piece(s)
709455	Tube Rack für Liquid Handling Station für 24 Sample Tubes 12 x 75mm, POM	1	piece(s)
709456	PCR Cooler Rack for Liquid Handling Station, 96-well	1	piece(s)

Robotic tips

Cat. No.	Description		From*	Unit
732146	Robotic tips, TipRack, PP, 1-50 µl, non-sterile, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732148	Robotic tips, TipRack, PP, 10-200 µl, non-sterile, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732150	Robotic tips, TipRack, PP, 10-300 µl, non-sterile, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732152	Robotic tips, TipRack, PP, 40-1000 µl, non-sterile, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack



Robotic filter tips

Cat. No.	Description		From*	Unit
732646	Robotic filter tips, TipRack, PP, 1-20 µl, non-sterile, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732650	Robotic filter tips, TipRack, PP, 10-100 µl, non-sterile, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732652	Robotic filter tips, TipRack, PP, 10-200 µl, non-sterile, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732654	Robotic filter tips, TipRack, PP, 40-1000 µl, non-sterile, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack



Robotic tips, BIO-CERT®

Cat. No.	Description		From*	Unit
732166	Robotic tips, TipRack, PP, 1-50 µl, BIO-CERT®, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732168	Robotic tips, TipRack, PP, 10-200 µl, BIO-CERT®, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732170	Robotic tips, TipRack, PP, 10-300 µl, BIO-CERT®, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732172	Robotic tips, TipRack, PP, 40-1000 µl, BIO-CERT®, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack



Robotic filter tips, BIO-CERT®

Cat. No.	Description		From*	Unit
732666	Robotic filter tips, TipRack, PP, 1-20 µl, BIO-CERT®, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732670	Robotic filter tips, TipRack, PP, 10-100 µl, BIO-CERT®, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732672	Robotic filter tips, TipRack, PP, 10-200 µl, BIO-CERT®, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack
732674	Robotic filter tips, TipRack, PP, 40-1000 µl, BIO-CERT®, DNA-, RNase-free	1 pack = 960 piece(s)	1	pack



Waste container for Liquid Handling Station

Cat. No.	Description		From*	Unit
709458	Waste container for Liquid Handling Station, PET		1	piece(s)



Temperature controlled shaker Teleshake 95



Cat. No.	Description	From*	Unit
709470	Teleshake 95 with fixed mounted basic plate for Liquid Handling Station	1	piece(s)

Controller for Teleshake 95



Cat. No.	Description	From*	Unit
709472	Single TEC Control + 1 black slot module, for Teleshake 95	1	piece(s)
709474	Multi TEC Control + 2 black slot modules, for Teleshake 95	1	piece(s)

Adapter for Teleshake 95



Cat. No.	Description	From*	Unit
709477	Universal adapter for Teleshake 95 for 96-well PCR-plates, low profile, Aluminum	1	piece(s)
709478	Universal adapter for Teleshake 95 for 96-well PCR-plates, standard, Aluminum	1	piece(s)
709479	Universal adapter for Teleshake 95 for 384-well PCR-plates, Aluminum	1	piece(s)
709480	Adapter for Teleshake 95 for 96-well BRAND F-, C- bottom plates, Aluminum	1	piece(s)
709481	Adapter for Teleshake 95 for 96-well BRAND V-bottom plates, Aluminum	1	piece(s)
709482	Adapter for Teleshake 95 for 96-well BRAND U-bottom plates, Aluminum	1	piece(s)

Reagent reservoirs, PP



Cat. No.	Description		From*	Unit
701450	Reservoir PP, 1-columns 220 ml, pyramidal bottom 96-well, DNase-, RNase free	1 pack = 10 piece(s)	1	pack
701452	Reservoir PP, 12-columns low profil, 6 ml, pyramidal bottom, DNase-, RNase free	1 pack = 10 piece(s)	1	pack
701454	Reservoir PP, 4-columns, 60 ml, pyramidal bottom, DNase-, RNase free	1 pack = 10 piece(s)	1	pack
701456	Reservoir PP, 6-columns, 40 ml, pyramidal bottom, DNase-, RNase free	1 pack = 10 piece(s)	1	pack
701458	Reservoir PP, 1 column low profile, 50 ml, flat bottom, DNase-, RNase free	1 pack = 10 piece(s)	1	pack

Reagent reservoirs with cover, PP



Cat. No.	Description		From*	Unit
701460	Reservoir with cover, PP, 40ml, graduated, non-sterile	1 pack = 24 piece(s)	1	pack
701462	Reservoir with cover, PP, 40ml, graduated, sterile	1 pack = 24 piece(s)	1	pack

Microcentrifuge tubes disposable, PP, 0,5 ml

Cat. No.	Description	From*	Unit
780507	Disposable microcentrifuge tube, PP, 0,5 ml, translucent, lid	1000	piece(s)
780536	Microcentrifuge tube PP 0,5 ml clear with lid locking	500	piece(s)



Microcentrifuge tubes disposable, PP, 1,5 ml

Cat. No.	Description	From*	Unit
780400	Disposable microcentrifuge tube, BIO-CERT®, PP 1,5 ml, neutral, lid	450	piece(s)
780500	Disposable microcentrifuge tube, PP, 1,5 ml, neutral, lid	500	piece(s)
780502	Disposable microcentrifuge tube, PP, 1,5 ml, neutral, lid	3000	piece(s)
780505	Disposable microcentrifuge tube, PP, 1,5 ml, neutral, w/o lid	2000	piece(s)
780521	Disposable microcentrifuge tube, PP, 1,5 ml, yellow, lid RCF max. 20000 at 20°C	500	piece(s)
780522	Disposable microcentrifuge tube, PP, 1,5 ml, blue, lid RCF max. 20000 at 20°C	500	piece(s)
780523	Disposable microcentrifuge tube, PP, 1,5 ml, green, lid RCF max. 20000 at 20°C	500	piece(s)
780524	Disposable microcentrifuge tube, PP, 1,5 ml, orange lid RCF max. 20000 at 20°C	500	piece(s)
780525	Disposable microcentrifuge tube, PP, 1,5 ml, amber, lid RCF max. 20000 at 20°C	500	piece(s)
780540	Disposable microcentrifuge tube, PP, 1,5 ml clear, with lid closure	1000	piece(s)



Microcentrifuge tubes disposable, PP, 2 ml

Cat. No.	Description	From*	Unit
780546	Disposable microcentrifuge tube, PP, 2,0 ml, transparent, with lid locking	500	piece(s)
780550	Disposable microcentrifuge tube, PP, 2,0 ml, transparent, with lid	500	piece(s)




Microcentrifuge tubes, PP, with screw cap, sterile

Cat. No.	Description	From*	Unit
780750	Micro tube (PP) with bulk screw cap (PP) with silicone sealing, 0,5 ml, self-standing, sterile	500	piece(s)
780751	Micro tube (PP) with bulk screw cap (PP) with silicone sealing, 1,5 ml, self-standing, sterile	500	piece(s)
780752	Micro tube (PP) with bulk screw cap (PP) with silicone sealing, 1,5 ml, round bottom, sterile	500	piece(s)
780753	Micro tube (PP) with bulk screw cap (PP) with silicone sealing, 2 ml, self-standing, sterile	500	piece(s)
780754	Micro tube (PP) with bulk screw cap (PP) with silicone sealing, 2 ml, round bottom, sterile	500	piece(s)
780755	Micro tubes with tamper-evident screw cap, PP, 0,5 ml, self-standing, sterile	500	piece(s)
780756	Micro tubes with tamper-evident screw cap, PP, 1,5 ml, self-standing, sterile	500	piece(s)
780757	Micro tubes with tamper-evident screw cap, PP, 1,5 ml, round bottom, sterile	500	piece(s)
780758	Micro tubes with tamper-evident screw cap, PP, 2 ml, self-standing, sterile	500	piece(s)
780759	Micro tubes with tamper-evident screw cap, PP, 2 ml, round bottom, sterile	500	piece(s)



Microcentrifuge tubes, PP, with screw cap, graduated



Cat. No.	Description	From*	Unit
780710	Micro tube (PP) with attached screw cap (PP) with silicone seal, 0,5 ml, self-standing, non graduated	1000	piece(s)
780711	Micro tube (PP) with attached screw cap (PP) with silicone seal, 1,5 ml, self-standing, graduated	1000	piece(s)
780712	Micro tube (PP) with attached screw cap (PP) with silicone seal, 1,5 ml, round bottom, graduated	1000	piece(s)
780713	Micro tube (PP) with attached screw cap (PP) with silicone seal, 2 ml, self-standing, graduated	1000	piece(s)
780714	Micro tube (PP) with attached screw cap (PP) with silicone seal, 2 ml, round bottom, graduated	1000	piece(s)

Microcentrifuge tubes, PP, with screw cap, ungraduated



Cat. No.	Description	From*	Unit
780700	Micro tube (PP) with attached screw cap (PE) with sealing cone, 0,5 ml, self-standing	1000	piece(s)
780701	Micro tube (PP) with attached screw cap (PE) with sealing cone, 1,5 ml, self-standing	1000	piece(s)
780702	Micro tube (PP) with attached screw cap (PE) with sealing cone, 1,5 ml, round bottom	1000	piece(s)
780703	Micro tube (PP) with attached screw cap (PE) with sealing cone, 2 ml, self-standing	1000	piece(s)
780704	Micro tube (PP) with attached screw cap (PE) with sealing cone, 2 ml, round bottom	1000	piece(s)

Microcentrifuge tubes, without cap, graduated, PP



Cat. No.	Description	From*	Unit
780730	Micro tube (PP) for bulk screw cap (PP) with silicone sealing, 0,5 ml, self-standing, non graduated	1000	piece(s)
780731	Micro tube (PP) for bulk screw cap (PP) with silicone sealing, 1,5 ml, self-standing	1000	piece(s)
780732	Micro tube (PP) for bulk screw cap (PP) with silicone sealing, 1,5 ml, round bottom	1000	piece(s)
780733	Micro tube (PP) for bulk screw cap (PP) with silicone sealing, 2 ml, self-standing	1000	piece(s)
780734	Micro tube (PP) for bulk screw cap (PP) with silicone sealing, 2 ml, round bottom	1000	piece(s)

Microcentrifuge tubes, without cap, ungraduated, PP



Cat. No.	Description	From*	Unit
780760	Micro tubes for screw caps with tamper-evident screw cap, PP, 0,5 ml, self-standing	1000	piece(s)
780761	Micro tubes for screw caps with tamper-evident screw cap, PP, 1,5 ml, self-standing	1000	piece(s)
780762	Micro tubes for screw caps with tamper-evident screw cap, PP, 1,5 ml, round bottom	1000	piece(s)
780763	Micro tubes for screw caps with tamper-evident screw cap, PP, 2 ml, self-standing	1000	piece(s)
780764	Micro tubes for screw caps with tamper-evident screw cap, PP, 2 ml, round bottom	1000	piece(s)

Microcentrifuge tube rack, PP

Cat. No.	Description	From*	Unit
780605	Microcentrifuge tube rack, PP, 210 x 70 x 37 mm grey for 20 tubes 1,5 ml	1	piece(s)
4341000	Microcentrifuge tube rack, PP, 265 x 126 x 38 mm for 84 tubes up to diameter 13 mm, white	5	piece(s)
4341001	Microcentrifuge tube rack, PP, 265 x 126 x 38 mm for 84 tubes up to diameter 13 mm, blue	5	piece(s)
4341002	Microcentrifuge tube rack, PP, 265 x 126 x 38 mm for 84 tubes up to diameter 13 mm, red	5	piece(s)
4341050	Microcentrifuge tube rack, PP, 265 x 126 x 38 mm for 128 tubes up to diameter 11 mm, white	5	piece(s)
4341051	Microcentrifuge tube rack, PP, 265 x 126 x 38 mm for 128 tubes up to diameter 11 mm, blue	5	piece(s)
4341052	Microcentrifuge tube rack, PP, 265 x 126 x 38 mm for 128 tubes up to diameter 11 mm, red	5	piece(s)



Adapter, PP, for microcentrifuge tube rack

Cat. No.	Description	From*	Unit
780608	Adapter for microcentrifuge tube rack 780605, 0,5 ml 1 pack = 20 piece(s)	1	pack



Mini cooler, PC

Cat. No.	Description	From*	Unit
114930	Mini cooler, PC, for 12 tubes 0,5-2,0 ml, bench temp. maintained 0°C, time 60 min.	1	piece(s)
114935	Mini cooler, PC, for 12 tubes 0,5-2,0 ml, bench temp. maintained -20°C, time 60 min.	1	piece(s)
114940	Mini cooler, PC, for 12 tubes 0,5-2,0 ml, bench temp. maintained -70°C, time 45 min.	1	piece(s)



Screw cap, PP, with silicone seal, non-sterile

Cat. No.	Description	From*	Unit
780740	Screw cap (PP) for micro tubes for bulk screw cap, white	1000	piece(s)
780741	Screw cap (PP) for micro tubes for bulk screw cap, blue	1000	piece(s)
780742	Screw cap (PP) for micro tubes for bulk screw cap, red	1000	piece(s)
780743	Screw cap (PP) for micro tubes for bulk screw cap, green	1000	piece(s)
780744	Screw cap (PP) for micro tubes for bulk screw cap, yellow	1000	piece(s)



Screw scap, PP, tamper-evident




Cat. No.	Description	From*	Unit
780770	Tamper-evident screw caps (PP) with silicone seal, for micro tubes 780760-64, transparent	1000	piece(s)
780771	Tamper-evident screw caps (PP) with silicone seal, for micro tubes 780760-64, blue	1000	piece(s)
780772	Tamper-evident screw caps (PP) with silicone seal, for micro tubes 780760-64, green	1000	piece(s)
780773	Tamper-evident screw caps (PP) with silicone seal, for micro tubes 780760-64, purple	1000	piece(s)
780774	Tamper-evident screw caps (PP) with silicone seal, for micro tubes 780760-64, red	1000	piece(s)
780775	Tamper-evident screw caps (PP) with silicone seal, for micro tubes 780760-64, yellow	1000	piece(s)

Cap inserts, PP



Cat. No.	Description	From*	Unit
780720	Cap insert, PP, for micro tubes, for identification, for all non-tamper evident screw caps, white	500	piece(s)
780721	Cap insert, PP, for micro tubes, for identification, for all non-tamper evident screw caps, blue	500	piece(s)
780722	Cap insert, PP, for micro tubes, for identification, for all non-tamper evident screw caps, red	500	piece(s)
780723	Cap insert, PP, for micro tubes, for identification, for all non-tamper evident screw caps, green	500	piece(s)
780724	Cap insert, PP, for micro tubes, for identification, for all non-tamper evident screw caps, yellow	500	piece(s)

PCR Single tubes, PP, 0,2 ml



Cat. No.	Description	From*	Unit
781300	PCR tubes, PP, colorless, thin-walled, 0,2 ml, indiv., attached domed cap	1000	piece(s)
781301	PCR tubes, PP, rose, thin-walled, 0,2 ml, indiv., attached flat cap	1000	piece(s)
781302	PCR tubes, PP, yellow, thin-walled, 0,2 ml, indiv., attached flat cap	1000	piece(s)
781303	PCR tubes, PP, green, thin-walled, 0,2 ml, indiv., attached flat cap	1000	piece(s)
781304	PCR tubes, PP, blue, thin-walled, 0,2 ml, indiv., attached flat cap	1000	piece(s)
781305	PCR tubes, PP, colorless, thin-walled, 0,2 ml, indiv., attached flat cap	1000	piece(s)

PCR Single tubes, PP, 0,5 ml



Cat. No.	Description	From*	Unit
781310	PCR tubes, PP, colorless, thin-walled, 0,5 ml, indiv., attached flat cap	1000	piece(s)
781311	PCR tubes, PP, rose, thin-walled, 0,5 ml, indiv., attached flat cap	1000	piece(s)
781312	PCR tubes, PP, yellow, thin-walled, 0,5 ml, indiv., attached flat cap	1000	piece(s)
781313	PCR tubes, PP, green, thin-walled, 0,5 ml, indiv., attached flat cap	1000	piece(s)
781314	PCR tubes, PP, blue, thin-walled, 0,5 ml, indiv., attached flat cap	1000	piece(s)

PCR strips, Strips of 8, PP, without cap

Cat. No.	Description		From*	Unit
781320	Strips of 8 PCR tubes, PP, colorless, thin-walled, 0,2 ml, w/o caps	1 pack = 125 piece(s)	1	pack
781321	Strips of 8 PCR tubes, PP, rose, thin-walled, 0,2 ml, w/o caps	1 pack = 125 piece(s)	1	pack
781322	Strips of 8 PCR tubes, PP, yellow, thin-walled, 0,2 ml, w/o caps	1 pack = 125 piece(s)	1	pack
781323	Strips of 8 PCR tubes, PP, green, thin-walled, 0,2 ml, w/o caps	1 pack = 125 piece(s)	1	pack
781324	Strips of 8 PCR tubes, PP, blue, thin-walled, 0,2 ml, w/o caps	1 pack = 125 piece(s)	1	pack
781325	Strips of 8 PCR tubes, PP, white, 0,2 ml, qPCR, separate lid strips	1 pack = 125 piece(s)	1	pack



PCR strips, Cap strips of 8

Cat. No.	Description		From*	Unit
781334	Strips of 8 caps, PP, clear, for PCR tubes, flat	1 pack = 125 piece(s)	1	pack
781340	Strips of 8 caps, PP, clear, for PCR tubes, domed	1 pack = 125 piece(s)	1	pack
781341	Strips of 8 caps, PP, rose, for PCR tubes, domed	1 pack = 125 piece(s)	1	pack
781342	Strips of 8 caps, PP, yellow, for PCR tubes, domed	1 pack = 125 piece(s)	1	pack
781343	Strips of 8 caps, PP, green, for PCR tubes, domed	1 pack = 125 piece(s)	1	pack
781344	Strips of 8 caps, PP, blue, for PCR tubes, domed	1 pack = 125 piece(s)	1	pack



PCR strips, Combi-pack of 8, detached cap strip

Cat. No.	Description		From*	Unit
781326	Strips of 8 PCR tubes, PP, with flat caps for qPCR, 250 tubes + 250 caps	1 pack = 250 piece(s)	1	pack
781327	Strips of 8 PCR tubes, PP, with domed caps, 250 tubes + 250 caps	1 pack = 250 piece(s)	1	pack



PCR strips, Strips of 8, PP, attached individual caps, 1 ridge

Cat. No.	Description		From*	Unit
781332	Strips of 8 PCR tubes, PP, attached single caps, 0,2 ml, for RealTime-PCR	1 pack = 120 piece(s)	1	pack
781333	Strips of 8 PCR tubes, PP, attached single caps, 0,15ml, Low Profile	1 pack = 120 piece(s)	1	pack



PCR strips, Strips of 8, PP, attached individual caps, 3 ridges




Cat. No.	Description		From*	Unit
781315	Strips of 8 PCR tubes with attached individual flat caps, 3 ridges between tubes, 0.2 ml, standard profile	1 pack = 120 piece(s)	1	pack
781316	Strips of 8 PCR tubes, PP, attached transparent single caps, 0,2 ml, Standard, white tubes	1 pack = 120 piece(s)	1	pack
781317	Strips of 8 PCR tubes, PP, attached transparent single caps, 0.15 ml, Low Profile, white tubes	1 pack = 120 piece(s)	1	pack
781318	Strips of 8 PCR tubes with attached individual flat caps, 3 ridges between tubes, 0.15 ml, low profile	1 pack = 120 piece(s)	1	pack

PCR strips, Strips of 8, PP, attached cap strip



Cat. No.	Description		From*	Unit
781330	Strips of 8 PCR tubes, PP, colorless, thin-walled, 0,2 ml, attached cap strips	1 pack = 125 piece(s)	1	pack

PCR strips, Strips of 12, PP, without cap



Cat. No.	Description		From*	Unit
781280	Strips of 12 PCR tubes, PP, colorless, thin-walled, 0,2 ml, w/o caps	1 pack = 125 piece(s)	1	pack

PCR strips, Cap strips of 12, PP, domed cap




Cat. No.	Description		From*	Unit
781290	Strips of 12 caps, PP, clear, for PCR tubes, domed	1 pack = 125 piece(s)	1	pack

PCR plate, 24-well, non-skirted



Cat. No.	Description		From*	Unit
781411	24-well PCR-plate, PP, 0,2 ml, thin-walled, flexible, no skirt	1 pack = 40 piece(s)	1	pack
781412	24-well PCR-plates, no skirt, PP, white, 0,2 ml	1 pack = 40 piece(s)	1	pack

PCR plate, 48-well, non-skirted



Cat. No.	Description		From*	Unit
781415	48-well PCR-plate, PP, 0,2 ml, thin-walled, flexible, no skirt	1 pack = 20 piece(s)	1	pack
781416	48-well PCR-plates, no skirt, PP, white, 0,2 ml	1 pack = 20 piece(s)	1	pack


PCR plate, 96-well, non-skirted

Cat. No.	Description		From*	Unit
781350	96-well PCR plate, non-skirted, standard profile, elev. rim, white, blue coding, cut corner H12, PP	1 pack = 50 piece(s)	1	pack
781354	96-well PCR plate, non-skirted, standard profile, elev. rim, white, blue coding, cut corner H12, PP	1 pack = 50 piece(s)	1	pack
781366	96-well PCR plate, non-skirted, low profile, clear, blue coding, cut corner H12, PP	1 pack = 50 piece(s)	1	pack
781367	96-well PCR plate, non-skirted, low profile, white, blue coding, cut corner H12, PP	1 pack = 50 piece(s)	1	pack
781368	96-well PCR plate, non-skirted, standard profile, clear, blue coding, cut corner A12, PP	1 pack = 50 piece(s)	1	pack
781369	96-well PCR plate, non-skirted, standard profile, white, blue coding, cut corner A12, PP	1 pack = 50 piece(s)	1	pack




PCR plate, 96-well, semi-skirted

Cat. No.	Description		From*	Unit
781357	96-well PCR plate, semi-skirted, standard profile, elev. rim, white, blue coding, cut corner A1, PP	1 pack = 50 piece(s)	1	pack
781371	96-well PCR plate, semi-skirted, low profile, clear, blue coding, cut corner A12, PP	1 pack = 50 piece(s)	1	pack
781372	96-well PCR plate, semi-skirted, low profile, white, blue coding, cut corner A12, PP	1 pack = 50 piece(s)	1	pack
781373	96-well PCR plate, semi-skirted, raised skirt, low profile, blue coding, clear, cut corner A1, PP	1 pack = 50 piece(s)	1	pack
781374	96-well PCR plate, semi-skirted, raised skirt, low profile, white, blue coding, cut corner A1, PP	1 pack = 50 piece(s)	1	pack
781375	96-well PCR plate, semi-skirted, standard profile, clear, blue coding, cut corner A12, PP	1 pack = 50 piece(s)	1	pack
781376	96-well PCR plate, semi-skirted, standard profile, white, blue coding, cut corner A12, PP	1 pack = 50 piece(s)	1	pack
781400	96-well PCR plate, semi-skirted, standard profile, elevated rim, clear, cut corner A1, PP	1 pack = 50 piece(s)	1	pack




PCR plate, 96-well, semi-skirted, for Roche® LightCycler® 480

Cat. No.	Description		From*	Unit
781364	96-well PCR plate, semi-skirted, white, 0,15 ml, low profile, embossed coding, cut corner H12, PP	1 pack = 50 piece(s)	1	pack
781365	96-well PCR plate, semi-skirted, white, 0,15 ml+ 50 films, black coding, cut corner H12, PP	1 pack = 50 piece(s)	1	pack



PCR plate, 96-well, skirted

Cat. No.	Description		From*	Unit
781377	96-well PCR plate, skirted, low profile, clear, black coding, cut corner A12, PP	1 pack = 50 piece(s)	1	pack
781378	96-well PCR plate, skirted, low profile, white, black coding, cut corner A12, PP	1 pack = 50 piece(s)	1	pack



PCR plate, 384-well, skirted



Cat. No.	Description		From*	Unit
781345	384-well PCR plate, 0,03ml, flexible, full skirt	1 pack = 50 piece(s)	1	pack
781347	384-well PCR plate, 0,03ml, flexible, full skirt, for RealTimePCR	1 pack = 50 piece(s)	1	pack

PCR plate, 384-well, skirted, for Roche® LightCycler® 480



Cat. No.	Description		From*	Unit
781358	384-well PCR-plates, PP, white, 0,03 ml	1 pack = 50 piece(s)	1	pack

PCR box/-rack




Cat. No.	Description		From*	Unit
781362	PCR-Box/-Rack, PP	1 pack = 5 piece(s)	1	pack

PCR Mini Cooler



Cat. No.	Description		From*	Unit
781260	PCR Mini Cooler, PP, with transparent lid		2	piece(s)

Life Science sealing films for PCR and qPCR



Cat. No.	Description		From*	Unit
781383	Filmstrips polyester, self adhesive, for PCR, qPCR and storage	1 pack = 50 piece(s) (each 8 strips)	1	pack
781390	Sealing film for microplates, self-adhesive, for PCR, PP	1 pack = 100 piece(s)	1	pack
781391	Sealing film for microplates, Polyester, pressure sensitive, for ELISA and Real-Time PCR	1 pack = 100 piece(s)	1	pack

Life Science sealing films for PCR and storage



Cat. No.	Description		From*	Unit
701367	Sealing film for storage, DMSO-resistant, allows visual inspection	1 pack = 100 piece(s)	1	pack

Life Science sealing films for storage and automation

Cat. No.	Description		From*	Unit
781381	Sealing film for microplates, self-adhesive, for cold storage, aluminium	1 pack = 100 piece(s)	1	pack
781382	Foilstrips Aluminum, self adhesive, for PCR and storage	1 pack = 50 piece(s) (each 6 strips)	1	pack



Life Science sealing films for automation processes

Cat. No.	Description		From*	Unit
701370	Sealing film for automation, PE/PP	1 pack = 50 piece(s)	1	pack
701374	Sealing film for plates, Vinyl, with pre-pierced zone, 100 pcs., non-sterile, 2 end tabs, alphanumeric code, -40 up to + 90°C	1 pack = 100 piece(s)	1	pack




Life Science sealing films for cell and tissue cultures

Cat. No.	Description		From*	Unit
701364	Sealing film for microplates, self-adhesive, for cell and tissue culture, gas-permeable, Rayon	1 pack = 100 piece(s)	1	pack
701365	Sealing film for microplates, self-adhesive, for cell and tissue culture, gas-permeable, Rayon, sterile	1 pack = 50 piece(s)	1	pack




Life Science sealing films for luminescence and fluorescence measuring

Cat. No.	Description		From*	Unit
701371	Sealing film for fluorescence measurement, black	1 pack = 50 piece(s)	1	pack
701372	Sealing film for luminescence measurement, white	1 pack = 50 piece(s)	1	pack



PCR plate and sealing film sets

Cat. No.	Description	From*	Unit
781230	PCR Set 1: 50x PCR plates, 96-well, non-skirted, low profile, transparent + 50x sealing films (qPCR), Polyester	1	piece(s)
781231	PCR Set 2: 50x PCR plates, 96-well, non-skirted, standard profile, transparent + 50x sealing films (qPCR), Polyester	1	piece(s)
781232	PCR Set 3: 50x PCR plates, 96-well, semi-skirted, low profile, transparent + 50x sealing films (qPCR), Polyester	1	piece(s)
781233	PCR Set 4: 50x PCR plates, 96-well, semi-skirted, raised skirt, low profile, transparent + 50x sealing films (qPCR), Polyester	1	piece(s)
781234	PCR Set 5: 50x PCR plates, 96-well, semi-skirted, standard profile, transparent + 50x sealing films (qPCR), Polyester	1	piece(s)



Sealing paddle for sealing films



Cat. No.	Description	From*	Unit
701381	Sealing paddle for plates, white, black BRAND-logo, PE	1	

PCR cap tool




Cat. No.	Description	From*	Unit
781419	Cap Tool for attaching and detaching caps on PCR plates, PA	1	piece(s)

Cap strip of 8 lids



Cat. No.	Description		From*	Unit
781413	PCR cap strips of 8, clear, for qPCR flat, for 781377-78, PP	1 pack = 300 piece(s)	1	pack
781414	PCR cap strips of 8, clear, domed for 781377-78	1 pack = 300 piece(s)	1	pack

Microplate, 96-well, U-bottom, PP



Cat. No.	Description		From*	Unit
701330	96-well Microplate, PP, U-bottom, 300 µl, non-sterile		100	piece(s)

Deep-well plate, 96-well, U-bottom, PP



Cat. No.	Description		From*	Unit
701340	Deep well plate, 96-well, stackable 1,2 ml PP non-sterile, Low Profile, alphan. identific. system	1 pack = 50 piece(s)	1	pack
701342	Deep well plate, 96-well, stackable 1,2 ml PP non-sterile, elevated wells, alphan. identific. system, DNA-, DNase- und RNase-free	1 pack = 32 piece(s)	1	pack
701346	Deep well plate, 96-well, stackable 0,5 ml, PP, non-sterile, alphan. identific. system, DNA-, DNase- und RNase-free	1 pack = 48 piece(s)	1	pack
701350	Deep well plate, 96-well, stackable 1,2 ml, PP, non-sterile, alphan. identific. system, DNA-, DNase- und RNase-free	1 pack = 24 piece(s)	1	pack
701354	Deep well plate, 96-well, stackable 2,2 ml, PP, non-sterile, alphan. identific. system, DNA-, DNase- und RNase-free	1 pack = 24 piece(s)	1	pack


Deep-well plate, 96-well, U-bottom, PS

Cat. No.	Description		From*	Unit
701352	Deep well plate, 96-well, stackable 1,1 ml, PS, non-sterile, alphan. identific. system	1 pack = 32 piece(s)	1	pack



Deep-well plate, 384-well, U-bottom, PP

Cat. No.	Description		From*	Unit
701355	Deep well plate, 384-well, stackable 300 µl, PP, non-sterile, alphan. identific. system	1 pack = 48 piece(s)	1	pack



96 Tube racks, PP

Cat. No.	Description		From*	Unit
781500	Tube, 1,2 ml, rack packed, PP, rack with 96 individual tubes non-sterile	1 pack = 10 piece(s)	1	pack
781510	Tube, 1,2 ml, rack packed, PP, rack with 12 strips of 8 tubes non-sterile	1 pack = 10 piece(s)	1	pack



96 Tube racks, PP, with lid

Cat. No.	Description		From*	Unit
781565	Rack with lid with 96 coded 0.65 ml tubes, for 96 tube rack. (Not available in the USA)	1 pack = 50 piece(s)	1	pack
781566	Rack with lid with 96 coded 1.2 ml tubes, for 96 tube rack. (Not available in the USA)	1 pack = 50 piece(s)	1	pack



Single tubes without coding

Cat. No.	Description		From*	Unit
781520	Tube, 1,2 ml, rack packed, PP, individual tubes, non-sterile, bulk	1 pack = 960 piece(s)	1	pack



Individual cap for tubes 1,2 ml

Cat. No.	Description		From*	Unit
781530	Caps, PE, for tubes 1,2 ml, non-sterile, bulk	1 pack = 960 piece(s)	1	pack



Tubes, strips of 8



Cat. No.	Description		From*	Unit
781525	Tube, 1,2 ml, rack packed, PP, 8 strip tubes, non-sterile, bulk	1 pack = 120 piece(s)	1	pack

Tube rack, empty



Cat. No.	Description		From*	Unit
781540	Tube 1,2ml PP 10 racks with grid, non-sterile (empty)	1 pack = 10 piece(s)	1	pack

Cap strips, strips of 8 for 96 tubes



Cat. No.	Description	From*	Unit	
781535	8 strip caps, PE, for tube strips 1,2 ml, PP, bulk, non-sterile	1 pack = 120 piece(s)	1	pack

Cryogenic tube, with external thread, screw-on cap with silicone seal



Cat. No.	Description	From*	Unit
114830	Cryogenic tube with screw cap, PP, 1,2 ml, with external thread, self-standing, y-sterile	1000	piece(s)
114831	Cryogenic tube with screw cap, PP, 2 ml, with external thread, round-bottom, y-sterile	1000	piece(s)
114832	Cryogenic tube with screw cap, PP, 2 ml, with external thread self-standing, y-sterile	1000	piece(s)
114833	Cryogenic tube with screw cap, PP, 3 ml, with external thread self-standing, y-sterile	1000	piece(s)
114834	Cryogenic tube with screw cap, PP, 4 ml, with external thread self-standing, y-sterile	1000	piece(s)
114835	Cryogenic tube with screw cap, PP, 5 ml, with external thread self-standing, y-sterile	1000	piece(s)

Cryogenic tube, with internal thread, screw-on cap with sealing cone



Cat. No.	Description	From*	Unit
114840	Cryogenic tube with screw cap, PP, 1,2 ml, with internal thread self-standing, y-sterile	1000	piece(s)
114841	Cryogenic tube with screw cap, PP, 2 ml, with internal thread self-standing, y-sterile	1000	piece(s)
114842	Cryogenic tube with screw cap, PP, 2 ml, with internal thread, round-bottom, y-sterile	1000	piece(s)
114843	Cryogenic tube with screw cap, PP, 4 ml, with internal thread, round-bottom, y-sterile	1000	piece(s)
114844	Cryogenic tube with screw cap, PP, 4 ml, with internal thread self-standing, y-sterile	1000	piece(s)
114845	Cryogenic tube with screw cap, PP, 5 ml, with internal thread, round-bottom, y-sterile	1000	piece(s)

Cap inserts for cryogenic tubes

Cat. No.	Description	From*	Unit
114850	Cap insert for cryogenic tubes, PP, white, for color coding	500	piece(s)
114851	Cap insert for cryogenic tubes, PP, blue, for color coding	500	piece(s)
114852	Cap insert for cryogenic tubes, PP, red, for color coding	500	piece(s)
114853	Cap insert for cryogenic tubes, PP, green, for color coding	500	piece(s)
114854	Cap insert for cryogenic tubes, PP, yellow, for color coding	500	piece(s)



Storage box for cryogenic tubes

Cat. No.	Description	From*	Unit
114862	Storage box, PC, for 1.2 ml and 2 ml cryogenic tubes, for 81 tubes	4	piece(s)
114864	Storage box, PC, for 3 ml, 4 ml and 5 ml cryogenic tubes (with external thread), for 81 tubes	5	piece(s)
114866	Storage box, PC, for 1.2 ml and 2 ml cryogenic tubes (with internal thread), for 100 tubes	4	piece(s)



Sealing mats for deep-well plates

Cat. No.	Description	From*	Unit
701357	Silicone cover for deep well plate 384-well 1 pack = 50 piece(s)	1	pack
701358	Mat cover for deep well plate 96-well 0,5 ml, non-sterile 1 pack = 50 piece(s)	1	pack
701360	Mat cover for deep-well plates 1.1 ml 96-well and 1.2 ml 96-well elevated rim, PE 1 pack = 24 piece(s)	1	pack
701362	Mat cover for deep well plate 96-well 2,2 ml, non-sterile 1 pack = 24 piece(s)	1	pack
701368	Mat cover for deep well plate 96-well 1,2 ml, non-sterile 1 pack = 50 piece(s)	1	pack



Cryogenic tube rack

Cat. No.	Description	From*	Unit
114860	Rack for cryogenic tubes, PP, stackable blue, for 50 cryogenic tubes	4	piece(s)



Cooling box

Cat. No.	Description	From*	Unit
156100	Ice bucket with lid, PU-foam, stackable, 4,5 l	1	piece(s)





Cell culture inserts, BRAND Insert 2in1, for 6-well

Cat. No.	Description		From*	Unit
782740	6-well Insert 2in1, PC-Membrane 0,4 µm	1 pack = 24 piece(s) (24 x 1)	1	pack
782741	6-well Insert 2in1, PC-membrane 0.4 µm multipack	1 pack = 24 piece(s) (4 x 6)	1	pack
782746	6-well Insert 2in1, PC-membrane 8.0 µm	1 pack = 24 piece(s) (24 x 1)	1	pack
782747	6-well Insert 2in1, PC-membrane 8.0 µm multipack	1 pack = 24 piece(s) (4 x 6)	1	pack
782750	6-well Insert 2in1, PET-membrane 0.4 µm	1 pack = 24 piece(s) (24 x 1)	1	pack
782751	6-well Insert 2in1, PET-membrane 0.4 µm multipack	1 pack = 24 piece(s) (4 x 6)	1	pack
782756	6-well Insert 2in1, PET-membrane 8.0 µm	1 pack = 24 piece(s) (24 x 1)	1	pack
782757	6-well Insert 2in1, PET-membrane 8.0 µm multipack	1 pack = 24 piece(s) (4 x 6)	1	pack

Cell culture inserts, BRAND Insert 2in1, for 12-well



Cat. No.	Description		From*	Unit
782720	12-well Insert 2in1, PC-membrane 0.4 µm	1 pack = 48 piece(s) (48 x 1)	1	pack
782721	12-well Insert 2in1, PC-membrane 0.4 µm, multipack	1 pack = 36 piece(s) (4 x 9)	1	pack
782726	12-well Insert 2in1, PC-membrane 8.0 µm	1 pack = 48 piece(s) (48 x 1)	1	pack
782727	12-well Insert 2in1, PC-membrane 8.0 µm, multipack	1 pack = 36 piece(s) (4 x 9)	1	pack
782730	12-well Insert 2in1, PET-membrane 0.4 µm	1 pack = 48 piece(s) (48 x 1)	1	pack
782731	12-well Insert 2in1, PET-membrane 0.4 µm, multipack	1 pack = 36 piece(s) (4 x 9)	1	pack
782736	12-well Insert 2in1, PET-membrane 8.0 µm	1 pack = 48 piece(s) (48 x 1)	1	pack
782737	12-well Insert 2in1, PET-membrane 8.0 µm, multipack	1 pack = 36 piece(s) (4 x 9)	1	pack

Cell culture inserts, BRAND Insert 2in1, for 24-well

Cat. No.	Description		From*	Unit
782700	24-well Insert 2in1, PC-membrane 0.4 µm	1 pack = 48 piece(s) (48 x 1)	1	pack
782701	24-well Insert 2in1, PC-membrane 0.4 µm, multipack	1 pack = 48 piece(s) (4 x 12)	1	pack
782706	24-well Insert 2in1, PC-membrane 8.0 µm	1 pack = 48 piece(s) (48 x 1)	1	pack
782707	24-well Insert 2in1, PC-membrane 8.0 µm, multipack	1 pack = 48 piece(s) (4 x 12)	1	pack
782710	24-well Insert 2in1, PET-membrane 0.4 µm	1 pack = 48 piece(s) (48 x 1)	1	pack
782711	24-well Insert 2in1, PET-membrane 0.4 µm, multipack	1 pack = 48 piece(s) (4 x 12)	1	pack
782716	24-well Insert 2in1, PET-membrane 8.0 µm	1 pack = 48 piece(s) (48 x 1)	1	pack
782717	24-well Insert 2in1, PET-membrane 8.0 µm, multipack	1 pack = 48 piece(s) (4 x 12)	1	pack



Insert Strips, smooth-walled, PC membrane

Cat. No.	Description		From*	Unit
782800	Insert strips 13mm for 24-/6-well BRANDplates®, smooth-walled, 4 Inserts/Strip, PC-membr. 0,4 µm	1 pack = 12 piece(s)	1	pack
782820	Insert strips 13mm for 24-/6-well BRANDplates®, smooth-walled, 4 Inserts/Strip, PC-membr. 1,0 µm	1 pack = 12 piece(s)	1	pack
782840	Insert strips 13mm for 24-/6-well BRANDplates®, smooth-walled, 4 Inserts/Strip, PC-membr. 3,0 µm	1 pack = 12 piece(s)	1	pack
782860	Insert strips 13mm for 24-/6-well BRANDplates®, smooth-walled, 4 Inserts/Strip, PC-membr. 8,0 µm	1 pack = 12 piece(s)	1	pack
782900	Insert Strips 9mm f.24/6-w.BRANDplates, PC-membr. 0,4µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack
782920	Insert Strips 9mm f.24/6-w.BRANDplates, PC-membr. 1,0µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack
782940	Insert Strips 9mm f.24/6-w.BRANDplates, PC-membr. 3,0µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack
782960	Insert Strips 9mm f.24/6-w.BRANDplates, PC-membr. 8,0µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack



Insert Strips, smooth-walled, PET membrane

Cat. No.	Description		From*	Unit
782810	Insert strips 13mm for 24-/6-well BRANDplates®, smooth-walled, 4 Inserts/Strip, PET-membr. 0,4 µm	1 pack = 12 piece(s)	1	pack
782870	Insert strips 13mm for 24-/6-well BRANDplates®, smooth-walled, 4 Inserts/Strip, PET-membr. 8,0 µm	1 pack = 12 piece(s)	1	pack
782910	Insert Strips 9mm f.24/6-w.BRANDplates, PET-membr. 0,4µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack
782930	Insert Strips 9mm f.24/6-w.BRANDplates, PET-membr. 1,0µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack
782950	Insert Strips 9mm f.24/6-w.BRANDplates, PET-membr. 3,0µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack
782970	Insert Strips 9mm f.24/6-w.BRANDplates, PET-membr. 8,0µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack



Insert Strips, with inlet channels (Inlet Opening System), PC membrane



Cat. No.	Description		From*	Unit
782801	Insert strips 13mm for 24-/6-well BRANDplates®, IOS, 4 Inserts/Strip, PC-membr. 0,4 µm	1 pack = 12 piece(s)	1	pack
782821	Insert strips 13mm for 24-/6-well BRANDplates®, IOS, 4 Inserts/Strip, PC-membr. 1,0 µm	1 pack = 12 piece(s)	1	pack
782841	Insert strips 13mm for 24-/6-well BRANDplates®, IOS, 4 Inserts/Strip, PC-membr. 3,0 µm	1 pack = 12 piece(s)	1	pack
782901	Insert Strips 9mm f.24/6-w.BRANDplates, IOS,PC-membr. 0,4µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack

Insert Strips, with inlet channels (Inlet Opening System), PET membrane



Cat. No.	Description		From*	Unit
782811	Insert strips 13mm for 24-/6-well BRANDplates®, IOS, 4 Inserts/Strip, PET-membr. 0,4 µm	1 pack = 12 piece(s)	1	pack
782861	Insert strips 13mm for 24-/6-well BRANDplates®, IOS, 4 Inserts/Strip, PC-membr. 8,0 µm	1 pack = 12 piece(s)	1	pack
782871	Insert strips 13mm for 24-/6-well BRANDplates®, IOS, 4 Inserts/Strip, PET-membr. 8,0 µm	1 pack = 12 piece(s)	1	pack
782911	Insert Strips 9mm f.24/6-w.BRANDplates, IOS,PET-membr. 0,4µm,4 Inserts/Strip	1 pack = 12 piece(s)	1	pack

Single Inserts, smooth-walled, PC membrane



Cat. No.	Description		From*	Unit
782806	Insert-Single 13mm for 24-/6-well BRANDplates®, smooth-walled, PC-membr., 0,4 µm	1 pack = 48 piece(s)	1	pack
782826	Insert-Single 13mm for 24-/6-well BRANDplates®, smooth-walled, PC-membr., 1,0 µm	1 pack = 48 piece(s)	1	pack
782846	Insert-Single 13mm for 24-/6-well BRANDplates®, smooth-walled, PC-membr., 3,0 µm	1 pack = 48 piece(s)	1	pack
782866	Insert-Single 13mm for 24-/6-well BRANDplates®, smooth-walled, PC-membr., 8,0 µm	1 pack = 48 piece(s)	1	pack

Single Inserts, smooth-walled, PET membrane



Cat. No.	Description		From*	Unit
782816	Insert-Single 13mm for 24-/6-well BRANDplates®, smooth-walled, PET-membr., 0,4 µm	1 pack = 48 piece(s)	1	pack
782876	Insert-Single 13mm for 24-/6-well BRANDplates®, smooth-walled, PET-membr., 8,0 µm	1 pack = 48 piece(s)	1	pack

Single Inserts, with inlet channels (Inlet Opening System), PET membrane



Cat. No.	Description		From*	Unit
782817	Insert-Single 13mm for 24-/6-well BRANDplates®, IOS, PET-membr., 0,4 µm	1 pack = 48 piece(s)	1	pack
782877	Insert-Single 13mm for 24-/6-well BRANDplates®, IOS, PET-membr., 8,0 µm	1 pack = 48 piece(s)	1	pack

Single Inserts, with inlet channels (Inlet Opening System), PC membrane

Cat. No.	Description		From*	Unit
782807	Insert-Single 13mm for 24-/6-well BRANDplates®, IOS, PC-membr., 0,4 µm	1 pack = 48 piece(s)	1	pack
782827	Insert-Single 13mm for 24-/6-well BRANDplates®, IOS, PC-membr., 1,0 µm	1 pack = 48 piece(s)	1	pack
782847	Insert-Single 13mm for 24-/6-well BRANDplates®, IOS, PC-membr., 3,0 µm	1 pack = 48 piece(s)	1	pack
782867	Insert-Single 13mm for 24-/6-well BRANDplates®, IOS, PC-membr., 8,0 µm	1 pack = 48 piece(s)	1	pack



Multiwell plates BRANDplates®, for inserts


Cat. No.	Description		From*	Unit
782880	BRANDplates for Inserts 24-well with lid pureGrade™ S, individually wrap.	1 pack = 10 piece(s)	1	pack
782881	BRANDplates for Inserts, 6-well with lid pureGrade™ S, individually wrap.	1 pack = 10 piece(s)	1	pack
782890	BRANDplates for Inserts 24-well with lid cellGrade™ plus, individually wrap.	1 pack = 10 piece(s)	1	pack
782891	BRANDplates for Inserts 6-well with lid cellGrade™ plus, individually wrap.	1 pack = 10 piece(s)	1	pack



Insert System, 6-well plates, filled with 6 insert strips

Cat. No.	Description		From*	Unit
782802	Insert System: 6-well plates BRANDplates® incl. 6 insert strips PC-membr., cellGrade™ plus, 0,4 µm	1 pack = 5 piece(s)	1	pack
782803	Insert System: 6-well plates BRANDplates® incl. 6 insert strips IOS, PC-membr., cellGrade™ plus, 0,4 µm	1 pack = 5 piece(s)	1	pack
782812	Insert System: 6-well plates BRANDplates® incl. 6 insert strips PET-membr., cellGrade™ plus, 0,4 µm	1 pack = 5 piece(s)	1	pack
782813	Insert System: 6-well plates BRANDplates® incl. 6 insert strips IOS, PET-membr., cellGrade™ plus, 0,4 µm	1 pack = 5 piece(s)	1	pack
782822	Insert System: 6-well plates BRANDplates® incl. 6 insert strips PC-membr., cellGrade™ plus, 1,0 µm	1 pack = 5 piece(s)	1	pack
782823	Insert System: 6-well plates BRANDplates® incl. 6 insert strips IOS, PC-membr., cellGrade™ plus, 1,0 µm	1 pack = 5 piece(s)	1	pack
782842	Insert System: 6-well plates BRANDplates® incl. 6 insert strips PC-membr., cellGrade™ plus, 3,0 µm	1 pack = 5 piece(s)	1	pack
782843	Insert System: 6-well plates BRANDplates® incl. 6 insert strips IOS, PC-membr., cellGrade™ plus, 3,0 µm	1 pack = 5 piece(s)	1	pack
782862	Insert System: 6-well plates BRANDplates® incl. 6 insert strips PC-membr., cellGrade™ plus, 8,0 µm	1 pack = 5 piece(s)	1	pack
782863	Insert System: 6-well plates BRANDplates® incl. 6 insert strips IOS, PC-membr., cellGrade™ plus, 8,0 µm	1 pack = 5 piece(s)	1	pack
782872	Insert System: 6-well plates BRANDplates® incl. 6 insert strips PET-membr., cellGrade™ plus, 8,0 µm	1 pack = 5 piece(s)	1	pack
782873	Insert System: 6-well plates BRANDplates® incl. 6 insert strips IOS, PET-membr., cellGrade™ plus, 8,0 µm	1 pack = 5 piece(s)	1	pack




BRANDplates[®] microplates, 96-well, pureGrade™

Cat. No.	Description		From*	Unit
781600	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS Standard transp., 'U'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781601	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS Standard transp., 'V'bottom 360 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781602	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS Standard transp., 'F'bottom 350 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781603	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS Standard transp., 'C'bottom 350 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781604	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS Standard, white, 'U'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781605	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS Standard, white, 'F'bottom 350 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781607	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS Standard black, 'U'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781608	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS Standard black, 'F'bottom 350 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781610	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS transp., b., white, 'F'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781611	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS transp., b., black, 'F'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781614	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS foil bottom UV transp., 350 µl, non-ster., no USA	1 pack = 50 piece(s)	1	pack
781615	Microplate, BRANDplates [®] , pureGrade™, 96-well, PS foil bottom UV transp., black, non-ster., no USA	1 pack = 50 piece(s)	1	pack

BRANDplates[®] microplates, 384-well, pureGrade™

Cat. No.	Description		From*	Unit
781620	Microplate, BRANDplates [®] , pureGrade™, 384-well, PS Standard transp., 'F'bottom 100 µl, non-sterile	1 pack = 50 piece(s)	1	pack
781621	Microplate, BRANDplates [®] , pureGrade™, 384-well, PS Standard, white, 'F'bottom 100 µl, non-sterile	1 pack = 50 piece(s)	1	pack
781622	Microplate, BRANDplates [®] , pureGrade™, 384-well, PS Standard, black, 'F'bottom 100 µl, non-sterile	1 pack = 50 piece(s)	1	pack
781626	Microplate, BRANDplates [®] , pureGrade™, 384-well, PS trans., b., white, 'F'bottom 120 µl, non-sterile	1 pack = 50 piece(s)	1	pack
781627	Microplate, BRANDplates [®] , pureGrade™, 384-well, PS trans., b., black, 'F'bottom 120 µl, non-sterile	1 pack = 50 piece(s)	1	pack

BRANDplates[®] microplates, 1536-well, pureGrade™

Cat. No.	Description		From*	Unit
781640	Microplate, BRANDplates [®] , pureGrade™, 1536-well, PS Standard trans., 'F'bottom 10 µl, non-sterile	1 pack = 50 piece(s)	1	pack
781641	Microplate, BRANDplates [®] , pureGrade™, 1536-well, PS Standard white, 'F'bottom 10 µl, non-sterile	1 pack = 50 piece(s)	1	pack
781642	Microplate, BRANDplates [®] , pureGrade™, 1536-well, PS Standard black, 'F'bottom 10 µl, non-sterile	1 pack = 50 piece(s)	1	pack

BRANDplates[®] microplates, 96-well, pureGrade[™] S

Cat. No.	Description		From*	Unit
781660	Microplate, BRANDplates [®] , pureGrade [™] S, 96-well, PS Standard transp., 'U'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
781661	Microplate, BRANDplates [®] , pureGrade [™] S, 96-well, PS Standard transp., 'V'bottom 360 µl, sterile	1 pack = 50 piece(s)	1	pack
781662	Microplate, BRANDplates [®] , pureGrade [™] S, 96-well, PS Standard transp., 'F'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
781663	Microplate, BRANDplates [®] , pureGrade [™] S, 96-well, PS Standard transp., 'C'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
781665	Microplate, BRANDplates [®] , pureGrade [™] S, 96-well, PS Standard white, 'F'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
781668	Microplate, BRANDplates [®] , pureGrade [™] S, 96-well, PS Standard black, 'F'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
781670	Microplate, BRANDplates [®] , pureGrade [™] S, 96-well, PS transp., b., white, 'F'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
781671	Microplate, BRANDplates [®] , pureGrade [™] S, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack



BRANDplates[®] microplates, 384-well, pureGrade[™] S

Cat. No.	Description		From*	Unit
781680	Microplate, BRANDplates [®] , pureGrade [™] S, 384-well, PS Standard, trans., 'F'bottom 100 µl, sterile	1 pack = 50 piece(s)	1	pack
781681	Microplate, BRANDplates [®] , pureGrade [™] S, 384-well, PS Standard, white, 'F'bottom 100 µl, sterile	1 pack = 50 piece(s)	1	pack
781682	Microplate, BRANDplates [®] , pureGrade [™] S, 384-well, PS Standard, black, 'F'bottom 100 µl, sterile	1 pack = 50 piece(s)	1	pack
781686	Microplate, BRANDplates [®] , pureGrade [™] S, 384-well, PS trans., b., white, 'F'bottom 120 µl, sterile	1 pack = 50 piece(s)	1	pack
781687	Microplate, BRANDplates [®] , pureGrade [™] S, 384-well, PS trans., b., black, 'F'bottom 120 µl, sterile	1 pack = 50 piece(s)	1	pack


BRANDplates[®] microplates, 1536-well, pureGrade[™] S

Cat. No.	Description		From*	Unit
781700	Microplate, BRANDplates [®] , pureGrade [™] S, 1536-well, PS Standard transp., 'F'bottom 10 µl, sterile	1 pack = 50 piece(s)	1	pack
781701	Microplate, BRANDplates [®] , pureGrade [™] S, 1536-well, PS Standard, white, 'F'bottom 10 µl, sterile	1 pack = 50 piece(s)	1	pack
781702	Microplate, BRANDplates [®] , pureGrade [™] S, 1536-well, PS Standard, black, 'F'bottom 10 µl, sterile	1 pack = 50 piece(s)	1	pack




BRANDplates® microplates, 96-well, immunoGrade™


Cat. No.	Description		From*	Unit
781720	Microplate, BRANDplates®, immunoGrade™, 96-well, PS Standard trans., 'U'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781721	Microplate, BRANDplates®, immunoGrade™, 96-well, PS Standard trans., 'V'bottom 360 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781722	Microplate, BRANDplates®, immunoGrade™, 96-well, PS Standard trans., 'F'bottom 350 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781724	Microplate, BRANDplates®, immunoGrade™, 96-well, PS Standard white, 'U'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781727	Microplate, BRANDplates®, immunoGrade™, 96-well, PS Standard black, 'U'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781731	Microplate, BRANDplates®, immunoGrade™, 96-well, PS trans., b., black, 'F'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781732	Microplate, BRANDplates®, immunoGrade™, 96-well, PS trans., b., black, 'F'bottom 330 µl, non-sterile	1 pack = 5 piece(s)	1	pack

BRANDplates® microplates, 384-well, immunoGrade™

Cat. No.	Description		From*	Unit
781740	Microplate, BRANDplates®, immunoGrade™, 384-well, PS Standard, trans., 'F'bottom 100 µl, non-sterile	1 pack = 50 piece(s)	1	pack
781741	Microplate, BRANDplates®, immunoGrade™, 384-well, PS Standard, white, 'F'bottom 100 µl, non-sterile	1 pack = 50 piece(s)	1	pack
781742	Microplate, BRANDplates®, immunoGrade™, 384-well, PS Standard, black, 'F'bottom 100 µl, non-sterile	1 pack = 50 piece(s)	1	pack

BRANDplates® microplates, 96-well, lipoGrade™


Cat. No.	Description		From*	Unit
781840	Microplate, BRANDplates®, lipoGrade™, 96-well, PS Standard trans., 'U'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781842	Microplate, BRANDplates®, lipoGrade™, 96-well, PS Standard trans., 'F'bottom 350 µl, non-sterile	1 pack = 100 piece(s)	1	pack

BRANDplates® microplates, 96-well, hydroGrade™

Cat. No.	Description		From*	Unit
781780	Microplate, BRANDplates®, hydroGrade™, 96-well, PS Standard trans., 'U'bottom 330 µl, non-sterile	1 pack = 100 piece(s)	1	pack
781782	Microplate, BRANDplates®, hydroGrade™, 96-well, PS Standard trans., 'F'bottom 350 µl, non-sterile	1 pack = 100 piece(s)	1	pack

BRANDplates[®] microplates, 96-well, cellGrade™

Cat. No.	Description		From*	Unit
781960	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS Standard transp., 'U'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
781961	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS Standard transp., 'V'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
781962	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS Standard transp., 'F'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
781963	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS Standard transp., 'C'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
781965	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS Standard, white, 'F'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
781968	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS Standard, black, 'F'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
781970	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS transp., b., white, 'U'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
781971	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
781974	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS transp., b., white, 'U'bottom 330 µl, sterile	1 pack = 5 piece(s)	1	pack
781975	Microplate, BRANDplates [®] , cellGrade™, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 5 piece(s)	1	pack


BRANDplates[®] microplates, 384-well, cellGrade™

Cat. No.	Description		From*	Unit
781980	Microplate, BRANDplates [®] , cellGrade™, 384-well, PS Standard, trans., 'F'bottom 100 µl, sterile	1 pack = 50 piece(s)	1	pack
781981	Microplate, BRANDplates [®] , cellGrade™, 384-well, PS Standard, white, 'F'bottom 100 µl, sterile	1 pack = 50 piece(s)	1	pack
781982	Microplate, BRANDplates [®] , cellGrade™, 384-well, PS Standard, black, 'F'bottom 100 µl, sterile	1 pack = 50 piece(s)	1	pack
781986	Microplate, BRANDplates [®] , cellGrade™, 384-well, PS trans., b., white, 'F'bottom 120 µl, sterile	1 pack = 50 piece(s)	1	pack
781987	Microplate, BRANDplates [®] , cellGrade™, 384-well, PS trans., b., black, 'F'bottom 120 µl, sterile	1 pack = 50 piece(s)	1	pack
781988	Microplate, BRANDplates [®] , cellGrade™, 384-well, PS trans., b., black, 'F'bottom 120 µl, sterile	1 pack = 5 piece(s)	1	pack
781989	Microplate, BRANDplates [®] , cellGrade™, 384-well, PS trans., b., white, 'F'bottom 120 µl, sterile	1 pack = 5 piece(s)	1	pack



BRANDplates[®] microplates, 1536-well, cellGrade™

Cat. No.	Description		From*	Unit
782000	Microplate, BRANDplates [®] , cellGrade™, 1536-well, PS Standard transp., 'F'bottom 10 µl, sterile	1 pack = 50 piece(s)	1	pack
782001	Microplate, BRANDplates [®] , cellGrade™, 1536-well, PS Standard, white, 'F'bottom 10 µl, sterile	1 pack = 50 piece(s)	1	pack
782002	Microplate, BRANDplates [®] , cellGrade™, 1536-well, PS Standard, black, 'F'bottom 10 µl, sterile	1 pack = 50 piece(s)	1	pack




BRANDplates® microplates, 96-well, cellGrade™ plus

Cat. No.	Description		From*	Unit
782022	Microplate, BRANDplates®, cellGrade™+, 96-well, PS Standard transp., 'F'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
782030	Microplate, BRANDplates®, cellGrade™+, 96-well, PS transp., b., white, 'F'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
782031	Microplate, BRANDplates®, cellGrade™+, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
782034	Microplate, BRANDplates®, cellGrade™, 96-well, PS transp., b., white, 'F'bottom 330 µl, sterile	1 pack = 5 piece(s)	1	pack
782035	Microplate, BRANDplates®, cellGrade™, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 5 piece(s)	1	pack

BRANDplates® microplates, 96-well, cellGrade™ premium


Cat. No.	Description		From*	Unit
782082	Microplate, BRANDplates®, cellGrade™ prem, 96-well, PS Standard transp., 'F'bottom 350 µl, sterile	1 pack = 50 piece(s)	1	pack
782090	Microplate, BRANDplates®, cellGrade™ prem, 96-well, PS transp., b., white, 'F'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
782091	Microplate, BRANDplates®, cellGrade™ prem, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 50 piece(s)	1	pack
782094	Microplate, BRANDplates®, cellGrade™ prem, 96-well, PS transp., b., white, 'F'bottom 330 µl, sterile	1 pack = 5 piece(s)	1	pack
782095	Microplate, BRANDplates®, cellGrade™ prem, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 5 piece(s)	1	pack

BRANDplates® microplates, 96-well, inertGrade™

Cat. No.	Description		From*	Unit
781900	Microplate, BRANDplates®, inertGrade™, 96-well, PS Standard transp., 'U'bottom 330 µl, sterile	1 pack = 40 piece(s)	1	pack
781902	Microplate, BRANDplates®, inertGrade™, 96-well, PS Standard transp., 'F'bottom 350 µl, sterile	1 pack = 40 piece(s)	1	pack
781910	Microplate, BRANDplates®, inertGrade™, 96-well, PS transp., b., white, 'F'bottom 330 µl, sterile	1 pack = 40 piece(s)	1	pack
781911	Microplate, BRANDplates®, inertGrade™, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 40 piece(s)	1	pack
781912	Microplate, BRANDplates®, inertGrade™, 96-well, PS transp., b., white, 'F'bottom 330 µl, sterile	1 pack = 5 piece(s)	1	pack
781913	Microplate, BRANDplates®, inertGrade™, 96-well, PS transp., b., black, 'F'bottom 330 µl, sterile	1 pack = 5 piece(s)	1	pack


Lid, for BRAND®plates

Cat. No.	Description		From*	Unit
782150	Lid with condensation rings for BRANDplates® 96-well, Standard	1 pack = 100 piece(s)	1	pack
782151	Lid w/o condensation rings for BRANDplates® 96-&384-well, trans., b. & Low Vol.	1 pack = 100 piece(s)	1	pack
782152	Lid w/o condensation rings for BRANDplates® 384-well, standard	1 pack = 50 piece(s)	1	pack
782153	Lid w/o condensation rings for BRANDplates® 1536-well, Standard	1 pack = 50 piece(s)	1	pack
782155	Lid with condens.rings for BRANDplates®, 96-well transp.bottom 100 pcs.bag of 5	1 pack = 100 piece(s)	1	pack




96-well strip plate (12 x F8), medium binding

Cat. No.	Description		From*	Unit
782300	Strip plates, 96-well, 12xF8, non dividable, pureGrade™, MB	1 pack = 100 piece(s)	1	pack
782301	Strip plates w.grid, 96-well, 12xF8, dividable, pureGrade™, MB	1 pack = 100 piece(s)	1	pack



96-well strip plate (12 x F8), high binding

Cat. No.	Description		From*	Unit
782305	Strip plates, 96-well, 12xF8, non dividable, immunoGrade™, HB	1 pack = 100 piece(s)	1	pack
782306	Strip plates with grid, 96-well, 12xF8, dividable, immunoGrade™, HB	1 pack = 100 piece(s)	1	pack



Media bottle



Cat. No.	Description	From*	Unit
122710	Media bottle 500 ml, GL 56, cap PP, low	4	piece(s)

Accessories for reagent bottle



Cat. No.	Description	From*	Unit
122750	Screw cap for media bottle 500 ml, GL 56, PP, white	10	piece(s)
122755	Thread adapter, PTFE, for media bottle, GL 45, autoclavable	1	piece(s)
122760	Vacuum filter funnel for media bottle, 0,2 µm sterile membrane PES GL 45 12 p.	12	piece(s)
122765	Vacuum filter funnel for media bottle, 0,45 µm sterile membrane PES GL 45 12 p.	12	piece(s)

Cuvettes, 4 clear sided



Cat. No.	Description	From*	Unit
759030	macro cuvettes, four clear sided, PS, minimum filling volume 2,5 ml, 100 pieces	100	piece(s)
759035	macro cuvettes, four clear sided, PS, minimum filling volume 2,5 ml, 500 pieces	500	piece(s)
759125	macro cuvettes, four clear sided, UV-transparent, minimum filling volume 2,5 ml, 100 pieces	100	piece(s)
759128	macro cuvettes, four clear sided, UV-transparent, minimum filling volume 2,5 ml, 500 pieces	500	piece(s)

Cuvettes, macro and semi-micro, PS



Cat. No.	Description	From*	Unit
759005	Macro-Cuvette, PS, minimum filling volume 2,5 ml	100	piece(s)
759015	Semi-micro cuvette, PS, minimum filling volume 1,5 ml	100	piece(s)

Cuvettes, macro and semi-micro, PMMA



Cat. No.	Description	From*	Unit
759105	Macro-Cuvette, PMMA, minimum filling volume 2,5 ml	100	piece(s)
759115	Semi-micro cuvette, PMMA, minimum filling volume 1,5 ml	100	piece(s)

UV-Cuvettes, macro

Cat. No.	Description	From*	Unit
759170	UV-Cuvette, macro, minimum filling volume 2,5 ml	100	piece(s)



UV-Cuvettes, semi-micro

Cat. No.	Description	From*	Unit
759150	UV-Cuvette, semi-micro, minimum filling volume 1,5 ml	100	piece(s)



UV-Cuvettes, micro

Cat. No.	Description	From*	Unit
759200	UV-Cuvette micro, center height 8.5 mm, filling volume 70 µl to 850 µl	100	piece(s)
759210	UV-Cuvette micro, center height 8.5 mm, filling volume 70 µl to 850 µl	500	piece(s)
759215	UV-Cuvette micro, center height 8.5 mm, filling volume 70 µl to 850 µl	100	piece(s)
759220	UV-Cuvette micro, center height 15 mm, filling volume 70 µl to 550 µl	100	piece(s)
759230	UV-Cuvette micro, center height 15 mm, filling volume 70 µl to 550 µl	500	piece(s)
759235	UV-Cuvette micro, center height 15 mm, filling volume 70 µl to 550 µl, individually wrapped	100	piece(s)



Cuvette rack, PP

Cat. No.	Description	From*	Unit
759500	Cuvette rack, PP, grey, for 16 cuvettes, 210 x 70 x 38 mm	1	piece(s)



Lid, for cuvettes

Cat. No.	Description	From*	Unit
759240	Cap for UV-Cuvette micro, blue	100	piece(s)
759241	Cap for UV-Cuvette micro, yellow	100	piece(s)
759242	Cap for UV-Cuvette micro, green	100	piece(s)
759243	Cap for UV-Cuvette micro, orange	100	piece(s)



Disposable stirring spatula

Cat. No.	Description	From*	Unit
759800	Disp. stirring spatula, PS, length 120 mm, rod diameter 3 mm	500	piece(s)

Bulb pipettes, BLAUBRAND®, class AS, 1 mark, AR-Glas®, DE-M

Cat. No.	Description	From*	Unit
29701	Bulb pipette, BLAUBRAND®, AS, DE-M, 0,5ml, one-mark	12	piece(s)
29702	Bulb pipette, BLAUBRAND®, AS, DE-M, 1ml, one-mark	12	piece(s)
29703	Bulb pipette, BLAUBRAND®, AS, DE-M, 2ml, one-mark	12	piece(s)
29704	Bulb pipette, BLAUBRAND®, AS, DE-M, 2,5ml, one-mark	6	piece(s)
29705	Bulb pipette, BLAUBRAND®, AS, DE-M, 3ml, one-mark	6	piece(s)
29706	Bulb pipette, BLAUBRAND®, AS, DE-M, 4ml, one-mark	6	piece(s)
29707	Bulb pipette, BLAUBRAND®, AS, DE-M, 5ml, one-mark	6	piece(s)
29708	Bulb pipette, BLAUBRAND®, AS, DE-M, 6ml, one-mark	6	piece(s)
29709	Bulb pipette, BLAUBRAND®, AS, DE-M, 7ml, one-mark	6	piece(s)
29710	Bulb pipette, BLAUBRAND®, AS, DE-M, 8ml, one-mark	6	piece(s)
29711	Bulb pipette, BLAUBRAND®, AS, DE-M, 9ml, one-mark	6	piece(s)
29712	Bulb pipette, BLAUBRAND®, AS, DE-M, 10ml, one-mark	6	piece(s)
29713	Bulb pipette, BLAUBRAND®, AS, DE-M, 15ml, one-mark	6	piece(s)
29714	Bulb pipette, BLAUBRAND®, AS, DE-M, 20ml, one-mark	6	piece(s)
29715	Bulb pipette, BLAUBRAND®, AS, DE-M, 25ml, one-mark	6	piece(s)
29716	Bulb pipette, BLAUBRAND®, AS, DE-M, 30ml, one-mark	6	piece(s)
29717	Bulb pipette, BLAUBRAND®, AS, DE-M, 40ml, one-mark	6	piece(s)
29718	Bulb pipette, BLAUBRAND®, AS, DE-M, 50ml, one-mark	6	piece(s)
29719	Bulb pipette, BLAUBRAND®, AS, DE-M, 100ml, one-mark	6	piece(s)


Bulb pipettes, BLAUBRAND®, class AS, 2 marks, AR-Glas®, DE-M

Cat. No.	Description	From*	Unit
29721	Bulb pipette BLAUBRAND®, AS, DE-M, 0,5 ml, two marks, AR-GLAS®	12	piece(s)
29722	Bulb pipette BLAUBRAND®, AS, DE-M, 1 ml, two marks, AR-GLAS®	12	piece(s)
29723	Bulb pipette BLAUBRAND®, AS, DE-M, 2 ml, two marks, AR-GLAS®	12	piece(s)
29725	Bulb pipette BLAUBRAND®, AS, DE-M, 3 ml, two marks, AR-GLAS®	6	piece(s)
29727	Bulb pipette BLAUBRAND®, AS, DE-M, 5 ml, two marks, AR-GLAS®	6	piece(s)
29732	Bulb pipette BLAUBRAND®, AS, DE-M, 10 ml, two marks, AR-GLAS®	6	piece(s)
29733	Bulb pipette BLAUBRAND®, AS, DE-M, 15 ml, two marks, AR-GLAS®	6	piece(s)
29734	Bulb pipette BLAUBRAND®, AS, DE-M, 20 ml, two marks, AR-GLAS®	6	piece(s)
29735	Bulb pipette BLAUBRAND®, AS, DE-M, 25 ml, two marks, AR-GLAS®	6	piece(s)
29738	Bulb pipette BLAUBRAND®, AS, DE-M, 50 ml, two marks, AR-GLAS®	6	piece(s)



Bulb pipettes, USP, BLAUBRAND®, class AS, 1 mark, AR-Glas®, DE-M


Cat. No.	Description	From*	Unit
30601	Bulb pipette USP BLAUBRAND®, AS, DE-M, 0,5 ml, one-mark, AR-GLAS®	12	piece(s)
30602	Bulb pipette USP BLAUBRAND®, AS, DE-M, 1 ml, one-mark, AR-GLAS®	12	piece(s)
30603	Bulb pipette USP BLAUBRAND®, AS, DE-M, 2 ml, one-mark, AR-GLAS®	12	piece(s)
30605	Bulb pipette USP BLAUBRAND®, AS, DE-M, 3 ml, one-mark, AR-GLAS®	6	piece(s)
30606	Bulb pipette USP BLAUBRAND®, AS, DE-M, 4 ml, one-mark, AR-GLAS®	6	piece(s)
30607	Bulb pipette USP BLAUBRAND®, AS, DE-M, 5 ml, one-mark, AR-GLAS®	6	piece(s)
30608	Bulb pipette USP BLAUBRAND®, AS, DE-M, 6 ml, one-mark, AR-GLAS®	6	piece(s)
30609	Bulb pipette USP BLAUBRAND®, AS, DE-M, 7 ml, one-mark, AR-GLAS®	6	piece(s)
30610	Bulb pipette USP BLAUBRAND®, AS, DE-M, 8 ml, one-mark, AR-GLAS®	6	piece(s)
30611	Bulb pipette USP BLAUBRAND®, AS, DE-M, 9 ml, one-mark, AR-GLAS®	6	piece(s)
30612	Bulb pipette USP BLAUBRAND®, AS, DE-M, 10 ml, one-mark, AR-GLAS®	6	piece(s)
30613	Bulb pipette USP BLAUBRAND®, AS, DE-M, 15 ml, one-mark, AR-GLAS®	6	piece(s)
30614	Bulb pipette USP BLAUBRAND®, AS, DE-M, 20 ml, one-mark, AR-GLAS®	6	piece(s)
30615	Bulb pipette USP BLAUBRAND®, AS, DE-M, 25 ml, one-mark, AR-GLAS®	6	piece(s)
30616	Bulb pipette USP BLAUBRAND®, AS, DE-M, 30 ml, one-mark, AR-GLAS®	6	piece(s)
30617	Bulb pipette USP BLAUBRAND®, AS, DE-M, 40 ml, one-mark, AR-GLAS®	6	piece(s)
30618	Bulb pipette USP BLAUBRAND®, AS, DE-M, 50 ml, one-mark, AR-GLAS®	6	piece(s)
30619	Bulb pipette USP BLAUBRAND®, AS, DE-M, 100 ml, one-mark, AR-GLAS®	6	piece(s)

Bulb pipettes, BLAUBRAND® ETERNA, class AS, 1 mark, AR-Glas®, DE-M


Cat. No.	Description	From*	Unit
30502	Bulb pipette BLAUBRAND® ETERNA, AS, DE-M, 1 ml, one-mark, AR-GLAS®	12	piece(s)
30503	Bulb pipette BLAUBRAND® ETERNA, AS, DE-M, 2 ml, one-mark, AR-GLAS®	12	piece(s)
30507	Bulb pipette BLAUBRAND® ETERNA, AS, DE-M, 5 ml, one-mark, AR-GLAS®	6	piece(s)
30512	Bulb pipette BLAUBRAND® ETERNA, AS, DE-M, 10 ml, one-mark, AR-GLAS®	6	piece(s)
30514	Bulb pipette BLAUBRAND® ETERNA, AS, DE-M, 20 ml, one-mark, AR-GLAS®	6	piece(s)
30515	Bulb pipette BLAUBRAND® ETERNA, AS, DE-M, 25 ml, one-mark, AR-GLAS®	6	piece(s)
30518	Bulb pipette BLAUBRAND® ETERNA, AS, DE-M, 50 ml, one-mark, AR-GLAS®	6	piece(s)

Bulb pipettes, SILBERBRAND ETERNA, class B, AR-GLAS®


Cat. No.	Description	From*	Unit
29501	Bulb pipette SILBERBRAND ETERNA, B, 0,5 ml, one-mark, AR-GLAS®	12	piece(s)
29502	Bulb pipette SILBERBRAND ETERNA, B, 1 ml, one-mark, AR-GLAS®	12	piece(s)
29503	Bulb pipette SILBERBRAND ETERNA, B, 2 ml, one-mark, AR-GLAS®	12	piece(s)
29507	Bulb pipette SILBERBRAND ETERNA, B, 5 ml, one-mark, AR-GLAS®	6	piece(s)
29512	Bulb pipette SILBERBRAND ETERNA, B, 10 ml, one-mark, AR-GLAS®	6	piece(s)
29514	Bulb pipette SILBERBRAND ETERNA, B, 20 ml, one-mark, AR-GLAS®	6	piece(s)
29515	Bulb pipette SILBERBRAND ETERNA, B, 25 ml, one-mark, AR-GLAS®	6	piece(s)
29518	Bulb pipette SILBERBRAND ETERNA, B, 50 ml, one-mark, AR-GLAS®	6	piece(s)
29519	Bulb pipette SILBERBRAND ETERNA, B, 100 ml, one-mark, AR-GLAS®	6	piece(s)


Bulb pipettes, PP

Cat. No.	Description	From*	Unit
30002	Bulb pipette, PP, 1 ml, one-mark	12	piece(s)
30003	Bulb pipette, PP, 2 ml, one-mark	12	piece(s)
30007	Bulb pipette, PP, 5 ml, one-mark	6	piece(s)
30012	Bulb pipette, PP, 10 ml, one-mark	6	piece(s)
30015	Bulb pipette, PP, 25 ml, one-mark	6	piece(s)
30018	Bulb pipette, PP, 50 ml, one-mark	6	piece(s)


Graduated pipettes, BLAUBRAND® ETERNA, type 3, AR-Glas®, DE-M

Cat. No.	Description	From*	Unit
27316	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 1:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27319	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 2:0,02 ml, total delivery, AR-GLAS®	12	piece(s)
28405	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 0,5:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
28406	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 1:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
28409	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 2:0,02 ml, total delivery, AR-GLAS®	12	piece(s)
28411	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 5:0,05 ml, total delivery, AR-GLAS®	12	piece(s)
28412	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 5:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
28413	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 10:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
28414	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 20:0,1 ml, total delivery, AR-GLAS®	6	piece(s)
28415	Graduated pipette, BLAUBRAND® ETERNA, AS, DE-M, type 3, 25:0,1 ml, total delivery, AR-GLAS®	6	piece(s)




Graduated pipettes, BLAUBRAND®, type graduated to contain, AR-Glas®, DE-M


Cat. No.	Description	From*	Unit
27702	Graduated pipette, BLAUBRAND®, A, DE-M, 0,1 ml:0,001 ml calibrated to contain, AR-GLAS®	12	piece(s)
27704	Graduated pipette, BLAUBRAND®, A, DE-M, 0,2 ml:0,002 ml calibrated to contain, AR-GLAS®	12	piece(s)

Graduated pipettes, SILBERBRAND ETERNA, type 2, class B



Cat. No.	Description	From*	Unit
27002	Graduated pipette, SILBERBRAND ETERNA, B, type 2, 1:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27004	Graduated pipette, SILBERBRAND ETERNA, B, type 2, 1:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27005	Graduated pipette, SILBERBRAND ETERNA, B, type 2, 2:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27006	Graduated pipette, SILBERBRAND ETERNA, B, type 2, 2:0,02 ml, total delivery, AR-GLAS®	12	piece(s)
27009	Graduated pipette, SILBERBRAND ETERNA, B, type 2, 5:0,05 ml, total delivery, AR-GLAS®	12	piece(s)
27010	Graduated pipette, SILBERBRAND ETERNA, B, type 2, 5:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27011	Graduated pipette, SILBERBRAND ETERNA, B, type 2, 10:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27013	Graduated pipette, SILBERBRAND ETERNA, B, type 2, 25:0,1 ml, total delivery, AR-GLAS®	6	piece(s)

Graduated pipettes, SILBERBRAND ETERNA, type 3, class B


Cat. No.	Description	From*	Unit
27070	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 1:0,01 ml, total delivery	12	piece(s)
27071	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 1:0,1 ml, total delivery	12	piece(s)
27072	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 2:0,01 ml, total delivery	12	piece(s)
27073	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 2:0,02 ml, total delivery	12	piece(s)
27075	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 5:0,05 ml, total delivery	12	piece(s)
27076	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 5:0,1 ml, total delivery	12	piece(s)
27077	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 10:0,1 ml, total delivery	12	piece(s)
27078	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 20:0,1 ml, total delivery	6	piece(s)
27079	Graduated pipette, SILBERBRAND ETERNA, B, type 3 (zero on top), 25:0,1 ml, total delivery	6	piece(s)

Graduated pipettes, BLAUBRAND®, class AS, type 1

Cat. No.	Description	From*	Unit
27721	Graduated pipette, BLAUBRAND®, AS, DE-M, type 1, 0,5:0,01 ml, partial delivery, AR-GLAS®	12	piece(s)
27722	Graduated pipette, BLAUBRAND®, AS, DE-M, type 1, 1:0,01 ml, partial delivery, AR-GLAS®	12	piece(s)
27724	Graduated pipette, BLAUBRAND®, AS, type 1, 2:0,01 ml, partial delivery, AR-GLAS®	12	piece(s)
27725	Graduated pipette, BLAUBRAND®, AS, DE-M, type 1, 2:0,02 ml, partial delivery, AR-GLAS®	12	piece(s)
27727	Graduated pipette, BLAUBRAND®, AS, DE-M, type 1, 5:0,05 ml, partial delivery, AR-GLAS®	12	piece(s)
27729	Graduated pipette, BLAUBRAND®, AS, DE-M, type 1, 10:0,1 ml, partial delivery, AR-GLAS®	12	piece(s)
27731	Graduated pipette, BLAUBRAND®, AS, DE-M, type 1, 25:0,1 ml, partial delivery, AR-GLAS®	6	piece(s)




Graduated pipettes, BLAUBRAND®, class AS, type 2, DE-M

Cat. No.	Description	From*	Unit
27816	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 0,5:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27817	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 1:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27818	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 1:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27819	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 2:0,02 ml, total delivery, AR-GLAS®	12	piece(s)
27820	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 2:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27821	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 5:0,05 ml, total delivery, AR-GLAS®	12	piece(s)
27822	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 5:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27823	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 10:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27824	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 20:0,1 ml, total delivery, AR-GLAS®	6	piece(s)
27825	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 25:0,1 ml, total delivery, AR-GLAS®	6	piece(s)
27826	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 50:0,5 ml, total delivery, AR-GLAS®	6	piece(s)
27827	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 1:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27828	Graduated pipette, BLAUBRAND®, AS, DE-M, type 2, 2:0,02 ml, total delivery, AR-GLAS®	12	piece(s)
27829	Graduated pipette, BLAUBRAND®, AS, type 2, 2:0,01 ml, total delivery, AR-GLAS®	12	piece(s)



Graduated pipettes, BLAUBRAND® type 2, class AS, USP, DE-M


Cat. No.	Description	From*	Unit
27506	Graduated pipette, USP, BLAUBRAND®, AS, DE-M, type 2, 1:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27509	Graduated pipette, USP, BLAUBRAND®, AS, DE-M, type 2, 2:0,02 ml, total delivery, AR-GLAS®	12	piece(s)
27511	Graduated pipette, USP, BLAUBRAND®, AS, DE-M, type 2, 5:0,05 ml, total delivery, AR-GLAS®	12	piece(s)
27513	Graduated pipette, USP, BLAUBRAND®, AS, DE-M, type 2, 10:0,1 ml, total delivery, AR-GLAS®	12	piece(s)

Graduated pipettes, BLAUBRAND®, class AS, type 3, DE-M


Cat. No.	Description	From*	Unit
27306	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 1:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27309	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 2:0,02 ml, total delivery, AR-GLAS®	12	piece(s)
27705	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 0,5:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27706	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 1:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27707	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 1:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27708	Graduated pipette, BLAUBRAND®, AS, type 3, 2:0,01 ml, total delivery, AR-GLAS®	12	piece(s)
27709	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 2:0,02 ml, total delivery, AR-GLAS®	12	piece(s)
27710	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 2:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27711	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 5:0,05 ml, total delivery, AR-GLAS®	12	piece(s)
27712	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 5:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27713	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 10:0,1 ml, total delivery, AR-GLAS®	12	piece(s)
27714	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 20:0,1 ml, total delivery, AR-GLAS®	6	piece(s)
27715	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 25:0,1 ml, total delivery, AR-GLAS®	6	piece(s)
27716	Graduated pipette, BLAUBRAND®, AS, DE-M, type 3, 50:0,5 ml, total delivery, AR-GLAS®	6	piece(s)

Graduated pipettes, PP, type 3

Cat. No.	Description	From*	Unit
27607	Graduated pipette, PP, 1 ml:0,1 ml, total delivery, diameter, 8 mm	12	piece(s)
27610	Graduated pipette, PP, 2 ml:0,1 ml, total delivery, diameter 8 mm	12	piece(s)
27612	Graduated pipette, PP, 5 ml:0,1 ml, total delivery, diameter 8 mm	12	piece(s)
27613	Graduated pipette, PP, 10 ml:0,1 ml, total delivery, diameter 10 mm	12	piece(s)
27614	Graduated pipette, PP, 10 ml:0,1 ml, total delivery, diameter 8 mm	12	piece(s)




Graduated pipettes, piston type, SILBERBRAND, AR-Glas®

Cat. No.	Description	From*	Unit
31106	Graduated pipette, piston type, SILBERBRAND, 1 ml:0,01 ml, calibrated 'Ex' to jet, AR-GLAS®	1	piece(s)
31109	Graduated pipette, piston type, SILBERBRAND, 2 ml:0,02 ml, calibrated 'Ex' to jet, AR-GLAS®	1	piece(s)
31111	Graduated pipette, piston type, SILBERBRAND, 5 ml:0,05 ml, calibrated 'Ex' to jet, AR-GLAS®	1	piece(s)
31113	Graduated pipette, piston type, SILBERBRAND, 10 ml:0,1 ml, calibrated 'Ex' to jet, AR-GLAS®	1	piece(s)
31115	Graduated pipette, piston type, SILBERBRAND, 25 ml:0,1 ml, calibrated 'Ex' to jet, AR-GLAS®	1	piece(s)



Graduated pipettes, serological, total delivery, SILBERBRAND ETERNA, AR-Glas®

Cat. No.	Description	From*	Unit
27107	Graduated pipette serological calibrated 'Ex' to jet, 1 ml:0,01 ml, tip dia. approx. 2 mm	12	piece(s)
27109	Graduated pipette serological calibrated 'Ex' to jet, 2 ml:0,02 ml, tip dia. approx. 2 mm	12	piece(s)
27112	Graduated pipette serological calibrated 'Ex' to jet, 5 ml:0,1 ml, tip dia. approx. 3 mm	12	piece(s)
27113	Graduated pipette serological calibrated 'Ex' to jet 10 ml:0,1 ml, tip dia. approx. 3 mm	12	piece(s)
27115	Graduated pipette serological calibrated 'Ex' to jet 25 ml:0,1 ml, tip dia. approx. 3 mm	6	piece(s)



Cotton cord for pipettes, 100 % cotton wool

Cat. No.	Description	From*	Unit
28205	Cotton roving for pipettes, 100% cotton wool, degreased, light weight, approx. 13g/10m	1	piece(s)



Pipette jar, aluminium, with inner silicone coating



Cat. No.	Description	From*	Unit
28920	Pipette box aluminum (150-200 mm)	1	piece(s)
28925	Pipette box aluminum (235-305 mm)	1	piece(s)
28930	Pipette box aluminum (315-385 mm)	1	piece(s)
28935	Pipette box aluminum (315-485 mm)	1	piece(s)

Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with PP/PE stopper



Cat. No.	Description	From*	Unit
37233	Volumetric flask, BLAUBRAND®, A, DE-M, 150 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
37234	Volumetric flask, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, W, NS 29/32, PP stopper	1	piece(s)
37236	Volumetric flask, BLAUBRAND®, A, DE-M, 10000 ml, Boro 3.3, NS 45/40, PE stopper	1	piece(s)
37238	Volumetric flask, BLAUBRAND®, A, DE-M, 5 ml, Boro 3.3, W, NS 10/19, PP stopper	2	piece(s)
37241	Volumetric flask, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, W, NS 10/19, PP stopper	2	piece(s)
37243	Volumetric flask, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, W, NS 10/19, PP stopper	2	piece(s)
37245	Volumetric flask, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, W, NS 12/21, PP stopper	2	piece(s)
37246	Volumetric flask, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, NS 10/19, PP stopper	2	piece(s)
37247	Volumetric flask, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, NS 10/19, PP stopper	2	piece(s)
37248	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, NS 12/21, PP stopper	2	piece(s)
37249	Volumetric flask, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 12/21, PP stopper	2	piece(s)
37250	Volumetric flask, BLAUBRAND®, A, DE-M, 200 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
37251	Volumetric flask, BLAUBRAND®, A, DE-M, 250 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
37252	Volumetric flask, BLAUBRAND®, A, DE-M, 500 ml, Boro 3.3, NS 19/26, PP stopper	2	piece(s)
37253	Volumetric flask, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, NS 24/29, PP stopper	1	piece(s)
37254	Volumetric flask, BLAUBRAND®, A, DE-M, 2000 ml, Boro 3.3, NS 29/32, PP stopper	1	piece(s)
37255	Volumetric flask, BLAUBRAND®, A, DE-M, 5000 ml, Boro 3.3, NS 34/35, PE stopper	1	piece(s)
37290	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, W, NS 14/23, PP stopper	2	piece(s)
37291	Volumetric flask, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
37293	Volumetric flask, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, W, NS 12/21, PP stopper	2	piece(s)

Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with PP stopper, PURprotect

Cat. No.	Description	From*	Unit
36548	Volumetric flask PUR-coated, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, W, NS 14/23, PP stopper	2	piece(s)
36549	Volumetric flask PUR-coated, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36550	Volumetric flask PUR-coated, BLAUBRAND®, A, DE-M, 200 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36551	Volumetric flask PUR-coated, BLAUBRAND®, A, DE-M, 250 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36552	Volumetric flask PUR-coated, BLAUBRAND®, A, DE-M, 500 ml, Boro 3.3, NS 19/26, PP stopper	2	piece(s)
36553	Volumetric flask PUR-coated, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, NS 24/29, PP stopper	1	piece(s)


Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with glass stopper

Cat. No.	Description	From*	Unit
37256	Volumetric flask, BLAUBRAND®, A, DE-M, 5 ml, Boro 3.3, W, NS 10/19, Boro 3.3, glass stopper	2	piece(s)
37257	Volumetric flask, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, NS 10/19, Boro 3.3, glass stopper	2	piece(s)
37258	Volumetric flask, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, NS 10/19, Boro 3.3, glass stopper	2	piece(s)
37259	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, NS 12/21, Boro 3.3, glass stopper	2	piece(s)
37260	Volumetric flask, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 12/21, Boro 3.3, glass stopper	2	piece(s)
37261	Volumetric flask, BLAUBRAND®, A, DE-M, 200 ml, Boro 3.3, NS 14/23, Boro 3.3, glass stopper	2	piece(s)
37262	Volumetric flask, BLAUBRAND®, A, DE-M, 250 ml, Boro 3.3, NS 14/23, Boro 3.3, glass stopper	2	piece(s)
37263	Volumetric flask, BLAUBRAND®, A, DE-M, 500 ml, Boro 3.3, NS 19/26, Boro 3.3, glass stopper	2	piece(s)
37264	Volumetric flask, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, NS 24/29, Boro 3.3, glass stopper	1	piece(s)
37265	Volumetric flask, BLAUBRAND®, A, DE-M, 2000 ml, Boro 3.3, NS 29/32, Boro 3.3, glass stopper	1	piece(s)
37266	Volumetric flask, BLAUBRAND®, A, DE-M, 5000 ml, Boro 3.3, NS 34/35, Boro 3.3, glass stopper	1	piece(s)
37267	Volumetric flask, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, W, NS 10/19, Boro 3.3, glass stopper	2	piece(s)
37268	Volumetric flask, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, W, NS 12/21, Boro 3.3, glass stopper	2	piece(s)
37288	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, W, NS 14/23, Boro 3.3, glass stopper	2	piece(s)
37289	Volumetric flask, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 14/23, Boro 3.3, glass stopper	2	piece(s)
37294	Volumetric flask, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, W, NS 12/21, Boro 3.3, glass stopper	2	piece(s)



Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, beaded rim


Cat. No.	Description	From*	Unit
37045	Volumetric flask, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, W, beaded rim	2	piece(s)
37046	Volumetric flask, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, beaded rim	2	piece(s)
37047	Volumetric flask, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, beaded rim	2	piece(s)
37048	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, beaded rim	2	piece(s)
37049	Volumetric flask, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, beaded rim	2	piece(s)
37050	Volumetric flask, BLAUBRAND®, A, DE-M, 200 ml, Boro 3.3, beaded rim	2	piece(s)
37051	Volumetric flask, BLAUBRAND®, A, DE-M, 250 ml, Boro 3.3, beaded rim	2	piece(s)
37052	Volumetric flask, BLAUBRAND®, A, DE-M, 500 ml, Boro 3.3, beaded rim	2	piece(s)
37053	Volumetric flask, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, beaded rim	1	piece(s)

Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with PP stopper, trapezoidal


Cat. No.	Description	From*	Unit
36401	Volumetric flask trapezoidal, BLAUBRAND®, A, DE-M, 1ml, Boro 3.3, NS 7/16 PP-stopper	2	piece(s)
36402	Volumetric flask trapezoidal, BLAUBRAND®, A, DE-M, 2ml, Boro 3.3, NS 7/16 PP-stopper	2	piece(s)
36403	Volumetric flask trapezoidal, BLAUBRAND®, A, DE-M, 5 ml, Boro 3.3, NS, 7/16, PP stopper	2	piece(s)
36404	Volumetric flask trapezoidal, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, W, NS 10/19, PP stopper	2	piece(s)
36406	Volumetric flask trapezoidal, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, NS 10/19, PP stopper	2	piece(s)
36407	Volumetric flask trapezoidal, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, NS 10/19, PP stopper	2	piece(s)
36408	Volumetric flask trapezoidal, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, NS 12/21, PP stopper	2	piece(s)

Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with PP stopper, amber

Cat. No.	Description	From*	Unit
37401	Volumetric flask, BLAUBRAND®, A, DE-M, 5 ml, Boro 3.3, W, NS 10/19, PP stopper, amber	2	piece(s)
37402	Volumetric flask, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, W, NS 10/19, PP stopper, amber	2	piece(s)
37403	Volumetric flask, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, NS 10/19, PP stopper, amber	2	piece(s)
37404	Volumetric flask, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, NS 10/19, PP stopper, amber	2	piece(s)
37405	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, NS 12/21, PP stopper, amber	2	piece(s)
37406	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, W, NS 14/23, PP stopper, amber	2	piece(s)
37407	Volumetric flask, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 14/23, PP stopper, amber	2	piece(s)
37408	Volumetric flask, BLAUBRAND®, A, DE-M, 200 ml, Boro 3.3, NS 14/23, PP stopper, amber	2	piece(s)
37409	Volumetric flask, BLAUBRAND®, A, DE-M, 250 ml, Boro 3.3, NS 14/23, PP stopper, amber	2	piece(s)
37410	Volumetric flask, BLAUBRAND®, A, DE-M, 500 ml, Boro 3.3, NS 19/26, PP stopper, amber	2	piece(s)
37411	Volumetric flask, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, NS 24/29, PP stopper, amber	1	piece(s)
37412	Volumetric flask, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, NS 12/21, PP stopper, amber	2	piece(s)
37413	Volumetric flask, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, NS 12/21, PP stopper, amber	2	piece(s)


Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with glass stopper, amber

Cat. No.	Description	From*	Unit
37438	Volumetric flask, BLAUBRAND®, A, DE-M, 5 ml, Boro 3.3, W, NS 10/19, glass stopper, amber	2	piece(s)
37443	Volumetric flask, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, W, NS 10/19, glass stopper, amber	2	piece(s)
37445	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, W, NS 14/23, glass stopper, amber	2	piece(s)
37446	Volumetric flask, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, NS 10/19, glass stopper, amber	2	piece(s)
37447	Volumetric flask, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, NS 10/19, glass stopper, amber	2	piece(s)
37448	Volumetric flask, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, NS 12/21, glass stopper, amber	2	piece(s)
37449	Volumetric flask, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 14/23, glass stopper, amber	2	piece(s)
37450	Volumetric flask, BLAUBRAND®, A, DE-M, 200 ml, Boro 3.3, NS 14/23, glass stopper, amber	2	piece(s)
37451	Volumetric flask, BLAUBRAND®, A, DE-M, 250 ml, Boro 3.3, NS 14/23, glass stopper, amber	2	piece(s)
37452	Volumetric flask, BLAUBRAND®, A, DE-M, 500 ml, Boro 3.3, NS 19/26, glass stopper, amber	2	piece(s)
37453	Volumetric flask, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, NS 24/29, glass stopper, amber	1	piece(s)



Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with PP stopper, USP


Cat. No.	Description	From*	Unit
36938	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 5ml, Boro 3.3, W, NS 10/19, PP stopper	2	piece(s)
36943	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 10ml, Boro 3.3, W, NS 10/19, PP stopper	2	piece(s)
36945	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 20ml, Boro 3.3, W NS 10/19, PP stopper	2	piece(s)
36947	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 25ml, Boro 3.3, NS 10/19, PP stopper	2	piece(s)
36948	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 50ml, Boro 3.3, NS 12/21, PP stopper	2	piece(s)
36949	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 100ml, Boro 3.3, W, NS 14/23, PP stopper	2	piece(s)
36950	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 200ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36951	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 250ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36952	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 500ml, Boro 3.3, NS 19/26, PP stopper	2	piece(s)
36953	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 1000ml, Boro 3.3, NS 24/29, PP stopper	1	piece(s)
36954	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 2000ml, Boro 3.3, NS 29/32, PP stopper	1	piece(s)

Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with glass stopper, USP


Cat. No.	Description	From*	Unit
36968	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 5ml, NS 10/19, Boro 3.3, glass stopper	2	piece(s)
36973	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 10ml, NS 10/19, Boro 3.3, glass stopper	2	piece(s)
36975	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 20ml, NS 10/19, Boro 3.3, glass stopper	2	piece(s)
36977	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 25ml, NS 10/19, Boro 3.3, glass stopper	2	piece(s)
36978	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 50ml, NS 12/21, Boro 3.3, glass stopper	2	piece(s)
36979	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 100ml, NS 14/23, Boro 3.3, glass stopper	2	piece(s)
36980	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 200ml, NS 14/23, Boro 3.3, glass stopper	2	piece(s)
36981	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 250ml, NS 14/23, Boro 3.3, glass stopper	2	piece(s)
36982	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 500ml, NS 19/26, Boro 3.3, glass stopper	2	piece(s)
36983	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 1000ml, NS 24/29, Boro 3.3, glass stopper	1	piece(s)
36984	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 2000ml, NS 29/32, Boro 3.3, glass stopper	1	piece(s)

Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with PP stopper, amber, USP

Cat. No.	Description	From*	Unit
37481	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 5 ml, Boro 3.3, NS 10/19, PP stopper, amber	2	piece(s)
37482	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, NS 10/19, PP stopper, amber	2	piece(s)
37483	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, NS 10/19, PP stopper, amber	2	piece(s)
37484	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, NS 10/19, PP stopper, amber	2	piece(s)
37485	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, NS 12/21, PP stopper, amber	2	piece(s)
37487	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 14/23, PP stopper, amber	2	piece(s)
37488	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 200 ml, Boro 3.3, NS 14/23, PP stopper, amber	2	piece(s)
37489	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 250 ml, Boro 3.3, NS 14/23, PP stopper, amber	2	piece(s)
37490	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 500 ml, Boro 3.3, NS 19/26, PP stopper, amber	2	piece(s)
37491	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, NS 24/29, PP stopper, amber	1	piece(s)


Measuring flask, BLAUBRAND®, class A, Boro 3.3, DE-M, with glass stopper, amber, USP

Cat. No.	Description	From*	Unit
37461	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 5 ml, Boro 3.3, NS 10/19, glass stopper, amber	2	piece(s)
37462	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 10 ml, Boro 3.3, NS 10/19, glass stopper, amber	2	piece(s)
37463	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 20 ml, Boro 3.3, NS 10/19, glass stopper, amber	2	piece(s)
37464	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 25 ml, Boro 3.3, NS 10/19, glass stopper, amber	2	piece(s)
37465	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 50 ml, Boro 3.3, NS 12/21, glass stopper, amber	2	piece(s)
37467	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 100 ml, Boro 3.3, NS 14/23, glass stopper, amber	2	piece(s)
37468	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 200 ml, Boro 3.3, NS 14/23, glass stopper, amber	2	piece(s)
37469	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 250 ml, Boro 3.3, NS 14/23, glass stopper, amber	2	piece(s)
37470	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 500 ml, Boro 3.3, NS 19/26, glass stopper, amber	2	piece(s)
37471	Volumetric flask, USP, BLAUBRAND®, A, DE-M, 1000 ml, Boro 3.3, NS 24/29, glass stopper, amber	1	piece(s)



Measuring flask, BLAUBRAND® ETERNA, class A, Boro 3.3, DE-M, with PP stopper


Cat. No.	Description	From*	Unit
36838	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 5ml, Boro 3.3, W, N10/19, PP stopper	2	piece(s)
36841	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 25ml, Boro 3.3, NS 10/19, PP stopper	2	piece(s)
36842	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 25ml, Boro 3.3, W, NS 12/21, PP stopper	2	piece(s)
36843	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 10ml, Boro 3.3, W, NS 10/19, PP stopper	2	piece(s)
36845	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 20ml, Boro 3.3, NS 10/19, PP stopper	2	piece(s)
36846	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 20ml, Boro 3.3, W, NS 12/21, PP stopper	2	piece(s)
36847	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 50ml, Boro 3.3, NS 12/21, PP stopper	2	piece(s)
36848	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 50ml, Boro 3.3, W, NS 14/23, PP stopper	2	piece(s)
36849	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 100ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36850	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 200ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36851	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 250ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36852	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 500ml, Boro 3.3, NS 19/26, PP stopper	2	piece(s)
36853	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 1000ml, Boro 3.3, NS 24/29, PP stopper	1	piece(s)
36854	Volumetric flask, BLAUBRAND® ETERNA, A, DE-M, 2000ml, Boro 3.3, NS 29/32, PP stopper	1	piece(s)

Measuring flask, SILBERBRAND, class B, Boro 3.3, with PP stopper


Cat. No.	Description	From*	Unit
36743	Volumetric flask, SILBERBRAND, B, 10 ml, Boro 3.3, W, NS 10/19, PP stopper	2	piece(s)
36747	Volumetric flask, SILBERBRAND, B, 25 ml, Boro 3.3, NS 10/19, PP stopper	2	piece(s)
36748	Volumetric flask, SILBERBRAND, B, 50 ml, Boro 3.3, NS 12/21, PP stopper	2	piece(s)
36749	Volumetric flask, SILBERBRAND, B, 100 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36750	Volumetric flask, SILBERBRAND, B, 200 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36751	Volumetric flask, SILBERBRAND, B, 250 ml, Boro 3.3, NS 14/23, PP stopper	2	piece(s)
36752	Volumetric flask, SILBERBRAND, B, 500 ml, Boro 3.3, NS 19/26, PP stopper	2	piece(s)
36753	Volumetric flask, SILBERBRAND, B, 1000 ml, Boro 3.3, NS 24/29, PP stopper	1	piece(s)
36754	Volumetric flask, SILBERBRAND, B, 2000 ml, Boro 3.3, NS 29/32, PP stopper	1	piece(s)
36755	Volumetric flask, SILBERBRAND, B, 5000 ml, Boro 3.3, NS 34/35, PP stopper	1	piece(s)

Measuring flask, PFA, class A, with screw cap

Cat. No.	Description	From*	Unit
36208	Volumetric flask, PFA, 10ml, A, with screw cap GL18	1	piece(s)
36220	Volumetric flask, PFA, 25ml, A, with screw cap, GL18	1	piece(s)
36228	Volumetric flask, PFA, 50ml, A, with screw cap, GL18	1	piece(s)
36238	Volumetric flask, PFA, 100ml, A, with screw cap, GL18	1	piece(s)
36248	Volumetric flask, PFA, 250ml, A, with screw cap, GL25	1	piece(s)
36254	Volumetric flask, PFA, 500ml, A, with screw cap, GL25	1	piece(s)



Measuring flask, PMP, class A, with PP stopper

Cat. No.	Description	From*	Unit
36170	Volumetric flask, PMP, transparent, 10 ml, NS 10/19, PP stopper	1	piece(s)
36172	Volumetric flask, PMP, transparent, 25 ml, NS 10/19, PP stopper	1	piece(s)
36174	Volumetric flask, PMP, transparent, 50 ml, NS 12/21, PP stopper	1	piece(s)
36176	Volumetric flask, PMP, transparent, 100 ml, NS 14/23, PP stopper	1	piece(s)
36178	Volumetric flask, PMP, transparent, 250 ml, NS 19/26, PP stopper	1	piece(s)
36180	Volumetric flask, PMP, transparent, 500 ml, NS 19/26, PP stopper	1	piece(s)
36182	Volumetric flask, PMP, transparent, 1000 ml, NS 24/29, PP stopper	1	piece(s)



Measuring flask, PMP, class B, with PP stopper

Cat. No.	Description	From*	Unit
36108	Volumetric flask, PMP (TPX) clear, 10 ml, NS 10/19, PP stopper	1	piece(s)
36120	Volumetric flask, PMP (TPX) clear, 25 ml, NS 10/19, PP stopper	1	piece(s)
36128	Volumetric flask, PMP (TPX) clear, 50 ml, NS 12/21, PP stopper	1	piece(s)
36138	Volumetric flask, PMP (TPX) clear, 100 ml, NS 14/23, PP stopper	1	piece(s)
36148	Volumetric flask, PMP (TPX) clear, 250 ml, NS 19/26, PP stopper	1	piece(s)
36154	Volumetric flask, PMP (TPX) clear, 500 ml, NS 19/26, PP stopper	1	piece(s)
36162	Volumetric flask, PMP (TPX) clear, 1000 ml, NS 24/29, PP stopper	1	piece(s)



Measuring flask, PP, class B, with PP stopper

Cat. No.	Description	From*	Unit
36008	Volumetric flask, PP translucent, 10 ml, NS 10/19, PP stopper	1	piece(s)
36020	Volumetric flask, PP translucent, 25 ml, NS 10/19, PP stopper	1	piece(s)
36028	Volumetric flask, PP translucent, 50 ml, NS 12/21, PP stopper	1	piece(s)
36038	Volumetric flask, PP translucent, 100 ml, NS 14/23, PP stopper	1	piece(s)
36048	Volumetric flask, PP translucent, 250 ml, NS 19/26, PP stopper	1	piece(s)
36054	Volumetric flask, PP translucent, 500 ml, NS 19/26, PP stopper	1	piece(s)
36062	Volumetric flask, PP translucent, 1000 ml, NS 24/29, PP stopper	1	piece(s)



Measuring flask, 3 marks, BLAUBRAND®, class A, Boro 3.3, DE-M


Cat. No.	Description	From*	Unit
38204	Volumetric flasks with 3 marks, DAkkS-calibrated, BLAUBRAND®, A, DE-M, 10 ml, NS 10/19, PP stopper	1	piece(s)
38206	Volumetric flasks with 3 marks, DAkkS-calibrated, BLAUBRAND®, A, DE-M, 25 ml, NS 10/19, PP stopper	1	piece(s)
38208	Volumetric flasks with 3 marks, DAkkS-calibrated, BLAUBRAND®, A, DE-M, 50 ml, NS 12/21, PP stopper	1	piece(s)
38210	Volumetric flasks with 3 marks, DAkkS-calibrated, BLAUBRAND®, A, DE-M, 100 ml, NS 14/23, PP stopper	1	piece(s)

Measuring flask for sugar analysis, SILBERBRAND, class B, Boro 3.3


Cat. No.	Description	From*	Unit
402038	Volumetric flask, SILBERBRAND, Boro 3.3, 100 ml, for sugar analysis acc. Kohlrausch	2	piece(s)
402046	Volumetric flask, SILBERBRAND, Boro 3.3, 200 ml, for sugar analysis acc. Kohlrausch	2	piece(s)

Measuring flask for oil content determination, SILBERBRAND, class B, Boro 3.3


Cat. No.	Description	From*	Unit
365538	Volumetric flask, SILBERBRAND, hollow glass stopper, NS 19/26, 100 ml, Boro 3.3, oil con. det.	1	piece(s)

Conical ground joint stoppers, amber, Boro 3.3


Cat. No.	Description	From*	Unit
144671	Conical joint stopper, amber, Boro 3.3, NS 10/19, honed	1	piece(s)
144672	Conical joint stopper, amber, Boro 3.3, NS 12/21, honed	1	piece(s)
144673	Conical joint stopper, amber, Boro 3.3, NS 14/23, honed	1	piece(s)
144674	Conical joint stopper, amber, Boro 3.3, NS 19/26, honed	1	piece(s)
144676	Conical joint stopper, amber, Boro 3.3, NS 24/29, honed	1	piece(s)

Conical ground joint stopper, hollow, Boro 3.3

Cat. No.	Description	From*	Unit
144605	Conical joint stopper, Boro 3.3, NS, 7/16, hollow, hexagonal grip	1	piece(s)
144610	Conical joint stopper, Boro 3.3, NS 10/19, hollow, hexagonal grip	1	piece(s)
144615	Conical joint stopper, Boro 3.3, NS 12/21, hollow, hexagonal grip	1	piece(s)
144620	Conical joint stopper, Boro 3.3, NS 14/23, hollow, hexagonal grip	1	piece(s)
144625	Conical joint stopper, Boro 3.3, NS 19/26, hollow, hexagonal grip	1	piece(s)
144630	Conical joint stopper, Boro 3.3, NS 24/29, hollow, hexagonal grip	1	piece(s)
144635	Conical joint stopper, Boro 3.3, NS 29/32, hollow, hexagonal grip	1	piece(s)
144640	Conical joint stopper, Boro 3.3, NS 34/35, hollow, hexagonal grip	1	piece(s)
144645	Conical joint stopper, Boro 3.3, NS 45/40, hollow, hexagonal grip	1	piece(s)



Conical ground joint stopper, Boro 3.3

Cat. No.	Description	From*	Unit
144515	Conical joint stopper, Boro 3.3, NS 12/21, solid, octagonal grip	1	piece(s)
144520	Conical joint stopper, Boro 3.3, NS 14/23, solid, octagonal grip	1	piece(s)
144525	Conical joint stopper, Boro 3.3, NS 19/26, solid, octagonal grip	1	piece(s)
144530	Conical joint stopper, Boro 3.3, NS 24/29, solid, octagonal grip	1	piece(s)
144535	Conical joint stopper, Boro 3.3, NS 29/32, semi-hollow, octagonal grip	1	piece(s)
144540	Conical joint stopper, Boro 3.3, NS 34/35, semi-hollow, octagonal grip	1	piece(s)
144545	Conical joint stopper, Boro 3.3, NS 45/40, semi-hollow, octagonal grip	1	piece(s)



Conical ground joint stopper, PP

Cat. No.	Description	From*	Unit
144405	Conical joint stopper, PP, NS, 7/16, square grip, hollow, sealing rings	10	piece(s)
144410	Conical joint stopper, PP, NS 10/19, square grip, hollow, sealing rings	10	piece(s)
144415	Conical joint stopper, PP, NS 12/21, square grip, hollow, sealing rings	10	piece(s)
144420	Conical joint stopper, PP, NS 14/23, square grip, hollow, sealing rings	10	piece(s)
144425	Conical joint stopper, PP, NS 19/26, square grip, hollow, sealing rings	10	piece(s)
144430	Conical joint stopper, PP, NS 24/29, square grip, hollow, sealing rings	10	piece(s)
144435	Conical joint stopper, PP, NS 29/32, square grip, hollow, sealing rings	10	piece(s)



Conical ground joint stopper, PE-LD

Cat. No.	Description	From*	Unit
144440	Conical joint stopper, PE-LD, NS 34/35, octagonal grip, hollow, sealing rings	1	piece(s)
144445	Conical joint stopper, PE-LD, NS 45/40, octagonal grip, hollow, sealing rings	1	piece(s)



Stopper, natural rubber



Cat. No.	Description	From*	Unit
144382	Stopper, natural rubber, red, height 20 mm, upper diameter 12 mm, lower diameter 8 mm	25	piece(s)
144384	Stopper, natural rubber, red, height 20 mm, upper diameter 14 mm, lower diameter 11 mm	25	piece(s)
144385	Stopper, natural rubber, red, height 20 mm, upper diameter 16 mm, lower diameter 12 mm	25	piece(s)
144386	Stopper, natural rubber, red, height 20 mm, upper diameter 18 mm, lower diameter 14 mm	25	piece(s)
144388	Stopper, natural rubber, red, height 25 mm, upper diameter 22 mm, lower diameter 17 mm	25	piece(s)
144390	Stopper, natural rubber, red, height 30 mm, upper diameter 27 mm, lower diameter 21 mm	25	piece(s)

Stopper, silicone




Cat. No.	Description	From*	Unit
144310	Stopper, silicone, height 20 mm, upper diameter 12 mm, lower diameter 8 mm	10	piece(s)
144315	Stopper, silicone, height 20 mm, upper diameter 14,5 mm, lower diameter 10,5 mm	10	piece(s)
144320	Stopper, silicone, height 20 mm, upper diameter 16,5 mm, lower diameter 12,5 mm	10	piece(s)
144325	Stopper, silicone, height 20 mm, upper diameter 18 mm, lower diameter 14 mm	10	piece(s)
144330	Stopper, silicone, height 25 mm, upper diameter 22 mm, lower diameter 17 mm	10	piece(s)
144340	Stopper, silicone, height 30 mm, upper diameter 27 mm, lower diameter 21 mm	10	piece(s)

Replacement screw caps, PFA, for volumetric flask PFA



Cat. No.	Description	From*	Unit
129250	Replacement screw-cap for vol.etric flask PFA, GL 18	1	piece(s)
129252	Replacement screw-cap for vol.etric flask PFA, GL 25	1	piece(s)

Graduated cylinder, tall form, BLAUBRAND®, class A, Boro 3.3, DE-M




Cat. No.	Description	From*	Unit
32105	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 5 ml: 0,1 ml, Boro 3.3	2	piece(s)
32108	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 10 ml: 0,2 ml, Boro 3.3,	2	piece(s)
32120	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 25 ml: 0,5 ml, Boro 3.3	2	piece(s)
32128	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 50 ml: 1 ml, Boro 3.3,	2	piece(s)
32138	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 100 ml: 1 ml, Boro 3.3	2	piece(s)
32148	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 250 ml: 2 ml, Boro 3.3	2	piece(s)
32154	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 500 ml: 5 ml, Boro 3.3	2	piece(s)
32162	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 1000 ml: 10 ml, Boro 3.3	1	piece(s)
32164	Graduated cylinder, tall form BLAUBRAND®, A, DE-M, 2000 ml: 20 ml, Boro 3.3	1	piece(s)

Graduated cylinder, tall form, BLAUBRAND®, class A, Boro 3.3, DE-M, USP

Cat. No.	Description	From*	Unit
32805	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 5 ml:0,1 ml, Boro 3.3	2	piece(s)
32808	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 10 ml:0,2 ml, Boro 3.3	2	piece(s)
32820	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 25 ml:0,5 ml, Boro 3.3	2	piece(s)
32828	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 50 ml:1 ml, Boro 3.3	2	piece(s)
32838	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 100 ml:1 ml, Boro 3.3	2	piece(s)
32848	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 250 ml:2 ml, Boro 3.3	2	piece(s)
32854	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 500 ml:5 ml, Boro 3.3	2	piece(s)
32862	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 1000 ml:10 ml, Boro 3.3	1	piece(s)
32864	Graduated cylinder, USP, tall form, BLAUBRAND®, A, DE-M, 2000 ml:20 ml, Boro 3.3	1	piece(s)



Graduated cylinder, tall form, BLAUBRAND® ETERNA, class A, Boro 3.3, DE-M


Cat. No.	Description	From*	Unit
32705	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 5 ml: 0,1 ml Boro 3.3	2	piece(s)
32708	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 10 ml: 0,2 ml Boro 3.3	2	piece(s)
32720	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 25 ml: 0,5 ml Boro 3.3	2	piece(s)
32728	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 50 ml: 1 ml Boro 3.3	2	piece(s)
32738	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 100 ml: 1 ml Boro 3.3	2	piece(s)
32748	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 250 ml: 2 ml Boro 3.3	2	piece(s)
32754	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 500 ml: 5 ml Boro 3.3	2	piece(s)
32762	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 1000 ml: 10 ml Boro 3.3	1	piece(s)
32764	Graduated cylinder, tall form BLAUBRAND® ETERNA, A, DE-M, 2000 ml: 20 ml Boro 3.3	1	piece(s)

Graduated cylinder, tall form, SILBERBRAND, class B, Boro 3.3


Cat. No.	Description	From*	Unit
31708	Graduated cylinder, tall form, SILBERBRAND, B, Boro 3.3, 10 ml: 0,2 ml	2	piece(s)
31720	Graduated cylinder, tall form, SILBERBRAND, B, Boro 3.3, 25 ml: 0,5 ml	2	piece(s)
31728	Graduated cylinder, tall form, SILBERBRAND, B, Boro 3.3, 50 ml: 1 ml	2	piece(s)
31738	Graduated cylinder, tall form, SILBERBRAND, B, Boro 3.3, 100 ml: 1 ml	2	piece(s)
31748	Graduated cylinder, tall form, SILBERBRAND, B, Boro 3.3, 250 ml: 2 ml	2	piece(s)
31754	Graduated cylinder, tall form, SILBERBRAND, B, Boro 3.3, 500 ml: 5 ml	2	piece(s)
31762	Graduated cylinder, tall form, SILBERBRAND, B, Boro 3.3, 1000 ml: 10 ml	1	piece(s)
31764	Graduated cylinder, tall form, SILBERBRAND, B, Boro 3.3, 2000 ml: 20 ml	1	piece(s)

Graduated cylinder, tall form, SILBERBRAND ETERNA, class B, Boro 3.3


Cat. No.	Description	From*	Unit
31905	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 5 ml: 0,1 ml, with glass base	2	piece(s)
31908	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 10 ml: 0,2 ml, with glass base	2	piece(s)
31920	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 25 ml: 0,5 ml, with glass base	2	piece(s)
31928	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 50 ml: 1 ml, with glass base	2	piece(s)
31938	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 100 ml: 1 ml, with glass base	2	piece(s)
31948	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 250 ml: 2 ml, with glass base	2	piece(s)
31954	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 500 ml: 5 ml, with glass base	2	piece(s)
31962	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 1000 ml: 10 ml, with glass base	1	piece(s)
31964	Graduated cylinder, tall form SILBERBRAND ETERNA, B, Boro 3.3, 2000 ml: 20 ml, with glass base	1	piece(s)

Graduated cylinder, short form, SILBERBRAND ETERNA, class B, Boro 3.3

Cat. No.	Description	From*	Unit
42008	Graduated cylinder, low form, SILBERBRAND ETERNA, 10 ml: 1 ml, Boro 3.3, graduated in amber	2	piece(s)
42020	Graduated cylinder, low form, SILBERBRAND ETERNA, 25 ml: 1 ml, Boro 3.3, graduated in amber	2	piece(s)
42028	Graduated cylinder, low form, SILBERBRAND ETERNA, 50 ml: 2 ml, Boro 3.3, graduated in amber	2	piece(s)
42038	Graduated cylinder, low form, SILBERBRAND ETERNA, 100 ml: 2 ml, Boro 3.3, graduated in amber	2	piece(s)
42048	Graduated cylinder, low form, SILBERBRAND ETERNA, 250 ml: 5 ml, Boro 3.3, graduated in amber	2	piece(s)
42054	Graduated cylinder, low form, SILBERBRAND ETERNA, 500 ml:10 ml, Boro 3.3, graduated in amber	2	piece(s)
42062	Graduated cylinder, low form, SILBERBRAND ETERNA, 1000 ml:20 ml, Boro 3.3, graduated in amber	1	piece(s)
42064	Graduated cylinder, low form, SILBERBRAND ETERNA, 2000 ml:50 ml, Boro 3.3, graduated in amber	1	piece(s)


Graduated cylinder, tall form, PMP, class A, DE-M

Cat. No.	Description	From*	Unit
35108	Graduated cylinder, tall form, A, PMP, 10ml:0,2ml, DE-M, with batch certificate	2	piece(s)
35120	Graduated cylinder, tall form, A, PMP, 25ml:0,5ml, DE-M, with batch certificate	2	piece(s)
35128	Graduated cylinder, tall form, A, PMP, 50ml:1ml, DE-M, with batch certificate	2	piece(s)
35138	Graduated cylinder, tall form, A, PMP, 100ml:1ml, DE-M, with batch certificate	2	piece(s)
35148	Graduated cylinder, tall form, A, PMP, 250ml:2ml, DE-M, with batch certificate	2	piece(s)
35154	Graduated cylinder, tall form, A, PMP, 500ml:5ml, DE-M, with batch certificate	2	piece(s)
35162	Graduated cylinder, tall form, A, PMP, 1000ml:10ml, DE-M, with batch certificate	1	piece(s)
35164	Graduated cylinder, tall form, A, PMP, 2000ml:20ml, DE-M, with batch certificate	1	piece(s)


Graduated cylinder, tall form, PMP, class B, blue graduation

Cat. No.	Description	From*	Unit
34708	Graduated cylinder, tall form, 10 ml: 0,2 ml PMP (TPX), graduated in blue	10	piece(s)
34720	Graduated cylinder, tall form, 25 ml: 0,5 ml PMP (TPX), graduated in blue	10	piece(s)
34728	Graduated cylinder, tall form, 50 ml: 1 ml PMP (TPX), graduated in blue	10	piece(s)
34738	Graduated cylinder, tall form, 100 ml: 1 ml PMP (TPX), graduated in blue	10	piece(s)
34748	Graduated cylinder, tall form, 250 ml: 2 ml PMP (TPX), graduated in blue	5	piece(s)
34754	Graduated cylinder, tall form, 500 ml: 5 ml PMP (TPX), graduated in blue	5	piece(s)
34762	Graduated cylinder, tall form, 1000 ml:10 ml PMP (TPX), graduated in blue	5	piece(s)
34764	Graduated cylinder, tall form, 2000 ml:20 ml PMP (TPX), graduated in blue	1	piece(s)



Graduated cylinder, tall form, PMP, class B, embossed scale



Cat. No.	Description	From*	Unit
34908	Graduated cylinder, tall form, 10 ml: 0,2 ml PMP (TPX), embossed scale	10	piece(s)
34920	Graduated cylinder, tall form, 25 ml: 0,5 ml PMP (TPX), embossed scale	10	piece(s)
34928	Graduated cylinder, tall form, 50 ml: 1 ml PMP (TPX), embossed scale	10	piece(s)
34938	Graduated cylinder, tall form, 100 ml: 1 ml PMP (TPX), embossed scale	10	piece(s)
34948	Graduated cylinder, tall form, 250 ml: 2 ml PMP (TPX), embossed scale	5	piece(s)
34954	Graduated cylinder, tall form, 500 ml: 5 ml PMP (TPX), embossed scale	5	piece(s)
34962	Graduated cylinder, tall form, 1000 ml:10 ml PMP (TPX), embossed scale	5	piece(s)
34964	Graduated cylinder, tall form, 2000 ml:20 ml PMP (TPX), embossed scale	1	piece(s)

Graduated cylinder, tall form, class B, PP, blue graduation



Cat. No.	Description	From*	Unit
34808	Graduated cylinder, tall form, 10 ml: 0,2 ml PP, graduated in blue	10	piece(s)
34820	Graduated cylinder, tall form, 25 ml: 0,5 ml PP, graduated in blue	10	piece(s)
34828	Graduated cylinder, tall form, 50 ml: 1 ml PP, graduated in blue	10	piece(s)
34838	Graduated cylinder, tall form, 100 ml: 1 ml PP, graduated in blue	10	piece(s)
34848	Graduated cylinder, tall form, 250 ml: 2 ml PP, graduated in blue	5	piece(s)
34854	Graduated cylinder, tall form, 500 ml: 5 ml PP, graduated in blue	5	piece(s)
34862	Graduated cylinder, tall form, 1000 ml:10 ml PP, graduated in blue	5	piece(s)
34864	Graduated cylinder, tall form, 2000 ml:20 ml PP, graduated in blue	1	piece(s)

Graduated cylinder, tall form, class B, PP, embossed scale



Cat. No.	Description	From*	Unit
35008	Graduated cylinder, tall form, 10 ml: 0,2 ml PP, embossed scale	10	piece(s)
35020	Graduated cylinder, tall form, 25 ml: 0,5 ml PP, embossed scale	10	piece(s)
35028	Graduated cylinder, tall form, 50 ml: 1 ml PP, embossed scale	10	piece(s)
35038	Graduated cylinder, tall form, 100 ml: 1 ml PP, embossed scale	10	piece(s)
35048	Graduated cylinder, tall form, 250 ml: 2 ml PP, embossed scale	5	piece(s)
35054	Graduated cylinder, tall form, 500 ml: 5 ml PP, embossed scale	5	piece(s)
35062	Graduated cylinder, tall form, 1000 ml:10 ml PP, embossed scale	5	piece(s)
35064	Graduated cylinder, tall form, 2000 ml:20 ml PP, embossed scale	1	piece(s)

Graduated cylinder, short form, PP

Cat. No.	Description	From*	Unit
41620	Graduated cylinder, low form embossed scale 25 ml: 0,5 ml, PP	10	piece(s)
41628	Graduated cylinder, low form embossed scale 50 ml: 1,0 ml, PP	10	piece(s)
41638	Graduated cylinder, low form embossed scale 100 ml: 2,0 ml, PP	10	piece(s)
41648	Graduated cylinder, low form embossed scale 250 ml: 5,0 ml, PP	5	piece(s)
41654	Graduated cylinder, low form embossed scale 500 ml: 10,0 ml, PP	5	piece(s)
41662	Graduated cylinder, low form embossed scale 1000 ml: 20,0 ml, PP	5	piece(s)



Dispenser, tilt measure with bottle

Cat. No.	Description	From*	Unit
43005	Automatic pipette, dispenser of Boro 3.3, 5 ml reserve bottle of soda-lime glass	1	piece(s)
43008	Automatic pipette, dispenser of Boro 3.3, 10 ml, reserve bottle of soda-lime glass	1	piece(s)
43016	Automatic pipette, dispenser of Boro 3.3, 20 ml, reserve bottle of soda-lime glass	1	piece(s)
43020	Automatic pipette, dispenser of Boro 3.3, 25 ml, reserve bottle of soda-lime glass	1	piece(s)
43028	Automatic pipette, dispenser of Boro 3.3, 50 ml, reserve bottle of soda-lime glass	1	piece(s)
43038	Automatic pipette, dispenser of Boro 3.3, 100 ml, reserve bottle of soda-lime glass	1	piece(s)



Dispenser, tilt measure without bottle

Cat. No.	Description	From*	Unit
43055	Automatic pipette, dosing attachment, Boro 3.3, w/o bottle, 5 ml.	1	piece(s)
43058	Automatic pipette, dosing attachment, Boro 3.3, w/o bottle, 10 ml.	1	piece(s)
43066	Automatic pipette, dosing attachment, Boro 3.3, w/o bottle, 20 ml.	1	piece(s)
43070	Automatic pipette, dosing attachment, Boro 3.3, w/o bottle, 25 ml.	1	piece(s)
43078	Automatic pipette, dosing attachment, Boro 3.3, w/o bottle, 50 ml.	1	piece(s)
43088	Automatic pipette, dosing attachment, Boro 3.3, w/o bottle, 100 ml.	1	piece(s)

Reservoir bottle

Cat. No.	Description	From*	Unit
126963	Reservoir bottle, soda-lime glass, clear 1000 ml, w/o dosing head, NS 29/32 DIN 12 463	1	piece(s)

Mixing cylinder, tall form, BLAUBRAND®, Boro 3.3, class A, DE-M


Cat. No.	Description	From*	Unit
32408	Mixing cylinder, BLAUBRAND®, A, DE-M, 10 ml: 0,2 ml Boro 3.3 NS10/19, PP stopper	2	piece(s)
32420	Mixing cylinder, BLAUBRAND®, A, DE-M, 25 ml: 0,5 ml Boro 3.3 NS14/23, PP stopper	2	piece(s)
32428	Mixing cylinder, BLAUBRAND®, A, DE-M, 50 ml: 1 ml Boro 3.3 NS19/26, PP stopper	2	piece(s)
32438	Mixing cylinder, BLAUBRAND®, A, DE-M, 100 ml: 1 ml Boro 3.3 NS24/29, PP stopper	2	piece(s)
32448	Mixing cylinder, BLAUBRAND®, A, DE-M, 250 ml: 2 ml Boro 3.3 NS29/32, PP stopper	2	piece(s)
32454	Mixing cylinder, BLAUBRAND®, A, DE-M, 500 ml: 5 ml Boro 3.3 NS34/35, PE stopper	2	piece(s)
32462	Mixing cylinder, BLAUBRAND®, A, DE-M, 1000 ml: 10 ml Boro 3.3 NS45/40, PE stopper	1	piece(s)

Mixing cylinder, tall form, SILBERBRAND ETERNA, Boro 3.3, class B


Cat. No.	Description	From*	Unit
33908	Mixing cylinder, SILBERBRAND ETERNA, B, 10 ml: 0,2 ml, Boro 3.3, NS10/19, PP stopper	2	piece(s)
33920	Mixing cylinder, SILBERBRAND ETERNA, B, 25 ml: 0,5 ml, Boro 3.3, NS14/23, PP stopper	2	piece(s)
33928	Mixing cylinder, SILBERBRAND ETERNA, B, 50 ml: 1 ml, Boro 3.3, NS19/26, PP stopper	2	piece(s)
33938	Mixing cylinder, SILBERBRAND ETERNA, B, 100 ml: 1 ml, Boro 3.3, NS24/29, PP stopper	2	piece(s)
33948	Mixing cylinder, SILBERBRAND ETERNA, B, 250 ml: 2 ml, Boro 3.3, NS29/32, PP stopper	2	piece(s)
33954	Mixing cylinder, SILBERBRAND ETERNA, B, 500 ml: 5 ml, Boro 3.3, NS34/35, PE stopper	2	piece(s)
33962	Mixing cylinder, SILBERBRAND ETERNA, B 1000 ml: 10 ml, Boro 3.3, NS45/40, PE stopper	1	piece(s)

Burettes, BLAUBRAND®, class AS, Boro 3.3, DE-M, lateral stopcock


Cat. No.	Description	From*	Unit
13844	Burette, BLAUBRAND®, AS, DE-M, 10:0,02 ml, Boro 3.3, lateral stopcock, PTFE key, blue Schellbach stripe	2	piece(s)
13847	Burette, BLAUBRAND®, AS, DE-M, 25:0,05 ml, Boro 3.3, lateral stopcock, PTFE key, blue Schellbach stripe	2	piece(s)
13848	Burette, BLAUBRAND®, AS, DE-M, 50:0,1 ml, Boro 3.3, lateral stopcock, PTFE key, blue Schellbach stripe	2	piece(s)

Burettes, BLAUBRAND®, class AS, Boro 3.3, DE-M, straight stopcock

Cat. No.	Description	From*	Unit
12484	Burette, BLAUBRAND®, AS, DE-M, 10:0,02 ml, Boro 3.3, straight stopcock, PTFE spindle, blue Schellbach stripe	2	piece(s)
12487	Burette, BLAUBRAND®, AS, DE-M, 25:0,05 ml, Boro 3.3, straight stopcock, PTFE spindle, blue Schellbach stripe	2	piece(s)
12488	Burette, BLAUBRAND®, AS, DE-M, 50:0,1 ml, Boro 3.3, straight stopcock, PTFE spindle, blue Schellbach stripe	2	piece(s)

**Micro burette Bang pattern, BLAUBRAND®, class AS, Boro 3.3, DE-M, lateral stopcock**

Cat. No.	Description	From*	Unit
24595	Microburette, BLAUBRAND®, AS, DE-M, 2:0,01 ml, Boro 3.3, lateral safety stopcock, blue Schellbach stripe	2	piece(s)
24597	Microburette, BLAUBRAND®, AS, DE-M, 5:0,01 ml, Boro 3.3, lateral safety stopcock, blue Schellbach stripe	2	piece(s)
24599	Microburette, BLAUBRAND®, AS, DE-M, 10:0,02 ml, Boro 3.3, lateral safety stopcock, blue Schellbach stripe	2	piece(s)

**Micro burette Bang pattern, BLAUBRAND®, class AS, Boro 3.3, DE-M, straight stopcock**

Cat. No.	Description	From*	Unit
24265	Microburette, BLAUBRAND®, AS, DE-M, 2:0,01 ml, Boro 3.3, straight safety stopcock, blue Schellbach stripe	2	piece(s)
24267	Microburette, BLAUBRAND®, AS, DE-M, 5:0,01 ml, Boro 3.3, straight safety stopcock, blue Schellbach stripe	2	piece(s)
24269	Microburette, BLAUBRAND®, AS, DE-M, 10:0,02 ml, Boro 3.3, straight safety stopcock, blue Schellbach stripe	2	piece(s)

**Burettes, SILBERBRAND, class B, Boro 3.3, lateral stopcock**

Cat. No.	Description	From*	Unit
13563	Burette, SILBERBRAND, B, lateral stopcock, 10 ml:0,02 ml, Boro 3.3, black Schellbach stripe, PTFE spindle, black Schellbach stripe	2	piece(s)
13566	Burette, SILBERBRAND, B, lateral stopcock, 25 ml:0,05 ml, Boro 3.3, black Schellbach stripe, PTFE spindle, black Schellbach stripe	2	piece(s)
13568	Burette, SILBERBRAND, B, lateral stopcock, 50 ml:0,1 ml, Boro 3.3, black Schellbach stripes, PTFE spindle, black Schellbach stripe	2	piece(s)



Burettes, SILBERBRAND, class B, Boro 3.3, without stopcock, 'Economy'



Cat. No.	Description	From*	Unit
10002	Burette length for burette 'Economy', SILBERBRAND, w/o stopcock, 25:0,1 ml, black graduations	2	piece(s)
10004	Burette length for burette 'Economy', SILBERBRAND, w/o stopcock, 25:0,1 ml, black graduations	2	piece(s)

Burettes, SILBERBRAND, class B, Boro 5.4, straight stopcock, amber


Cat. No.	Description	From*	Unit
13536	Burette, SILBERBRAND, B, straight stopcock, 25 ml:0,1, ml, amber glass, Boro 5.4, PTFE key	2	piece(s)
13538	Burette, SILBERBRAND, B, straight stopcock, 50 ml:0,1, ml, amber glass, Boro 5.4, PTFE key	2	piece(s)

Compact-Burettes, BLAUBRAND®, class AS, Boro 3.3, DE-M


Cat. No.	Description	From*	Unit
13913	Compact burette, BLAUBRAND®, cl.AS, DE-M, 10:0,02ml, Boro 3.3,PTFE stopcock, blue Schellbach stripe	1	piece(s)
13916	Compact burette, BLAUBRAND®, cl.AS, DE-M, 25:0,05ml, Boro 3.3,PTFE stopcock, blue Schellbach stripe	1	piece(s)
13918	Compact burette, BLAUBRAND®, cl.AS, DE-M, 50:0,1ml, Boro 3.3,PTFE stopcock, blue Schellbach stripe	1	piece(s)

Compact-Burettes, BLAUBRAND®, class AS, Boro 3.3, DE-M, amber


Cat. No.	Description	From*	Unit
13963	Compact burette, BLAUBRAND®, cl.AS, DE-M, 10:0,02ml, Boro 3.3,PTFE stopcock, Schellbach	1	piece(s)
13966	Compact burette, BLAUBRAND®, cl.AS, DE-M, 25:0,05ml, Boro 3.3,PTFE stopcock, Schellbach	1	piece(s)
13968	Compact burette, BLAUBRAND®, cl.AS, DE-M, 50:0,1ml, Boro 3.3,PTFE stopcock, Schellbach	1	piece(s)

Compact-Burettes, SILBERBRAND, class B, Boro 3.3



Cat. No.	Description	From*	Unit
13903	Compact burette, SILBERBRAND, 10:0,05ml, Boro 3.3,PTFE stopcock, black Schellbach stripe	1	piece(s)
13906	Compact burette, SILBERBRAND, 25:0,1ml, Boro 3.3, PTFE stopcock, black Schellbach stripe	1	piece(s)
13908	Compact burette, SILBERBRAND, 50:0,1ml, Boro 3.3, PTFE stopcock, black Schellbach stripe	1	piece(s)

Compact-Burettes, SILBERBRAND, class B, Boro 3.3, 'Economy'

Cat. No.	Description	From*	Unit
10012	Compact-Burette 'Economy', SILBERBRAND, 25 ml, with PTFE stopcock and PP tip	2	piece(s)
10014	Compact-Burette 'Economy', SILBERBRAND, 50 ml, with PTFE stopcock and PP tip	2	piece(s)


Compact-Burettes, SILBERBRAND, class B, Boro 5.4, amber

Cat. No.	Description	From*	Unit
13926	Compact burette, SILBERBRAND, 25:0,1 ml, amber, PTFE stopcock	1	piece(s)
13928	Compact burette, SILBERBRAND, 50:0,1 ml, amber, PTFE stopcock	1	piece(s)


Spare burette tube, for compact burette

Cat. No.	Description	From*	Unit
13933	Spare burette length for compact burette 10 ml, SILBERBRAND, Boro 3.3	1	piece(s)
13936	Spare burette length for compact burette 25 ml, SILBERBRAND, Boro 3.3	1	piece(s)
13938	Spare burette length for compact burette 50 ml, SILBERBRAND, Boro 3.3	1	piece(s)
13943	Spare burette length for compact burette 10 ml, BLAUBRAND®, Boro 3.3	1	piece(s)
13946	Spare burette length for compact burette 25 ml, BLAUBRAND®, Boro 3.3	1	piece(s)
13948	Spare burette length for compact burette 50 ml, BLAUBRAND®, Boro 3.3	1	piece(s)
13956	Spare burette length for compact burette 25 ml, SILBERBRAND, Boro 5.4, amber glas	1	piece(s)
13958	Spare burette length for compact burette 50 ml, SILBERBRAND, Boro 5.4, amber glas	1	piece(s)


Burette clamps

Cat. No.	Description	From*	Unit
16515	Burette clamp, Zinc alloy, nickel-plated for 1 burette, clamp PVC-coated	1	piece(s)
16520	Burette clamp, Zinc alloy, nickel-plated for 2 burettes, clamp PVC-coated	1	piece(s)
578000	Burette clamp, aluminum die casting for 2 burettes, with roll mounting	1	piece(s)


Burette holder, PP, white

Cat. No.	Description	From*	Unit
16505	Burette clamp, PP, white, for 1 burette	1	piece(s)
16510	Burette clamp, PP, white, for 2 burettes	1	piece(s)



Burette cover, PP

Cat. No.	Description	From*	Unit
16400	Burette cover, PP, inner dia. 20 mm, height 30 mm	10	piece(s)

Burette clamp for micro burettes, Bang pattern

Cat. No.	Description	From*	Unit
915	Burette clamp for micro burettes, Bang pattern	1	piece(s)

Needle-valve stopcock for burettes, Boro 3.3

Cat. No.	Description	From*	Unit
82120	Needle-valve stopcock/burette lateral stopcock for burette cap. 2-10 ml, valve open 0-2,5 mm	1	piece(s)
82121	Needle-valve stopcock/burette lateral stopcock for burette cap. 25 ml, valve open 0-2,5 mm	1	piece(s)
82122	Needle-valve stopcock/burette lateral stopcock for burette cap. 50 ml, valve open 0-2,5 mm	1	piece(s)
82220	Straight needle-valve stopcock for Burette 2-10ml, Boro 3.3, 0-2,5mm	1	piece(s)
82221	Straight needle-valve stopcock for Burette 25ml, Boro 3.3, 0-2,5mm	1	piece(s)
82222	Straight needle-valve stopcock for Burette 50ml, Boro 3.3, 0-2,5 mm	1	piece(s)

Spare valve spindle, PTFE, with retention device

Cat. No.	Description	From*	Unit
82170	Needle-valve stopcock/spare spindle reten. dev./snap ring for valve open. 0-2,5 mm	1	piece(s)

Burette stopcock, Boro 3.3, straight stopcock

Cat. No.	Description	From*	Unit
81805	Burette stopcock, lateral stopcock for burette cap. 2-10 ml, glass key, Boro 3.3	1	piece(s)
81807	Burette stopcock, lateral stopcock for burette cap. 25 ml, glass key, Boro 3.3	1	piece(s)
81809	Burette stopcock, lateral stopcock for burette cap. 50 ml, glass key, Boro 3.3	1	piece(s)
81815	Burette stopcock, straight stopcock for burette cap. 2-10 ml, glass key, Boro 3.3	1	piece(s)
81817	Burette stopcock, straight stopcock for burette cap. 25 ml, glass key, Boro 3.3	1	piece(s)
81819	Burette stopcock, straight stopcock for burette cap. 50 ml, glass key, Boro 3.3	1	piece(s)

Burette stopcock, Boro, brown

Cat. No.	Description	From*	Unit
81146	Lateral burette stopcock/spare key for burette cap. 25 ml, securing device, Boro 3.3, amber	1	piece(s)
81147	Lateral burette stopcock/spare key for burette cap. 50 ml, securing device, Boro 3.3, amber	1	piece(s)
81907	Burette stopcock, lateral stopcock for burette cap. 25 ml, glass key, amber	1	piece(s)
81909	Burette stopcock, lateral stopcock for burette cap. 50 ml, glass key, amber	1	piece(s)
81917	Burette stopcock, straight stopcock for burette cap. 25 ml, glass key, amber	1	piece(s)
81919	Burette stopcock, straight stopcock for burette cap. 50 ml, glass key, amber	1	piece(s)
81927	Burette stopcock, lateral stopcock for burette cap. 25+50 ml, PTFE key, amber	1	piece(s)
81937	Burette stopcock, straight stopcock for burette cap. 25 ml, PTFE key, amber	1	piece(s)
81939	Burette stopcock, straight stopcock for burette cap. 50 ml, PTFE key, amber	1	piece(s)



Burette stopcock, PTFE

Cat. No.	Description	From*	Unit
11800	PTFE-stopcock, for Burette length, with PP-tip, complete	1	piece(s)
11805	Stopcock for compact/automatic burettes PTFE, w/o burette tip with screwings and sealings	1	piece(s)



Spare key with retention

Cat. No.	Description	From*	Unit
81053	Straight-bore stopcock/spare key, solid, bore 2,5 mm stopc.3 NS/12 securing dev.	1	piece(s)
81140	Lateral burette stopcock/spare key for burette cap., 2-10 ml, securing device, Boro 3.3	1	piece(s)
81141	Lateral burette stopcock/spare key for burette cap. 25 ml, securing device, Boro 3.3	1	piece(s)
81142	Lateral burette stopcock/spare key for burette cap. 50 ml, securing device, Boro 3.3	1	piece(s)
81167	Lateral burette stopcock/spare key for burette cap. 25+50 ml, securing device, PTFE	1	piece(s)
81247	Straight burette stopcock/spare key for burette cap. 25+50 ml, securing device, Boro 3.3, amber	1	piece(s)
81248	Straight burette stopcock/spare key for burette cap. 25+50 ml, securing device, PTFE	1	piece(s)



Discharge tip for compact burettes 'Economy'

Cat. No.	Description	From*	Unit
11600	Burette tip, for burette length, PP, with screw cap, for PTFE stopcock	1	piece(s)



Spare burette tip, for compact burette and compact automatic burette, Boro 3.3



Cat. No.	Description	From*	Unit
11510	Spare burette tip for compact- and automatic burette, 10 ml, Boro 3.3, clear glass	1	piece(s)
11511	Spare burette tip for compact- and automatic burette, 25 ml, Boro 3.3, clear glass	1	piece(s)
11512	Spare burette tip for compact- and automatic burette, 50 ml, Boro 3.3, clear glass	1	piece(s)
11515	Spare burette tip for compact- and automatic burette 25 + 50 ml, Boro 5.4, amber glass	1	piece(s)

Automatic burette Pellet pattern, BLAUBRAND®, class AS, DE-M, with intermediate stopcock



Cat. No.	Description	From*	Unit
22761	Auto. burette, BLAUBRAND®, AS, DE-M, 10:0,02 ml, i-stopcock, lateral stopcock, PTFE stopcock, blue Schellbach stripe	1	piece(s)
22762	Auto. burette, BLAUBRAND®, AS, DE-M, 25:0,05 ml, i-stopcock, lateral stopcock, PTFE stopcock, blue Schellbach stripe	1	piece(s)
22763	Auto. burette, BLAUBRAND®, AS, DE-M, 50:0,1 ml, i-stopcock, lateral stopcock, PTFE stopcock, blue Schellbach stripe	1	piece(s)
22764	Auto. burette, BLAUBRAND®, AS, DE-M, 10:0,02 ml, i-stopcock, lateral stopcock, PTFE stopcock, bottle, blue Schellbach stripe	1	piece(s)
22767	Auto. burette, BLAUBRAND®, AS, DE-M, 25:0,05 ml, i-stopcock, lateral stopcock, PTFE stopcock, bottle, blue Schellbach stripe	1	piece(s)
22768	Auto. burette, BLAUBRAND®, AS, DE-M, 50:0,1 ml, i-stopcock, lateral stopcock, PTFE stopcock, bottle, blue Schellbach stripe	1	piece(s)

Automatic burette Pellet pattern, BLAUBRAND®, class AS, DE-M, without intermediate stopcock



Cat. No.	Description	From*	Unit
22521	Auto. burette, BLAUBRAND®, AS, DE-M, 10:0,02 ml, lateral safety stopcock, blue Schellbach stripe	1	piece(s)
22522	Auto. burette, BLAUBRAND®, AS, DE-M, 25:0,05 ml, lateral safety stopcock, blue Schellbach stripe	1	piece(s)
22523	Auto. burette, BLAUBRAND®, AS, DE-M, 50:0,1 ml, lateral safety stopcock, blue Schellbach stripe	1	piece(s)
22524	Auto. burette, BLAUBRAND®, AS, DE-M, 10:0,02 ml, lateral safety stopcock, w.bottle, blue Schellbach stripe	1	piece(s)
22528	Auto. burette, BLAUBRAND®, AS, DE-M, 25:0,05 ml, lateral safety stopcock, w.bottle, blue Schellbach stripe	1	piece(s)
22530	Auto. burette, BLAUBRAND®, AS, DE-M, 50:0,1 ml, lateral safety stopcock, w.bottle, blue Schellbach stripe	1	piece(s)

Automatic burette Pellet pattern, SILBERBRAND, class B, without intermediate stopcock

Cat. No.	Description	From*	Unit
21911	Auto. burette, SILBERBRAND, 10:0,02 ml, lateral safety stopcock, black Schellbach stripe	1	piece(s)
21912	Auto. burette, SILBERBRAND, 25:0,05 ml, lateral safety stopcock, black Schellbach stripe	1	piece(s)
21913	Auto. burette, SILBERBRAND, 50:0,1 ml, lateral safety stopcock, black Schellbach stripe	1	piece(s)
21914	Auto. burette, SILBERBRAND, 10:0,02 ml, lateral safety stopcock, w.bottle, black Schellbach stripe	1	piece(s)
21917	Auto. burette, SILBERBRAND, 25:0,05 ml, lateral safety stopcock, w.bottle, black Schellbach stripe	1	piece(s)
21918	Auto. burette, SILBERBRAND, 50:0,1 ml, lateral safety stopcock, w. bottle, black Schellbach stripe	1	piece(s)


Automatic burette Pellet pattern, SILBERBRAND, class B, with intermediate stopcock

Cat. No.	Description	From*	Unit
22332	Auto. burette, SILBERBRAND, 25:0,1 ml, i-stopcock, lateral stopcock PTFE amber, bottle	1	piece(s)
22334	Auto. burette, SILBERBRAND, 50:0,1 ml, i-stopcock, lateral stopcock PTFE amber, bottle	1	piece(s)
22336	Auto. burette, SILBERBRAND, 25:0,1 ml, i-stopcock, lateral stopcock PTFE k. amber	1	piece(s)
22338	Auto. burette, SILBERBRAND, 50:0,1 ml, i-stopcock, lateral stopcock PTFE k. amber	1	piece(s)


Automatic burette Pellet pattern, SILBERBRAND, class B, without intermediate stopcock

Cat. No.	Description	From*	Unit
22312	Auto. burette, SILBERBRAND, 25:0,1 ml, lateral stopcock PTFE amber, w.bottle	1	piece(s)
22314	Auto. burette, SILBERBRAND, 50:0,1 ml, lateral stopcock PTFE amber, w.bottle	1	piece(s)
22316	Auto. burette, SILBERBRAND, 25:0,1 ml, amber lateral stopcock PTFE key	1	piece(s)
22318	Auto. burette, SILBERBRAND, 50:0,1 ml, amber lateral stopcock PTFE key	1	piece(s)


Compact-Automatic burette, BLAUBRAND®, class AS, DE-M, with glass bottle

Cat. No.	Description	From*	Unit
23919	Compact auto. burette, BLAUBRAND®, AS, DE-M, 10:0.02 ml, 2000ml glass-bottle, blue Schellbach stripe	1	piece(s)
23920	Compact auto. burette, BLAUBRAND®, AS, DE-M, 25:0.05 ml, 2000 ml glass bottle, blue Schellbach stripe	1	piece(s)
23921	Compact auto. burette, BLAUBRAND®, AS, DE-M, 50:0.1 ml, 2000 ml glass bottle, blue Schellbach stripe	1	piece(s)



Compact-Automatic burette, BLAUBRAND®, class AS, DE-M, with PE bottle

Cat. No.	Description	From*	Unit
23819	Compact auto. burette, BLAUBRAND®, DE-M, Schellbach, Boro 3.3, w. PE-LD bottle 1000ml, blue Schellbach stripe	1	piece(s)
23820	Compact auto. burette, BLAUBRAND®, DE-M, Schellbach, Boro 3.3, w. PE-LD bottle 1000ml, blue Schellbach stripe	1	piece(s)
23821	Compact auto. burette, BLAUBRAND®, DE-M, Schellbach, Boro 3.3, w. PE-LD bottle 1000ml, blue Schellbach stripe	1	piece(s)

Compact-Automatic burette, BLAUBRAND®, amber, class AS, DE-M, with glass bottle


Cat. No.	Description	From*	Unit
23939	Compact auto. burette, BLAUBRAND®, AS, DE-M, 10:0.02 ml, 2000ml glass-bottle	1	piece(s)
23940	Compact auto. burette, BLAUBRAND®, AS, DE-M, 25:0.05 ml, 2000 ml glass bottle	1	piece(s)
23941	Compact auto. burette, BLAUBRAND®, AS, DE-M, 50:0.1 ml, 2000 ml glass bottle	1	piece(s)

Compact-Automatic burette, BLAUBRAND®, amber, class AS, DE-M, with PE bottle

Cat. No.	Description	From*	Unit
23870	Compact auto. burette, BLAUBRAND®, DE-M, Schellbach, Boro 3.3, w. PE-LD bottle 1000ml	1	piece(s)
23871	Compact auto. burette, BLAUBRAND®, DE-M, Schellbach, Boro 3.3, w. PE-LD bottle 1000ml	1	piece(s)
23872	Compact auto. burette, BLAUBRAND®, DE-M, Schellbach, Boro 3.3, w. PE-LD bottle 1000ml	1	piece(s)

Compact-Automatic burette, SILBERBRAND, class B, with PE bottle

Cat. No.	Description	From*	Unit
23809	Compact auto. burette, SILBERBRAND, Schellbach, Boro 3.3, w. PE-LD bottle 1000 ml, black Schellbach stripe	1	piece(s)
23810	Compact auto. burette, SILBERBRAND, Schellbach, Boro 3.3, w. PE-LD bottle 1000 ml, black Schellbach stripe	1	piece(s)
23811	Compact auto. burette, SILBERBRAND, Schellbach, Boro 3.3, w. PE-LD bottle 1000 ml, black Schellbach stripe	1	piece(s)

Compact-Automatic burette, SILBERBRAND, amber, class B, with PE bottle

Cat. No.	Description	From*	Unit
23829	Compact auto. burette, SILBERBRAND, amber, 25 ml, Boro 5.4, w. PE-LD bottle 1000 ml	1	piece(s)
23830	Compact auto. burette, SILBERBRAND, amber, 50 ml, Boro 5.4, w. PE-LD bottle 1000 ml	1	piece(s)

Automatic burette Schilling pattern, SILBERBRAND, class B, with reservoir bottle

Cat. No.	Description	From*	Unit
23753	Auto. burette Schilling SILBERBRAND 10 ml, Schellbach, Boro 3.3, w. PE-LD bottle 500 ml, black Schellbach stripe	1	piece(s)
23755	Auto. burette Schilling SILBERBRAND 15 ml, Schellbach, Boro 3.3, w. PE-LD bottle 500 ml, black Schellbach stripe	1	piece(s)
23756	Auto. burette Schilling SILBERBRAND 25 ml, Schellbach, Boro 3.3, w. PE-LD bottle 1000 ml, black Schellbach stripe	1	piece(s)
23758	Auto. burette Schilling SILBERBRAND 50 ml, Schellbach, Boro 3.3, w. PE-LD bottle 1000 ml, black Schellbach stripe	1	piece(s)



Automatic burette Schilling pattern, SILBERBRAND, amber, class B, with reservoir bottle (brown)

Cat. No.	Description	From*	Unit
23766	Auto. burette Schilling SILBERBRAND 25 ml, Schellbach, AR-GLAS®, w. PE-LD bottle 1000 ml	1	piece(s)
23768	Auto. burette Schilling SILBERBRAND 50 ml, Schellbach, AR-GLAS®, w. PE-LD bottle 1000 ml	1	piece(s)



Spare burette for automatic burette Dr. Schilling pattern, SILBERBRAND

Cat. No.	Description	From*	Unit
23713	Spare burette f. aut.burette Schilling 10 ml, SILBERBRAND, Schellbach, Boro 3.3, length 330 mm	1	piece(s)
23715	Spare burette f. aut.burette Schilling 15 ml, SILBERBRAND, Schellbach, Boro 3.3, length 330 mm	1	piece(s)
23716	Spare burette f. aut.burette Schilling 25 ml, SILBERBRAND, Schellbach, Boro 3.3, length 380 mm	1	piece(s)
23718	Spare burette f. aut.burette Schilling 50 ml, SILBERBRAND, Schellbach, Boro 3.3, length 580 mm	1	piece(s)
23723	Spare burette f. aut.burette Schilling 25 ml SILBERBRAND Schellb., AR-GLAS®, amber glass, l. 390 mm	1	piece(s)
23724	Spare burette f. aut.burette Schilling 50 ml SILBERBRAND Schellb., AR-GLAS®, amber glass, l. 580 mm	1	piece(s)




Bottles for automatic burettes

Cat. No.	Description	From*	Unit
23310	Spare reservoir bottle for aut. burettes 2000 ml, clear glass, Boro 3.3, NS 29/32	1	piece(s)
23320	Spare reservoir bottle for aut. burettes 2000 ml, amber glass, Boro 3.3, NS 29/32	1	piece(s)
126965	Washbottle, BISTABIL, soda-lime glass, clear, 2000 ml, w/o head, NS 29/32	1	piece(s)
127065	Washbottle, BISTABIL, soda-lime glass, amber, 2000 ml, w/o head, NS 29/32	1	piece(s)
129055	Bottle for Schilling burettes, PE-LD, 500 ml	1	piece(s)
129060	Bottle for Schilling burettes, PE-LD, 1000 ml	1	piece(s)
130260	Bottle for Schilling burette, PE-LD, amber, 1000 ml, with screw cap and tube bushing	1	piece(s)



Rubber bellow for automatic burettes



Cat. No.	Description	From*	Unit
23400	Rubber bellows for auto. burettes single bulb bellow with valve/tubing	1	piece(s)

Joint clip for dispenser, tilt measure



Cat. No.	Description	From*	Unit
55618	Joint clip for dispensers, tilt measure Stainless steel, neck size 29/32	1	piece(s)

Spare filling tube for automatic burettes

Cat. No.	Description	From*	Unit
11525	Spare silicone tubing for automatic burettes and compact automatic burettes, PVC, transparent	1	piece(s)

Spare key for automatic burettes

Cat. No.	Description	From*	Unit
81255	Spare key for automatic burettes Pellet pattern, glass key, Boro 5.4	1	piece(s)
81256	Spare key for automatic burettes Pellet pattern, glass key, Boro 5.4, amber glass	1	piece(s)
81265	Spare key for automatic burettes Pellet pattern, PTFE key	1	piece(s)

Spare burette tube for compact automatic burettes



Cat. No.	Description	From*	Unit
23833	Spare burette length for compact autom. burette SILBERBRAND, 10 ml, Boro 3.3	1	piece(s)
23836	Spare burette length for compact autom. burette SILBERBRAND, 25 ml, Boro 3.3	1	piece(s)
23838	Spare burette length for compact autom. burette SILBERBRAND, 50 ml, Boro 3.3	1	piece(s)
23843	Spare burette length for compact autom. burette BLAUBRAND®, 10 ml, Boro 3.3	1	piece(s)
23846	Spare burette length for compact autom. burette BLAUBRAND®, 25 ml, Boro 3.3	1	piece(s)
23848	Spare burette length for compact autom. burette BLAUBRAND®, 50 ml, Boro 3.3	1	piece(s)
23866	Spare burette length for compact autom. burette SILBERBRAND, 25 ml, Boro 5.4, amber glass	1	piece(s)
23868	Spare burette length for compact autom. burette SILBERBRAND, 50 ml, Boro 5.4, amber glass	1	piece(s)

Base for plastic bottles

Cat. No.	Description	From*	Unit
23725	Base for auto. burette Schilling for plastic bottle 500 ml, diameter 74-76 mm	1	piece(s)
23728	Base for auto. burette Schilling for 1000 ml plastic bottle, diameter 94-96 mm	1	piece(s)



Burette support

Cat. No.	Description	From*	Unit
23882	Support for compact automatic burette rectangular base plate	1	piece(s)



Pumphead for compact automatic burettes, Boro 3.3

Cat. No.	Description	From*	Unit
23875	Pump head for glass reservoir bottle, for compact auto. burette, clear glass	1	piece(s)
23876	Pump head for glass reservoir bottle, for compact auto. burette, amber glass	1	piece(s)



Discharge tip for automatic burettes Dr. Schilling pattern, with silicone tubing

Cat. No.	Description	From*	Unit
11500	Glass discharge tip, clear glass, for automatic burette Schilling, with silicone tubing	10	piece(s)
11505	Glass discharge tip, amber glass, for automatic burette Schilling, with silicone tubing	10	piece(s)



Fitting for automatic burette Dr. Schilling pattern

Cat. No.	Description	From*	Unit
23745	Fitting for aut. burette Dr. Schilling f. 10 ml micrometer screw + push-button incl. support clamp	1	piece(s)
23746	Fitting for aut. burette Dr. Schilling f. 15 ml micrometer screw + push-button incl. support clamp	1	piece(s)
23748	Fitting for aut. burette Dr. Schilling f. 25-50 ml micrometer screw+push button incl. support clamp	1	piece(s)



Sockets, PTFE



Cat. No.	Description	From*	Unit
51467	Ground socket, PTFE NS 29/32, grip collar and sealing rings	1	piece(s)

Ground joint sleeves, PTFE



Cat. No.	Description	From*	Unit
51422	Manguito p. esmerilados PTFE para tamaño, PTFE for NS 29/32, thin-walled	10	piece(s)

Conical joint clips, POM



Cat. No.	Description	From*	Unit
55646	Conical joint clip, POM, for NS 29/32, red	1	piece(s)

Sealing film Seal-R-film™

Cat. No.	Description	From*	Unit
701620	BRAND Seal-R-film™ sealing film, length 38m, width 100mm	12	piece(s)
701622	BRAND Seal-R-film™ sealing film, length 75m, width 50mm	24	piece(s)
701624	BRAND Seal-R-film™ sealing film, length 75m, width 100mm	12	piece(s)
701626	BRAND Seal-R-film™ sealing film, length 15m, width 500mm	6	piece(s)



Cutter, for Sealing film

Cat. No.	Description	From*	Unit
701650	Sealing film cutter for rolls 50/100 mm width	1	piece(s)



Erlenmeyer flask, wide neck, Boro 3.3

Cat. No.	Description	From*	Unit
92814	Erlenmeyer flask, wide neck, Boro 3.3, 25 ml, with beaded rim and graduation	10	piece(s)
92817	Erlenmeyer flask, wide neck, Boro 3.3, 50 ml, with beaded rim and graduation	10	piece(s)
92824	Erlenmeyer flask, wide neck, Boro 3.3, 100 ml, with beaded rim and graduation	10	piece(s)
92836	Erlenmeyer flask, wide neck, Boro 3.3, 250 ml, with beaded rim and graduation	10	piece(s)
92839	Erlenmeyer flask, wide neck, Boro 3.3, 300 ml, with beaded rim and graduation	10	piece(s)
92844	Erlenmeyer flask, wide neck, Boro 3.3, 500 ml, with beaded rim and graduation	10	piece(s)
92854	Erlenmeyer flask, wide neck, Boro 3.3, 1000 ml, with beaded rim and graduation	10	piece(s)
92863	Erlenmeyer flask, wide neck, Boro 3.3, 2000 ml, with beaded rim and graduation	10	piece(s)



Erlenmeyer flask, narrow neck, Boro 3.3

Cat. No.	Description	From*	Unit
92714	Erlenmeyer flask, narrow neck, Boro 3.3, 25 ml, with beaded rim and graduation	10	piece(s)
92717	Erlenmeyer flask, narrow neck, Boro 3.3, 50 ml, with beaded rim and graduation	10	piece(s)
92724	Erlenmeyer flask, narrow neck, Boro 3.3, 100 ml, with beaded rim and graduation	10	piece(s)
92736	Erlenmeyer flask, narrow neck, Boro 3.3, 250 ml, with beaded rim and graduation	10	piece(s)
92739	Erlenmeyer flask, narrow neck, Boro 3.3, 300 ml, with beaded rim and graduation	10	piece(s)
92744	Erlenmeyer flask, narrow neck, Boro 3.3, 500 ml, with beaded rim and graduation	10	piece(s)
92754	Erlenmeyer flask, narrow neck, Boro 3.3, 1000 ml, with beaded rim and graduation	10	piece(s)
92763	Erlenmeyer flask, narrow neck, Boro 3.3, 2000 ml, with beaded rim and graduation	10	piece(s)
92768	Erlenmeyer flask, narrow neck, Boro 3.3, 3000 ml, with beaded rim and graduation	2	piece(s)
92773	Erlenmeyer flask, narrow neck, Boro 3.3, 5000 ml, with beaded rim and graduation	1	piece(s)



Erlenmeyer flask, wide neck, PP



Cat. No.	Description	From*	Unit
93170	Erlenmeyer flask, wide neck, PP, 50 ml, NS 34/35 screw cap/graduated	6	piece(s)
93172	Erlenmeyer flask, wide neck, PP, 100 ml, NS 34/35 screw cap/graduated	6	piece(s)
93175	Erlenmeyer flask, wide neck, PP, 250 ml, NS 45/40 screw cap/graduated	6	piece(s)
93180	Erlenmeyer flask, wide neck, PP, 500 ml, NS 45/40 screw cap/graduated	6	piece(s)
93185	Erlenmeyer flask, wide neck, PP, 1000 ml, NS 45/40 screw cap/graduated	4	piece(s)

Beaker, tall form, Boro 3.3



Cat. No.	Description	From*	Unit
90617	Beaker, tall form, Boro 3.3, 50 ml, with graduation and spout	10	piece(s)
90624	Beaker, tall form, Boro 3.3, 100 ml, with graduation and spout	10	piece(s)
90629	Beaker, tall form, Boro 3.3, 150 ml, with graduation and spout	10	piece(s)
90636	Beaker, tall form, Boro 3.3, 250 ml, with graduation and spout	10	piece(s)
90641	Beaker, tall form, Boro 3.3, 400 ml, with graduation and spout	10	piece(s)
90648	Beaker, tall form, Boro 3.3, 600 ml, with graduation and spout	10	piece(s)
90654	Beaker, tall form, Boro 3.3, 1000 ml, with graduation and spout	10	piece(s)
90663	Beaker, tall form, Boro 3.3, 2000 ml, with graduation and spout	10	piece(s)
90668	Beaker, tall form, Boro 3.3, 3000 ml, with graduation and spout	2	piece(s)

Beaker, low form, Boro 3.3



Cat. No.	Description	From*	Unit
91214	Beaker, low form, Boro 3.3, 25 ml, with graduation and spout	10	piece(s)
91217	Beaker, low form, Boro 3.3, 50 ml, with graduation and spout	10	piece(s)
91224	Beaker, low form, Boro 3.3, 100 ml, with graduation and spout	10	piece(s)
91229	Beaker, low form, Boro 3.3, 150 ml, with graduation and spout	10	piece(s)
91236	Beaker, low form, Boro 3.3, 250 ml, with graduation and spout	10	piece(s)
91241	Beaker, low form, Boro 3.3, 400 ml, with graduation and spout	10	piece(s)
91248	Beaker, low form, Boro 3.3, 600 ml, with graduation and spout	10	piece(s)
91253	Beaker, low form, Boro 3.3, 800 ml, with graduation and spout	10	piece(s)
91254	Beaker, low form, Boro 3.3, 1000 ml, with graduation and spout	10	piece(s)
91263	Beaker, low form, Boro 3.3, 2000 ml, with graduation and spout	10	piece(s)
91268	Beaker, low form, Boro 3.3, 3000 ml, with graduation and spout	4	piece(s)
91273	Beaker, low form, Boro 3.3, 5000 ml, with graduation and spout	1	piece(s)

Beaker, low form, Boro 3.3



Cat. No.	Description	From*	Unit
91210	Beaker, low form, Boro 3.3, 5 ml, w/o graduation, with spout	10	piece(s)
91212	Beaker, low form, Boro 3.3, 10 ml, w/o graduation, with spout	10	piece(s)

Beaker, short form, PMP, blue printed scale

Cat. No.	Description	From*	Unit
87502	Beaker, low form, PMP, 10 ml: 2 ml, blue graduation, beaded rim and spout	12	piece(s)
87504	Beaker, low form, PMP, 25 ml: 5 ml, blue graduation, beaded rim and spout	12	piece(s)
87506	Beaker, low form, PMP, 50 ml: 10 ml, blue graduation, beaded rim and spout	12	piece(s)
87508	Beaker, low form, PMP, 100 ml: 10 ml, blue graduation, beaded rim and spout	12	piece(s)
87510	Beaker, low form, PMP, 150 ml: 20 ml, blue graduation, beaded rim and spout	12	piece(s)
87512	Beaker, low form, PMP, 250 ml: 25 ml, blue graduation, beaded rim and spout	6	piece(s)
87514	Beaker, low form, PMP, 400 ml: 50 ml, blue graduation, beaded rim and spout	6	piece(s)
87516	Beaker, low form, PMP, 500 ml: 50 ml, blue graduation, beaded rim and spout	6	piece(s)
87518	Beaker, low form, PMP, 600 ml: 50 ml, blue graduation, beaded rim and spout	6	piece(s)
87520	Beaker, low form, PMP, 1000 ml:100 ml, blue graduation, beaded rim and spout	6	piece(s)
87522	Beaker, low form, PMP, 2000 ml:200 ml, blue graduation, beaded rim and spout	6	piece(s)
87524	Beaker, low form, PMP, 3000 ml:250 ml, blue graduation, beaded rim and spout	4	piece(s)
87526	Beaker, low form, PMP, 5000 ml:500 ml, blue graduation, beaded rim and spout	4	piece(s)



Beaker, low form, PMP, embossed scale

Cat. No.	Description	From*	Unit
87704	Beaker, low form, PMP, 25 ml: 5 ml, embossed graduation, beaded rim and spout	12	piece(s)
87706	Beaker, low form, PMP, 50 ml: 10 ml, embossed graduation, beaded rim and spout	12	piece(s)
87708	Beaker, low form, PMP, 100 ml: 10 ml, embossed graduation, beaded rim and spout	12	piece(s)
87710	Beaker, low form, PMP, 150 ml: 20 ml, embossed graduation, beaded rim and spout	12	piece(s)
87712	Beaker, low form, PMP, 250 ml: 25 ml, embossed graduation, beaded rim and spout	6	piece(s)
87714	Beaker, low form, PMP, 400 ml: 50 ml, embossed graduation, beaded rim and spout	6	piece(s)
87716	Beaker, low form, PMP, 500 ml: 50 ml, embossed graduation, beaded rim and spout	6	piece(s)
87718	Beaker, low form, PMP, 600 ml: 50 ml, embossed graduation, beaded rim and spout	6	piece(s)
87720	Beaker, low form, PMP, 1000 ml:100 ml, embossed graduation, beaded rim and spout	6	piece(s)
87722	Beaker, low form, PMP, 2000 ml:200 ml, embossed graduation, beaded rim and spout	6	piece(s)
87724	Beaker, low form, PMP, 3000 ml:250 ml, embossed graduation, beaded rim and spout	4	piece(s)
87726	Beaker, low form, PMP, 5000 ml:500 ml, embossed graduation, beaded rim and spout	4	piece(s)



Beaker, low form, PP, blue printed scale



Cat. No.	Description	From*	Unit
87602	Beaker, low form, PP, 10 ml: 2 ml, blue graduation, beaded rim and spout	12	piece(s)
87604	Beaker, low form, PP, 25 ml: 5 ml, blue graduation, beaded rim and spout	12	piece(s)
87606	Beaker, low form, PP, 50 ml: 10 ml, blue graduation, beaded rim and spout	12	piece(s)
87608	Beaker, low form, PP, 100 ml: 10 ml, blue graduation, beaded rim and spout	12	piece(s)
87610	Beaker, low form, PP, 150 ml: 20 ml, blue graduation, beaded rim and spout	12	piece(s)
87612	Beaker, low form, PP, 250 ml: 25 ml, blue graduation, beaded rim and spout	6	piece(s)
87614	Beaker, low form, PP, 400 ml: 50 ml, blue graduation, beaded rim and spout	6	piece(s)
87616	Beaker, low form, PP, 500 ml: 50 ml, blue graduation, beaded rim and spout	6	piece(s)
87618	Beaker, low form, PP, 600 ml: 50 ml, blue graduation, beaded rim and spout	6	piece(s)
87620	Beaker, low form, PP, 1000 ml:100 ml, blue graduation, beaded rim and spout	6	piece(s)
87622	Beaker, low form, PP, 2000 ml:200 ml, blue graduation, beaded rim and spout	6	piece(s)
87624	Beaker, low form, PP, 3000 ml:250 ml, blue graduation, beaded rim and spout	4	piece(s)
87626	Beaker, low form, PP, 5000 ml:500 ml, blue graduation, beaded rim and spout	4	piece(s)

Beaker, low form, PP, embossed scale



Cat. No.	Description	From*	Unit
87804	Beaker, low form, PP, 25 ml: 5 ml, embossed graduation, beaded rim and spout	12	piece(s)
87806	Beaker, low form, PP, 50 ml: 10 ml, embossed graduation, beaded rim and spout	12	piece(s)
87808	Beaker, low form, PP, 100 ml: 10 ml, embossed graduation, beaded rim and spout	12	piece(s)
87810	Beaker, low form, PP, 150 ml: 20 ml, embossed graduation, beaded rim and spout	12	piece(s)
87812	Beaker, low form, PP, 250 ml: 25 ml, embossed graduation, beaded rim and spout	6	piece(s)
87814	Beaker, low form, PP, 400 ml: 50 ml, embossed graduation, beaded rim and spout	6	piece(s)
87816	Beaker, low form, PP, 500 ml: 50 ml, embossed graduation, beaded rim and spout	6	piece(s)
87818	Beaker, low form, PP, 600 ml: 50 ml, embossed graduation, beaded rim and spout	6	piece(s)
87820	Beaker, low form, PP, 1000 ml:100 ml, embossed graduation, beaded rim and spout	6	piece(s)
87822	Beaker, low form, PP, 2000 ml:200 ml, embossed graduation, beaded rim and spout	6	piece(s)
87824	Beaker, low form, PP, 3000 ml:250 ml, embossed graduation, beaded rim and spout	4	piece(s)
87826	Beaker, low form, PP, 5000 ml:500 ml, embossed graduation, beaded rim and spout	4	piece(s)

Beaker, low form, ETFE, blue printed scale



Cat. No.	Description	From*	Unit
90224	Beaker, low form, 25 ml:5 ml, with graduation and spout, ETFE	1	piece(s)
90228	Beaker, low form, 50 ml:10 ml, with graduation and spout, ETFE	1	piece(s)
90238	Beaker, low form, 100 ml:20 ml, with graduation and spout, ETFE	1	piece(s)
90248	Beaker, low form, 250 ml:50 ml, with graduation and spout, ETFE	1	piece(s)
90252	Beaker, low form, 400 ml:50 ml, with graduation and spout, ETFE	1	piece(s)
90254	Beaker, low form, 500 ml:100 ml, with graduation and spout, ETFE	1	piece(s)
90256	Beaker, low form, 600 ml:100 ml, with graduation and spout, ETFE	1	piece(s)
90262	Beaker, low form, 1000 ml:100 ml, with graduation and spout, ETFE	1	piece(s)

Beaker, low form, PTFE, without graduation

Cat. No.	Description	From*	Unit
90320	Beaker, low form, 25 ml, without graduation, with spout, PTFE	1	piece(s)
90328	Beaker, low form, 50 ml, without graduation, with spout, PTFE	1	piece(s)
90338	Beaker, low form, 100 ml, without graduation, with spout, PTFE	1	piece(s)
90344	Beaker, low form, 150 ml, without graduation, with spout, PTFE	1	piece(s)
90348	Beaker, low form, 250 ml, without graduation, with spout, PTFE	1	piece(s)
90352	Beaker, low form, 400 ml, without graduation, with spout, PTFE	1	piece(s)
90354	Beaker, low form, 500 ml, without graduation, with spout, PTFE	1	piece(s)
90362	Beaker, low form, 1000 ml, without graduation, with spout, PTFE	1	piece(s)
90364	Beaker, low form, 2000 ml, without graduation, with spout, PTFE	1	piece(s)



Beaker, low form, PFA

Cat. No.	Description	From*	Unit
90120	Griffin beaker with spout, low form, PFA, 25:5 ml	1	piece(s)
90128	Griffin beaker with spout, low form, PFA, 50:10 ml	1	piece(s)
90138	Griffin beaker with spout, low form, PFA, 100:20 ml	1	piece(s)
90148	Griffin beaker with spout, low form, PFA, 250:50 ml	1	piece(s)
90154	Griffin beaker with spout, low form, PFA, 500:100 ml	1	piece(s)
90162	Griffin beaker with spout, low form, PFA, 1000:100 ml	1	piece(s)



Graduated beaker, PP, blue printed scale

Cat. No.	Description	From*	Unit
40802	Graduated beaker, PP, blue graduation, 50 ml: 2 ml, with handle and spout	24	piece(s)
40804	Graduated beaker, PP, blue graduation, 100 ml: 2 ml, with handle and spout	24	piece(s)
40806	Graduated beaker, PP, blue graduation, 250 ml: 5 ml, with handle and spout	12	piece(s)
40808	Graduated beaker, PP, blue graduation, 500 ml: 10 ml, with handle and spout	12	piece(s)
40810	Graduated beaker, PP, blue graduation, 1000 ml: 20 ml, with handle and spout	6	piece(s)
40812	Graduated beaker, PP, blue graduation, 2000 ml: 50 ml, with handle and spout	6	piece(s)
40814	Graduated beaker, PP, blue graduation, 3000 ml: 50 ml, with handle and spout	6	piece(s)
40816	Graduated beaker, PP, blue graduation, 5000 ml: 100 ml, with handle and spout	6	piece(s)



Graduated beaker, PP, embossed scale



Cat. No.	Description	From*	Unit
41002	Graduated beaker, PP, embossed graduation, 50 ml: 2 ml, with handle and spout	24	piece(s)
41004	Graduated beaker, PP, embossed graduation, 100 ml: 2 ml, with handle and spout	24	piece(s)
41006	Graduated beaker, PP, embossed graduation, 250 ml: 5 ml, with handle and spout	12	piece(s)
41008	Graduated beaker, PP, embossed graduation, 500 ml: 10 ml, with handle and spout	12	piece(s)
41010	Graduated beaker, PP, embossed graduation, 1000 ml: 20 ml, with handle and spout	6	piece(s)
41012	Graduated beaker, PP, embossed graduation, 2000 ml: 50 ml, with handle and spout	6	piece(s)
41014	Graduated beaker, PP, embossed graduation, 3000 ml: 50 ml, with handle and spout	6	piece(s)
41016	Graduated beaker, PP, embossed graduation, 5000 ml: 100 ml, with handle and spout	6	piece(s)

Funnels, short stem, Boro 3.3, fluted



Cat. No.	Description	From*	Unit
145611	Funnel, short stem, Boro 3.3, fluted out. diameter 70 mm, stem diameter 8 mm, length 70 mm	10	piece(s)
145616	Funnel, short stem, Boro 3.3, fluted out. diameter 80 mm, stem diameter 10 mm, length 80 mm	10	piece(s)
145619	Funnel, short stem, Boro 3.3, fluted out. diameter 100 mm, stem diameter 10 mm, length 100 mm	10	piece(s)
145631	Funnel, short stem, Boro 3.3, fluted out. diameter 150 mm, stem diameter 16 mm, length 150 mm	10	piece(s)
145638	Funnel, short stem, Boro 3.3, fluted out. diameter 200 mm, stem diameter 26 mm, length 170 mm	1	piece(s)

Funnels, short stem, Boro 3.3



Cat. No.	Description	From*	Unit
145505	Funnel, short stem, Boro 3.3, out. diameter 35 mm, stem diameter 6 mm, length 35 mm	1	piece(s)
145507	Funnel, short stem, Boro 3.3, out. diameter 40 mm, stem diameter 6 mm, length 40 mm	1	piece(s)
145515	Funnel, short stem, Boro 3.3, out. diameter 50 mm, stem diameter 7 mm, length 50 mm	1	piece(s)
145520	Funnel, short stem, Boro 3.3, out. diameter 55 mm, stem diameter 8 mm, length 55 mm	1	piece(s)
145525	Funnel, short stem, Boro 3.3, out. diameter 60 mm, stem diameter 8 mm, length 60 mm	1	piece(s)
145530	Funnel, short stem, Boro 3.3, out. diameter 70 mm, stem diameter 8 mm, length 70 mm	1	piece(s)
145535	Funnel, short stem, Boro 3.3, out. diameter 80 mm, stem diameter 9 mm, length 80 mm	1	piece(s)
145540	Funnel, short stem, Boro 3.3, out. diameter 100 mm, stem diameter 10 mm, length 100 mm	1	piece(s)
145545	Funnel, short stem, Boro 3.3, out. diameter 120 mm, stem diameter 16 mm, length 120 mm	1	piece(s)
145550	Funnel, short stem, Boro 3.3, out. diameter 150 mm, stem diameter 16 mm, length 150 mm	1	piece(s)
145560	Funnel, short stem, Boro 3.3, out. diameter 200 mm, stem diameter 26 mm, length 175 mm	1	piece(s)

Powder funnel, Boro 3.3

Cat. No.	Description	From*	Unit
146506	Powder funnel, Boro 3.3 glass, top diameter 50 mm, bottom diameter 18 mm	1	piece(s)
146508	Powder funnel, Boro 3.3 glass, top diameter 60 mm, bottom diameter 18 mm	1	piece(s)
146511	Powder funnel, Boro 3.3 glass, top diameter 70 mm, bottom diameter 20 mm	1	piece(s)
146516	Powder funnel, Boro 3.3 glass, top diameter 80 mm, bottom diameter 22 mm	1	piece(s)
146517	Powder funnel, Boro 3.3 glass, top diameter 90 mm, bottom diameter 25 mm	1	piece(s)
146519	Powder funnel, Boro 3.3 glass, top diameter 100 mm, bottom diameter 25 mm	1	piece(s)
146523	Powder funnel short, wide stem, Boro 3.3, top diameter 120 mm, bottom diameter 30 mm	1	piece(s)



Standard ground joint funnel, PP

Cat. No.	Description	From*	Unit
94905	Standard ground joint funnel, PP, funnel out. dia. 40 mm, stem for NS 14/23	10	piece(s)
94910	Standard ground joint funnel, PP, funnel out. dia. 50 mm, stem for NS 19/26	10	piece(s)
94915	Standard ground joint funnel, PP, funnel out. dia. 80 mm, stem for NS 29/32	5	piece(s)



Separatory funnel, PP

Cat. No.	Description	From*	Unit
140004	Separating funnel, PP, 125 ml, screw cap PP, con. PTFE stopcock	1	piece(s)
140005	Separating funnel, PP, 250 ml, screw cap PP, con. PTFE stopcock	1	piece(s)
140006	Separating funnel, PP, 500 ml, screw cap PP, con. PTFE stopcock	1	piece(s)
140007	Separating funnel, PP, 1000 ml, screw cap PP, con. PTFE stopcock	1	piece(s)




Funnels, PP

Cat. No.	Description	From*	Unit
147010	Funnel PP, outer diameter 40 mm, stem diameter 5 mm, length 35 mm	24	piece(s)
147015	Funnel PP, outer diameter 50 mm, stem diameter 9 mm, length 43 mm	24	piece(s)
147020	Funnel PP, outer diameter 75 mm, stem diameter 10 mm, length 55 mm	12	piece(s)
147025	Funnel PP, outer diameter 100 mm, stem diameter 13 mm, length 77 mm	12	piece(s)
147030	Funnel PP, outer diameter 120 mm, stem diameter 14 mm, length 90 mm	12	piece(s)
147035	Funnel PP, outer diameter 150 mm, stem diameter 17 mm, length 95 mm	12	piece(s)



Funnels, PP, fluted interior



Cat. No.	Description	From*	Unit
148215	Funnel with internal fluting, PP, out. diameter 45 mm, stem diameter 5 mm length 44 mm	20	piece(s)
148220	Funnel with internal fluting, PP, out. diameter 65 mm, stem diameter 9 mm length 63 mm	20	piece(s)
148225	Funnel with internal fluting, PP, out. diameter 80 mm, stem diameter 10 mm, length 69 mm	20	piece(s)
148235	Funnel with internal fluting, PP, out. diameter 100 mm, stem diameter 11 mm, length 82 mm	1	piece(s)
148237	Funnel with internal fluting, PP, out. diameter 120 mm, stem diameter 11 mm, length 85 mm	1	piece(s)
148240	Funnel with internal fluting, PP, out. diameter 150 mm, stem diameter 14 mm, length 115 mm	1	piece(s)
148245	Funnel with internal fluting, PP, out. diameter 180 mm, stem diameter 14 mm, length 143 mm	1	piece(s)

Powder funnel, PP



Cat. No.	Description	From*	Unit
148020	Powder funnel short, wide stem, PP out. diameter 60 mm stem diameter 15 mm, length 19 mm	20	piece(s)
148025	Powder funnel short, wide stem, PP out. diameter 80 mm stem diameter 15 mm, length 23 mm	20	piece(s)
148030	Powder funnel short, wide stem, PP out. diameter 100 mm stem diameter 25 mm, length 23 mm	20	piece(s)
148035	Powder funnel short, wide stem, PP out. diameter 120 mm stem diameter 30 mm, length 28 mm	1	piece(s)
148040	Powder funnel short, wide stem, PP out. diameter 150 mm stem diameter 35 mm, length 42 mm	1	piece(s)
148045	Powder funnel short, wide stem, PP out. diameter 180 mm stem diameter 40 mm, length 48 mm	1	piece(s)

Funnel holder, PP



Cat. No.	Description	From*	Unit
140300	Holder for separating funnel, PP, for rods of 8 - 14 mm diameter	1	piece(s)
148500	Funnel holder, PP for 1 funnel	5	piece(s)
148502	Funnel holder, PP for 2 funnels	5	piece(s)

Filter funnel support



Cat. No.	Description	From*	Unit
74506	Filter funnel support for 2 funnels, base plate, PP, 250 x 140 mm	1	piece(s)
74511	Filter funnel support for 4 funnels, base plate, PP, 450 x 140 mm	1	piece(s)

Filter funnels, Boro 3.3

Cat. No.	Description	From*	Unit
464003	Filter funnel, Boro 3.3, 50 ml, 3 D 3	1	piece(s)
464004	Filter funnel, Boro 3.3, 50 ml, 3 D 4	1	piece(s)
464713	Filter funnel, Boro 3.3, 75 ml, 11 D 3	1	piece(s)
464714	Filter funnel, Boro 3.3, 75 ml, 11 D 4	1	piece(s)
465423	Filter funnel, Boro 3.3, 125 ml, 17 D 3	1	piece(s)
465424	Filter funnel, Boro 3.3, 125 ml, 17 D 4	1	piece(s)
466133	Filter funnel, Boro 3.3, 500 ml, 25 D 3	1	piece(s)
466134	Filter funnel, Boro 3.3, 500 ml, 25 D 4	1	piece(s)



Filter crucible, Boro 3.3

Cat. No.	Description	From*	Unit
458021	Filter crucible, Boro 3.3, 30 ml, 1 D 1	10	piece(s)
458022	Filter crucible, Boro 3.3, 30 ml, 1 D 2	10	piece(s)
458023	Filter crucible, Boro 3.3, 30 ml, 1 D 3	10	piece(s)
458024	Filter crucible, Boro 3.3, 30 ml, 1 D 4	10	piece(s)
458731	Filter crucible, Boro 3.3, 50 ml, 2 D 1	10	piece(s)
458732	Filter crucible, Boro 3.3, 50 ml, 2 D 2	10	piece(s)
458733	Filter crucible, Boro 3.3, 50 ml, 2 D 3	10	piece(s)
458734	Filter crucible, Boro 3.3, 50 ml, 2 D 4	10	piece(s)



Filter adapter, Boro 3.3

Cat. No.	Description	From*	Unit
463616	Filter adapter Boro 3.3, suitable filter crucibles 1 D, upper diameter 41 mm, length 125 mm	10	piece(s)
463732	Filter adapter Boro 3.3, suitable filter crucibles 2 D, upper diameter 50 mm, length 132 mm	10	piece(s)



Rubber gasket, conical, for filter funnels

Cat. No.	Description	From*	Unit
462500	Rubber gasket con. EPDM for filter funnels and filter flasks size 22-84 mm	1 pack = 8 piece(s)	1 pack
462512	Rubber gasket con. EPDM for filter funnels and filter flasks, size 22	10	piece(s)
462617	Rubber gasket con. EPDM for filter funnels and filter flasks, size 28	10	piece(s)
462723	Rubber gasket con. EPDM for filter funnels and filter flasks, size 35	10	piece(s)
462827	Rubber gasket con. EPDM for filter funnels and filter flasks, size 42	10	piece(s)
462932	Rubber gasket con. EPDM for filter funnels and filter flasks, size 50	10	piece(s)
463036	Rubber gasket con. EPDM for filter funnels and filter flasks, size 63	10	piece(s)
463139	Rubber gasket con. EPDM for filter funnels and filter flasks, size 71	10	piece(s)
463243	Rubber gasket con. EPDM for filter funnels and filter flasks, size 84	10	piece(s)



Rubber sleeve, for filter crucibles



Cat. No.	Description	From*	Unit
462026	Rubber sleeve for filter crucibles 1 D, outer diameter 41 mm	10	piece(s)
462231	Rubber sleeve for filter crucibles 2 D, outer diameter 49 mm	10	piece(s)

Filter flask, Boro 3.3, with lateral socket



Cat. No.	Description	From*	Unit
73714	Filter flask, with lateral socket, 250 ml, Boro 3.3	1	piece(s)
73719	Filter flask, with lateral socket, 500 ml, Boro 3.3	1	piece(s)
73724	Filter flask, with lateral socket, 1000 ml, Boro 3.3	1	piece(s)
73729	Filter flask, with lateral socket, 2000 ml, Boro 3.3	1	piece(s)

Desiccators, Boro 3.3, with knob lid



Cat. No.	Description	From*	Unit
65031	Desiccator with knob lid, nominal size 150 mm, diameter 215 mm, Boro 3.3	1	piece(s)
65038	Desiccator with knob lid, nominal size 200 mm, diameter 270 mm, Boro 3.3	1	piece(s)
65043	Desiccator with knob lid, nominal size 250 mm, diameter 320 mm, Boro 3.3	1	piece(s)

Desiccators, Boro 3.3, with socket in the lid



Cat. No.	Description	From*	Unit
65231	Desiccator socket, NS 24/29, nominal size 150 mm, diameter 215 mm, Boro 3.3	1	piece(s)
65238	Desiccator socket, NS 24/29, nominal size 200 mm, diameter 270 mm, Boro 3.3	1	piece(s)
65243	Desiccator socket, NS 24/29, nominal size 250 mm, diameter 320 mm, Boro 3.3	1	piece(s)

Desiccators, PP, Lid PC



Cat. No.	Description	From*	Unit
65805	Desiccator/lid PC, base/desiccant tray PP, nominal size 150 mm, with venting stopper	1	piece(s)
65810	Desiccator/lid PC, base/desiccant tray PP, nominal size 200 mm, with venting stopper	1	piece(s)
65815	Desiccator/lid PC, base/desiccant tray PP, nominal size 250 mm, with venting stopper	1	piece(s)

Sealing ring, CR

Cat. No.	Description	From*	Unit
65820	Desiccator-sealing ring, CR, nominal size 150 mm	1	piece(s)
65822	Desiccator-sealing ring, CR, nominal size 200 mm	1	piece(s)
65824	Desiccator-sealing ring, CR, nominal size 250 mm	1	piece(s)



Desiccator plate

Cat. No.	Description	From*	Unit
65965	Dessicator/dessicator plate, porcelain, nominal size 150 mm, diameter 140 mm	1	piece(s)
65975	Dessicator/dessicator plate, porcelain, nominal size 200 mm, diameter 190 mm	1	piece(s)
65980	Dessicator/dessicator plate, porcelain, nominal size 250 mm, diameter 240 mm	1	piece(s)
66015	Dessicator/dessicator plate, PP, nominal size 150 mm, diameter 140 mm	1	piece(s)
66025	Dessicator/dessicator plate, PP, nominal size 200 mm, diameter 190 mm	1	piece(s)
66030	Dessicator/dessicator plate, PP, nominal size 250 mm, diameter 240 mm	1	piece(s)



Venting stopper, PC, for desiccators

Cat. No.	Description	From*	Unit
65804	Dessicator/venting stopper, PC, with non-return valve	1	piece(s)



Novus knob lid, Boro 3.3, for desiccators

Cat. No.	Description	From*	Unit
65631	Lid with knob, nominal size 150 mm, diameter 215 mm, Boro 3.3	1	piece(s)
65638	Lid with knob, nominal size 200 mm, diameter 270 mm, Boro 3.3	1	piece(s)
65643	Lid with knob, nominal size 250 mm, diameter 320 mm, Boro 3.3	1	piece(s)



Novus lid with socket in the lid for desiccators, Boro 3.3

Cat. No.	Description	From*	Unit
65731	Lid with socket, NS 24/29, nominal size 150 mm, diameter 215 mm, Boro 3.3	1	piece(s)
65738	Lid with socket, NS 24/29, nominal size 200 mm, diameter 270 mm, Boro 3.3	1	piece(s)
65743	Lid with socket, NS 24/29, nominal size 250 mm, diameter 320 mm, Boro 3.3	1	piece(s)



Needle-valve stopcock, Boro 3.3, for desiccators



Cat. No.	Description	From*	Unit
82415	Dessicator/needle-valve stopc. NS 24/29, for lid/lateral socket, Boro 3.3	1	piece(s)

Water jet pumps, PP



Cat. No.	Description	From*	Unit
159600	Water jet filter pump, PP, connection R 3/4", reducing adapter R 1/2", nozzle 10-12 mm	1	piece(s)

Tubing adaper, PP, for water jet pump



Cat. No.	Description	From*	Unit
159665	Reducing adapter for art.-no. 1596.00, R 3/8", PP	1	piece(s)
159670	Reducing adapter for art.-no. 1596.00, M 22x1 thread, PP	1	piece(s)

Tubings, PVC



Cat. No.	Description	From*	Unit
143255	Special laboratory tubing, PVC, inner diameter 2,0 mm, outer diameter 4,0, width 1,0 mm	20	meter
143260	Special laboratory tubing, PVC, inner diameter 4,0 mm, outer diameter 6,0, width 1,0 mm	20	meter
143265	Special laboratory tubing, PVC, inner diameter 5,0 mm, outer diameter 8,0, width 1,5 mm	20	meter
143270	Special laboratory tubing, PVC, inner diameter 6,0 mm, outer diameter 9,0, width 1,5 mm	20	meter
143275	Special laboratory tubing, PVC, inner diameter 7,0 mm, outer diameter 11,0, width 2,0 mm	20	meter
143280	Special laboratory tubing, PVC, inner diameter 8,0 mm, outer diameter 12,0, width 2,0 mm	20	meter
143285	Special laboratory tubing, PVC, inner diameter 10,0 mm, outer diameter 14,0, width 2,0 mm	20	meter

Tubings, silicone

Cat. No.	Description	From*	Unit
143352	Tubing Silicone, inner diameter 2,0, outer diameter 4,0, width 1,0 mm	25	meter
143355	Tubing Silicone, inner diameter 3,0, outer diameter 5,0, width 1,0 mm	25	meter
143356	Tubing Silicone, inner diameter 4,0, outer diameter 6,0, width 1,0 mm	25	meter
143357	Tubing Silicone, inner diameter 4,0, outer diameter 7,0, width 1,5 mm	25	meter
143358	Tubing Silicone, inner diameter 5,0, outer diameter 8,0, width 1,5 mm	25	meter
143359	Tubing Silicone, inner diameter 6,0, outer diameter 9,0, width 1,5 mm	25	meter
143360	Tubing Silicone, inner diameter 6,0, outer diameter 10,0, width 2,0 mm	25	meter
143361	Tubing Silicone, inner diameter 7,0, outer diameter 10,0, width 1,5 mm	25	meter
143362	Tubing Silicone, inner diameter 8,0, outer diameter 12,0, width 2,0 mm	25	meter
143364	Tubing Silicone, inner diameter 10,0, outer diameter 15,0, width 2,5 mm	25	meter



Flow indicator, SAN

Cat. No.	Description	From*	Unit
134080	Flow indicator, SAN, for tubing inner dia. 6-11 mm	5	piece(s)



Spring clip

Cat. No.	Description	From*	Unit
76015	Spring clip, spring steel, zinc/PE-coated grip, diameter 13, mm, bore 3,0 mm approx.	1	piece(s)
76020	Spring clip, spring steel, zinc/PE-coated grip, diameter 16, mm, bore 3,0 mm approx.	1	piece(s)
76025	Spring clip, spring steel, zinc/PE-coated grip, diameter 19, mm, bore 3,5 mm approx.	1	piece(s)
76035	Spring clip, spring steel, zinc/PE-coated grip, diameter 25, mm, bore 3,5 mm approx.	1	piece(s)
76040	Spring clip, spring steel, zinc/PE-coated grip, diameter 29, mm, bore 4,0 mm approx.	1	piece(s)
76045	Spring clip, spring steel, zinc/PE-coated grip, diameter 32, mm, bore 4,0 mm approx.	1	piece(s)
76050	Spring clip, spring steel, zinc/PE-coated grip, diameter 42, mm, bore 4,0 mm approx.	1	piece(s)
76055	Spring clip, spring steel, zinc/PE-coated grip, diameter 51, mm, 2 bores 4,0 mm approx.	1	piece(s)



Tubing connector, straight, PP

Cat. No.	Description	From*	Unit
152400	Universal tubing, connector, PP, straight for tubing, inner diameter 5 - 17 mm, length 110 mm	10	piece(s)



Tubing adapters, straight



Cat. No.	Description	From*	Unit
152505	Tubing adapter, PE-HD, for tubing, inner diameter 8-10/ 11-14 mm, length 66 mm	10	piece(s)
152510	Tubing adapter, PE-HD, for tubing, inner diameter 12-14/18-20,5 mm, length 73 mm	10	piece(s)
152605	Tubing adapter, PP, for tubing, inner diameter 3- 5/ 6-10 mm, length 55 mm	10	piece(s)
152615	Tubing adapter, PP, for tubing, inner diameter 6-10/ 9-13 mm, length 59 mm	10	piece(s)
152705	Tubing connector, PP, straight, for tubing, inner diameter 3 - 5 mm, length 47 mm	10	piece(s)
152710	Tubing connector, PP, straight, for tubing, inner diameter 6 - 10 mm, length 55 mm	10	piece(s)
152715	Tubing connector, PP, straight, for tubing, inner diameter 9 - 13 mm, length 57 mm	10	piece(s)

Non-return valve, PE-HD



Cat. No.	Description	From*	Unit
134002	Non-return valve, PE-HD, for tubing inner dia. 6-9 mm, length 78 mm	10	piece(s)

Tubing connector, PP, T-shape



Cat. No.	Description	From*	Unit
152800	Tubing connector, PP, T-shape, for tubing, inner diameter 3 - 4 mm, total length 20 mm	20	piece(s)
152803	Tubing connector, PP, T-shape, for tubing, inner diameter 4 - 5 mm, total length 30 mm	20	piece(s)
152805	Tubing connector, PP, T-shape, for tubing, inner diameter 6 - 7 mm, total length 40 mm	20	piece(s)
152810	Tubing connector, PP, T-shape, for tubing, inner diameter 8 - 9 mm, total length 55 mm	20	piece(s)
152815	Tubing connector, PP, T-shape, for tubing, inner diameter 10 - 11 mm, total length 69 mm	20	piece(s)

Tubing connector, PP, Y-shape



Cat. No.	Description	From*	Unit
152900	Tubing connector, PP, Y-shape, for tubing, inner diameter 3 mm, total length 24 mm	20	piece(s)
152903	Tubing connector, PP, Y-shape, for tubing, inner diameter 4 - 5 mm, total length 35 mm	20	piece(s)
152905	Tubing connector, PP, Y-shape, for tubing, inner diameter 6 - 7 mm, total length 54 mm	20	piece(s)
152907	Tubing connector, PP, Y-shape, for tubing, inner diameter 8 - 9 mm, total length 64 mm	20	piece(s)
152910	Tubing connector, PP, Y-shape, for tubing, inner diameter 10 - 11 mm, total length 74 mm	20	piece(s)

Stopcocks, PTFE, straight-bore

Cat. No.	Description	From*	Unit
88610	Stopcock one-way with nozzles, for tub. inner dia. 3 mm bore, 2 mm, PTFE	1	piece(s)
88615	Stopcock one-way with nozzles, for tub. inner dia. 4 mm bore, 2 mm, PTFE	1	piece(s)
88620	Stopcock one-way with nozzles, for tub. inner dia. 5 mm bore, 3 mm, PTFE	1	piece(s)
88625	Stopcock one-way with nozzles, for tub. inner dia. 6 mm bore, 3 mm, PTFE	1	piece(s)
88630	Stopcock one-way with nozzles, for tub. inner dia. 7 mm bore, 4 mm, PTFE	1	piece(s)
88635	Stopcock one-way with nozzles, for tub. inner dia. 8 mm bore, 4 mm, PTFE	1	piece(s)
88640	Stopcock one-way with nozzles, for tub. inner dia. 9 mm bore, 5 mm, PTFE	1	piece(s)



Stopcocks, PTFE, T-bore

Cat. No.	Description	From*	Unit
88650	Stopcock T-bore with nozzles, for tub. inner dia. 3 mm bore, 2 mm, PTFE	1	piece(s)
88655	Stopcock T-bore with nozzles, for tub. inner dia. 4 mm bore, 2 mm, PTFE	1	piece(s)
88660	Stopcock T-bore with nozzles, for tub. inner dia. 5 mm bore, 3 mm, PTFE	1	piece(s)
88665	Stopcock T-bore with nozzles, for tub. inner dia. 6 mm bore, 3 mm, PTFE	1	piece(s)
88670	Stopcock T-bore with nozzles, for tub. inner dia. 7 mm bore, 4 mm, PTFE	1	piece(s)
88675	Stopcock T-bore with nozzles, for tub. inner dia. 8 mm bore, 4 mm, PTFE	1	piece(s)
88680	Stopcock T-bore with nozzles, for tub. inner dia. 9 mm bore, 5 mm, PTFE	1	piece(s)



Stopcock grease, silicone free

Cat. No.	Description	From*	Unit
61610	Stopcock grease, high vacuum quality, silicone-free, operating temp. -40/+320°C, 60 g	1	piece(s)



Silicone paste, UNISILKON TK M 1011

Cat. No.	Description	From*	Unit
61605	Silicone paste UNISILKON TK M 1011, operating temperature -40/+160°C, 50 g	1	piece(s)



Magnetic stirring bars, cylindrical, PTFE



Cat. No.	Description	From*	Unit
137100	Magnetic stirring bar PTFE, 2,5 x 2,5 mm, cylindrical	10	piece(s)
137101	Magnetic stirring bar PTFE, 3,5 x 3,5 mm, cylindrical	10	piece(s)
137102	Magnetic stirring bar PTFE, 5 x 2 mm, cylindrical	10	piece(s)
137103	Magnetic stirring bar PTFE, 6 x 3 mm, cylindrical	10	piece(s)
137104	Magnetic stirring bar PTFE, 8 x 3 mm, cylindrical	10	piece(s)
137105	Magnetic stirring bar PTFE, 12 x 4,5 mm, cylindrical	10	piece(s)
137106	Magnetic stirring bar PTFE, 7 x 2 mm, cylindrical	10	piece(s)
137107	Magnetic stirring bar PTFE, 8 x 1,7 mm, cylindrical	10	piece(s)
137108	Magnetic stirring bar PTFE, 10 x 3 mm, cylindrical	10	piece(s)
137109	Magnetic stirring bar PTFE, 13 x 3 mm, cylindrical	10	piece(s)
137110	Magnetic stirring bar PTFE, 15 x 4,5 mm, cylindrical	10	piece(s)
137111	Magnetic stirring bar PTFE, 15 x 1,7 mm, cylindrical	10	piece(s)
137113	Magnetic stirring bar PTFE, 10 x 6 mm, cylindrical	10	piece(s)
137114	Magnetic stirring bar PTFE, 15 x 6 mm, cylindrical	10	piece(s)
137115	Magnetic stirring bar PTFE, 20 x 6 mm, cylindrical	10	piece(s)
137120	Magnetic stirring bar PTFE, 25 x 6 mm, cylindrical	10	piece(s)
137125	Magnetic stirring bar PTFE, 30 x 6 mm, cylindrical	10	piece(s)
137127	Magnetic stirring bar PTFE, 35 x 6 mm, cylindrical	10	piece(s)
137130	Magnetic stirring bar PTFE, 40 x 8 mm, cylindrical	10	piece(s)
137132	Magnetic stirring bar PTFE, 45 x 8 mm, cylindrical	10	piece(s)
137135	Magnetic stirring bar PTFE, 50 x 8 mm, cylindrical	10	piece(s)
137140	Magnetic stirring bar PTFE, 60 x 9 mm, cylindrical	10	piece(s)
137145	Magnetic stirring bar PTFE, 70 x 9 mm, cylindrical	10	piece(s)
137150	Magnetic stirring bar PTFE, 80 x 9 mm, cylindrical	10	piece(s)
137155	Magnetic stirring bar PTFE, 108 x 27 mm, cylindrical	1	piece(s)
137160	Magnetic stirring bar PTFE, 159 x 27 mm, cylindrical	1	piece(s)

Magnetic stirring bars, cylindrical with pivot ring, PTFE



Cat. No.	Description	From*	Unit
137404	Magnetic stirring bar PTFE, 8 x 3 mm, cylindrical, small pivot ring	10	piece(s)
137405	Magnetic stirring bar PTFE, 12 x 4,5 mm, cylindrical, small pivot ring	10	piece(s)
137410	Magnetic stirring bar PTFE, 15 x 4,5 mm, cylindrical, small pivot ring	10	piece(s)
137415	Magnetic stirring bar PTFE, 20 x 6 mm, cylindrical, small pivot ring	10	piece(s)
137420	Magnetic stirring bar PTFE, 25 x 6 mm, cylindrical, small pivot ring	10	piece(s)
137425	Magnetic stirring bar PTFE, 30 x 6 mm, cylindrical, small pivot ring	10	piece(s)
137427	Magnetic stirring bar PTFE, 35 x 6 mm, cylindrical, small pivot ring	10	piece(s)
137430	Magnetic stirring bar PTFE, 40 x 8 mm, cylindrical, small pivot ring	10	piece(s)
137432	Magnetic stirring bar PTFE, 45 x 8 mm, cylindrical, small pivot ring	10	piece(s)
137435	Magnetic stirring bar PTFE, 50 x 8,5 mm, cylindrical, small pivot ring	10	piece(s)
137445	Magnetic stirring bar PTFE, 70 x 8,5 mm, cylindrical, small pivot ring	10	piece(s)

Magnetic stirring bars, triangular, PTFE

Cat. No.	Description	From*	Unit
137810	Magnetic stirring bar PTFE, 25 x 8 mm, triangular	10	piece(s)
137812	Magnetic stirring bar PTFE, 25 x 14 mm, triangular	10	piece(s)
137815	Magnetic stirring bar PTFE, 35 x 10 mm, triangular	10	piece(s)
137817	Magnetic stirring bar PTFE, 40 x 14 mm, triangular	10	piece(s)
137820	Magnetic stirring bar PTFE, 50 x 12 mm, triangular	10	piece(s)
137822	Magnetic stirring bar PTFE, 55 x 14 mm, triangular	10	piece(s)
137826	Magnetic stirring bar PTFE, 80 x 16 mm, triangular	10	piece(s)



Magnetic stirring bars, octagonal with pivot ring, PTFE

Cat. No.	Description	From*	Unit
137607	Magnetic stirring bar PTFE, 13 x 8 mm, octagonal, pivot ring	10	piece(s)
137608	Magnetic stirring bar PTFE, 15 x 8 mm, octagonal, pivot ring	10	piece(s)
137610	Magnetic stirring bar PTFE, 25 x 8 mm, octagonal, pivot ring	10	piece(s)
137612	Magnetic stirring bar PTFE, 25 x 10 mm, octagonal, pivot ring	10	piece(s)
137615	Magnetic stirring bar PTFE, 38 x 8 mm, octagonal, pivot ring	10	piece(s)
137617	Magnetic stirring bar PTFE, 38 x 10 mm, octagonal, pivot ring	10	piece(s)
137618	Magnetic stirring bar PTFE, 51 x 8 mm, octagonal, pivot ring	10	piece(s)
137620	Magnetic stirring bar PTFE, 51 x 10 mm, octagonal, pivot ring	10	piece(s)
137622	Magnetic stirring bar PTFE, 64 x 10 mm, octagonal, pivot ring	10	piece(s)
137625	Magnetic stirring bar PTFE, 75 x 12 mm, octagonal, pivot ring	10	piece(s)



Magnetic stirring bars, double-ended, PTFE

Cat. No.	Description	From*	Unit
137930	Magnetic stirring bar PTFE length 37 mm, double-ended	10	piece(s)
137935	Magnetic stirring bar PTFE length 54 mm, double-ended	10	piece(s)



Magnetic stirring bars, oval, PTFE



Cat. No.	Description	From*	Unit
137300	Magnetic stirring bar PTFE, 10 x 5 mm, oval	10	piece(s)
137301	Magnetic stirring bar PTFE, 15 x 6 mm, oval	10	piece(s)
137302	Magnetic stirring bar PTFE, 20 x 10 mm, oval	10	piece(s)
137305	Magnetic stirring bar PTFE, 25 x 12 mm, oval	10	piece(s)
137307	Magnetic stirring bar PTFE, 32 x 15 mm, oval	10	piece(s)
137310	Magnetic stirring bar PTFE, 34 x 15 mm, oval	10	piece(s)
137312	Magnetic stirring bar PTFE, 40 x 20 mm, oval	10	piece(s)
137315	Magnetic stirring bar PTFE, 50 x 20 mm, oval	10	piece(s)
137320	Magnetic stirring bar PTFE, 65 x 20 mm, oval	10	piece(s)
137325	Magnetic stirring bar PTFE, 70 x 20 mm, oval	10	piece(s)

Magnetic stirring bars, crosshead, PTFE



Cat. No.	Description	From*	Unit
137630	Magnetic stirring bar PTFE, dia. 10 mm, h. 5 mm, crosshead	10	piece(s)
137632	Magnetic stirring bar PTFE, dia. 20 mm, h. 9 mm, crosshead	10	piece(s)
137634	Magnetic stirring bar PTFE, dia. 25 mm, h. 10 mm, crosshead	10	piece(s)
137636	Magnetic stirring bar PTFE, dia. 30 mm, h. 12 mm, crosshead	10	piece(s)
137638	Magnetic stirring bar PTFE, dia. 38 mm, h. 15 mm, crosshead	10	piece(s)

Magnetic stirring bar retriever, PP



Cat. No.	Description	From*	Unit
137750	Magnetic stirring bar retriever 300 x 10 mm, PP, with ring	1	piece(s)
137755	Magnetic stirring bar retriever 450 x 10 mm, PP, with ring	1	piece(s)

Magnetic stirring bar retriever, PTFE



Cat. No.	Description	From*	Unit
137700	Magnetic stirring bar retriever 150 x 8 mm, PTFE	1	piece(s)
137710	Magnetic stirring bar retriever 250 x 8 mm, PTFE	1	piece(s)
137720	Magnetic stirring bar retriever 350 x 8 mm, PTFE	1	piece(s)

Stirring rods, AR-Glas®

Cat. No.	Description	From*	Unit
135805	Stirring rod, AR-GLAS®, 200 mm x 5 mm fused closely at both ends	50	piece(s)
135810	Stirring rod, AR-GLAS®, 200 mm x 6 mm fused closely at both ends	50	piece(s)
135815	Stirring rod, AR-GLAS®, 250 mm x 6 mm fused closely at both ends	50	piece(s)
135820	Stirring rod, AR-GLAS®, 250 mm x 8 mm fused closely at both ends	50	piece(s)
135825	Stirring rod, AR-GLAS®, 300 mm x 8 mm fused closely at both ends	50	piece(s)

Stirring rods, spatulate, PTFE

Cat. No.	Description	From*	Unit
136215	Stirring rod, PTFE 200 mm x 8 mm, with spatulate ends	1	piece(s)
136220	Stirring rod, PTFE 250 mm x 8 mm, with spatulate ends	1	piece(s)
136225	Stirring rod, PTFE 300 mm x 8 mm, with spatulate ends	1	piece(s)

Stirring rods, with steel core, PTFE

Cat. No.	Description	From*	Unit
136315	Stirring rod, PTFE 200 mm x 6 mm, with steel core	1	piece(s)
136320	Stirring rod, PTFE 250 mm x 6 mm, with steel core	1	piece(s)
136325	Stirring rod, PTFE 300 mm x 6 mm, with steel core	1	piece(s)
136335	Stirring rod, PTFE 400 mm x 6 mm, with steel core	1	piece(s)

Sample cups, PS, for Technicon-Analyzer, CE-IVD



Cat. No.	Description	From*	Unit
115015	Sample vial, PS, glass clear, IVD 1,5 ml f. Technicon-Analyzer	1000	piece(s)
115016	Sample vial, PS, glass clear, IVD 2,0 ml f. Technicon-Analyzer	1000	piece(s)
115017	Sample vial, PS, glass clear, IVD 4,0 ml f. Technicon-Analyzer	1000	piece(s)

Press-on lid, PE, for Technicon-Analyzer, CE-IVD



Cat. No.	Description	From*	Unit
115020	Push-on cap for Technicon Analyzer, IVD PE suitable for 1,5/2,0 ml sample tubes	1000	piece(s)

Sample cup, PS, COULTER COUNTER®, with lid PE, CE-IVD



Cat. No.	Description	From*	Unit
722055	Sample vial, PS, with lid, PE, for Coulter-Counter, 20 ml, CE-IVD	1000	piece(s)

Sample cups, PS



Cat. No.	Description	From*	Unit
722060	Sample vial, PS, with lid, PE, for Leuko-Ery Counter, 12 ml	1000	piece(s)

Jar with snap-on lid, PS



Cat. No.	Description	From*	Unit
62315	Jar, PS, with snap-on lid, PE-LD, approx., 25 ml, conical shape	1000	piece(s)

Jar with push-on lid, PP



Cat. No.	Description	From*	Unit
61815	Jar with push-on lid, PP, approx. 50 ml, max. dia. 56 mm, h. 25 mm	10	piece(s)
61820	Jar with push-on lid, PP, approx. 80 ml max. dia. 75 mm, h. 30 mm	10	piece(s)

Jar with screw cap, PP

Cat. No.	Description	From*	Unit
62310	Jar with screw cap, PP, conical shape, approx., 30 ml, max.diameter 57 mm, height 32 mm	500	piece(s)

**Urine beaker, PP, CE-IVD**

Cat. No.	Description	From*	Unit
758900	Urine beaker, PP, grad. to 125 ml, red press-on lid, PE	1000	piece(s)

**Press-on lid, PP, for urine beaker**

Cat. No.	Description	From*	Unit
758902	Press-on lid for urine beaker 7589 01, PE, red	1000	piece(s)

**Urine beaker, PP, Screw cap, PE, CE-IVD**

Cat. No.	Description	From*	Unit
758905	Urine beaker, PP, screw cap grad. to 125 ml, non-sterile, green cap	1000	piece(s)

**Urine beaker, PP, without lid, CE-IVD**

Cat. No.	Description	From*	Unit
758901	Urine beaker w/o lid PP grad. to 125ml	1000	piece(s)

Universal container, PP, Screw cap, PE, y-sterile

Cat. No.	Description	From*	Unit
758910	Universal container, PP, screw cap PE (yellow), 125 ml, y-sterile	240	piece(s)



Feces container, PS, CE-IVD



Cat. No.	Description	From*	Unit
62305	Faeces container, PS, approx. 30 ml, non-sterile (blue cap), CE-IVD	400	piece(s)

Petri dishes, PS



Cat. No.	Description	From*	Unit
452000	Disposable petri dish, PS, with lid, lid dia. 94 mm overall, h. 16 mm, w/o vent	480	piece(s)
452005	Disposable petri dish, PS, with lid, lid dia. 94 mm overall, h. 16 mm, with vent	480	piece(s)
452010	Disposable petri dish, PS, with lid, lid dia. 55 mm overall, h. 16 mm, with vent	1620	piece(s)
452015	Disposable petri dish, PS, with lid, lid dia. 55 mm overall, h. 16 mm, w/o vent	1620	piece(s)

Petri dishes, soda-lime glass



Cat. No.	Description	From*	Unit
455701	Petri dish, soda-lime glass, lid diameter 40 mm, height of dish 12 mm	10	piece(s)
455717	Petri dish, soda-lime glass, lid diameter 60 mm, height of dish 15 mm	10	piece(s)
455732	Petri dish, soda-lime glass, lid diameter 80 mm, height of dish 15 mm	10	piece(s)
455742	Petri dish, soda-lime glass, lid diameter 100 mm, height of dish 15 mm	10	piece(s)
455743	Petri dish, soda-lime glass, lid diameter 100 mm, height of dish 20 mm	10	piece(s)
455751	Petri dish, soda-lime glass, lid diameter 150 mm, height of dish 25 mm	10	piece(s)


Sterilization indicator tape



Cat. No.	Description	From*	Unit
61750	Sterilization indicator tape for temp. ind. 50 m x 19 mm, crepe paper, self-adhesive	1	piece(s)

Inoculation loops, PS, γ -sterilized

Cat. No.	Description	From*	Unit
452201	Inoculation loop with needle, PS γ -sterilized 1, μ l white for inocul. nutrient media	1000	piece(s)
452210	Inoculation loop with needle, PS γ -sterilized 10, μ l blue for inocul. nutrient media	1000	piece(s)
452215	Inoculation loop with needle, PS, γ -ray sterilized, 1+10 μ l, yellow, for inocul. of nutrient media	1000	piece(s)



Disposal bags, PP

Cat. No.	Description	From*	Unit
759705	Disposal bag, PP, l. 300 mm, width 200 mm	100	piece(s)



Disposal bags, PA

Cat. No.	Description	From*	Unit
759710	Disposal bag, PA, l. 300 mm, width 200 mm	50	piece(s)



Rack for disposal bag

Cat. No.	Description	From*	Unit
759700	Stand for disposal bag, 100 bags of PP	1	piece(s)
759701	Stand for disposal bag, 50 bags of PA	1	piece(s)
759703	Stand for disposal bags, white	1	piece(s)



Pasteur Pipettes glass, soda-lime glass

Cat. No.	Description	From*	Unit
747715	Pasteur pipette, soda-lime glass, total length approx. 145 mm cap. approx. 2 ml	1000	piece(s)
747720	Pasteur pipette, soda-lime glass, total length approx. 225 mm cap. approx. 2 ml	1000	piece(s)
747725	Pasteur pipette, soda-lime glass, total length approx. 270 mm cap. approx. 1,5 ml	1000	piece(s)
747730	Pasteur pipette, soda-lime glass, total length approx. 310 mm cap. approx. 1,5 ml	1000	piece(s)



Pasteur pipettes, PE-LD



Cat. No.	Description	From*	Unit
747750	Pasteur pipette, PE-LD, ungrad., suction volume max. 3,0 ml	500	piece(s)
747755	Pasteur pipette, PE-LD, 1/0,25 ml, suction volume max. 3,5 ml	500	piece(s)
747760	Pasteur pipette, PE-LD, 3/0,5 ml, suction volume max. 3,5 ml	500	piece(s)
747765	Pasteur pipette, PE-LD, 2/0,5 ml, suction volume max. 2,0 ml	500	piece(s)
747770	Pasteur pipette, PE-LD, ungrad., suction volume max. 4,0 ml	500	piece(s)
747775	Pasteur pipette, PE-LD, ungrad., suction volume max. 1,0 ml	500	piece(s)

Micro haematocrit capillaries, CE-IVD



Cat. No.	Description	From*	Unit
749311	Micro haematocrit capillary, sodium hep., w/o cal. mark inner dia. 1,1 - 1,2 mm, code red	1000	piece(s)
749321	Micro haematocrit capillary, non-hep., w/o cal. mark inner dia. 1,1 - 1,2 mm, code blue	1000	piece(s)

Micropipettes, BLAUBRAND®, intraMark, DE-M, CE



Cat. No.	Description		From*	Unit
708707	Micropipettes, BLAUBRAND®, DE-M, intraMark, 1-2-3-4-5 µl, white	1 jar = 250 piece(s)	4	jar
708709	Micropipettes, BLAUBRAND®, DE-M, intraMark, 10 µl, orange	1 jar = 250 piece(s)	4	jar
708718	Micropipettes, BLAUBRAND®, DE-M, intraMark, 20 µl, black	1 jar = 250 piece(s)	4	jar
708722	Micropipettes, BLAUBRAND®, DE-M, intraMark, 25 µl, 2 x white	1 jar = 250 piece(s)	4	jar
708727	Micropipettes, BLAUBRAND®, DE-M, intraMark, 40 µl, 2 x red	1 jar = 250 piece(s)	4	jar
708728	Micropipettes, BLAUBRAND®, DE-M, intraMark, 20 and 40 µl, 2 x red	1 jar = 250 piece(s)	4	jar
708733	Micropipettes, BLAUBRAND®, DE-M, intraMark, 50 µl, green	1 jar = 250 piece(s)	4	jar
708744	Micropipettes, BLAUBRAND®, DE-M, intraMark, 100 µl, blue	1 jar = 250 piece(s)	4	jar
708745	Micropipettes, BLAUBRAND®, DE-M, intraMark, 50 and 100 µl, blue	1 jar = 250 piece(s)	4	jar
708757	Micropipettes, BLAUBRAND®, DE-M, intraMark, 200 µl, red	1 jar = 200 piece(s)	4	jar

Micropipettes, BLAUBRAND®, intraEND, DE-M, CE

Cat. No.	Description		From*	Unit
709101	Micropipettes, BLAUBRAND®, DE-M, intraEND, 1 µl	4 jar = 250 piece(s)	4	jar
709103	Micropipettes, BLAUBRAND®, DE-M, intraEND, 2 µl	4 jar = 250 piece(s)	4	jar
709105	Micropipettes, BLAUBRAND®, DE-M, intraEND, 3 µl	4 jar = 250 piece(s)	4	jar
709107	Micropipettes, BLAUBRAND®, DE-M, intraEND, 5 µl	4 jar = 250 piece(s)	4	jar
709109	Micropipettes, BLAUBRAND®, DE-M, intraEND, 10 µl	4 jar = 250 piece(s)	4	jar
709118	Micropipettes, BLAUBRAND®, DE-M, intraEND, 20 µl	4 jar = 250 piece(s)	4	jar
709122	Micropipettes, BLAUBRAND®, DE-M, intraEND, 25 µl	4 jar = 250 piece(s)	4	jar
709133	Micropipettes, BLAUBRAND®, DE-M, intraEND, 50 µl	4 jar = 250 piece(s)	4	jar
709144	Micropipettes, BLAUBRAND®, DE-M, intraEND, 100 µl	4 jar = 250 piece(s)	4	jar



Pipette holder, for micropipettes intraEND

Cat. No.	Description		From*	Unit
708605	Pipette holder for intraEND and Delbrück micropipettes		10	piece(s)



Pipetting aid, for micropipettes intraEND, PET tube, SI adapter, TPE suction bellow

Cat. No.	Description		From*	Unit
709110	Pipetting aid BLAUBRAND®, intraEND micropip. glass tube, rub. adapter, rub. teat with venting hole		10	piece(s)



Haematocrit sealing compound, 24 positions, CE-IVD

Cat. No.	Description		From*	Unit
749510	Haematocrit sealing compound, 24 numbered pos. for two capillaries each	1 pack = 10 piece(s)	1	pack



Dropping pipette, PE-LD, with graduation

Cat. No.	Description		From*	Unit
125410	Dropping pipette with bellows, PE-LD, 1,5 ml, length 133 mm, graduated ml: 0.5; 1; 1.5		100	piece(s)
125420	Dropping pipette with bellows, PE-LD, 5 ml, length 194 mm, graduated ml: 0.5; 1; 2; 3; 4; 5		100	piece(s)



Dropping pipette, PE-LD



Cat. No.	Description	From*	Unit
125400	Dropping pipette with pipetting bulb, PE-LD, approx. 1,8 ml, length approx. 98 mm	250	piece(s)

Staining troughs with tray, soda-lime glass




Cat. No.	Description	From*	Unit
472000	Staining trough, soda-lime glass, tray for 10 slides, 88 x 70 x 40 mm	10	piece(s)
472200	Staining trough, soda-lime glass, with lid 105 x 85 x 70 mm	10	piece(s)

Staining troughs with tray, PMP



Cat. No.	Description	From*	Unit
474300	Staining trough, PMP, glass clear, 101 x 83 x 70 mm without tray, with two lids	4	piece(s)

Staining troughs with tray, POM



Cat. No.	Description	From*	Unit
471800	Staining trough, complete, with tray for 25 slides	5	piece(s)

Staining troughs, Schiefferdecker, soda-lime glass



Cat. No.	Description	From*	Unit
472500	Staining trough, Schiefferdecker pattern, for 10 slides 76 x 26 mm, soda-lime glass	10	piece(s)

Staining troughs, Schiefferdecker, PMP



Cat. No.	Description	From*	Unit
474410	Staining trough, PMP, Schiefferdecker pattern, for 20 slides 76x26 mm, with lid	4	piece(s)

Staining troughs, Hellendahl, with lid

Cat. No.	Description	From*	Unit
472600	Staining trough, Hellendahl pattern, for 16 slides 76 x 26 mm, soda-lime glass	10	piece(s)



Staining troughs, Hellendahl, soda-lime glass, extended

Cat. No.	Description	From*	Unit
472700	Staining trough, 'Hellendahl extended', for 16 slides 76 x 26 mm, soda-lime glass	10	piece(s)




Staining troughs, Hellendahl, PMP, extended

Cat. No.	Description	From*	Unit
474400	Staining trough, PMP, glass clear Hellendahl, extended, for 16 slides 76x26 mm, with lid	4	piece(s)



Staining troughs, Coplin, with lid

Cat. No.	Description	From*	Unit
472800	Staining trough, Coplin pattern, round, for 10 slides 76 x 26 mm, soda-lime glass	10	piece(s)



Cover glasses, pure white

Cat. No.	Description	From*	Unit
470045	Cover glass, boros. glass, thickn. no 1, 18 x 18 mm, square, water-white	2000	piece(s)
470050	Cover glass, boros. glass, thickn. no 1, 20 x 20 mm, square, water-white	2000	piece(s)
470055	Cover glass, boros. glass, thickn. no 1, 22 x 22 mm, square, water-white	2000	piece(s)
470060	Cover glass, boros. glass, thickn. no 1, 24 x 24 mm, square, water-white	2000	piece(s)
470816	Cover glass, boros. glass, thickn. no 1, 24 x 40 mm, rectangular, water-white	1000	piece(s)
470819	Cover glass, boros. glass, thickn. no 1, 24 x 50 mm, rectangular, water-white	1000	piece(s)
470820	Cover glass, boros. glass, thickn. no 1, 24 x 60 mm, rectangular, water-white	1000	piece(s)




Haemocytometer cover glasses for counting chambers, CE-IVD




Cat. No.	Description	From*	Unit
723014	Haemocytometer cover glass for counting chambers, borosilicate glass, 24 x 24 x 0,4 mm, IVD	100	piece(s)
723015	Haemocytometer cover glass for counting chambers, borosilicate glass, 20 x 26 x 0,4 mm, IVD	100	piece(s)
723016	Haemocytometer cover glass for counting chambers, borosilicate glass, 22 x 30 x 0,4 mm, IVD	100	piece(s)

Counting chambers, Neubauer, BLAUBRAND®, CE




Cat. No.	Description	From*	Unit
718605	Counting chamber BLAUBRAND® Neubauer pattern, w/o clips double rul.	1	piece(s)
718620	Counting chamber BLAUBRAND® Neubauer pattern, with clips double rul.	1	piece(s)

Counting chambers, Neubauer, Improved, BLAUBRAND®, CE



Cat. No.	Description	From*	Unit
717805	Counting chamber BLAUBRAND® Neubauer improved, w/o clips double rul.	1	piece(s)
717820	Counting chamber BLAUBRAND® Neubauer improved, with clips double rul.	1	piece(s)

Counting chambers, Neubauer, Improved, bright-line, BLAUBRAND®, CE




Cat. No.	Description	From*	Unit
717810	Counting chamber BLAUBRAND® Neubauer improved, w/o clips bright-line	1	piece(s)

Counting chambers, Thoma pattern, BLAUBRAND®, CE



Cat. No.	Description	From*	Unit
718005	Counting chamber BLAUBRAND® Thoma pattern, w/o clips double rul.	1	piece(s)
718020	Counting chamber BLAUBRAND® Thoma pattern, with clips double rul.	1	piece(s)


Counting chambers, Bürker pattern, BLAUBRAND®, CE



Cat. No.	Description	From*	Unit
718905	Counting chamber BLAUBRAND® Bürker pattern, w/o clips double rul.	1	piece(s)
718920	Counting chamber BLAUBRAND® Bürker pattern, with clips double rul.	1	piece(s)


Counting chambers, Malassez pattern, BLAUBRAND®, CE

Cat. No.	Description	From*	Unit
719005	Counting chamber BLAUBRAND® Malassez pattern, w/o clips double rul.	1	piece(s)




Counting chambers, Bürker-Türk pattern, BLAUBRAND®, CE

Cat. No.	Description	From*	Unit
719505	Counting chamber BLAUBRAND® Bürker-Türk, w/o clips double rul.	1	piece(s)
719520	Counting chamber BLAUBRAND® Bürker-Türk, with clips double rul.	1	piece(s)



Counting chambers, Fuchs-Rosenthal pattern, BLAUBRAND®, CE

Cat. No.	Description	From*	Unit
719805	Counting chamber BLAUBRAND® Fuchs-Rosenthal, w/o clips double rul.	1	piece(s)
719820	Counting chamber BLAUBRAND® Fuchs-Rosenthal, with clips double rul.	1	piece(s)



Counting chambers, Nageotte pattern, BLAUBRAND®, CE

Cat. No.	Description	From*	Unit
721305	Counting chamber BLAUBRAND® Nageotte pattern, w/o clips, double rul.	1	piece(s)



Microscope slides, 1 mm

Cat. No.	Description	From*	Unit
474701	Microscope slide, half-white, cut edges, standard pack	2500	piece(s)
474702	Microscope slide, half-white, cut edges, twin frosted end, standard pack	2500	piece(s)
474743	Microscope slide, half-white, ground edges, standard pack	2500	piece(s)
474744	Microscope slide, half-white, ground edges, twin frosted end, standard pack	2500	piece(s)



Microscope slides, white, 1,2 - 1,5 mm

Cat. No.	Description	From*	Unit
475505	Cavity slides, white, thick. 1,2 - 1,5 mm, 76x26 mm ground edges with 1 concavity	50	piece(s)
475535	Cavity slides, white, thick. 1,2 - 1,5 mm, 76x26 mm ground edges with 2 concavities	50	piece(s)
475565	Cavity slides, white, thick. 1,2 - 1,5 mm, 76x26 mm ground edges with 3 concavities	50	piece(s)



Spotting tile, soda-lime glass




Cat. No.	Description	From*	Unit
473500	Spotting tile, soda-lime glass, 130 x 100 x 6 mm, 12 cavities diameter 20 - 22 mm	1	piece(s)

Transport box for microscope slides, PP



Cat. No.	Description	From*	Unit
476900	Slide box for storage/transport PP round, for 5 thick or 10 thin slides	10	piece(s)

Slide box, angled, PS



Cat. No.	Description	From*	Unit
475800	Slide box, PS, for, 25 slides 76 x 26 mm	1	piece(s)
475900	Slide box, PS, for, 50 slides 76 x 26 mm	1	piece(s)
476000	Slide box, PS, for 100 slides 76 x 26 mm	1	piece(s)

Trays, POM, for staining troughs

Cat. No.	Description	From*	Unit
471400	Trays for staining troughs, POM, 190 x 166 x 40 mm, for 25 slides 76 x 26 mm	5	piece(s)
471500	Staining trough/through for slides, trough separate, POM	5	piece(s)

Trays, PP, for staining troughs

Cat. No.	Description	From*	Unit
474305	Tray, PP, for staining through, PMP, for 20 slides	2	piece(s)

Wire, for staining trough tray



Cat. No.	Description	From*	Unit
473100	Wire handle for stain. trough tray 4720 00 made of stainless steel	10	piece(s)

Storage/dispenser box

Cat. No.	Description	From*	Unit
131900	Storage/dispenser box, PMMA, large, 165 x 152 x 355 mm	1	piece(s)
131902	Storage/dispenser box, PMMA, small, 165 x 152 x 178 mm	1	piece(s)



Test tubes, Boro 3.3, with PP stopper

Cat. No.	Description	From*	Unit
114508	Test tube, Boro 3.3, with PP stopper, 10 ml: 0,1 ml, NS 12/21, 15 x 165 mm	10	piece(s)
114516	Test tube, Boro 3.3, with PP stopper, 20 ml: 0,2 ml, NS 14/23, 17 x 200 mm	10	piece(s)
114521	Test tube, Boro 3.3, with PP stopper, 25 ml: 0,2 ml, NS 14/23, 17 x 220 mm	10	piece(s)
114524	Test tube, Boro 3.3, with PP stopper, 30 ml: 0,5 ml, NS 14/23, 19 x 220 mm	10	piece(s)



Centrifuge tubes, PC, rimless

Cat. No.	Description	From*	Unit
781029	Centrifuge tube, PC, round, 50 ml, diameter 35 mm, height 98 mm, graduation, subd. 1 ml	10	piece(s)
781039	Centrifuge tube, PC, round, 100 ml, diameter 41 mm, height 115 mm, graduation, subd. 2 ml	10	piece(s)
781040	Centrifuge tube, PC, round, 100 ml, diameter 45 mm, height 98 mm, graduation, subd. 1 ml	10	piece(s)



Centrifuge tubes, PP, with screw cap

Cat. No.	Description	From*	Unit
114817	Centrifuge tube, PP, graduated, 15 ml (graduated up to 13 ml), non-sterile, w/o base, screw cap	750	piece(s)
114818	Centrifuge tube, PP, graduated, 15 ml (graduated up to 13 ml), y-sterilized, w/o base, screw cap	750	piece(s)
114820	Centrifuge tube, PP, graduated, 50 ml, non-sterile, w/o base, screw cap	300	piece(s)
114821	Centrifuge tube, PP, graduated, 50 ml, y-sterilized, w/o base, screw cap	300	piece(s)
114822	Centrifuge tube, PP, graduated, 50 ml, non-sterile, with base, screw cap	250	piece(s)
114823	Centrifuge tube, PP, graduated, 50 ml, y-sterilized, with base, screw cap	250	piece(s)




Centrifuge tubes, PP, with beaded rim, without stopper


Cat. No.	Description	From*	Unit
115342	Centrifuge tube, PP, 10 ml, 16x100 mm, cylindrical, with rim, without stopper	250	piece(s)
115346	Centrifuge tube, PP, 26 ml, 24x90 mm, cylindrical, with rim, without stopper	50	piece(s)
115348	Centrifuge tube, PP, 30 ml, 20x100 mm, cylindrical, with rim, without stopper	50	piece(s)
115350	Centrifuge tube, PP, 48 ml, 30x100 mm, cylindrique, with rim, without stopper	25	piece(s)
115352	Centrifuge tube, PP, 75 ml, 35x100 mm, cylindrical, with rim, without stopper	20	piece(s)
115354	Centrifuge tube, PP, 110 ml, 40x120 mm, cylindrical, with rim, without stopper	20	piece(s)
115356	Centrifuge tube, PP, 160 ml, 45x120 mm, cylindrical, with rim, without stopper	10	piece(s)

Centrifuge tubes, ungraduated

Cat. No.	Description	From*	Unit
778012	Centrifuge tube, AR-GLAS®, approx. 15 ml, conical, beaded rim, ungrad., length 113 mm	100	piece(s)
778013	Centrifuge tube, Boro 3.3, approx. 15 ml, conical, beaded rim, ungrad., length 113 mm	100	piece(s)
779012	Centrifuge tube, AR-GLAS®, approx. 15 ml, cone-point bottom rimlength ungrad. length 98 mm	100	piece(s)
779022	Centrifuge tube, AR-GLAS®, approx. 15 ml, round bottom, rimless, ungrad., length 98 mm	100	piece(s)



Centrifuge tubes, graduated

Cat. No.	Description	From*	Unit
778312	Centrifuge tube, AR-GLAS®, approx. 15 ml, conical, beaded rim, grad. white, 0-10 ml	10	piece(s)
778313	Centrifuge tube, Boro 3.3, approx. 15 ml, conical, beaded rim, grad. white, 0-10 ml	10	piece(s)
778712	Centrifuge tube, AR-GLAS®, approx. 15 ml, conical, beaded rim, grad. white, 0-15 ml	10	piece(s)
778713	Centrifuge tube, Boro 3.3, approx. 15 ml, conical, beaded rim, grad. white, 0-15 ml	10	piece(s)



Stopper for centrifuge tubes

Cat. No.	Description	From*	Unit
115360	Stopper, PE for centrifuge tubes art.-no. 1153 42	100	piece(s)
115366	Stopper, PE for centrifuge tubes art.-no. 1153 48	100	piece(s)
115368	Stopper, PE for centrifuge tubes art.-no. 1153 46	100	piece(s)
115370	Stopper, PE for centrifuge tubes art.-no. 1153 50	100	piece(s)
115372	Stopper, PE for centrifuge tubes art.-no. 1153 52	100	piece(s)
115374	Stopper, PE for centrifuge tubes art.-no. 1153 54	100	piece(s)
115376	Stopper, PE for centrifuge tubes art.-no. 1153 56	100	piece(s)



Culture tubes, soda-lime glass, rimless

Cat. No.	Description	From*	Unit
114105	Culture tube, soda-lime glass, rimless, 10 x 75 mm, wall thickness 0,6 mm	250	piece(s)
114106	Culture tube, soda-lime glass, rimless, 12 x 75 mm, wall thickness 0,6 mm	250	piece(s)
114110	Culture tube, soda-lime glass, rimless, 12 x 100 mm, wall thickness 0,6 mm	144	piece(s)
114115	Culture tube, soda-lime-glass, rimless, 16 x 100 mm, wall thickness 0,7 mm	78	piece(s)
114120	Culture tube, soda-lime glass, rimless, 16 x 125 mm, wall thickness 0,8 mm	105	piece(s)
114125	Culture tube, soda-lime glass, rimless, 16 x 160 mm, wall thickness 0,8 mm	100	piece(s)
114130	Culture tube, soda-lime glass, rimless, 18 x 180 mm, wall thickness 0,8 mm	121	piece(s)



Culture tubes, soda-lime glass, with screw cap, PP


Cat. No.	Description	From*	Unit
113931	Culture tube, soda-lime glass, screw cap PP, 12 x 100 mm wall thickn. 1 mm, GL 14	100	piece(s)
113935	Culture tube, soda-lime glass, screw cap PP, 16 x 100 mm, wall thickness 1 mm, GL 18	100	piece(s)
113941	Culture tube, soda-lime glass, screw cap PP, 16 x 160 mm, wall thickness 1 mm, GL 18	100	piece(s)
113943	Culture tube, soda-lime glass, screw cap PP, 18 x 180 mm, wall thickness 1 mm, GL 18	100	piece(s)

Sample tube, PS


Cat. No.	Description	From*	Unit
114715	Sample tube, PS, glass clear, 16 x 100 mm	2000	piece(s)
114750	Sample tube (coagulometer), PS, glass clear, 11 x 55 mm	5000	piece(s)
114760	Sample tube, PS, glass clear, 11,5 x 75 mm	4000	piece(s)

Grip stopper, PE-LD


Cat. No.	Description	From*	Unit
114720	Grip stopper, PE-LD, for sample tubes cat. no. 114715	1000	piece(s)
114730	Grip stopper, PE-LD, for sample tubes cat. no. 114750/114760	1000	piece(s)

Sedimentation tube, PS


Cat. No.	Description	From*	Unit
114740	Sedimentation tube, PS, 16 x 105 mm, conical bottom	2000	piece(s)

Test tube rack, PP

Cat. No.	Description	From*	Unit
4340000	Test tube rack, PP, white, 265 x 126 x 75 mm for 84 tubes up to diameter 13 mm	5	piece(s)
4340001	Test tube rack, PP, blue, 265 x 126 x 75 mm for 84 tubes up to diameter 13 mm	5	piece(s)
4340002	Test tube rack, PP, red, 265 x 126 x 75 mm for 84 tubes up to diameter 13 mm	5	piece(s)
4340010	Test tube rack, PP, white, 265 x 126 x 75 mm for 55 tubes up to diameter 18 mm	5	piece(s)
4340011	Test tube rack, PP, blue, 265 x 126 x 75 mm for 55 tubes up to diameter 18 mm	5	piece(s)
4340012	Test tube rack, PP, red, 265 x 126 x 75 mm for 55 tubes up to diameter 18 mm	5	piece(s)
4340020	Test tube rack, PP, white, 265 x 126 x 75 mm for 40 tubes up to diameter 20 mm	5	piece(s)
4340021	Test tube rack, PP, blue, 265 x 126 x 75 mm for 40 tubes up to diameter 20 mm	5	piece(s)
4340022	Test tube rack, PP, red, 265 x 126 x 75 mm for 40 tubes up to diameter 20 mm	5	piece(s)
4340030	Test tube rack, PP, white, 265 x 126 x 88 mm for 32 tubes up to diameter 25 mm	5	piece(s)
4340031	Test tube rack, PP, blue, 265 x 126 x 88 mm for 32 tubes up to diameter 25 mm	5	piece(s)
4340032	Test tube rack, PP, red, 265 x 126 x 88 mm for 32 tubes up to diameter 25 mm	5	piece(s)
4340040	Test tube rack, PP, white, 265 x 126 x 88 mm for 21 tubes up to diameter 30 mm	5	piece(s)
4340041	Test tube rack, PP, blue, 265 x 126 x 88 mm for 21 tubes up to diameter 30 mm	5	piece(s)
4340042	Test tube rack, PP, red, 265 x 126 x 88 mm for 21 tubes up to diameter 30 mm	5	piece(s)
4340060	Test tube rack, PP, white, 265 x 126 x 75 mm for 55 tubes up to diameter 16 mm	5	piece(s)
4340061	Test tube rack, PP, blue, 265 x 126 x 75 mm for 55 tubes up to diameter 16 mm	5	piece(s)
4340062	Test tube rack, PP, red, 265 x 126 x 75 mm for 55 tubes up to diameter 16 mm	5	piece(s)



Test tube rack, PTFE

Cat. No.	Description	From*	Unit
115510	Test tube rack, PTFE f. 21 tub., up to dia.13 mm, 180x60x60 mm	1	piece(s)
115515	Test tube rack, PTFE f. 10 tub., up to dia.19 mm, 180x60x70 mm	1	piece(s)
115520	Test tube rack, PTFE f. 4 tub., up to dia.30 mm, 180x60x80 mm	1	piece(s)



Pipette tray



Cat. No.	Description	From*	Unit
28300	Pipette tray, PVC, l x w x h: 280 x 220 x 39 mm	1	piece(s)
28305	Pipette tray, PVC, l x w x h: 420 x 300 x 30 mm	1	piece(s)
768515	Pipette tray PVC f. small-vol. pipettes, 355 x 300 x 45 mm, nine oblong dividers	1	piece(s)

Pipette stand



Cat. No.	Description	From*	Unit
30900	Pipette stand, PP for 94 pipettes, rotary, 230 x 450 mm	1	piece(s)

Trays (photographic trays), PP



Cat. No.	Description	From*	Unit
156640	Tray (photogr. tray), PP, white stackable, 225 x 180 x 45 mm	1	piece(s)
156642	Tray (photogr. tray), PP, white stackable, 300 x 240 x 70 mm	1	piece(s)
156644	Tray (photogr. tray), PP, white stackable, 370 x 310 x 75 mm	1	piece(s)
156646	Tray (photogr. tray), PP, white stackable, 510 x 410 x 120 mm	1	piece(s)
156648	Tray (photogr. tray), PP, white stackable, 625 x 530 x 140 mm	1	piece(s)
156650	Tray (photogr. tray), PP, white stackable, 840 x 645 x 160 mm	1	piece(s)

Tidy tray, PVC



Cat. No.	Description	From*	Unit
768505	Tidy tray with compartments, PVC, 5 compartments, 400 x 300 x 65 mm	1	piece(s)
768510	Tidy tray with compartments, PVC, 12 compartments, 400 x 300 x 65 mm	1	piece(s)


Scoop, PP, white



Cat. No.	Description	From*	Unit
139115	Scoop, PP, white, length 160 mm, 50 ml	12	piece(s)
139120	Scoop, PP, white, length 200 mm, 100 ml	12	piece(s)
139125	Scoop, PP, white, length 260 mm, 250 ml	6	piece(s)
139130	Scoop, PP, white, length 315 mm, 500 ml	6	piece(s)
139135	Scoop, PP, white, length 385 mm, 1000 ml	6	piece(s)

Scoop, PE

Cat. No.	Description	From*	Unit
139035	Scoop, PE, natural colored, length 350 mm, approx. 500 ml	6	piece(s)
139040	Scoop, PE, natural colored, length 400 mm, approx. 1000 ml	6	piece(s)



Spatula, impact-resistant PS

Cat. No.	Description	From*	Unit
139810	Laboratory spatula, impact PS, 150 mm, double spatula	10	piece(s)
139815	Laboratory spatula, impact PS, 180 mm, double spatula	10	piece(s)
139817	Laboratory spatula, impact PS, 180 mm, spatula / spoon	10	piece(s)
139820	Laboratory spatula, impact PS, 210 mm, spatula / spoon	10	piece(s)



Spatula, PE-HD, coated steel core

Cat. No.	Description	From*	Unit
136100	Stirring spatula, PE-HD, with steel core l. approx. 260 mm, dia. approx. 7 mm	1	piece(s)



Spatula, PP

Cat. No.	Description	From*	Unit
136110	Stirring spatula, PP, l. 245 mm, dia. 6 mm	10	piece(s)



Sampling dipper, PTFE

Cat. No.	Description	From*	Unit
90438	Sampling dipper with spout 100 ml, dia. 55 mm, h. 65 mm, PTFE	1	piece(s)
90448	Sampling dipper with spout 250 ml, dia. 70 mm, h. 95 mm, PTFE	1	piece(s)
90454	Sampling dipper with spout 500 ml, dia. 80 mm, h. 125 mm, PTFE	1	piece(s)
90462	Sampling dipper with spout 1000 ml, dia. 105 mm, h. 155 mm, PTFE	1	piece(s)



Forceps, PMP

Cat. No.	Description	From*	Unit
112805	Forcep, PMP, pointed ends, length 115 mm, autoclavable at 121 °C	10	piece(s)
112810	Forcep, PMP, pointed ends, length 145 mm, autoclavable at 121 °C	10	piece(s)



Extension rod, PTFE, for sampling dipper



Cat. No.	Description	From*	Unit
90470	Extension rod for sampling dippers, length 600 mm, PTFE	1	piece(s)

Forceps, POM



Cat. No.	Description	From*	Unit
113020	Forcep, POM, glass-fibre reinforced length 250 mm, rounded ends	5	piece(s)

Forceps, PTFE



Cat. No.	Description	From*	Unit
113125	Forcep, PTFE, rectangular ends length 100 mm	1	piece(s)
113130	Forcep, PTFE, rectangular ends length 150 mm	1	piece(s)
113135	Forcep, PTFE, rectangular ends length 200 mm	1	piece(s)

Weighing dishes, PS



Cat. No.	Description	From*	Unit
155540	Weighing dish, PS, square shape, 7 ml, 45 x 45 x 7 mm	500	piece(s)
155542	Weighing dish, PS, square shape, 100 ml, 80 x 80 x 23 mm	500	piece(s)
155544	Weighing dish, PS, square shape, 250 ml, 140 x 140 x 21 mm	500	piece(s)

Weighing dishes, PS



Cat. No.	Description	From*	Unit
155546	Weighing dish, PS, diamond shape, 25 ml, 80 x 50 x 12 mm	500	piece(s)
155548	Weighing dish, PS, diamond shape, 30 ml, 80 x 60 x 14 mm	500	piece(s)
155550	Weighing dish, PS, diamond shape, 100 ml, 125 x 100 x 19 mm	500	piece(s)

Bucket, PE-HD



Cat. No.	Description	From*	Unit
71772	Bucket, PE-HD, w/o lid and spout, 5 l, h. 240 mm with graduation and handle	1	piece(s)
71774	Bucket, PE-HD, w/o lid and spout, 10 l, h. 300 mm with graduation and handle	1	piece(s)

Bucket, PP

Cat. No.	Description	From*	Unit
72376	Bucket, PP, w/o lid, with spout, 12 l, h. 300 mm reinforced rim and handle	1	piece(s)
72378	Bucket, PP, w/o lid, with spout, 15 l, h. 340 mm reinforced rim and handle	1	piece(s)

**Press-on lid, PE-LD, for buckets**

Cat. No.	Description	From*	Unit
71771	Push-on lid for buckets, PE-LD, for, 5 l buckets	1	piece(s)
71776	Push-on lid for buckets, PE-LD for 10 l buckets	1	piece(s)

**Mortars with pestle, MF, white**

Cat. No.	Description	From*	Unit
110705	Mortar with pestle, MF, white, with foot/spout, 125 x 80 mm, autoclavable at 121 °C	1	piece(s)
110710	Mortar with pestle, MF, white, with foot/spout, 150 x 90 mm, autoclavable at 121 °C	1	piece(s)

**Dropping bottle, soda-lime glass, amber glass**

Cat. No.	Description	From*	Unit
124529	Dropping bottle, soda-lime glass, amber, 50 ml, with dropping pipette and rubber cap	1	piece(s)
124539	Dropping bottle, soda-lime glass, amber, 100 ml, with dropping pipette and rubber cap	1	piece(s)

**Dropping bottle, soda-lime glass**

Cat. No.	Description	From*	Unit
124633	Dropping bottle, soda-lime glass, clear, 50 ml, with dropping pipette and rubber cap	1	piece(s)
124639	Dropping bottle, soda-lime glass, clear, 100 ml, with dropping pipette and rubber cap	1	piece(s)



Laboratory bottle, Boro 3.3



Cat. No.	Description	From*	Unit
122520	Laboratory bottle, Boro 3.3, 25 ml, graduated, GL 25, screw cap (PP)	10	piece(s)
122528	Laboratory bottle, Boro 3.3, 50 ml, graduated, GL 32, screw cap (PP)	10	piece(s)
122538	Laboratory bottle, Boro 3.3, 100 ml, graduated, GL 45, screw cap (PP)	10	piece(s)
122548	Laboratory bottle, Boro 3.3, 250 ml, graduated, GL 45, screw cap (PP)	10	piece(s)
122554	Laboratory bottle, Boro 3.3, 500 ml, graduated, GL 45, screw cap (PP)	10	piece(s)
122562	Laboratory bottle, Boro 3.3, 1000 ml, graduated, GL 45, screw cap (PP)	10	piece(s)
122564	Laboratory bottle, Boro 3.3, 2000 ml, graduated, GL 45, screw cap (PP)	10	piece(s)
122570	Laboratory bottle, Boro 3.3, 5000 ml, graduated, GL 45, screw cap (PP)	1	piece(s)
122574	Laboratory bottle, Boro 3.3, 10000 ml, graduated, GL 45, screw cap (PP)	1	piece(s)
122578	Laboratory bottle, Boro 3.3, 20000 ml, graduated, GL 45, screw cap (PP)	1	piece(s)

Threaded bottle, amber glass



Cat. No.	Description	From*	Unit
704012	Threaded bottle amber, without coating, 100 ml thread 32 with screw cap PP	1	piece(s)
704014	Threaded bottle amber, without coating, 250 ml thread 32 with screw cap PP	1	piece(s)
704016	Threaded bottle amber, without coating, 500 ml thread 32 with screw cap PP	1	piece(s)
704018	Threaded bottle amber, without coating, 1000 ml thread 45 with screw cap PP	1	piece(s)
704020	Threaded bottle amber, without coating, 2500 ml thread 45 with screw cap PP	1	piece(s)

Threaded bottle, amber glass, ethylene-acrylate coating



Cat. No.	Description	From*	Unit
704002	Threaded bottle amber Ethyl-Acrylate coat. 100ml thread 32 screw cap PP temp. 80 C.	1	piece(s)
704004	Threaded bottle amber Acrylate coat. 250 ml thread 32 screw cap PP	1	piece(s)
704006	Threaded bottle amber Acrylate coat. 500 ml thread 32 screw cap PP	1	piece(s)
704008	Threaded bottle amber Acrylate coat. 1000 ml thread 45 screw cap PP	1	piece(s)
704010	Threaded bottle amber Acrylate coat. 2500 ml thread 45 screw cap PP	1	piece(s)

Wash bottle, PFA, technical quality, narrow neck



Cat. No.	Description	From*	Unit
143848	Washing bottle, PFA, narrow-mouth, 250 ml, ETFE screw caps, GL 25, FEP wash elements	1	piece(s)
143854	Washing bottle, PFA, narrow-mouth, 500 ml, ETFE screw caps, GL 25, FEP wash elements	1	piece(s)
143862	Washing bottle, PFA, narrow-mouth, 1000 ml, ETFE screw caps, GL 32, FEP wash elements	1	piece(s)

Wash bottle, PE-LD, narrow neck

Cat. No.	Description	From*	Unit
144148	Wash bottle, PE-LD, 250 ml, GL 25 interchangeable wash head	50	piece(s)
144154	Wash bottle, PE-LD, 500 ml, GL 25 interchangeable wash head	50	piece(s)
144162	Wash bottle, PE-LD, 1000 ml, GL 28 interchangeable wash head	25	piece(s)



Wash bottle, PE-LD, wide neck

Cat. No.	Description	From*	Unit
144018	Wash bottle, PE-LD, wide-neck, 250 ml, Dist. water, white, with valve	5	piece(s)
144019	Wash bottle, PE-LD, wide-neck, 250 ml, Acetone, red, with valve	5	piece(s)
144020	Wash bottle, PE-LD, wide-neck, 250 ml, Methanol, green, with valve	5	piece(s)
144021	Wash bottle, PE-LD, wide-neck, 250 ml, Isopropanol, blue, with valve	5	piece(s)
144022	Wash bottle, PE-LD, wide-neck, 250 ml, Ethanol, orange, with valve	5	piece(s)
144024	Wash bottle, PE-LD, 500 ml, Dist. water, white, with valve	5	piece(s)
144025	Wash bottle, PE-LD, 500 ml, Acetone, red, with valve	5	piece(s)
144026	Wash bottle, PE-LD, 500 ml, Methanol, green, with valve	5	piece(s)
144027	Wash bottle, PE-LD, 500 ml, Isopropanol, blue, with valve	5	piece(s)
144028	Wash bottle, PE-LD, 500 ml, Ethanol, orange, with valve	5	piece(s)
144078	Wash bottle, PE-LD, wide-neck, 250 ml, Dist. water, white	5	piece(s)
144079	Wash bottle, PE-LD, wide-neck, 250 ml, Acetone, red	5	piece(s)
144080	Wash bottle, PE-LD, wide-neck, 250 ml, Methanol, green	5	piece(s)
144081	Wash bottle, PE-LD, wide-neck, 250 ml, Isopropanol, blue	5	piece(s)
144082	Wash bottle, PE-LD, wide-neck, 250 ml, Ethanol, orange	5	piece(s)
144083	Wash bottle, PE-LD, wide-neck, 250 ml, neutral, yellow	5	piece(s)
144084	Wash bottle, PE-LD, wide-neck, 500 ml, Dist. water, white	5	piece(s)
144085	Wash bottle, PE-LD, wide-neck, 500 ml, Acetone, red	5	piece(s)
144086	Wash bottle, PE-LD, wide-neck, 500 ml, Methanol, green	5	piece(s)
144087	Wash bottle, PE-LD, wide-neck, 500 ml, Isopropanol, blue	5	piece(s)
144088	Wash bottle, PE-LD, wide-neck, 500 ml, Ethanol, orange	5	piece(s)
144089	Wash bottle, PE-LD, wide-neck, 500 ml, neutral, yellow	5	piece(s)



Aspirator bottle, PE-HD, narrow neck

Cat. No.	Description	From*	Unit
131170	Aspirator bottle, PE-HD, narrow-mouth, 5 l, with stopcock, screw cap and carrying handle	1	piece(s)
131174	Aspirator bottle, PE-HD, narrow-mouth, 10 l, with stopcock, screw cap and carrying handle	1	piece(s)



Storage bottle, PE-HD, narrow neck

Cat. No.	Description	From*	Unit
130870	Storage bottle, PE-HD, narrow-mouth, 5 l, with screw cap and carrying handle	1	piece(s)
130874	Storage bottle, PE-HD, narrow-mouth, 10 l, with screw cap and carrying handle	1	piece(s)

Storage bottle, PE-HD, wide neck

Cat. No.	Description	From*	Unit
130970	Storage bottle, PE-HD, wide-neck, 5 l, with screw cap and carrying handle	1	piece(s)
130974	Storage bottle, PE-HD, wide-neck, 10 l, with screw cap and carrying handle	1	piece(s)

Bottle, PE-LD, narrow neck

Cat. No.	Description	From*	Unit
129428	Bottle, PE-LD, narrow-mouth, 50 ml, GL 18, with cap	100	piece(s)
129438	Bottle, PE-LD, narrow-mouth, 100 ml, GL 18, with cap	50	piece(s)
129448	Bottle, PE-LD, narrow-mouth, 250 ml, GL 25, with cap	50	piece(s)
129454	Bottle, PE-LD, narrow-mouth, 500 ml, GL 25, with cap	50	piece(s)
129462	Bottle, PE-LD, narrow-mouth, 1000 ml, GL 28, with cap	25	piece(s)
129464	Bottle, PE-LD, narrow-mouth, 2000 ml, GL 28, with cap	25	piece(s)

Bottle, PE-LD, wide neck

Cat. No.	Description	From*	Unit
129628	Bottle, PE-LD, wide-neck, 50 ml, GL 32, with cap	100	piece(s)
129638	Bottle, PE-LD, wide-neck, 100 ml, GL 32, with cap	50	piece(s)
129648	Bottle, PE-LD, wide-neck, 250 ml, GL 40, with cap	50	piece(s)
129654	Bottle, PE-LD, wide-neck, 500 ml, GL 50, with cap	50	piece(s)
129662	Bottle, PE-LD, wide-neck, 1000 ml, GL 65, with cap	25	piece(s)
129664	Bottle, PE-LD, wide-neck, 2000 ml, GL 65, with cap	25	piece(s)

Bottle, amber, PE-LD, wide neck

Cat. No.	Description	From*	Unit
129728	Bottle, PE-LD, wide neck, 50 ml, GL 32, amber, opaque, screw cap	1	piece(s)
129738	Bottle, PE-LD, wide neck, 100 ml, GL 32, amber, opaque, screw cap	1	piece(s)
129748	Bottle, PE-LD, wide neck, 250 ml, GL 32, amber, opaque, screw cap	1	piece(s)
129754	Bottle, PE-LD, wide neck, 500 ml, GL 50, amber, opaque, screw cap	1	piece(s)
129762	Bottle, PE-LD, wide neck, 1000 ml, GL 65, amber, opaque, screw cap	1	piece(s)
129764	Bottle, PE-LD, wide neck, 2000 ml, GL 65, amber, opaque, screw cap	1	piece(s)



Bottle, amber, PE-LD, narrow neck

Cat. No.	Description	From*	Unit
130228	Bottle, PE-LD, narrow-mouth, 50 ml, GL 18, amber, opaque, screw cap	1	piece(s)
130238	Bottle, PE-LD, narrow-mouth, 100 ml, GL 18, amber, opaque, screw cap	1	piece(s)
130248	Bottle, PE-LD, narrow-mouth, 250 ml, GL 25, amber, opaque, screw cap	1	piece(s)
130254	Bottle, PE-LD, narrow-mouth, 500 ml, GL 25, amber, opaque, screw cap	1	piece(s)
130262	Bottle, PE-LD, narrow-mouth, 1000 ml, GL 28, amber, opaque, screw cap	1	piece(s)
130264	Bottle, PE-LD, narrow-mouth, 2000 ml, GL 28, amber, opaque, screw cap	1	piece(s)



Square bottle, PE-HD, wide neck

Cat. No.	Description	From*	Unit
130412	Square bottle, PE-HD, wide-neck, 100 ml, GL 32, screw cap PP	1	piece(s)
130414	Square bottle, PE-HD, wide-neck, 250 ml, GL 45, screw cap PP	1	piece(s)
130418	Square bottle, PE-HD, wide-neck, 500 ml, GL 54, screw cap PP	1	piece(s)
130422	Square bottle, PE-HD, wide-neck, 1000 ml, GL 54, screw cap PP	1	piece(s)
130424	Square bottle, PE-HD, wide-neck, 1500 ml, GL 80, screw cap PP	1	piece(s)
130426	Square bottle, PE-HD, wide-neck, 2500 ml, GL 80, screw cap PP	1	piece(s)
130428	Square bottle, PE-HD, wide-neck, 4000 ml, GL W 80, screw cap PP	1	piece(s)



Bottle, PTFE, wide neck

Cat. No.	Description	From*	Unit
130538	Bottle, PTFE, wide-neck, 100 ml, mouth inner dia. 35 mm, screw cap	1	piece(s)
130548	Bottle, PTFE, wide-neck, 250 ml, mouth inner dia. 42 mm, screw cap	1	piece(s)
130554	Bottle, PTFE, wide-neck, 500 ml, mouth inner dia. 52 mm, screw cap	1	piece(s)
130562	Bottle, PTFE, wide-neck, 1000 ml, mouth inner dia. 57 mm, screw cap	1	piece(s)
130564	Bottle, PTFE, wide-neck, 2000 ml, mouth inner dia. 60 mm, screw cap	1	piece(s)



Bottle, PFA, wide neck



Cat. No.	Description	From*	Unit
128748	Bottle, PFA, wide-neck, 250 ml, with thread S40 and screw cap	1	piece(s)
128754	Bottle, PFA, wide-neck, 500 ml, with thread S40 and screw cap	1	piece(s)
128762	Bottle, PFA, wide-neck, 1000 ml, with thread S40 and screw cap	1	piece(s)
128764	Bottle, PFA, wide-neck, 2000 ml, with thread S40 and screw cap	1	piece(s)

Bottle, PFA, narrow neck



Cat. No.	Description	From*	Unit
128928	Bottle, PFA, narrow-mouth, 50 ml, with thread S28 and screw cap	1	piece(s)
128938	Bottle, PFA, narrow-mouth, 100 ml, with thread S28 and screw cap	1	piece(s)
128948	Bottle, PFA, narrow-mouth, 250 ml, with thread S28 and screw cap	1	piece(s)
128954	Bottle, PFA, narrow-mouth, 500 ml, with thread S28 and screw cap	1	piece(s)
128962	Bottle, PFA, narrow-mouth, 1000 ml, with thread S28 and screw cap	1	piece(s)

Dropping bottle, PE-LD, with dropper nozzle



Cat. No.	Description	From*	Unit
125216	Dropping bottle, PE-LD, 20 ml, GL 14, with dropper nozzle	100	piece(s)
125224	Dropping bottle, PE-LD, 30 ml, GL 14, with dropper nozzle	100	piece(s)
125228	Dropping bottle, PE-LD, 50 ml, GL 18, with dropper nozzle	100	piece(s)
125238	Dropping bottle, PE-LD, 100 ml, GL 18, with dropper nozzle	50	piece(s)
125248	Dropping bottle, PE-LD, 250 ml, GL 25, with dropper nozzle	50	piece(s)
125254	Dropping bottle, PE-LD, 500 ml, GL 25, with dropper nozzle	50	piece(s)

Dropping bottle, PTFE, with dropper nozzle and threaded cap



Cat. No.	Description	From*	Unit
125316	Dropping bottle, PTFE, 25 ml, with dropper nozzle and screw cap	1	piece(s)
125328	Dropping bottle, PTFE, 50 ml, with dropper nozzle and screw cap	1	piece(s)

Atomizer, PP



Cat. No.	Description	From*	Unit
144180	Atomizer, white, 400 ml, bottle including lid made from PP, atomizer head made from PE-LD	5	piece(s)
144190	Atomizer, white, 850 ml, bottle including lid made from PP, atomizer head made from PE-LD	5	piece(s)

Stopcock for aspirator bottle

Cat. No.	Description	From*	Unit
131100	Stopcock, PP, 3/4" connection, for storage- / aspirator bottles and drums	1	piece(s)




Rubber cap

Cat. No.	Description	From*	Unit
124700	Rubber teat, pipetting aid for Pasteur pipettes and dropping bottles of glass	100	piece(s)



Dropper nozzle, PE-LD, for threaded bottles

Cat. No.	Description	From*	Unit
129230	Closure with dropper nozzle, PE-LD for bottles with thread, GL 14	100	piece(s)
129232	Closure with dropper nozzle, PE-LD for bottles with thread, GL 18	50	piece(s)
129234	Closure with dropper nozzle, PE-LD for bottles with thread, GL 25	50	piece(s)



Spare wash head, ETFE, with delivery tube FEP

Cat. No.	Description	From*	Unit
129273	Wash head: screw cap ETFE + del. tube FEP for wash bottles of PFA, GL 25	1	piece(s)
129275	Replacement wash heads ETFE with FEP delivery tube, GL 32	1	piece(s)



Wash head, PE-LD, for threaded bottles

Cat. No.	Description	From*	Unit
129222	Wash head, PE-LD for bottles with thread, GL 25	50	piece(s)
129224	Wash head, PE-LD for bottles with thread, GL 28	25	piece(s)



Screw cap, PFA, for sample jars PFA

Cat. No.	Description	From*	Unit
129260	Replacement screw cap for narrow-mouth PFA bottles, S 28	1	piece(s)
129262	Replacement screw cap for wide-neck PFA bottles, S 40	1	piece(s)



Screw cap, PP

Cat. No.	Description	From*	Unit
130320	Screw cap, PP, thread GL W 32, for wide-neck bottles	1	piece(s)
130322	Screw cap, PP, thread GL W 45, for wide-neck bottles	1	piece(s)
130324	Screw cap, PP, thread GL W 54, for wide-neck bottles	1	piece(s)
130326	Screw cap, PP, thread GL W 80, for wide-neck bottles	1	piece(s)

Screw cap, PE-LD

Cat. No.	Description	From*	Unit
129202	Screw cap, PE-LD for bottles with thread, GL 18	50	piece(s)
129204	Screw cap, PE-LD for bottles with thread, GL 25	50	piece(s)
129206	Screw cap, PE-LD for bottles with thread, GL 28	25	piece(s)
129208	Screw cap, PE-LD for bottles with thread, GL 32	25	piece(s)
129210	Screw cap, PE-LD for bottles with thread, GL 40	25	piece(s)
129212	Screw cap, PE-LD for bottles with thread, GL 50	25	piece(s)
129214	Screw cap, PE-LD for bottles with thread, GL 65	25	piece(s)

Pouring ring, PP, for laboratory bottles

Cat. No.	Description	From*	Unit
122662	Pouring ring, PP, for laboratory bottles GL 32	10	piece(s)
122664	Pouring ring, PP, for laboratory bottles GL 45	10	piece(s)

Screw cap, PP, for laboratory bottles

Cat. No.	Description	From*	Unit
122620	Screw cap, PP, for laboratory bottles GL 25, integral glass ring	10	piece(s)
122628	Screw cap, PP, for laboratory bottles GL 32	10	piece(s)
122638	Screw cap, PP, for laboratory bottles GL 45	10	piece(s)

Watch glass, soda-lime glass

Cat. No.	Description	From*	Unit
150010	Watch glass, soda-lime glass, dia. 40 mm ground edges stress-reduced	10	piece(s)
150015	Watch glass, soda-lime glass, dia. 50 mm ground edges stress-reduced	10	piece(s)
150020	Watch glass, soda-lime glass, dia. 60 mm ground edges stress-reduced	10	piece(s)
150030	Watch glass, soda-lime glass, dia. 70 mm ground edges stress-reduced	10	piece(s)
150040	Watch glass, soda-lime glass, dia. 80 mm ground edges stress-reduced	10	piece(s)
150045	Watch glass, soda-lime glass, dia. 90 mm ground edges stress-reduced	10	piece(s)
150050	Watch glass, soda-lime glass, dia. 100 mm ground edges stress-reduced	10	piece(s)
150060	Watch glass, soda-lime glass, dia. 120 mm ground edges stress-reduced	10	piece(s)
150065	Watch glass, soda-lime glass, dia. 125 mm ground edges stress-reduced	10	piece(s)
150070	Watch glass, soda-lime glass, dia. 150 mm ground edges stress-reduced	10	piece(s)
150080	Watch glass, soda-lime glass, dia. 200 mm ground edges stress-reduced	10	piece(s)



Watch glass, PTFE

Cat. No.	Description	From*	Unit
150900	Watch glass, PTFE, diameter 20 mm	1	piece(s)
150902	Watch glass, PTFE, diameter 30 mm	1	piece(s)
150904	Watch glass, PTFE, diameter 40 mm	1	piece(s)
150905	Watch glass, PTFE, diameter 50 mm	1	piece(s)
150907	Watch glass, PTFE, diameter 65 mm	1	piece(s)
150910	Watch glass, PTFE, diameter 75 mm	1	piece(s)
150912	Watch glass, PTFE, diameter 80 mm	1	piece(s)
150915	Watch glass, PTFE, diameter 100 mm	1	piece(s)
150920	Watch glass, PTFE, diameter 125 mm	1	piece(s)
150922	Watch glass, PTFE, diameter 150 mm	1	piece(s)



Disinfection tissues mikrozid®

Cat. No.	Description	From*	Unit
44860	mikrozid® universal wipes premium disinfection wipes 1 pack = 100 piece(s)	6	pack



Universal detergent, Mucaso® liquid concentrate

Cat. No.	Description	From*	Unit
44003	Mucaso®, universal detergent, 2 l, liquid concentrate	5	piece(s)
44005	Mucaso®, universal detergent, 5 l, liquid concentrate	1	piece(s)



Edisonite® CLASSIC, Universal detergent, powdered

Cat. No.	Description	From*	Unit
44205	EDISONITE® CLASSIC, universal detergent, 5 kg, powder compound	1	piece(s)
44301	EDISONITE® CLASSIC, universal detergent, 1 kg, powder compound	1	piece(s)

Edisonite® SUPER, Universal detergent, powdered

Cat. No.	Description	From*	Unit
44361	EDISONITE® Super, universal detergent, 1 kg, powder compound	1	piece(s)
44365	EDISONITE® Super, universal detergent, 5 kg, powder compound	1	piece(s)

Disinfecting detergent, Mucocit®-T liquid concentrate

Cat. No.	Description	From*	Unit
44822	Mucocit® -T, instrument disinfecting detergent, 2 l, liquid concentrate	5	piece(s)
44825	Mucocit® -T, instrument disinfecting detergent, 5 l, liquid concentrate	1	piece(s)

Pursept® A Xpress, Disinfecting detergent

Cat. No.	Description	From*	Unit
44830	Pursept® -A Xpress, surface disinfecting detergent, 1 l spray bottle	10	piece(s)
44835	Pursept® -A Xpress, surface disinfecting detergent, 5 l, liquid	1	piece(s)

Atomizer head for 1 l Pursept® A Xpress bottle

Cat. No.	Description	From*	Unit
44832	Spray head for disinfection spray 1 l Pursept® -A Xpress	1	piece(s)

Disinfecting detergent, Pursept®-AF liquid concentrate

Cat. No.	Description	From*	Unit
44850	Pursept® -AF, surface disinfecting detergent, 2 l, liquid concentrate	5	piece(s)
44855	Pursept® -AF, surface disinfecting detergent, 5 l, liquid concentrate	1	piece(s)

schülke wipes Disinfecting detergent

Cat. No.	Description	From*	Unit
44870	schülke wipes 70, fleece wipes, 230x280mm, Roll with 70 wipes, dry	6	piece(s)



Dispenser box for schülke wipes 70

Cat. No.	Description	From*	Unit
44875	Dispenser box for schülke wipes 70, fleece wipes 230x280mm, empty	1	piece(s)



Pipette rinsing apparatus, PE-HD

Cat. No.	Description	From*	Unit
29120	Pipette rinser, PE-HD, for pipette length 460 mm, 13, l	1	piece(s)
29125	Pipette rinser, PE-HD, for pipette length 500 mm, 17, l	1	piece(s)



Draining rack, PS

Cat. No.	Description	From*	Unit
95105	Draining rack, PS, light grey, 450 x 650 mm, 72 pegs, dia. 16 mm	1	piece(s)



Soaking jar, PE-HD

Cat. No.	Description	From*	Unit
29210	Soaking jar, PE-HD, for pipette length 460 mm, 10, l	1	piece(s)
29215	Soaking jar, PE-HD, for pipette length 600 mm, 12,5 l	1	piece(s)



Pipette basket, PE-HD

Cat. No.	Description	From*	Unit
29010	Pipette basket, PE-HD, with handle, for pipette rinsers, height with handle 650 mm	1	piece(s)
29015	Pipette basket, PE-HD, with handle, for pipette rinsers, height with handle 850 mm	1	piece(s)



Special pegs, draining rack



Cat. No.	Description	From*	Unit
95111	Special pegs for drain. rack PS, out. dia. 6 mm, l. 100 mm 1 pack = 11 piece(s)	1	pack

Accessories for water still



Cat. No.	Description	From*	Unit
584700	Hosecock clamp, nickel plated brass for tubes up to 14 mm outer diameter	1	piece(s)

Spare cap sedimentation cones

Cat. No.	Description	From*	Unit
388002	Rack for two Imhoff sedimentation cones, made of glass or plastics 300x130x400 mm	1	

Rack for sedimentation cones



Cat. No.	Description	From*	Unit
388050	Rack for two Imhoff sedimentation cones, made of plastics, 300 x 130 x 315 mm	1	piece(s)
388060	Rack for two Imhoff sedimentation cones, made of glass or plastics 300x130x400 mm	1	piece(s)

Spring clip for oxygen flasks Winkler pattern



Cat. No.	Description	From*	Unit
386138	Fastening clip for oxygen flask, Winkler for capacity 100 - 150 ml	1	piece(s)
386148	Fastening clip for oxygen flask, Winkler for capacity 250 - 300 ml	1	piece(s)

ASTM centrifuge tubes, BLAUBRAND®, Boro 3.3



Cat. No.	Description	From*	Unit
362038	ASTM centrifuge tube, BLAUBRAND®, Boro 3.3, 100 ml, cylindrical, conical bottom, norm ASTM D91	2	piece(s)
362138	ASTM centrifuge tube, BLAUBRAND®, Boro 3.3, 100 ml, pear-shaped, cylindrical bottom, standard D96	2	piece(s)
362338	ASTM centrifuge tube, BLAUBRAND®, Boro 3.3, 100 ml, cylindrical, conical bottom, standard D96	2	piece(s)

Oxygen flasks Winkler pattern

Cat. No.	Description	From*	Unit
386038	Oxygen flask, Winkler, soda-lime glass, 100-150 ml, with glass stopper, NS 14/23	2	piece(s)
386048	Oxygen flask, Winkler, soda-lime glass, 250-300 ml, with glass stopper, NS 19/26	2	piece(s)



Imhoff sedimentation cones Imhoff pattern, SAN

Cat. No.	Description	From*	Unit
388000	Imhoff sedimentation cone, grad. 0 - 1000 ml, SAN, with screw cap	1	piece(s)



Imhoff Sedimentation Cones Imhoff pattern, SILBERBRAND

Cat. No.	Description	From*	Unit
387362	Imhoff sedimentation cone, SILBERBRAND, grad. 0 - 100 ml, mark 1000 ml, w/o stopcock	4	piece(s)
387462	Imhoff sedimentation cone, SILBERBRAND, grad. 0 - 1000 ml, mark 1000 ml, w/o stopcock	4	piece(s)
387662	Imhoff sedimentation cone, SILBERBRAND, grad. 0 - 100 ml, mark 1000 ml, with stopcock	1	piece(s)



Density bottle Gay-Lussac pattern, uncalibrated

Cat. No.	Description	From*	Unit
43205	Density bottle Gay-Lussac, 5 ml, uncalibrated, Boro 3.3	2	piece(s)
43208	Density bottle Gay-Lussac, 10 ml, uncalibrated, Boro 3.3	2	piece(s)
43220	Density bottle Gay-Lussac, 25 ml, uncalibrated, Boro 3.3	2	piece(s)
43228	Density bottle Gay-Lussac, 50 ml, uncalibrated, Boro 3.3	2	piece(s)
43238	Density bottle Gay-Lussac, 100 ml, uncalibrated, Boro 3.3	2	piece(s)



Density bottle Gay-Lussac pattern, BLAUBRAND®, calibrated

Cat. No.	Description	From*	Unit
43305	Density bottle Gay-Lussac, BLAUBRAND®, 5 ml, calibrated, Boro 3.3	1	piece(s)
43308	Density bottle Gay-Lussac, BLAUBRAND®, 10 ml, calibrated, Boro 3.3	1	piece(s)
43320	Density bottle Gay-Lussac, BLAUBRAND®, 25 ml, calibrated, Boro 3.3	1	piece(s)
43328	Density bottle Gay-Lussac, BLAUBRAND®, 50 ml, calibrated, Boro 3.3	1	piece(s)
43338	Density bottle Gay-Lussac, BLAUBRAND®, 100 ml, calibrated, Boro 3.3	1	piece(s)



BRAND GMBH + CO KG General Terms and Conditions

for orders from outside shop.brand.de

1. General

- 1.1. These General Terms and Conditions (GT&C) are intended for use in commercial transactions between businesses.
- 1.2. These GT&C shall apply for all, including future, contracts with the customer, except for online business via shop.brand.de/en/. Other terms and conditions shall not become part of the contract, even if BRAND does not expressly object to them. Amendments to and changes of the contract must be in written form. The waiver of the requirement for written form shall only be possible in writing. This shall not apply to individual contractual agreements. The language of the contract shall be German and/ or English. In the event of a discrepancy between the German language version of these GT&C and a version in any language, the German language version shall prevail.
- 1.3. BRAND offers are subject to change and non-binding. BRAND reserves the right to make technical improvements to BRAND products.
- 1.4. BRAND may electronically store and process data necessary for the purpose of processing the contract.
- 1.5. A set-off by the customer shall not be permitted unless the counterclaims are undisputed or legally established, or pecuniary counterclaims arising from the right to refuse payment pursuant to Section 320 Bürgerliches Gesetzbuch (BGB) (German Civil Code).
- 1.6. Orders with a goods value of below € 250 shall be subject to a minimum order surcharge of € 50. Delivery shall be undertaken generally in packaging units (PU) according to the currently valid price list. For deliveries within five (5) working days or for order values up to € 500, BRAND reserves the right to waive an order confirmation.
- 1.7. For commercial transactions with customers having no general place of jurisdiction in Germany and between businesses, public law legal persons or special funds under public law the place of jurisdiction shall be the court responsible in Frankfurt am Main, Germany. BRAND shall also be entitled to appeal to the court responsible for the head office of the customer. BRAND shall, furthermore, as plaintiff have the right to invoke the Arbitration Court at the Chamber of Commerce and Industry in Frankfurt am Main, Germany. The Arbitration Court shall, in this case, make the final judgment in accordance with the Rules of Arbitration of the Chamber of Commerce and Industry in Frankfurt am Main without recourse to the ordinary courts of law. The instigation of legal dunning proceedings by BRAND shall not signify the exertion of its right of choice; it shall be admissible in all cases.
- 1.8. German law shall apply exclusively under the exclusion of the conflict of laws principles of Private International Law and the UN Convention on Contracts for the International Sale of Goods (CISG).

2. Delivery

- 2.1 The place of performance shall be the factory of BRAND in Wertheim, Germany. The risk shall transfer to the customer when the goods for delivery are packed and ready for pick-up (EXW (Incoterms® 2010 ex works)). This shall also apply to partial deliveries or where BRAND has performed additional services, such as shipping; costs for transport, packaging or insurance; exportation and installation. This shall also apply in case of delivery to a consignment warehouse of the customer.
- 2.2 Insofar as BRAND has agreed to orders on call, the customer must take delivery prior to the agreed date, otherwise within six (6) months.

3. Delivery Period, Force Majeure, Delay

- 3.1 Delivery times shall be ex works. Delivery periods or delivery dates are subject to the customer providing relevant information and documents such as drawings, permits or approvals, opening letters of credit as agreed, making agreed advanced payments and complying with any obligation on time. In particular, the customer is obliged to immediately provide all information required for export, import or shipment (e.g. final recipient, final destination and intended use), documents, approvals and certificates which may be required for the fulfillment of BRAND's obligations. If delays occur due to required approval, examination or information procedures, delivery periods and delivery dates shall be extended accordingly if BRAND is not solely responsible for the delay.
- 3.2 Delivery shall be subject to BRAND receiving its own supplies punctually and in good order.
- 3.3 If delays of delivery periods are due to Force Majeure, e.g. natural disasters, mobilization, war, riot or similar events beyond control of BRAND, e.g. strike or lock-out, delivery periods shall be extended by the time during which the aforementioned event or its effects persist.
- 3.4 BRAND shall immediately inform the customer about delivery barriers or delays and their expected duration.
- 3.5 BRAND shall be considered to be in default of delivery only if the customer has issued BRAND with a reminder, has set a reasonable extension period which has elapsed. The customer shall be obliged to immediately inform BRAND in writing of any likely consequences of delay.
- 3.6 In the case of delay damages, BRAND's liability for compensation shall be limited to 10% of the value of the delayed delivery/ service. The limitation shall not apply in cases of wilful intent, gross negligence and/ or injury to life, limb or health.

4. Prices, Terms of Payment

- 4.1 Prices shall be EXW (Incoterms® 2010 ex works), Wertheim and exclusive of statutory VAT, if applicable. Costs of packaging, transportation, freight and insurance shall be borne by the customer. Prices shall also be exclusive of the cost of returning and recycling/ disposing of old equipment.
- 4.2 Invoices shall be payable to BRAND account in EUROs (€) without deductions and free of charges and expenses. Payment shall be made immediately or by the date stated. The determinant factor shall be the receipt of payment. Cheques and bills of exchange shall only be accepted on account of performance and at the cost of the customer.
- 4.3 In the case of customers, with whom BRAND is working for the first time or with whom BRAND does not work regularly, after delay in payment or in the case of reasonable doubt as to the creditworthiness of the customer, BRAND shall reserve the right to make individual deliveries dependent on a pre-payment or a security deposit to the value of the invoice amount.
- 4.4 Should the period between conclusion of contract and agreed delivery exceed four (4) months, so may BRAND, at its discretion, demand a reasonable additional charge equivalent to the increase in its costs up until delivery.
- 4.5 In the case of an agreed return of goods that are free of defects, the customer shall be charged a checking and processing fee of 20% of the invoice amount (minimum € 50).
- 4.6 Should the customer be in arrears with payment, BRAND debt claims against him shall be due immediately, and BRAND shall not be obliged to make any further deliveries based on current delivery contracts.
- 4.7 If payment is delayed, BRAND shall charge - notwithstanding further damage compensation claims – interest on arrears at the statutory rate.
- 4.8 BRAND may offset amounts payable to the customer (e.g. from credit notes) against BRAND claims against the customer.

5. Retention of Title, Assignment of Future Claims

- 5.1 The goods delivered shall remain property of BRAND until the complete and unconditional payment. Should BRAND still have further claims against the customer arising from the business relationship, BRAND shall then retain its property rights until payment of such claims has been effected.
- 5.2 The customer may neither use goods subject to retention of title nor combine them with other objects to which a third party may have rights. Should, however, goods subject to retention of title become, through their combination with other objects, part of a new (complete) item, BRAND shall be a proportional co-owner of this new item directly, even if this latter component is regarded as the main component. BRAND's proportion of co-ownership shall be determined by the ratio of the invoice value of the goods to the value of the new item at the time of combination.
- 5.3 The customer may resell goods subject to retention of title in the course of his normal business as long as he has not assigned, pledged or otherwise encumbered his claims from the resale.
- 5.4 The customer shall assign to BRAND in advance as collateral any claims against his customers from the resale of the goods subject to retention of title (see Clause 5.3) and/or newly formed items (see Clause 5.2) to the value of BRAND's invoice for the goods subject to retention of title. As long as the customer is not in default of payment for the goods subject to retention of title, he may collect the assigned claims in the ordinary course of business. He may, however, only use the proportional proceeds for the payment to BRAND for the goods subject to retention of title.
- 5.5 At the customer's request, BRAND shall release collateral at its discretion if and to the extent that the nominal value of the collateral exceeds 120% of the nominal value of its outstanding debt claims against the customer.
- 5.6 The customer shall be required to inform BRAND immediately of any attachments, seizures or any other third-party dispositions relating to the goods that are reserved or co-owned by BRAND.

- 5.7 In the event of delay in payment, failure to pay bills of exchange or cheques, or failure or recall of a payment made via SEPA Direct Debit Scheme, suspension of payments or insolvency of the customer or of the end buyer, the rights of the customer under Clause 5.3 shall no longer be valid. The customer must then immediately inform the buyer of BRAND's extended retention of title; he may use the proportional proceeds relating to the assignment only to pay for the delivered goods. BRAND shall be entitled to collect the assigned receivables itself.

- 5.8 In the event of customer's culpable breach of contractual obligations, in particular for the cases covered in Clause 5.7, BRAND shall be entitled to withdraw from the contract and/ or, without withdrawing from the contract, demand the return of any goods subject to retention of title still in possession of the customer and to collect the assigned receivables itself. In order to ascertain the rights of BRAND, BRAND shall be entitled to have all of customer's documents/ books concerning the reserved rights of BRAND examined by a person who is subject to the professional duty of confidentiality.

6. Warranty, Limitation of Liability

- 6.1 BRAND warrants that its delivered goods (including any agreed installation) are free of defects at time of risk transfer. The required quality, durability and use of BRAND's delivered goods are based solely on the agreed written specification, product description and/ or operating manuals. Any information beyond this, in particular in preliminary discussions, advertising and/or referencing industrial standards shall only become part of the contract if they are expressly referenced in writing.
- 6.2 Should the customer require the delivered goods for purposes other than those agreed, he must take responsibility himself for examining their special suitability for this - also in terms of product safety – and ensure their compliance with all relevant technical, legal or regulatory provisions before the intended use. BRAND shall not be liable for any usability that was not expressly confirmed by BRAND in writing. In the case of material or design requirements of the customer, BRAND shall accept no liability for the suitability or permissibility of the desired materials or designs. Compliance with safety-related and occupational health regulations depends on the location and operating conditions of which BRAND has no prior knowledge. Action for ensuring compliance shall therefore be the responsibility of the customer or his buyer.
- 6.3 BRAND shall not be liable for the consequences of improper handling, use, maintenance and operation of the delivered goods; the consequences of normal wear and tear, in particular of wearing parts, such as pistons, seals, valves; the breakage of glass, plastic or ceramic parts; for the consequences of chemical, electrochemical or electrical influences; or non-observance of the operating instructions.
- 6.4 If a notice of defect is justified, BRAND shall initially only be required to provide cure. Any such cure shall be, at the discretion of BRAND, either rectification of the defect or delivery of goods free of defects. Further warranty claims shall only apply in the event of rejection, impossibility or failure of the cure. If cure is possible only at disproportionate expense pursuant to Section 439 Para. 4 BGB (German Civil Code), BRAND will bear the cost necessary for the cure according to Section 439 Para. 2, 3 BGB (German Civil Code) up to a limit of 150 % of the value of the delivered goods free from defects.
- 6.5 The customer must, immediately upon receipt of the goods, inspect them carefully, also in terms of product safety, and notify obvious defects immediately in writing; any hidden defects must be immediately notified upon discovery. The customer must notify the carrier immediately of any transport damage. Failure to observe the testing and notification obligation shall void any customer claims for defects.
- 6.6 BRAND's liability for slight negligence shall be limited to claims for injury to life, limb or health, to claims under the Produkthaftungsgesetz (German Product Liability Act) or to claims of culpable breach of fundamental contractual obligations through which the purpose of the contract is endangered. Otherwise, its liability for slightly negligent breach of fundamental contractual obligations is limited to the typically occurring damages which BRAND could have foreseen when the contract was concluded.
- 6.7 Should the customer use the delivered goods in conjunction with environmentally harmful, toxic, radioactive or otherwise hazardous materials, he must notify BRAND about any such materials and decontaminate the delivered goods prior to returning them to BRAND. If applicable, BRAND may charge any necessary costs for decontamination/cleaning and disposal to the customer's account.

7. Limitation Period

The warranty period shall be one year and starts from the date of delivery of the goods to the customer. The same shall apply for claims for damages, irrespective of their legal basis. The limitation periods of Section 438 Para. 1 Nos. 1 and 2 and Section 634a Para. 1 No. 2 of the BGB (German Civil Code) shall remain unaffected. The warranty period for any rights of recourse pursuant to Section 445a BGB (German Civil Code) shall be one year and starts from the date of delivery of the goods to the customer. The suspension of expiry pursuant to Section 445b BGB (German Civil Code) shall end no later than three years after BRAND has delivered the goods to the customer. The restrictions of the limitation periods shall not apply to claims based on fraudulent concealment of a defect, for claims under the Produkthaftungsgesetz (German Product Liability Act) or for damages resulting from injury to life, limb or health and other damages based on intent or gross negligence. The limitation period in respect of replaced or repaired goods shall only commence anew if BRAND admitted the defectiveness of the replaced or repaired goods.

8. Software Use

- 8.1 If software is included in the scope of a delivery, the customer shall be granted a non-exclusive right to use the software and its associated documentation. It is provided for use on the designated delivery item. The use of the software on more than one system shall be prohibited.
- 8.2 The customer shall only be entitled to copy, transfer or translate the software or to convert it from object code to source code to the extent permitted by law (Sections 69a et seq. Urheberrechtsgesetz – German Copyright Act). The customer undertakes to refrain from removing manufacturer information, in particular copyright notices, or from changing these without BRAND's prior express consent or the prior express consent of the software supplier.
- 8.3 All other rights to the software and the documentation including copies thereof shall remain with BRAND and/ or the software supplier. The issue of sublicenses is not permitted.

9. Installation

- 9.1 Installation costs may be invoiced on a monthly basis. Fixed installation prices shall only cover the work that has been agreed upon.
- 9.2 The customer shall be responsible for providing the following at his own expense: lighting, motive power, if necessary, compressed air; water; electrical power for welding and heating, including the necessary connections; electrical installations to connect the products supplied by BRAND; the devices required (such as lifting equipment); a lockable room that can be used for storing materials; tools and clothing during the installation.

10. Spare Parts, Maintenance/Repair and Calibration

- 10.1 For spare parts and maintenance, repair and calibration services, the current repair and exchange price list shall apply.
- 10.2 Insofar as there is an obligation on the part of BRAND to maintain/supply spare parts, then this obligation shall be limited to a period of five (5) years from the date of delivery. If spare parts are not manufactured by BRAND, or are no longer available on the market, for example electronic components, or if the raw material for their production is no longer available, the obligation of BRAND to deliver spare parts shall lapse.
- 10.3 For calibration and maintenance, expendable items from BRAND production are normally used.
- 10.4 Maintenance and calibration services can only be provided if the customer has declared the devices sent to be safe to work on from a health hazard perspective.
- 10.5 For repair/service values of up to € 50, BRAND reserves the right not to provide a separate cost estimate.

11. Legal Reservation, Industrial Property Rights, Confidentiality

- 11.1 BRAND reserves ownership and all industrial property rights and copyrights to all moulds, tools or other devices, samples, pictures, and business and technical documents produced or provided by BRAND. This also applies where the customer has wholly or in part taken on the costs hereof. The customer may use these only in the manner agreed with BRAND. Without BRAND's written consent, the customer may not himself manufacture contractual objects delivered or have the same manufactured by third parties.
- 11.2 Insofar as BRAND delivers goods according to the designs or other requirements specified by the customer (models, patterns etc.), the customer shall be liable to BRAND by default for ensuring that, through the manufacture and delivery of these goods, the industrial property rights or other rights of third parties are not infringed. If the customer is at fault he shall reimburse BRAND all damage resulting from any such infringement of rights.
- 11.3 Any information acquired from this business relationship and not deemed to be public knowledge must not be disclosed by the customer to third parties.

Status as of: May 2019

