

Features	Technical data
• KERN PEJ: Automatic internal adjustment, guarantees high degree of accuracy and makes the balance independent of its location of use • KERN PES: Adjusting program CAL for quick setting of the balance accuracy, external test weights at an additional price, see <i>test weights</i> • Metal housing: robust and sturdy • Weighing with tolerance range: a visual and audible signal helps with portioning, dispensing or grading • Draught shield standard for models with weighing plate size A, weighing space W×D×H 170×150×100 mm • A, B Underfloor weighing: Possibility of load support with built-in loop on the underside of the balance • A, B Protective working cover included with delivery	• Fluorescent display, bright with high contrast, digit height 14 mm • Dimensions weighing surface, stainless steel A W×D 140×120 mm B W×D 200×200 mm, see larger picture C W×D 250×220 mm • Overall dimensions (without draught shield) W×D×H A, B 220×333×93 mm C 260×330×113 mm • Permissible ambient temperature 10 °C/30 °C

www.kern-sohn.com · Order Hotline +49 7433 9933-0