

# KERN 313-020-100

**KERN**

Wooden box for weight sets (Cylindrical/compact, 1 mg - 50 g)



## Pictograms

### STANDARD



## Category

|                  |               |
|------------------|---------------|
| Brand            | KERN          |
| Product category | Cases & Boxes |
| Product group    | Wooden box    |

## Measuring System

|                                         |                     |
|-----------------------------------------|---------------------|
| Case suitable for OIML classes          | E1+E2               |
| Suitable weight shape for case          | Cylindrical/compact |
| Case suitable for nominal weight values | 1 mg - 50 g         |

## Construction

|                                 |               |
|---------------------------------|---------------|
| Dimension (W×D×H)               | 152×152×66 mm |
| Material                        | beech         |
| Case - Smallest possible weight | 1 mg          |
| Case - Highest possible weight  | 50 g          |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Delivery time                | 1 d            |
| Dimensions packaging (W×D×H) | 152×152×66 mm  |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,70 kg        |
| Gross weight approx.         | 0,70 kg        |
| Shipping weight              | 0,68 kg        |

## Product information

|                      |                          |
|----------------------|--------------------------|
| GTIN/EAN number      | 4045761003071            |
| Trade number         | 84239090                 |
| REACH classification | no prohibited substances |