

# KERN 317-020-400

**KERN**

Plastic box for individual weights (Cylindrical/compact, 1 g - 2 g)



## Pictograms

### STANDARD



## Category

|                  |               |
|------------------|---------------|
| Brand            | KERN          |
| Product category | Cases & Boxes |
| Product group    | Plastic box   |

## Measuring System

|                                         |                     |
|-----------------------------------------|---------------------|
| Case suitable for OIML classes          | E1, E2              |
| Suitable weight shape for case          | Cylindrical/compact |
| Case suitable for nominal weight values | 1 g - 2 g           |

## Construction

|                                 |          |
|---------------------------------|----------|
| Dimension (Ø×H)                 | 24×51 mm |
| Material                        | plastic  |
| Case - Smallest possible weight | 1 g      |
| Case - Highest possible weight  | 2 g      |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Delivery time                | 1 d            |
| Dimensions packaging (W×D×H) | 24×24×51 mm    |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,05 kg        |
| Gross weight approx.         | 0,05 kg        |
| Shipping weight              | 0,1 kg         |

## Product information

|                      |                          |
|----------------------|--------------------------|
| GTIN/EAN number      | 4045761081116            |
| Trade number         | 84239090                 |
| REACH classification | no prohibited substances |