

KERN 317-050-400

KERN

Plastic box for individual weights (Cylindrical/compact, 20 g)



Pictograms

STANDARD



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

Measuring System

Case suitable for OIML classes	E1, E2
Suitable weight shape for case	Cylindrical/compact
Case suitable for nominal weight values	20 g

Construction

Dimension (Ø×H)	24×51 mm
Material	plastic
Case - Smallest possible weight	20 g
Case - Highest possible weight	20 g

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	24×24×51 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Product information

GTIN/EAN number	4045761081352
Trade number	84239090
REACH classification	no prohibited substances