

KERN TGC 150-2S05

KERN

Slim pocket balance with large stainless steel weighing plate and practical tare pan



Category

Brand	KERN
Product category	Laboratory balance
Product group	Pocket balance
Product family	TGC

Measuring System

Weighing system	Strain gauge
Weighing capacity [Max]	150 g
Readability [d]	0,01 g
Reproducibility	0,2 g
Linearity	± 0,3 g
Resolution	15.000
Tare range	150 g

Adjustment options Adjusting with external weight

Recommended adjusting weight 50 g (M2); 100 g (M2)

Possible calibration points 150 g

Stabilization time 2 s

Eccentric loading at 1/3 [Max] 0,3 g

Construction type of the scale Single-range balance

Units g
gn
dwt
ozt
ct
oz

Weighing capacity [Max] (ct) 750 ct

Readability carat [d] (ct) 0,05 ct

Tolerance (% of [Max]) 0,2%

Maximum creep (15 minutes) 40 mg

Maximum creep (30 minutes) 80 mg

Default unit g

Display

Display type	LCD
Display backlight	✓
Display digit height	12 mm
Languages of the user interface	Symbol language

Construction

Dimension housing (W×D×H)	100×130×18 mm
Dimensions weighing surface (Ø)	81 mm
Dimensions weighing surface	81 mm
Material housing	plastic
Material weighing plate	stainless steel
Material display housing	plastic

Functions

IP protection - complete device	-
Auto-off interval(s) in battery mode/rechargeable battery mode	4 min
Tare function	manual (multi)
Number of keys for operation	4

Power Supply

Supplied power supply	Battery
Battery / accumulator type	Li-Metal
Battery	2×CR2032
Battery connection	Pad + Spring
Battery operating time - backlight on	8,75 h
Battery operating time - backlight off	22 h
Battery operating time	8,75 h
Battery capacity	210 mAh
Battery voltage	30 V
Battery - Number of cells	10
Battery - Individual weight	17400 g
Battery - Total weight	300 g
Battery - Construction	Individual cells
Battery capacity (Wh)	0,63 Wh
Battery - Lithium content	550 mg
Battery - Tests passed UN 38.3	✓
Battery/Accu connected	not connected, enclosed

Environmental conditions

Ambient temperature [Min]	0 °C
Ambient temperature [Max]	40 °C
Humidity of environment [Min]	80 %
Humidity of environment [Max]	90 %

KERN TGC 150-2S05



Slim pocket balance with large stainless steel weighing plate and practical tare pan

Storage temperature [Min]	-20 °C
---------------------------	--------

Storage temperature [Max]	60 °C
---------------------------	-------

Approval

CE mark	✓
---------	---

Services (optional)

Article number for adjustment at the location of installation	961-247S05
---	------------

Article number for DAkkS calibration	963-127S05.
--------------------------------------	-------------

Article number for certificate of conformity	969-517
--	---------

Packing & Shipping

Delivery time	1 d
---------------	-----

Dimensions packaging (W×D×H)	100×130×16,5 mm
------------------------------	-----------------

Shipping method	Parcel service
-----------------	----------------

Net weight approx.	0,75 kg
--------------------	---------

Gross weight approx.	0,80 kg
----------------------	---------

Shipping weight	0,8 kg
-----------------	--------

Product information

GTIN/EAN number	4045761217553
-----------------	---------------

Trade number	90160010
--------------	----------

REACH classification	no prohibited substances
----------------------	--------------------------

Pictograms

STANDARD



OPTION

