

Data Sheet

Page 1/1



30M X 0.25MM ID - BP5MS 0.25UM GC Column

Part Number: TJ-054310
Manufacturer: Trajan Scientific
Quantity: 1 pc

Description

30M X 0.25MM ID - BP5MS 0.25UM GC-Säule

Technical Information

Brand:	BP	Base Material:	5% Phenyl 95% polysilphenylene siloxane
Product Type:	GC Column	Temperature Minimum [°C]:	-40
Length:	30 m	Temperature Maximum [°C]:	330
Internal Diameter:	0.25 mm	Mode of Separation:	Gas Chromatography
Film Thickness [µm]:	0.25 µm		
Chemistry:	5 (5HT, 5MS)		
Chemistry Producers Term:	BP5MS		
USP Class:	G27, G36		
Country of Origin:	Australia		