

Data Sheet

Page 1/1



25M X 0.22MM ID - HT5 0.1UM GC Column

Part Number: TJ-054632
Manufacturer: Trajan Scientific
Quantity: 1 pc

Description

25M X 0.22MM ID - HT5 0.1UM GC-Säule

Technical Information

Product Type:	GC Column	Base Material:	5% Phenylpolycaboranesiloxane
Length:	25 m	Mode of Separation:	Gas Chromatography
Internal Diameter:	0.22 mm		
Film Thickness [µm]:	0.1 µm		
Chemistry:	5 (5HT, 5MS)		
Chemistry Producers Term:	HT5		
Country of Origin:	Australia		